



Our kids are our passion.™

Freedom Designs, Inc.  
New Product Press Release  
P.R.O. CG Center of Gravity Tilt-in-Space Mobility  
October 23, 2015

Freedom Designs, Inc. of Simi Valley, California is recognized as the market leader in the development, manufacture, sales and distribution of customized seating systems and specialized mobility products worldwide. Since 1981 Freedom Designs has made available products that have helped to change the lives of people with disabilities. Freedom's latest "NEW" product introduction the **P.R.D. CG™** (precise rotational orientation) tilt-in-space wheelchair, is yet another innovative breakthrough that's designed to help make the activities of daily living just a bit easier.

The **P.R.D. CG™** is not "complex" in fact, it's just the opposite! Developed around an appreciation for simplicity, it's so easy to adjust to achieve the "perfect" fit. When fit, it's almost effortless to tilt and un-tilt. The **P.R.D. CG™** telescoping tilt arcs, simple CG adjustment and lightweight compact frame all work in concert to enhance its maneuverability making the challenge of pushing, turning, controlling and navigating small spaces, less challenging.

The **P.R.D. CG™** tilts from 0° to 50° in 2° increments. The tilt position is easily identified from the care givers position and can be restricted to meet specific tilt angle requirements. The frame is fully adjustable in width, depth, height and CG. It is WC19 approved for transit and available in pediatric and adult sizes. With a weight capacity of 250lbs., the **P.R.D. CG™** robust yet lightweight design is durable, dependable and transportable.

For more detailed information on the "NEW" **P.R.D. CG™** or any of Freedom's other products contact your Freedom Designs sales representative, product provider or referral source professional. Or call Freedom customer service at 1-800-331-8551 or visit us at [www.freedomdesigns.com](http://www.freedomdesigns.com).

