

# Easy Conforming Back

Clear Form

Email Form to Customer Service



## Freedom Designs, Incorporated

2241 N Madera Road, Simi Valley, CA 93065  
 phone: 800-331-8551 fax: 888-582-1509  
 e: customerservice@freedomdesigns.com



Yes, you can.®

Quote Order

Buyer: \_\_\_\_\_  
 Co./Account: \_\_\_\_\_  
 Location: \_\_\_\_\_  
 Contact #: \_\_\_\_\_  
 P.O. #: \_\_\_\_\_  
 Marked For: \_\_\_\_\_  
 Fax: \_\_\_\_\_  
 Email: \_\_\_\_\_

### Wheelchair Specifications

Model: \_\_\_\_\_  
 Width: \_\_\_\_\_  
 Depth: \_\_\_\_\_

### Upholstery Selections Black lycra standard. Stretch fabrics recommended.

| Body Contact Surfaces |         | Non-Contact Surfaces |         |
|-----------------------|---------|----------------------|---------|
| Number:               | Choose: | Number:              | Choose: |
| Monogram Option       |         | Embroidery Option    |         |
| Name:                 |         | Design:              |         |
| Script                | Block   |                      |         |
| Thread Color:         |         |                      |         |

## Easy Conforming Standard

|           |  |
|-----------|--|
| 102-12000 | 4 J's and L's and 4 Glide Locks                      |
| 102-17000 | 4 J's and L's and 2 Glide Locks with 2 Split Collars |
| 102-20000 | 4 Easy Clamps  |
| EZ CONF2  | 2 J's and L's and 2 Glide Locks and Hinges           |



E2617

|        |
|--------|
| \$1262 |
| \$1262 |
| \$1466 |
| \$1262 |

## Easy Conforming Extended

|               |  |
|---------------|--|
| 102-13000     | 4 J's and L's and 4 Glide Locks                      |
| 102-18000     | 4 J's and L's and 2 Glide Locks with 2 Split Collars |
| 102-30000     | 4 Easy Clamps  |
| EZ CONF-EXT 2 | 2 J's and L's and 2 Glide Locks and Hinges           |



E2617

|        |
|--------|
| \$1262 |
| \$1262 |
| \$1466 |
| \$1262 |

## Back Modifications

|              |   |        |             |  |        |
|--------------|---|--------|-------------|--|--------|
| MOD-2QQ      | Split Outer Wings                         | \$ 280 | MOD-2SS     | Laterals Upholstered                                 | \$ 234 |
|              | <b>Standard Easy Conforming Back Only</b> |        |             | Separate from center panel                           |        |
| MOD-2UU      | Wings Molded Together                     | \$ 280 | MOD-2S      | Free Form (box style)                                | \$ 303 |
| MOD-2Econo4  | Hardware Substitution 4-Point             | \$ 202 | MOD-2Econo2 | Hardware Substitution 2-Point                        | \$ 68  |
|              | Econo-Eze                                 |        |             | Econo-Eze  |        |
| MOD-2OVER18  | Oversize Back Over 18" wide or 18" high   | \$ 81  | MOD-2OVER24 | Oversize Back Over 24" wide or 24" high              | \$ 136 |
| MOD-2OVER 22 | Oversize Back Over 22" wide or 22" high   | \$ 110 |             |  |        |
| 77PLATEMATE  | Headrest Adapter - Universal              | \$ 52  | 7733        | Adaptaplate - Headrest Adapter - Universal & Freedom | \$ 76  |

## Back Foam


|                |              |               |              |                        |                |
|----------------|--------------|---------------|--------------|------------------------|----------------|
| 3# Foam        | \$ 55 per 1" | Firm Sunmate  | \$127 per 1" | Pudgee                 | \$273 per 1"   |
| 4# Foam        | \$ 55 per 1" | X-Soft T-Foam | \$137 per 1" | Waffle                 | \$409 per 1/2" |
| 5# Foam        | \$ 75 per 1" | Soft T-Foam   | \$137 per 1" | Waffle (over 18x16)    | \$815 per 1/2" |
| Soft Sunmate   | \$127 per 1" | Medium T-Foam | \$137 per 1" | Free Float             | \$409 per 1/2" |
| Medium Sunmate | \$127 per 1" | Firm T-Foam   | \$137 per 1" | FreeFloat (over 18x16) | \$815 per 1/2" |



# Easy Conforming Back

## Dimensions - Standard

| Width | A    | B    | C  | D | E   | F  | G |
|-------|------|------|----|---|-----|----|---|
| 11/12 | 7    | 8    | 13 | 3 | 2   | 7  |   |
| 13/14 | 9.5  | 10.5 | 15 | 4 | 3   | 9  |   |
| 15/16 | 10.5 | 11.5 | 17 | 4 | 3   | 11 |   |
| 17/18 | 11   | 12.5 | 19 | 4 | 3   | 13 |   |
| 19/20 | 12   | 13.5 | 21 | 4 | 3.5 | 15 |   |
| 21/22 | 13   | 14.5 | 23 | 4 | 3.5 | 17 |   |
| 23/24 | 14   | 15.5 | 25 | 4 | 3.5 | 19 |   |
| 25/26 | 14   | 15.5 | 27 | 4 | 3.5 | 21 |   |



Comments:

## Custom Dimensions

|   | Left      | Right |
|---|-----------|-------|
| A |           |       |
| B |           |       |
| C |           |       |
| D |           |       |
| E |           |       |
| F | See Chart |       |
| G |           |       |

G Measurement on Standard is for Split-Wing measurement

## Dimensions - Extended

| Width | A   | B    | C  | D | E   | F  | G    |
|-------|-----|------|----|---|-----|----|------|
| 11/12 | 5   | 6    | 13 | 3 | 2   | 7  | 8    |
| 13/14 | 6   | 7    | 15 | 4 | 3   | 9  | 10.5 |
| 15/16 | 7   | 7.5  | 17 | 4 | 3   | 11 | 11.5 |
| 17/18 | 7.5 | 8    | 19 | 4 | 3   | 13 | 12.5 |
| 19/20 | 8   | 9    | 21 | 4 | 3.5 | 15 | 13.5 |
| 21/22 | 9   | 10.5 | 23 | 4 | 3.5 | 17 | 14.5 |
| 23/24 | 9   | 10.5 | 25 | 4 | 3.5 | 19 | 14.5 |
| 25/26 | 9   | 10.5 | 27 | 4 | 3.5 | 21 | 14.5 |

