

NEW!

Simple...

Freedom[®] Simply Seating

It's never been easier

**A single part number for the seat...
and a single part number for the back.**



Our kids are our passion.[™]


New from Freedom®
A simple solution
that offers so much.



Order the EB1-Kit
and the seat includes...

- ✓ **NEW Boltaron seat base platform**
 - fire retardant
 - incontinence proof
 - super strong, 1/2" thick
- ✓ **"T" nuts for tracks, hinge and abductor**
- ✓ **Tracks for hardware and accessories**
- ✓ **Seat attachment hardware**
- ✓ **Seat-to-back hinge**
- ✓ **4 lb. support foam**
- ✓ **NEW "Bravo Silver" anti-microbial cover***
 - 3 sides zippered
 - helps control bacterial growth
 - latex free
 - highly incontinence resistant
 - fire retardant
 - FDA approved

* Freedom's **NEW "Bravo Silver"** fabric is another important step we have taken to help keep kids healthy. An antibacterial fabric, "Bravo Silver" utilizes silver technology to create a barrier to help control bacteria growth.



Order the EB2-Kit
and the back includes...

- ✓ **NEW Boltaron back base platform**
 - fire retardant
 - incontinence proof
 - super strong, 1/2" thick
- ✓ **"T" nuts for tracks, hinge and headrest**
- ✓ **Tracks for hardware and accessories**
- ✓ **Back attachment hardware**
- ✓ **Seat-to-back hinge**
- ✓ **5 lb. support foam**
- ✓ **Height adjustable shoulder strap risers**
- ✓ **Universal headrest adapt-a-plate**
- ✓ **"Air Mesh" comfort cover**
 - 3 sides zippered
 - breathable in warm conditions
 - incontinence resistance liner option
 - latex free
 - fire retardant

NOTE: Non- "standard" options shown throughout include:

- **NEW 7711 Adjust-Eze** multi-positional removeable headrest
- Swing away lateral trunk supports
- Freedom Select anterior trunk support
- Fixed lateral hip supports
- Flip down medial thigh support



Our kids are our passion.™

Freedom Designs Incorporated
 2241 N. Madera Road, Simi Valley, CA 93065

phone: 800-331-8551

fax: 888-582-1509

e: CustomerService@freedomdesigns.com

www.freedomdesigns.com

Freedom Designs Incorporated (FDI) maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without notification. Check with FDI for the latest information.

Copyright © 2012 Freedom Designs, Inc. All rights reserved. Printed in the USA. FSSS 22012.



Manufactured in the U.S.A.