

Toddler Seating & Positioning Order Form

"New"

email it!



Freedom Designs, Incorporated
 2241 N Madera Road, Simi Valley, CA 93065
 phone: 800-331-8551 fax: 888-582-1509
 e: customerservice@freedomdesigns.com
 www.freedomdesigns.com

Buyer: _____
 Co./Account: _____
 Location: _____
 Contact #: _____
 P.O. #: _____
 Marked For: _____
 Fax: _____
 Email: _____

Quote Order

Size 1: Seat: 10" Wide x 10" Deep Back: 10" Wide x 10" Tall
 Size 2: Seat: 11" Wide x 11" Deep Back: 11" Wide x 11" Tall
 Size 3: Seat: 12" Wide x 12" Deep Back: 12" Wide x 12" Tall
 Custom Size: Seat: _____ W x _____ D Back: _____ W x _____ T


Seating Specifications

- A Back Height: _____
- B Lateral Height: _____
- C Lateral Width: _____
- D Hip Width (rear): _____
- E Set Up Seat Depth: _____
- F Hip Width (front): _____
- G "E" + 3" for Growth: _____
- H Seat Width: _____
- I Back Width: _____

Important Note: To ensure quote accuracy, seat and back width, depth and height dimensions must be provided.

Finished Dimensions

Standard cut is Boltaron with foam rolled over edge.



Wheelchair Specifications

Model: _____
 Width: _____
 Depth: _____

| Freedom Designs Upholstery Selections | Seat Base | Back Base | Headrest Pad | Lateral Pads | Hip & Knee Pads |
|---|---|--------------------------|--------------------------|--------------------------|--------------------------|
| 1. Choose Contact / Non-Contact 2. Choose Upholstery 3. Choose Component | | | | | |
| Body Contact Surface | | | | | |
| Bravo "Silver" Anti-Microbial Upholstery | | | | | |
| Choose Standard in Black Only | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Choose Reversed in Black Only | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Cushioned AirMesh Breathable Upholstery | | | | | |
| Choose Standard | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Choose Incontinence Resistance Liner | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| "Other" Upholstery Options | | | | | |
| Number: _____ Choose | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Body Non-Contact Surface | | | | | |
| Bravo Standard Upholstery | | | | | |
| Choose Standard in Black Only | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Choose Reversed in Black Only | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Cushioned AirMesh Breathable Upholstery | | | | | |
| Choose Standard | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Choose Incontinence Resistance Liner | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| "Other" Upholstery Options | | | | | |
| Number: _____ Choose | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| Monogram Option | | | | | |
| Name: _____ | When ordering choose either Monogram or Embroidery | | | | |
| Script <input type="checkbox"/> | | | | | |
| Block <input type="checkbox"/> | | | | | |
| Thread Color: _____ | | | | | |
| Embroidery Option | | | | | |
| Design: _____ | | | | | |
| <input type="checkbox"/> Freedom Install Program Purchase Order #: _____ of frame for install Freedom Install Program does not include power wheelchairs. | | | | | |



Toddler Seating & Positioning Order Form






Purchase Order #

Marked for:

Toddler Seating Construction Includes:





- 1/4" thick Boltaron base material up to 14" x 14" size • Cut full width unless otherwise specified • 2" x 2" back post notches • Full standard "T" nut pattern • 1/2" 4 lb. base foam with 1/2" soft Sunmate foam top layer is standard on Seat & Back • 3 sided zipper closed cover • Backs feature slightly rounded top and bottom corners • Backs are available with flat base material only
- **MAXIMUM WEIGHT CAPACITY 75 lbs.** •

Toddler Solid Seat Insert Packages (E2609)

- | | | | |
|---|---|---|---|
| <p><input type="checkbox"/> A EZ LITE1-TDLR 4-Point Toddler Ez-Lite Adjustable Mounting Hdwe With Seat to Back Hinge \$492</p> <p>4 - 504C Toddler EZ-Lite Adjustable Hardware</p> <p>2 - SEZFC Select EZ Frame Clamps</p> <p>1 - 430D Seat to Back Hinge Toddler 2" Offset</p> |  | <p><input type="checkbox"/> B FASTCLAMP1-TDLR 2-Point Toddler Fast-Clamp Adjustable Mounting Hdwe With Seat to Back Hinge \$484</p> <p>2 - 590A Toddler Fast Clamp Adjustable hardware</p> <p>2 - SEZFC Select EZ Frame Clamps</p> <p>2 - 504C Toddler EZ-Lite Adjustable Hardware</p> <p>1 - 430D Seat to Back Hinge Toddler 2" Offset</p> |  |
| <p><input type="checkbox"/> C EE1-TDLR 2-Point Mini Econo-Eze Adjustable Mounting Hdwe With Seat to Back Hinge \$633</p> <p>2 - 462G Mini Econo Eze Adjustable Knob Hardware</p> <p>2 - 439A Seat Mount Bracket Mini Econo Eze</p> <p>2 - 504C Toddler Ez-Lite Adjustable Hardware</p> <p>1 - 430D Seat to Back Hinge Toddler 2" Offset</p> |  | | |
| <p><input type="checkbox"/> MOD-1K Anti-Thrust Foam Modification \$187</p> <p>• Height of Anti-Thrust shelf is 1/2"</p> <p>• C= First 1/3 of set-up seat depth (see "E" page 1)</p> | | <p><input type="checkbox"/> MOD-1WW-TDLR Abduction Flared Toddler Seat Modification N/C</p> | |
|  | | <p><input type="checkbox"/> MOD-1Y Contour Foam Modification \$320</p> <p>• Adductor height is 1/2"</p> <p>• Height of Anti-Thrust shelf is 1/2"</p> <p>• C= First 1/3 of set-up seat depth (see "E" page 1)</p> | |
| | |  | |

Notes:

Toddler Solid Back Insert Packages (E2617)

- | | | | |
|--|---|--|---|
| <p><input type="checkbox"/> A FASTCLAMP2-TDLR 2-Point Toddler Fast-Clamp Adjustable Mounting Hdwe \$434</p> <p>2 - 590A Toddler Fast Clamp Adjustable Hardware</p> <p>2 - SEZFC Select EZ Frame Clamps</p> |  | <p><input type="checkbox"/> B EE2-TDLR 2-Point Mini Econo-Eze Adjustable Mounting Hardware \$633</p> <p>2 - 462G Mini Econo Eze Adjustable Knob Hardware</p> <p>2 - 437E Back Mount Bracket Mini Econo Eze</p> |  |
| <p><input type="checkbox"/> MOD-2H I-Cut Modification \$114</p> | | <p><input type="checkbox"/> MOD-2IA Curved Foam Modification \$103</p> | |
| <p><input type="checkbox"/> MOD-2J T-Cut Modification \$114</p> | | | |
|  | |  | |
| <p>A _____</p> <p>B _____</p> <p>C _____</p> <p>D _____</p> <p>E _____</p> | | <p>E _____</p> <p>F _____</p> | |

Notes:




Toddler Seating & Positioning Order Form

Purchase Order # _____


Marked for: _____

Toddler Medial Pelvic Support Package – Each

| Part Number | Description | HCCPS | Price | Shape | Standard Pad Sizing | Custom Pad Size |
|------------------------------------|---|---------|-------|---|--|--|
| <input type="checkbox"/> EE3A-TDLR | Depth Adjustable - Removable Hdwe.– Choice of Pad | (E0957) | \$232 | <input type="checkbox"/> Wedge <input type="checkbox"/> Oval | <input type="checkbox"/> Small – 3L x 3H x 3-2 Wedge |  A _____ B _____ C _____ D _____ |
| <input type="checkbox"/> MOD-3E | Medial Pelvic Support Pad Only | (E0957) | \$167 | <input type="checkbox"/> Round | | |
| <input type="checkbox"/> 438T | Removable, Adjustable Toddler Medial Hdwe | (E0128) | \$ 65 | | | |


Notes:

Toddler Lateral Pelvic Support Package – Pair

| Part Number | Description | HCCPS | Price | Location | Standard Pad Sizing | Custom Pad Size |
|------------------------------------|-----------------------------------|---------|-------|---|--|--|
| <input type="checkbox"/> EE4A-TDLR | Fixed, Flush Toddler Hdwe & Pads | (E0956) | \$232 | <input type="checkbox"/> Hip <input type="checkbox"/> Knee | <input type="checkbox"/> Set Up Seat Depth x 3H |  A _____ B _____ C _____ D _____ |
| <input type="checkbox"/> MOD-4B | Lateral Pelvic Pads Only | (E0956) | \$160 | | Pad Profile <input type="checkbox"/> Standard (1" thick) | |
| <input type="checkbox"/> MOD-4A | Contour Pad Modification | (K0108) | \$176 | | <input type="checkbox"/> Thin (1/2" thick) | |
| <input type="checkbox"/> 450T-2 | Fixed, Flush Toddler Lateral Hdwe | (K0108) | \$ 72 | | | |

Notes:

Toddler Lateral Thoracic Support Package – Pair

| Part Number | Description | HCCPS | Price | Shape | Pad Profile | Std Pad Size | Custom Pad Size |
|------------------------------------|--|-----------------|-------|---|--|----------------------------------|---|
| <input type="checkbox"/> EE5A-TDLR | Fixed, Flush Toddler Hdwe & Pads | (E0956) | \$232 | <input type="checkbox"/> Flat | <input type="checkbox"/> Standard (1" thick) | <input type="checkbox"/> 4L x 3H | A _____ B _____ |
| <input type="checkbox"/> EE5E-TDLR | Lift & Swing-Away, Flush Toddler Hdwe & Pads | (E0956 & E1028) | \$465 | <input type="checkbox"/> Curved <input type="checkbox"/> Reverse Curve | <input type="checkbox"/> Thin (1/2" thick) | <input type="checkbox"/> 5L x 4H | |
| <input type="checkbox"/> MOD-5G | Lateral Thoracic Pads Only | (E0956) | \$168 | | | |  |
| <input type="checkbox"/> MOD-5F | Protective Boots for Hardware | (K0108) | \$ 53 | | | | |
| <input type="checkbox"/> MOD-5J | Extended Hardware | (K0108) | \$ 82 | | | | |
| <input type="checkbox"/> 450T-2 | Fixed, Flush Toddler Lateral Hdwe | (K0108) | \$ 72 | | | | |
| <input type="checkbox"/> 455T-2 | Lift, Swing-Away Toddler Lateral Hdwe | (E1028) | \$298 | | | | |

Notes:

Toddler Headrest Support – Each

| Part Number | Headrest Hardware Assemblies | HCCPS | Price | Headrest Pads | | |
|------------------------------------|--|---------|-------|---|--|---|
| <input type="checkbox"/> 7700-TDLR | Height, Depth Adjustable - Removable Toddler Hdwe | (E1028) | \$141 | | | |
| <input type="checkbox"/> 7714-TDLR | Height, Depth, & Angle Adjustable - Removable Toddler Hdwe | (E1028) | \$180 | | | |
| <input type="checkbox"/> 7756 | Chubbie Headrest Pad Only | (E0955) | \$184 | <input type="checkbox"/> X-Sm – 8W x 3H | <input type="checkbox"/> Small – 10W x 4H | <input type="checkbox"/> Med – 12W x 5H |
| <input type="checkbox"/> 7760 | Skinny Headrest Pad Only | (E0955) | \$184 | <input type="checkbox"/> X-Sm – 6W x 3.5H | <input type="checkbox"/> Small – 8W x 4.5H | <input type="checkbox"/> Med – 10W x 4.5H |

Notes:




Toddler Seating & Positioning Order Form

Purchase Order #

Freedom Select Anterior Pelvic Supports

Marked for:




| FREEDOM SELECT 2-POINT PELVIC POSITIONING BELTS - PADDED | | | | | | | | |
|--|------|----------|-----|----|------|---|---------|--|
| PUSH BUTTON BUCKLE | | | | | | | | |
| Part Number | Size | Buckle | A | B | C | D | Price | |
| <input type="checkbox"/> SFDI-3A2MP-XS | XS | MICRO 1" | 1 | 50 | 1.75 | 6 | \$75.00 | |
| <input type="checkbox"/> SFDI-3A2MP-S | S | MICRO 1" | 1 | 55 | 2.25 | 7 | \$75.00 | |
| <input type="checkbox"/> SFDI-3A2P-S | S | STD 1.5" | 1.5 | 55 | 2.25 | 7 | \$75.00 | |
| <input type="checkbox"/> SFDI-3A2P-CUSTOM | | QUOTE | | | | | \$95.00 | |

| FREEDOM SELECT 2-POINT PELVIC POSITIONING BELTS - PADDED | | | | | | | | |
|--|------|----------|-----|----|------|---|---------|--|
| SIDE RELEASE BUTTON BUCKLE | | | | | | | | |
| Part Number | Size | Buckle | A | B | C | D | Price | |
| <input type="checkbox"/> SFDI-3A2MS-XS | XS | MICRO 1" | 1 | 50 | 1.75 | 6 | \$55.00 | |
| <input type="checkbox"/> SFDI-3A2MS-S | S | MICRO 1" | 1 | 55 | 2.25 | 7 | \$55.00 | |
| <input type="checkbox"/> SFDI-3A2S-S | S | STD 1.5" | 1.5 | 55 | 2.25 | 7 | \$55.00 | |
| <input type="checkbox"/> SFDI-3A2S-CUSTOM | | QUOTE | | | | | \$75.00 | |

Notes:

Freedom Select Anterior Trunk Supports



| Part Number | Description | A | B | C | D | E | F | G | Price |
|---|--|-----|-----|------|---|-----|------|-------|--------|
| <input type="checkbox"/> SFDI-6SNZ-S | Select <i>Standard</i> Small | 7.5 | 4.5 | 3.25 | 6 | 11 | 5.25 | 13.75 | \$ 130 |
| <input type="checkbox"/> SFDI-6SNZ-1-S | Select <i>Standard</i> Dual Rear-Pull Small | 7.5 | 4.5 | 3.25 | 6 | 11 | 5.25 | 13.75 | \$ 130 |
| <input type="checkbox"/> SFDI-6NNZ-XS | Select <i>Narrow</i> X-Small | 6 | 3 | 2 | 2 | 6.5 | 4.25 | 10.5 | \$ 130 |
| <input type="checkbox"/> SFDI-6NNZ-1-XS | Select <i>Narrow</i> Dual Rear-Pull X-Small | 6 | 3 | 2 | 2 | 6.5 | 4.25 | 10.5 | \$ 130 |
| <input type="checkbox"/> SFDI-6NNZ-S | Select <i>Narrow</i> Small | 7.5 | 4 | 2.75 | 3 | 11 | 5.25 | 13.75 | \$ 130 |
| <input type="checkbox"/> SFDI-6NNZ-1-S | Select <i>Narrow</i> Dual Rear-Pull Small | 7.5 | 4 | 2.75 | 3 | 11 | 5.25 | 13.75 | \$ 130 |
| <input type="checkbox"/> SFDI-6SZ-S | Select <i>Standard</i> Zipper Small | 7.5 | 4.5 | 3.25 | 6 | 11 | 5.25 | 13.75 | \$ 135 |
| <input type="checkbox"/> SFDI-6SZ-1-S | Select <i>Standard</i> Zipper Dual Rear-Pull Small | 7.5 | 4.5 | 3.25 | 6 | 11 | 5.25 | 13.75 | \$ 135 |
| <input type="checkbox"/> SFDI-6SZ-S | Select <i>Narrow</i> Zipper Small | 7.5 | 4 | 2.75 | 3 | 11 | 5.25 | 13.75 | \$ 135 |
| <input type="checkbox"/> SFDI-6SZ-1-S | Select <i>Narrow</i> Zipper Dual Rear-Pull Small | 7.5 | 4 | 2.75 | 3 | 11 | 5.25 | 13.75 | \$ 135 |
| <input type="checkbox"/> SFDI-6CUSTOM | Select Anterior Support Custom Sizing | | | | | | | | \$ 190 |

Notes:



*** FREE * FREE * FREE * FREE * FREE * FREE ***

Freedom Designs & Pindot

Are Pleased To Announce

“A Sizzling Summer Sensation”

And For a Limited Time It's...

FREE * FREE * FREE

When ordering any

Freedom Wheelchair and

a complete Freedom or Pindot

seating system including seat,

back, and attachment hardware,

simply provide the frame Quote

or frame PO number

and receive a

Freedom Canopy and

Freedom Backpack

*** FREE * FREE * FREE * FREE * FREE * FREE ***