



Our kids are our passion.™

Freedom 2Kids

Manual Wheelchair



Yes, you can.

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: _____ Buyer: _____
 Company: _____
 Marked For: _____
 Email: _____

Quote Order PO# _____
 Location: _____
 Phone: _____
 Fax: _____

| Description | Part No. | Retail | Description | Part No. | Retail |
|-------------|----------|--------|-------------|----------|--------|
|-------------|----------|--------|-------------|----------|--------|

Frame Options

FD2K Folding (E1236) **FD2K-30270 \$2085**
Transport Optional or (K0005)

Note: Transportation Option selected separately below.
Note: Includes Freedom Future Fit Program.

FD2K Rigid (E1235) **FD2K-30270R \$2085**
Transport Optional or (K0005)

Note: Transportation Option selected separately below.
Note: Includes Freedom Future Fit Program.

Transportation Option

Note: Used when Transport Optional frame is selected above.

Transportation Package **FD2K-99998 \$ 250**

Note: Must be installed by Freedom Designs manufacturing.
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.

Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

Frame Colors

| | | | |
|-----------------|-----|----------------|-----|
| Sugar Cookie | n/c | Pacific Blue | n/c |
| Neon Yellow | n/c | Blue Blast | n/c |
| Yellow | n/c | Electric Blue | n/c |
| Bright Pink | n/c | Galaxy Blue | n/c |
| Pink Parfait | n/c | Passion Purple | n/c |
| Sugar Plum | n/c | Octane Orange | n/c |
| Candy Violet | n/c | Red Sunset | n/c |
| Candy Raspberry | n/c | Really Red | n/c |
| Extreme Green | n/c | Fiery Red | n/c |
| Sparkle Green | n/c | Twilight | n/c |
| Mystic Teal | n/c | Stormy Gray | n/c |
| | | Black | n/c |

Frame Dimensions

Specify Frame Width Setting

| | | |
|----------|------------|-----|
| 11" Wide | FD2K-30210 | n/c |
| 12" Wide | FD2K-30220 | n/c |
| 13" Wide | FD2K-30230 | n/c |
| 14" Wide | FD2K-30240 | n/c |
| 15" Wide | FD2K-30290 | n/c |
| 16" Wide | FD2K-30260 | n/c |
| 17" Wide | FD2K-30330 | n/c |
| 18" Wide | FD2K-30280 | n/c |

Note: Frame widths measured across and from outside of seat frame tubes.

Back Post Locations

Specify Frame Type and Frame Length

| | | |
|---|-------------------|------------|
| Regular Frame Length - Adjusts 10"-14" | FD2K-30544 | n/c |
| 10" | FD2K-10 | n/c |
| 11" | FD2K-11 | n/c |
| 12" | FD2K-12 | n/c |
| 13" | FD2K-13 | n/c |
| 14" | FD2K-14 | n/c |

Specify Back Post Location (if other than specified above) _____"

Note: 24" wheels not available with Regular frame.

Long Frame Length - Adjusts 13"-17" **FD2K-30545** **n/c**

| | | |
|-----|---------|-----|
| 13" | FD2K-13 | n/c |
| 15" | FD2K-15 | n/c |
| 16" | FD2K-16 | n/c |
| 17" | FD2K-17 | n/c |

Specify Back Post Location (if other than specified above) _____"

Note: All 24" wheels only with 3", 4" or 5" casters and Long frame.

Note: FD2K has 5" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating/hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back/hardware thickness of 2" = back post location set at 14".

Sling Seat Upholstery

| | | |
|------------------------|-------------------|------------|
| Omit Seat Sling | FD2K-39004 | |
| Seat Sling | FD2K-39011 | n/c |

Note: Not Available with ANY Single Post armrests.

Note: See Page 4 for Seat Pan options.

Specify Seat Sling Length

| | | |
|-----|------------|-------|
| 10" | FD2K-30500 | \$ 56 |
| 11" | FD2K-30501 | \$ 56 |
| 12" | FD2K-30502 | n/c |
| 13" | FD2K-30503 | \$ 56 |
| 14" | FD2K-30504 | n/c |
| 15" | FD2K-30505 | \$ 56 |
| 16" | FD2K-30506 | n/c |
| 17" | FD2K-30507 | \$ 56 |

Seat Cushions

| | | |
|------------------------|-------------------|--------------|
| 1" Seat Cushion | FD2K-39003 | \$ 87 |
| 2" Seat Cushion | FD2K-39007 | \$ 97 |

Note: Includes 4# foam & black knit cover with zipper closure.

Back Post Tubes

Specify Push Handle, Back Post and Extensions

Straight Back Post Tubes - No Push Handles **FD2K-30600** **n/c**

Specify Back Post Height

| | | |
|----------|------------|-----|
| 10" High | FD2K-30610 | n/c |
| 11" High | FD2K-30611 | n/c |
| 12" High | FD2K-30612 | n/c |
| 13" High | FD2K-30613 | n/c |
| 14" High | FD2K-30614 | n/c |
| 15" High | FD2K-30615 | n/c |





FD2K Rider Weight Capacity is 250 lbs.

Page 1 of 5

Freedom 2 Kids

€€

Date: _____ PO# _____ Marked For: _____

| Description | Part No. | Retail | Description | Part No. | Retail |
|---|------------|--------|--|---------------------------|--------|
| Back Post Tubes (continued) | | | Back Post Tubes (continued) | | |
| Standard 90° Bend Push Handles and Back Post Tubes | FD2K-30518 | n/c | Sling Back Upholstery | | |
| Ergo Bend Push Handles and Back Post Tubes | FD2K-30519 | n/c | Omit Sling Back Upholstery | FD2K-39005 | |
| Note: All push/stroller handle extensions not recommended with Ergo Bend. | | | Standard Sling Back Upholstery | FD2K-39006 | n/c |
| 8° Bend Back Canes | FD2K-30511 | \$ 26 | Note: Not available with Flip-Down Rigidizing Bar. Note: Sling back height must be selected below. | | |
| Back Post Height | | | Tension Adjustable Sling Back Uph. | FD2K-39008 | \$ 236 |
| Note: Back Posts are measured from top of seat rail to top of handle bend. | | | Note: Not available with Flip-Down Rigidizing Bar. Note: Sling back height must be selected below. | | |
| Short Push Handle - 20" high | FD2K-30515 | n/c | 10" High | FD2K-30910 | n/c |
| Medium Push Handle - 23" high | FD2K-30514 | n/c | 11" High | FD2K-30911 | n/c |
| Tall Push Handle - 26" high | FD2K-30516 | n/c | 12" High | FD2K-30912 | n/c |
| | | | 13" High | FD2K-30913 | n/c |
| | | | 14" High | FD2K-30914 | n/c |
| | | | 15" High | FD2K-30915 | n/c |
| | | | 16" High | FD2K-30916 | n/c |
| | | | 17" High | FD2K-30917 | n/c |
| | | | 18" High | FD2K-30918 | n/c |
| | | | Note: Sling back height must be selected when ordering Standard or Adjustable Tension Sling Back Upholstery. | | |
| Stroller Handle Extensions | | | Legrests, Footrests and Extensions | | |
| Note: All stroller handle ext. add approximately 6.5" to push handle height. | | | Legrests | | |
| Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad. | | | Omit Legrest Hangers & Footrests | FD2K-33000 | |
|  | | | 75° Pop-Up Swing Away Hangers | FD2K-30284 | n/c |
| -30560 | | | Note: Not available with high mount bracket. | | |
|  | | | 90° Pop-Up Swing Away Hangers | FD2K-30281 | n/c |
| -30570/5 | | | Note: Not available with 8" adjustable extensions. | | |
|  | | | 90° Pop-Up Swing Away Elevating Legrests (ELRs) | FD2K-30279 | \$ 255 |
| -30580/5 | | | Note: Not available with all length adjustable extensions. Note: Not available with 6.5" heavy duty extensions. | | |
|  | | | 75° Latch Style Swing Away Hangers | FD2K-30282 | n/c |
| -30517 | | | Note: Not available with high mount bracket. | | |
| Push Handle Extension - Width Adjustable | FD2K-30560 | n/c | 90° Latch Style Swing Away Hangers | FD2K-30283 | n/c |
| Note: Not removable, but allows frame to fold. Note: Mounts direct to back post push handle. | | | Note: Not available with 8" adjustable extensions. | | |
| Push Button, Angle-Adjustable Folding Stroller Handle Extension | FD2K-30570 | \$ 112 | Footplates & Footboards | | |
| Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below. Note: Not removable, but allows frame to fold. Note: Mounts direct to back post push handle. | | | Omit Footrest Plates | FD2K-36404 | |
| Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty | FD2K-30575 | \$ 179 | Composite Flip-Up | FD2K-17070 | n/c |
| Note: Available 10W-22W wheelchair frames. Note: Not removable, but allows frame to fold. Note: Same as above but with metal hinge mechanism. Note: Mounts direct to back post push handle. | | | Kids Aluminum Angle Adjustable, Flip-up | (K0040) FD2K-36020 | \$ 160 |
| Push Button, Angle-Adjustable, Folding Push Handle Extension | FD2K-30580 | \$ 112 | Note: For frames 14" wide and narrower. | | |
| Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below. Note: Not removable, but allows frame to fold. Note: Mounts direct to straight back post tubes. | | | Heel Loops | (E0951) FD2K-36022 | \$ 41 |
| Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty | FD2K-30585 | \$ 179 | Calf Strap - Webbing | (K0038) FD2K-36024 | \$ 41 |
| Note: Available 10W-22W wheelchair frames. Note: Not removable, but allows frame to fold. Note: Same as above but with metal hinge mechanism. Note: Mounts direct to straight back post tubes. | | | Adult Aluminum Angle Adj. Flip-up w/ Heel Loops | (K0040/ E0951) FD2K-36021 | \$ 180 |
| Individual Stroller Handle Extensions Removable - Pair | FD2K-30517 | \$ 102 | Note: For frames 15" wide and greater. | | |
| Note: Not available with Straight Back Posts. | | | Omit Heel Loops | FD2K-36023 | n/c |
| | | | Aluminum Angle Adjustable Footboard Removable | (K0040) FD2K-36029 | \$ 160 |

Continue To Next Page

Date: _____ PO# _____ Marked For: _____

| Description | Part No. | Retail | Description | Part No. | Retail |
|--|--------------|--------|---|------------|--------|
| Legrests, Footrests and Extensions(continued) | | | Armrests and Armrest Pads (continued) | | |
| MAPS Footplate Mounting Hardware | | | Frame Mounted Armrests | | |
| MAPS (Multi Axis Positioning System) - Pair | FD2K-36040 | \$ 306 | Single Post Fixed Height, Removable - Short (K0015) | FD2K-38222 | \$ 225 |
| MAPS (Multi Axis Positioning System) - Left | FD2K-36040-L | \$ 153 | Single Post Fixed Height, Removable - Tall (K0015) | FD2K-38224 | \$ 225 |
| MAPS (Multi Axis Positioning System) - Right | FD2K-36040-R | \$ 153 | Note: Single Post Fixed Height arms include a removable side guard. | | |
| <p>Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.</p> <p>Note: Must choose aluminum footplates.</p> <p>Note: Not available on 10" & 11" wide frames</p> <p>Note: MAPS not available with composite footplates or aluminum footboard.</p> | | | Note: Adjustable in height with tools only. See matrix for height ranges. | | |
| Legrest Extensions | | | Note: Not available with ANY clothing guard accessory. | | |
| 4" Short Length Adjustable Extensions | FD2K-30285 | n/c | Single Post Height Adjustable, Removable - Short (E0973) | FD2K-38217 | \$ 265 |
| 6" Medium Length Adjustable Extensions | FD2K-30286 | n/c | Single Post Height Adjustable, Removable - Tall (E0973) | FD2K-38203 | \$ 265 |
| 8" Long Length Adjustable Extensions | FD2K-30287 | n/c | Note: Not available with ANY clothing guard accessory. | | |
| <p>Note: Not available with all 90° hangers.</p> | | | Desk Length Pads | FD2K-38213 | n/c |
| 4.5" Heavy Duty ELR Extensions | FD2K-36401 | \$ 56 | Full Length Pads | FD2K-38214 | n/c |
| 4.5" Heavy Duty Extensions | FD2K-36402 | \$ 56 | Specify Seat Rail to Top of Armrest Measurement _____ " | | |
| <p>Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p> | | | Front Wheels - Casters | | |
| 6.5" Heavy Duty Extensions | FD2K-36403 | \$ 66 | Omit Caster Wheels | FD2K-34333 | |
| <p>Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p> | | | 3" Polyurethane | FD2K-34888 | \$ 56 |
| High Mount Brackets | FD2K-36407 | \$ 102 | 3" Polyurethane Light-Up | FD2K-34111 | \$ 56 |
| <p>Note: Includes flip-up mounting hardware for future growth adjustment.</p> <p>Note: Not available with all 75° hangers.</p> | | | <p>Note: 3" casters recommended for low seat to floor heights.</p> | | |
| Specify Knee to Heel Measurement _____ " (see matrix page) | | | 4" Polyurethane | FD2K-40000 | \$ 66 |
| Armrests and Armrest Pads | | | 4" Polyurethane Light-Up | FD2K-41000 | \$ 101 |
| Omit Armrests | FD2K-38209 | | 4" x 1.4" Poly, Aluminum Mag | FD2K-43000 | \$ 92 |
| Back Post Mounted Armrests | | | 5" x 1" Poly, Plastic Mag | FD2K-34555 | n/c |
| <p>Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.</p> | | | 5" Polyurethane Light-Up | FD2K-35050 | \$ 101 |
| 10.5" Long Flip-Back | FD2K-38204 | n/c | 5" x 1.5" Poly, Plastic Mag | FD2K-34222 | \$ 87 |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | 5" x 1.4" Poly, Aluminum Mag | FD2K-42000 | \$ 92 |
| 13.5" Long Flip-Back | FD2K-38202 | n/c | 6" x 1" Poly, Plastic Mag | FD2K-34666 | n/c |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | 6" x 1.5" Poly, Plastic Mag | FD2K-34999 | \$ 87 |
| 10.5" Long Flip-Back Locking (E0973) | FD2K-38201 | \$ 245 | 6" x 1.4" Poly, Aluminum Mag | FD2K-44000 | \$ 92 |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | 6" x 2" Poly with Tread | FD2K-34444 | \$ 92 |
| 13.5" Long Flip-Back Locking (E0973) | FD2K-38200 | \$ 245 | 6" x 1" Pneumatic Plastic Mag (E2214/E2215) | FD2K-34667 | \$ 77 |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | 7" x 1" Poly, Plastic Mag | FD2K-34777 | \$ 77 |
| 10.5" to 13.5" Long Flip-Back, Telescoping (E0973) | FD2K-38208 | \$ 240 | 7" x 2" Polyurethane | FD2K-34200 | \$ 118 |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | <p>Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock.</p> | | |
| 10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973) | FD2K-38210 | \$ 265 | Caster Stem Spacer - 0.75" | FD2K-34778 | n/c |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | Rear Wheels - Tires - Accessories | | |
| 17" to 20" Long Flip-Back Locking, Telescoping (E0973) | FD2K-38215 | \$ 265 | <p>Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.</p> | | |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | Omit Rear Wheels and Axles | FD2K-31524 | |
| Foam Wraps - Standard | FD2K-38216 | n/c | Mag Wheels | | |
| Desk Length Pads | FD2K-38213 | n/c | 12" Low Profile Poly Mag (E2218) | FD2K-31526 | \$ 36 |
| Full Length Pads | FD2K-38214 | n/c | <p>Note: Not available with Hand Activated Drum Brake/Wheel Lock.</p> <p>Note: Not available with airless inserts.</p> <p>Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock.</p> | | |
| <p>Note: Full length pad on 13.5" or longer armrests only.</p> | | | 12" Semi Pneumatic Plastic Mag (K0069) | FD2K-31523 | \$ 168 |
| | | | <p>Note: Not available with airless inserts.</p> <p>Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock.</p> | | |
| | | | 12" Pneumatic Plastic Mag | FD2K-31525 | n/c |
| | | | 12" Airless Inserts (E2213) | FD2K-31019 | \$ 77 |
| | | | <p>Note: All 12" wheels not available with all handrims.</p> <p>Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock.</p> | | |

Date: _____ PO# _____

Marked For: _____

| Description | Part No. | Retail | Description | Part No. | Retail |
|--|------------|--------|---|------------|-------------------------|
| Rear Wheels - Tires - Accessories (continued) | | | Rear Wheels - Tires - Accessories (continued) | | |
| 16" Pneumatic Plastic Mag | FD2K-31530 | \$ 41 | <i>Note: 16" & 18" spoke wheels do not include handrims.</i> | | |
| 16" Airless Inserts (E2213) | FD2K-31020 | \$ 77 | 20" Pneumatic Spoke Wheel with Aluminum Handrim | FD2K-31540 | \$ 204 |
| 16" Aluminum Handrims | FD2K-31014 | \$ 169 | 20" Airless Inserts (E2213) | FD2K-31022 | \$ 77 |
| <i>Note: 16" mag wheels not available with plastic coated handrims.</i> | | | <i>Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.</i> | | |
| 18" Pneumatic Plastic Mag | FD2K-31533 | \$ 41 | 22" Pneumatic Spoke Wheel with Aluminum Handrim | FD2K-31541 | \$ 204 |
| 18" Airless Inserts (E2213) | FD2K-31021 | \$ 77 | 22" Airless Inserts (E2213) | FD2K-31023 | \$ 77 |
| 18" Aluminum Handrims | FD2K-31015 | \$ 169 | <i>Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.</i> | | |
| <i>Note: 18" mag wheels not available with plastic coated handrims.</i> | | | 24" Pneumatic Spoke Wheel with Aluminum Handrim | FD2K-31542 | \$ 204 |
| 20" Pneumatic Plastic Mag with Aluminum Handrim | FD2K-31527 | n/c | 24" Airless Inserts (E2213) | FD2K-31024 | \$ 77 |
| 20" Airless Inserts (E2213) | FD2K-31022 | \$ 77 | <i>Note: All 24" wheels not available with Regular frame.</i> | | |
| 20" Plastic Coated Handrims | FD2K-31016 | \$ 119 | <i>Note: All 24" wheels only with 3", 4" or 5" casters and Long frame.</i> | | |
| 22" Pneumatic Plastic Mag with Aluminum Handrim | FD2K-31528 | n/c | <i>Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.</i> | | |
| 22" Airless Inserts (E2213) | FD2K-31023 | \$ 77 | Omit Handrims from 20", 22" or 24" Wheels | FD2K-31006 | |
| 22" Plastic Coated Handrims | FD2K-31017 | \$ 119 | Natural Fit Handrims, Super Grip - 22" Spoke wheel | FD2K-31009 | \$ 347 |
| 24" Pneumatic Plastic Mag with Aluminum Handrim | FD2K-31529 | n/c | Natural Fit Handrims, Super Grip - 24" Spoke wheel | FD2K-31010 | \$ 347 |
| 24" Airless Inserts (E2213) | FD2K-31024 | \$ 77 | Camber Spacer Kit | FD2K-31003 | \$ 56 |
| 24" Plastic Coated Handrims | FD2K-31018 | \$ 119 | 3 Degrees | FD2K-31011 | n/c |
| <i>Note: All 24" wheels not available with Regular frame.</i> | | | 6 Degrees | FD2K-31012 | n/c |
| <i>Note: All 24" wheels only with 3", 4" or 5" casters and Long frame.</i> | | | | | |
| Omit Handrims from 20", 22" or 24" Wheels | FD2K-31006 | | Wheel Locks | | |
| Natural Fit Handrims, Super Grip - 24" Mag wheel | FD2K-31008 | \$ 347 | Push to Lock - Standard | FD2K-31501 | n/c |
| | | | Pull to Lock | FD2K-31506 | n/c |
| | | | Extension Handles (E0961) | FD2K-31801 | \$ 61 |
| | | | Hand Activated Brake and Wheel Lock System (E2228) | FD2K-31521 | \$ 700 |
| | | | <i>Note: Includes manual wheel locks to be used with public transportation.</i> | | |
| | | | <i>Note: Not available with 12" mag wheels with low pro poly tires.</i> | | |
| | | | <i>Note: Used to control direction, speed, brake and lock wheels.</i> | | |
| | | | <i>Note: Activation handles are mounted to the back post push handles.</i> | | |
| | | | Rear Anti-Tippers | | |
| | | | Anti-Tip Assembly - Fixed (E0971) | FD2K-31600 | \$ 97 |
| | | | Anti-Tip Assembly - Swing-Away (E0971) | FD2K-31602 | \$ 97 |
| | | | Seat Pan Options | | |
| | | | Fixed Solid Seat Pan - Bolteron (E2231) | FD2K-39012 | \$ 150 |
| | | | <i>Note: Not available with folding frame.</i> | | |
| | | | Fixed Solid Seat Pan - Aluminum (E2231) | FD2K-39030 | \$ 170 |
| | | | <i>Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".</i> | | |
| | | | <i>Note: Not available with folding frame.</i> | | |
| | | | Removable Solid Seat Pan - Bolteron (E2231) | FD2K-39013 | \$ 170 |
| | | | Removable Solid Seat Pan - Aluminum (E2231) | FD2K-39032 | \$ 255 |
| | | | <i>Note: Removable Seat Pan raises the effective seat to floor height by 5/8".</i> | | |
| | | | <i>Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.</i> | | |
| | | | <i>Note: All seat pans include 2" x 2" notches for growth adjustment.</i> | | |
| | | | <i>Note: All seat pans made to fit frame dimensions unless specified below.</i> | | |
| | | | Seat Width _____ | X | Seat Depth _____ |

Date: _____ PO# _____

Marked For: _____

| Description | Part No. | Retail |
|-------------|----------|--------|
|-------------|----------|--------|

| Description | Part No. | Retail |
|-------------|----------|--------|
|-------------|----------|--------|

Accessories

- Clothing Guards – Fabric** FD2K-CGF \$ 71
Note: Not available with ANY single post armrests.
- Clothing Guards – Quick-Release, Rigid** FD2K-CGR \$ 122
Note: Not available with ANY single post armrests.
- Medical Necessities Pouch** FD2K-23207 \$ 68
Note: Black nylon.
- Footrest Hanger Pads - Pair** SFDI-9HP \$ 92
Note: Wrap-around pads made from black neoprene.
- Frame Pads - Pair** 10911-2 \$ 82
Note: Wrap-around pads made from black neoprene.
Note: Typically used to pad exposed back post tubes.
- Flip-Down Rigidizing Bar** FD2K-23210 \$ 129
Note: Mounts to back post tubes to add rigidity to frame.
Note: Not available with ANY Sling Back Upholstery.
Note: Does NOT adjust in width to grow with frame.
- Freedom Backpack** FD2K-10900 \$ 56
Note: Black Nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom Backpack.
- Freedom Canopy - Flip-Back** FD2K-10905 \$ 128
Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.

Anterior Support Surfaces – (Trays)

- 3/8" Lexan Tray - Channel Bracket Hardware** (E0950) 08-2000 \$ 271
- 3/8" Lexan Tray - Adjustable Knob Hardware & Ribs** (E0950) 08-6000 \$ 219
- 1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs** (E0950) 08-13000 \$ 183
- 1/4" Lexan Tray - Hook & Loop Straps** (E0950) 08-16000 \$ 106

Footrest / Extension Lengths

| Extension/Hanger | 75° Pop | 75° Lch | 90° Pop | 90° Lch | ELR (90°) |
|------------------|---------|-------------|---------|---------|-----------|
| Hanger Only | 6"-8" | 5.5"-8.5" | 4"-8" | 4"-8" | 6"-10" |
| High Mount | N/A | N/A | 1.5"-3" | N/A | N/A |
| 4.5" H D Ext. | 9"-12" | 9.5"-12.5" | 9"-12" | 9"-12" | 13"-15" |
| 6.5" H D Ext. | 9"-15" | 9.5"-15.5" | 9"-15" | 9"-15" | N/A |
| 4" Adjustable | 10"-13" | 9.5"-12.5" | 11"-13" | 10"-13" | N/A |
| 6" Adjustable | 12"-15" | 10.5"-14.5" | 13"-15" | 12"-15" | N/A |
| 8" Adjustable | 14"-17" | 12.5"-16.5" | N/A | N/A | N/A |

MAPS Footrest / Extension Lengths

| Extension/Hanger | 75° Pop | 75° Lch | 90° Pop | 90° Lch | ELR (90°) |
|------------------|-------------|---------|-------------|------------|--------------------|
| Hanger Only | 5.5"-7.5" | 5"-8" | 5.5"-7.5" | 4.5"-7.5" | 4.5"- 8.5" & 11.5" |
| High Mount | N/A | N/A | N/A | N/A | N/A |
| 4.5" H D Ext. | 8.5"-11.5" | 9"-12" | 8.5"-11.5" | 8.5"-11.5" | 12.5"-15.5" |
| 6.5" H D Ext. | 8.5"-13.5" | 9"-14" | 8.5"-13.5" | 8.5"-13.5" | N/A |
| 4" Adjustable | 8.5"-11.5" | 9"-12" | 9.5"-11.5" | 8.5"-11.5" | N/A |
| 6" Adjustable | 10.5"-13.5" | 10"-14" | 11.5"-13.5" | 8.5"-13.5" | N/A |
| 8" Adjustable | 12.5"-15.5" | 12"-16" | N/A | N/A | N/A |

Note: Length from top of footrest plate to top of seat rail.
Note: Without extensions bottom hanger hole not accessible for hardware.
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

FD2K Rider Weight Capacity is 250 lbs.

Freedom Designs Armrest Heights

| Armrest and Pad | Foam | Desk | Full |
|---------------------------------|--------------|--------------|--------------|
| Flip Back | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Flip Back Locking | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Telescoping Flip Back | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Flip Back Locking Telescoping | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Single Post Height Adj. Short | N/A | 8" - 14" | 8" - 14" |
| Single Post Height Adj. Tall | N/A | 10" - 16" | 10" - 16" |
| Dual Post Fixed Height | N/A | 9.5" | 9.5" |
| Dual Post Flip Back Height Adj. | N/A | 9" - 13" | 9" - 13" |

Note: All heights from top of seat frame to top of armrest pad.

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Choose Seat to Floor Height

Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, and front caster size.
Note: Caster sizes vary slightly in height by type.
Note: Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".
Note: Seat pan choice increases seat rail to floor height 3/8" to 3/4".

FD2K Wheel / Caster Seat to Floor Heights

| Wheel | Caster | Height | Arc Over Rail |
|-------|--------|-------------|---------------|
| 12" | 3" | 14"-14.5" | 0" |
| 12" | 4" | 15" | 0" |
| 12" | 5" | 17" | 0" |
| 12" | 6" | 17.5" | 0" |
| 16" | 3" | 13"-14.5" | 3"-1.5" |
| 16" | 4" | 13.5"-14" | 2.5"-2" |
| 16" | 5" | 14"-14.5" | 2"-1.5" |
| 16" | 6" | 14.5"-15" | 1.5"-1" |
| 16" | 7" | 15"-15.5" | 1" |
| 18" | 3" | 14"-15.5" | 4"-2.5" |
| 18" | 4" | 15.5" | 3" |
| 18" | 5" | 16" | 2" |
| 18" | 6" | 16.5" | 1.5" |
| 18" | 7" | 17" | 1" |
| 20" | 3" | 16"-17.5" | 4"-3.5" |
| 20" | 4" | 15" | 3.5" |
| 20" | 5" | 15.5"-17.5" | 4.5"-2.5" |
| 20" | 6" | 17"-17.5" | 3"-2.5" |
| 20" | 7" | 18" | 2" |
| 22" | 3" | 15.5" | 6.5" |
| 22" | 4" | 16"-17" | 6"-5" |
| 22" | 5" | 16.5"-17.5" | 5.5"-4.5" |
| 22" | 6" | 17"-18.5" | 5"-3.5" |
| 22" | 7" | 18"-18.5" | 4"-3.5" |
| 24" | 3" | 16" | 8" |
| 24" | 4" | 16.5"-17.5" | 7.5"-6.5" |
| 24" | 5" | 17"-18" | 7"-6" |
| 24" | 6" | N/A | N/A |
| 24" | 7" | N/A | N/A |