



Our kids are our passion.™

# Freedom 2

# Manual Wheelchair



Yes, you can.™

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: \_\_\_\_\_ Buyer: \_\_\_\_\_  
 Company: \_\_\_\_\_  
 Marked For: \_\_\_\_\_  
 Email: \_\_\_\_\_

Quote Order PO# \_\_\_\_\_  
 Location: \_\_\_\_\_  
 Phone: \_\_\_\_\_  
 Fax: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Frame Options

**FD2 Folding** (E1236) FD2-30270 **\$2085**  
**Transport Optional** or (K0005)  
 Note: Transportation Option selected separately below.  
 Note: Includes Freedom Future Fit Program.

**FD2 Rigid** (E1235) FD2-30270R **\$2085**  
**Transport Optional** or (K0005)  
 Note: Transportation Option selected separately below.  
 Note: Includes Freedom Future Fit Program.

### Transportation Option

Note: Used when Transport Optional frame is selected above.  
**Transportation Package** FD2-99998 **\$ 250**  
 Note: Must be installed by Freedom Designs manufacturing.  
 Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.  
 Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

### Frame Colors

Sugar Cookie	n/c	Pacific Blue	n/c
Neon Yellow	n/c	Blue Blast	n/c
Yellow	n/c	Electric Blue	n/c
Bright Pink	n/c	Galaxy Blue	n/c
Pink Parfait	n/c	Passion Purple	n/c
Sugar Plum	n/c	Octane Orange	n/c
Candy Violet	n/c	Red Sunset	n/c
Candy Raspberry	n/c	Really Red	n/c
Extreme Green	n/c	Fiery Red	n/c
Sparkle Green	n/c	Twilight	n/c
Mystic Teal	n/c	Stormy Gray	n/c
		Black	n/c

### Frame Dimensions

#### Specify Frame Width Setting

11" Wide	FD2-30210	n/c
12" Wide	FD2-30220	n/c
13" Wide	FD2-30230	n/c
14" Wide	FD2-30240	n/c
15" Wide	FD2-30290	n/c
16" Wide	FD2-30260	n/c
17" Wide	FD2-30330	n/c
18" Wide	FD2-30280	n/c

Note: Frame widths measured across and from outside of seat frame tubes.

### Back Post Locations

#### Specify Frame Type and Frame Length

<b>Regular Frame Length - Adjusts 15"-19"</b>	FD2-30558	n/c
15"	FD2-15	n/c
16"	FD2-16	n/c
17"	FD2-17	n/c
18"	FD2-18	n/c
19"	FD2-19	n/c

**Specify Back Post Location** (if other than specified above) \_\_\_\_\_ "

**Long Frame Length - Adjusts 17"-21"** FD2-30510 n/c

17"	FD2-17	n/c
18"	FD2-18	n/c
19"	FD2-19	n/c
20"	(E2203) FD2-20	\$ 500
21"	(E2203) FD2-21	\$ 500

**Specify Back Post Location** (if other than specified above) \_\_\_\_\_ "

Note: FD2 has 5" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating/hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back/hardware thickness of 2" = back post location set at 14".

### Sling Seat Upholstery

<b>Omit Seat Sling</b>	FD2-39004	
<b>Seat Sling</b>	FD2-39011	n/c

Note: Not Available with ANY Single Post armrests.

Note: See Page 4 for Seat Pan options.

#### Specify Seat Sling Length

15"	FD2-39015	\$ 56
16"	FD2-39016	n/c
17"	FD2-39017	\$ 56
18"	FD2-39018	n/c
19"	FD2-39021	\$ 77
20"	FD2-39022	\$ 77
21"	FD2-39023	\$ 77

#### Seat Cushions

<b>1" Seat Cushion</b>	FD2-39003	\$ 87
<b>2" Seat Cushion</b>	FD2-39007	\$ 97

Note: Includes 4# foam & black knit cover with zipper closure.

### Back Post Tubes

#### Specify Push Handle, Back Post and Extensions

<b>Straight Back Post Tubes - No Push Handles</b>	FD2-30600	n/c
---	-----------	-----

#### Specify Back Post Height

10" High	FD2-30610	n/c
11" High	FD2-30611	n/c
12" High	FD2-30612	n/c
13" High	FD2-30613	n/c
14" High	FD2-30614	n/c
15" High	FD2-30615	n/c

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail
<b>Back Post Tubes (continued)</b>		
<b>Standard 90° Bend Push Handles and Back Post Tubes</b>	FD2-30518	n/c
<b>Ergo Bend Push Handles and Back Post Tubes</b>	FD2-30519	n/c
<b>Note:</b> All push/stroller handle extensions not recommended with Ergo Bend.		
<b>8° Bend Back Canes</b>	FD2-30511	\$ 26
<b>Seating Dynamics Rocker Back</b>	(K0108) FD2-30532	\$1734
<b>Note:</b> Angle adjustable from 72° to 108° in 6° increments.		
<b>Note:</b> Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.		
<b>Note:</b> Locks into a non-dynamic position when necessary.		
<b>Note:</b> Does not allow the back post tubes to fold forward.		
<b>Back Post Height</b>		
<b>Note:</b> Back Posts are measured from top of seat rail to top of handle bend.		
<b>Short Push Handle - 20" high</b>	FD2-30515	n/c
<b>Medium Push Handle - 23" high</b>	FD2-30514	n/c
<b>Tall Push Handle - 26" high</b>	FD2-30516	n/c
<b>Stroller Handle Extensions</b>		
<b>Note:</b> All stroller handle ext. add approximately 6.5" to push handle height.		
<b>Note:</b> All stroller/push handle extensions except Individual include a black neoprene wrap around pad.		
	<b>-30560</b>	
	<b>-30570/5</b>	
	<b>-30580/5</b>	
	<b>-30517</b>	
<b>Push Handle Extension - Width Adjustable</b>	FD2-30560	n/c
<b>Note:</b> Not removable, but allows frame to fold.		
<b>Note:</b> Mounts direct to back post push handle.		
<b>Push Button, Angle-Adjustable Folding Stroller Handle Extension</b>	FD2-30570	\$ 112
<b>Note:</b> Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.		
<b>Note:</b> Not removable, but allows frame to fold.		
<b>Note:</b> Mounts direct to back post push handle.		
<b>Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty</b>	FD2-30575	\$ 179
<b>Note:</b> Available 10W-22W wheelchair frames.		
<b>Note:</b> Not removable, but allows frame to fold.		
<b>Note:</b> Same as above but with metal hinge mechanism.		
<b>Note:</b> Mounts direct to back post push handle.		
<b>Push Button, Angle-Adjustable, Folding Push Handle Extension</b>	FD2-30580	\$ 112
<b>Note:</b> Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.		
<b>Note:</b> Not removable, but allows frame to fold.		
<b>Note:</b> Mounts direct to straight back post tubes.		
<b>Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty</b>	FD2-30585	\$ 179
<b>Note:</b> Available 10W-22W wheelchair frames.		
<b>Note:</b> Not removable, but allows frame to fold.		
<b>Note:</b> Same as above but with metal hinge mechanism.		
<b>Note:</b> Mounts direct to straight back post tubes.		
<b>Individual Stroller Handle Extensions Removable - Pair</b>	FD2-30517	\$ 102
<b>Note:</b> Not available with Straight Back Posts.		

Description	Part No.	Retail
<b>Back Post Tubes (continued)</b>		
<b>Sling Back Upholstery</b>		
<b>Omit Sling Back Upholstery</b>	FD2-39005	
<b>Standard Sling Back Upholstery</b>	FD2-39006	n/c
<b>Note:</b> Not available with Flip-Down Rigidizing Bar.		
<b>Note:</b> Sling back height must be selected below.		
<b>Tension Adjustable Sling Back Uph.</b>	FD2-39008	\$ 236
<b>Note:</b> Not available with Flip-Down Rigidizing Bar.		
<b>Note:</b> Sling back height must be selected below.		
10" High	FD2-30910	n/c
11" High	FD2-30911	n/c
12" High	FD2-30912	n/c
13" High	FD2-30913	n/c
14" High	FD2-30914	n/c
15" High	FD2-30915	n/c
16" High	FD2-30916	n/c
17" High	FD2-30917	n/c
18" High	FD2-30918	n/c
<b>Note:</b> Sling back height must be selected when ordering Standard or Adjustable Tension Sling Back Upholstery.		
<b>Legrests, Footrests and Extensions</b>		
<b>Legrests</b>		
<b>Omit Legrest Hangers &amp; Footrests</b>	FD2-33000	
<b>75° Pop-Up Swing Away Hangers</b>	FD2-30284	n/c
<b>Note:</b> Not available with high mount bracket.		
<b>90° Pop-Up Swing Away Hangers</b>	FD2-30281	n/c
<b>Note:</b> Not available with 8" adjustable extensions.		
<b>90° Pop-Up Swing Away Elevating Legrests (ELRs)</b>	FD2-30279	\$ 255
<b>Note:</b> Not available with all length adjustable extensions.		
<b>Note:</b> Not available with 6.5" heavy duty extensions.		
<b>75° Latch Style Swing Away Hangers</b>	FD2-30282	n/c
<b>Note:</b> Not available with high mount bracket.		
<b>90° Latch Style Swing Away Hangers</b>	FD2-30283	n/c
<b>Note:</b> Not available with 8" adjustable extensions.		
<b>Footplates &amp; Footboards</b>		
<b>Omit Footrest Plates</b>	FD2-36404	
<b>Composite Flip-Up</b>	FD2-17070	n/c
<b>Kids Aluminum Angle Adjustable, Flip-up</b>	(K0040) FD2-36020	\$ 160
<b>Note:</b> For frames 14" wide and narrower.		
<b>Heel Loops</b>	(E0951) FD2-36022	\$ 41
<b>Calf Strap - Webbing</b>	(K0038) FD2-36024	\$ 41
<b>Adult Aluminum Angle Adj. Flip-up w/ Heel Loops</b>	(K0040/ E0951) FD2-36021	\$ 180
<b>Note:</b> For frames 15" wide and greater.		
<b>Omit Heel Loops</b>	FD2-36023	n/c
<b>Aluminum Angle Adjustable Footboard Removable</b>	(K0040) FD2-36029	\$ 160

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Legrests, Footrests and Extensions(continued)</b> <b>MAPS Footplate Mounting Hardware</b> <b>MAPS</b> (Multi Axis Positioning System) - Pair FD2-36040 \$ 306 <b>MAPS</b> (Multi Axis Positioning System) - Left FD2-36040-L \$ 153 <b>MAPS</b> (Multi Axis Positioning System) - Right FD2-36040-R \$ 153 Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning. Note: Must choose aluminum footplates. Note: Not available on 10" & 11" wide frames Note: MAPS not available with composite footplates or aluminum footboard.			<b>Armrests and Armrest Pads (continued)</b> <b>Frame Mounted Armrests</b> <b>Single Post Fixed Height, Removable - Short</b> (K0015) FD2-38222 \$ 225 <b>Single Post Fixed Height, Removable - Tall</b> (K0015) FD2-38224 \$ 225 Note: Single Post Fixed Height arms include a removable side guard. Note: Adjustable in height with tools only. See matrix for height ranges. Note: Not available with ANY clothing guard accessory. <b>Single Post Height Adjustable, Removable - Short</b> (E0973) FD2-38217 \$ 265 <b>Single Post Height Adjustable, Removable - Tall</b> (E0973) FD2-38203 \$ 265 Note: Not available with ANY clothing guard accessory. <b>Desk Length Pads</b> FD2-38213 n/c <b>Full Length Pads</b> FD2-38214 n/c Specify Seat Rail to Top of Armrest Measurement _____ "		
<b>Legrest Extensions</b> <b>4" Short Length Adjustable Extensions</b> FD2-30285 n/c <b>6" Medium Length Adjustable Extensions</b> FD2-30286 n/c <b>8" Long Length Adjustable Extensions</b> FD2-30287 n/c Note: Not available with all 90° hangers. <b>4.5" Heavy Duty ELR Extensions</b> FD2-36401 \$ 56 <b>4.5" Heavy Duty Extensions</b> FD2-36402 \$ 56 Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes. <b>6.5" Heavy Duty Extensions</b> FD2-36403 \$ 66 Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes. <b>High Mount Brackets</b> FD2-36407 \$ 102 Note: Includes flip-up mounting hardware for future growth adjustment. Note: Not available with all 75° hangers. Specify Knee to Heel Measurement _____ " (see matrix page)			<b>Front Wheels - Casters</b> <b>Omit Caster Wheels</b> FD2-34333 <b>3" Polyurethane</b> FD2-34888 \$ 56 <b>3" Polyurethane Light-Up</b> FD2-34111 \$ 56 Note: 3" casters recommended for low seat to floor heights. <b>4" Polyurethane</b> FD2-40000 \$ 66 <b>4" Polyurethane Light-Up</b> FD2-41000 \$ 101 <b>4" x 1.4" Poly, Aluminum Mag</b> FD2-43000 \$ 92 <b>5" x 1" Poly, Plastic Mag</b> FD2-34555 n/c <b>5" Polyurethane Light-Up</b> FD2-35050 \$ 101 <b>5" x 1.5" Poly, Plastic Mag</b> FD2-34222 \$ 87 <b>5" x 1.4" Poly, Aluminum Mag</b> FD2-42000 \$ 92 <b>6" x 1" Poly, Plastic Mag</b> FD2-34666 n/c <b>6" x 1.5" Poly, Plastic Mag</b> FD2-34999 \$ 87 <b>6" x 1.4" Poly, Aluminum Mag</b> FD2-44000 \$ 92 <b>6" x 2" Poly with Tread</b> FD2-34444 \$ 92 <b>6" x 1" Pneumatic Plastic Mag</b> (E2214/E2215) FD2-34667 \$ 77 <b>7" x 1" Poly, Plastic Mag</b> FD2-34777 \$ 77 <b>7" x 2" Polyurethane</b> FD2-34200 \$ 118 Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock. <b>Caster Stem Spacer - 0.75"</b> FD2-34778 n/c		
<b>Armrests and Armrest Pads</b> <b>Omit Armrests</b> FD2-38209 <b>Back Post Mounted Armrests</b> Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed. <b>10.5" Long Flip-Back</b> FD2-38204 n/c Note: Full length pad on 13.5" or longer armrests only. <b>13.5" Long Flip-Back</b> FD2-38202 n/c Note: Full length pad on 13.5" or longer armrests only. <b>10.5" Long Flip-Back Locking</b> (E0973) FD2-38201 \$ 245 Note: Full length pad on 13.5" or longer armrests only. <b>13.5" Long Flip-Back Locking</b> (E0973) FD2-38200 \$ 245 Note: Full length pad on 13.5" or longer armrests only. <b>10.5" to 13.5" Long Flip-Back, Telescoping</b> (E0973) FD2-38208 \$ 240 Note: Full length pad on 13.5" or longer armrests only. <b>10.5" to 13.5" Long Flip-Back Locking, Telescoping</b> (E0973) FD2-38210 \$ 265 Note: Full length pad on 13.5" or longer armrests only. <b>17" to 20" Long Flip-Back Locking, Telescoping</b> (E0973) FD2-38215 \$ 265 <b>Foam Wraps - Standard</b> FD2-30216 n/c <b>Desk Length Pads</b> FD2-38213 n/c <b>Full Length Pads</b> FD2-38214 n/c Note: Full length pad on 13.5" or longer armrests only.			<b>Rear Wheels - Tires - Accessories</b> Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles. <b>Omit Rear Wheels and Axles</b> FD2-31524 <b>Mag Wheels</b> <b>12" Low Profile Poly Mag</b> (E2218) FD2-31526 \$ 36 Note: Not available with Hand Activated Drum Brake/Wheel Lock. Note: Not available with airless inserts. Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock. <b>12" Semi Pneumatic Plastic Mag</b> (K0069) FD2-31523 \$ 168 Note: Not available with airless inserts. Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock. <b>12" Pneumatic Plastic Mag</b> FD2-31525 n/c <b>12" Airless Inserts</b> (E2213) FD2-31019 \$ 77 Note: All 12" wheels not available with all handrims. Note: All 7" casters not available with all 12" rear wheels and Hand Activated Drum Brake / Wheel Lock.		

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Rear Wheels - Tires - Accessories (continued)</b>			<b>Rear Wheels - Tires - Accessories (continued)</b>		
16" Pneumatic Plastic Mag	FD2-31530	\$ 41	20" Pneumatic Spoke Wheel with Aluminum Handrim	FD2-31540	\$ 204
16" Airless Inserts (E2213)	FD2-31020	\$ 77	20" Airless Inserts (E2213)	FD2-31022	\$ 77
16" Aluminum Handrims	FD2-31014	\$ 169	Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
Note: 16" mag wheels not available with plastic coated handrims.			22" Pneumatic Spoke Wheel with Aluminum Handrim	FD2-31541	\$ 204
18" Pneumatic Plastic Mag	FD2-31533	\$ 41	22" Airless Inserts (E2213)	FD2-31023	\$ 77
18" Airless Inserts (E2213)	FD2-31021	\$ 77	Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
18" Aluminum Handrims	FD2-31015	\$ 169	24" Pneumatic Spoke Wheel with Aluminum Handrim	FD2-31542	\$ 204
Note: 18" mag wheels not available with plastic coated handrims.			24" Airless Inserts (E2213)	FD2-31024	\$ 77
20" Pneumatic Plastic Mag with Aluminum Handrim	FD2-31527	n/c	Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
20" Airless Inserts (E2213)	FD2-31022	\$ 77	Omit Handrims from 20", 22" or 24" Wheels	FD2-31006	
20" Plastic Coated Handrims	FD2-31016	\$ 119	Natural Fit Handrims, Super Grip - 22" Spoke wheel	FD2-31009	\$ 347
22" Pneumatic Plastic Mag with Aluminum Handrim	FD2-31528	n/c	Natural Fit Handrims, Super Grip - 24" Spoke wheel	FD2-31010	\$ 347
22" Airless Inserts (E2213)	FD2-31023	\$ 77	Camber Spacer Kit	FD2-31003	\$ 56
22" Plastic Coated Handrims	FD2-31017	\$ 119	3 Degrees	FD2-31011	n/c
24" Pneumatic Plastic Mag with Aluminum Handrim	FD2-31529	n/c	6 Degrees	FD2-31012	n/c
24" Airless Inserts (E2213)	FD2-31024	\$ 77	<b>Wheel Locks</b>		
24" Plastic Coated Handrims	FD2-31018	\$ 119	Push to Lock - Standard	FD2-31501	n/c
Omit Handrims from 20", 22" or 24" Wheels	FD2-31006		Pull to Lock	FD2-31506	n/c
Natural Fit Handrims, Super Grip - 24" Mag wheel	FD2-31008	\$ 347	Extension Handles (E0961)	FD2-31801	\$ 61
<b>Spoke Wheels</b>			Hand Activated Brake and Wheel Lock System (E2228)	FD2-31521	\$ 700
16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	FD2-31534	n/c	Note: Includes manual wheel locks to be used with public transportation.		
Note: Rim, hub & spokes are anodized silver.			Note: Not available with 12" mag wheels with low pro poly tires.		
16" Airless Inserts (E2213)	FD2-31020	\$ 77	Note: Used to control direction, speed, brake and lock wheels.		
16" Silver Aluminum Handrim	FD2-31030	\$ 77	Note: Activation handles are mounted to the back post push handles.		
16" Black Plastic Coated Handrim	FD2-31032	\$ 77	<b>Rear Anti-Tippers</b>		
Note: 16" & 18" spoke wheels do not include handrims.			Anti-Tip Assembly - Fixed (E0971)	FD2-31600	\$ 97
16" Black Pneumatic Tire & Black Spoke Aluminum Wheel	FD2-31535	n/c	Anti-Tip Assembly - Swing-Away (E0971)	FD2-31602	\$ 97
Note: Rim, hub & spokes are anodized black.			<b>Seat Pan Options</b>		
16" Airless Inserts (E2213)	FD2-31020	\$ 77	Fixed Solid Seat Pan - Bolteron (E2231)	FD2-39010	\$ 150
Note: 16" & 18" spoke wheels do not include handrims.			Note: Not available with folding frame.		
16" Black Aluminum Handrim	FD2-31031	\$ 77	Fixed Solid Seat Pan - Aluminum (E2231)	FD2-39030	\$ 170
16" Black Plastic Coated Handrim	FD2-31032	\$ 77	Note: Not available with folding frame.		
Note: 16" & 18" spoke wheels do not include handrims.			Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".		
18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	FD2-31537	n/c	Removable Solid Seat Pan - Bolteron (E2231)	FD2-39012	\$ 170
Note: Rim, hub & spokes are anodized silver.			Removable Solid Seat Pan - Aluminum (E2231)	FD2-39032	\$ 255
18" Airless Inserts (E2213)	FD2-31021	\$ 77	Note: Removable Seat Pan raises the effective seat to floor height by 5/8".		
18" Silver Aluminum Handrim	FD2-31034	\$ 77	Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.		
18" Black Plastic Coated Handrim	FD2-31036	\$ 77	Note: All seat pans include 2" x 2" notches for growth adjustment.		
Note: 16" & 18" spoke wheels do not include handrims.			Note: All seat pans made to fit frame dimensions unless specified below.		
18" Black Pneumatic Tire & Black Spoke Aluminum Wheel	FD2-31538	n/c	Seat Width _____" X Seat Depth _____"		
Note: Rim, hub & spokes are anodized black.					
18" Airless Inserts (E2213)	FD2-31021	\$ 77			
18" Black Aluminum Handrim	FD2-31035	\$ 77			
18" Black Plastic Coated Handrim	FD2-31036	\$ 77			
Note: 16" & 18" spoke wheels do not include handrims.					

FD2 Rider Weight Capacity is 250 lbs.

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

**Accessories**

<b>Clothing Guards – Fabric</b> <b>Note:</b> Not available with ANY single post armrests.	FD2-CGF	\$ 71
<b>Clothing Guards – Quick-Release, Rigid</b> <b>Note:</b> Not available with ANY single post armrests.	FD2-CGR	\$ 122
<b>Medical Necessities Pouch</b> <b>Note:</b> Black nylon.	FD2-23207	\$ 68
<b>Footrest Hanger Pads - Pair</b> <b>Note:</b> Wrap-around pads made from black neoprene.	SFDI-9HP	\$ 92
<b>Frame Pads - Pair</b> <b>Note:</b> Wrap-around pads made from black neoprene. <b>Note:</b> Typically used to pad exposed back post tubes.	10911-2	\$ 82
<b>Flip-Down Rigidizing Bar</b> <b>Note:</b> Mounts to back post tubes to add rigidity to frame. <b>Note:</b> Not available with ANY Sling Back Upholstery. <b>Note:</b> Does NOT adjust in width to grow with frame.	FD2-23210	\$ 129
<b>Freedom Backpack</b> <b>Note:</b> Black Nylon only with Freedom Designs logo. <b>Note:</b> Monogramming is not available on Freedom Backpack.	FD2-10900	\$ 56
<b>Freedom Canopy - Flip-Back</b> <b>Note:</b> Canopy assembly includes removable mounting hardware. <b>Note:</b> Canopy is available in black only with Freedom Designs logo. <b>Note:</b> Canopy includes rear flap with zippered storage pocket. <b>Note:</b> Canopy includes clear-view window.	FD2-10905	\$ 128

**Anterior Support Surfaces – (Trays)**

<b>3/8" Lexan Tray - Channel Bracket Hardware</b>	(E0950) 08-2000	\$ 271
<b>3/8" Lexan Tray - Adjustable Knob Hardware &amp; Ribs</b>	(E0950) 08-6000	\$ 219
<b>1/4" Lexan Tray - EZ Lite Channel Hardware &amp; Ribs</b>	(E0950) 08-13000	\$ 183
<b>1/4" Lexan Tray - Hook &amp; Loop Straps</b>	(E0950) 08-16000	\$ 106

**Footrest / Extension Lengths**

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

**MAPS Footrest / Extension Lengths**

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

**Note:** Length from top of footrest plate to top of seat rail.  
**Note:** Without extensions bottom hanger hole not accessible for hardware.  
**Note:** Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

**Freedom Designs Armrest Heights**

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

**Note:** All heights from top of seat frame to top of armrest pad.

**HCPCS CODES** are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Choose Seat to Floor Height**

**Note:** Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, and front caster size.  
**Note:** Caster sizes vary slightly in height by type.  
**Note:** Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".  
**Note:** Seat pan choice increases seat rail to floor height 3/8" to 3/4".

**FD2 Wheel / Caster Seat to Floor Heights**

Wheel	Caster	Height	Arc Over Rail
12"	3"	14"-14.5"	0"
12"	4"	15"	0"
12"	5"	17"	0"
12"	6"	17.5"	0"
16"	3"	13"-14.5"	3"-1.5"
16"	4"	13.5"-14"	2.5"-2"
16"	5"	14"-14.5"	2"-1.5"
16"	6"	14.5"-15"	1.5"-1"
16"	7"	15"-15.5"	1"
18"	3"	14"-15.5"	4"-2.5"
18"	4"	15.5"	3"
18"	5"	16"	2"
18"	6"	16.5"	1.5"
18"	7"	17"	1"
20"	3"	16"-17.5"	4"-3.5"
20"	4"	15"	3.5"
20"	5"	15.5"-17.5"	4.5"-2.5"
20"	6"	17"-17.5"	3"-2.5"
20"	7"	18"	2"
22"	3"	15.5"	6.5"
22"	4"	16"-17"	6"-5"
22"	5"	16.5"-17.5"	5.5"-4.5"
22"	6"	17"-18.5"	5"-3.5"
22"	7"	18"-18.5"	4"-3.5"
24"	3"	16"	8"
24"	4"	16.5"-17.5"	7.5"-6.5"
24"	5"	17"-18"	7"-6"
24"	6"	17.5"-18.5"	6.5"-5.5"
24"	7"	18"-19"	6"-5.5"

**FD2 Rider Weight Capacity is 250 lbs.**