



Freedom SP3

2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

Bill to:
 Date: _____ PO#: _____
 Company: _____
 Buyer: _____
 Marked For: _____


Ship to:
 Company: _____
 Address: _____
 City: _____ State: _____ Zip: _____
 Phone: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------


Frame Type

- SP3—Folding** SP3-30270 **\$2158**
 (E1238) or (K0005)
 HCPCS codes for frames 10"W—18"W
- SP3—Non-Folding (Rigid)** SP3-30270R **\$2158**
 (E1237) or (K0005)
 HCPCS codes for frames 10"W—18"W

Folding Frame



Rigid Frame




Transport Option

- Transport Package** SP3-99998 **\$150**
 Available with FDI WC19 approved seating or on frame only.
- No Transport Package** SP3-99999-1 **n/c**


Configuration Set-Up

- Standard—Rear Wheel Drive** SP3-STD **n/c**
- Reverse Configuration—Front Wheel Drive** SP3-REV **n/c**

Standard Configuration




Reverse Configuration



Frame Width

Specify Frame Width:

- 10" Wide SP3-30250 **n/c**
- 11" Wide SP3-30210 **n/c**
- 12" Wide SP3-30220 **n/c**
- 13" Wide SP3-30230 **n/c**
- 14" Wide SP3-30240 **n/c**
- 15" Wide SP3-30290 **n/c**
- 16" Wide SP3-30260 **n/c**
- 17" Wide SP3-30330 **\$ 75**
- 18" Wide SP3-30280 **\$ 75**



← Width →

Frame Length

Specify Frame Length:

- X-Small (11" - 14")** SP3-30507 **n/c**
Set-up Depth:
 - 11" Deep SP3-11 **n/c**
 - 12" Deep SP3-12 **n/c**
 - 13" Deep SP3-13 **n/c**
 - 14" Deep SP3-14 **n/c**
- Small (14" - 17")** SP3-30509 **n/c**
Set-up Depth:
 - 14" Deep SP3-14 **n/c**
 - 15" Deep SP3-15 **n/c**
 - 16" Deep SP3-16 **n/c**
 - 17" Deep SP3-17 **n/c**
- Medium (16" - 19")** SP3-30508 **n/c**
Set-up Depth:
 - 16" Deep SP3-16 **n/c**
 - 17" Deep SP3-17 **n/c**
 - 18" Deep SP3-18 **n/c**
 - 19" Deep SP3-19 **n/c**
- Long (18" - 21")** SP3-30510 **\$ 75**
Set-up Depth:
 - 18" Deep SP3-18 **n/c**
 - 19" Deep SP3-19 **n/c**
 - 20" Deep SP3-20 **n/c**
 - 21" Deep SP3-21 **n/c**

Frame Colors

<input type="radio"/> Sugar Cookie	<input type="radio"/> Pacific Blue
<input type="radio"/> Neon Yellow	<input type="radio"/> Blue Blast
<input type="radio"/> Yellow	<input type="radio"/> Electric Blue
<input type="radio"/> Bright Pink	<input type="radio"/> Galaxy Blue
<input type="radio"/> Pink Parfait	<input type="radio"/> Passion Purple
<input type="radio"/> Sugar Plum	<input type="radio"/> Octane Orange
<input type="radio"/> Vigorous Violet	<input type="radio"/> Red Sunset
<input type="radio"/> Candy Violet	<input type="radio"/> Really Red
<input type="radio"/> Candy Raspberry	<input type="radio"/> Fiery Red
<input type="radio"/> Extreme Green	<input type="radio"/> Twilight
<input type="radio"/> Sparkle Green	<input type="radio"/> Stormy Gray
<input type="radio"/> Mystic Teal	<input type="radio"/> Black
<input type="radio"/> Blizzard Blue	

Date: _____ PO#: _____

Marked For _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles
Hand Grip Bend Options**

- Traditional Bend Push Handles** SP3-30518 n/c
- Ergo Bend Push Handles** SP3-30519 n/c



Specify height from rear seat rail to top of bend on cane:

- Short (19" High)** SP3-30515 n/c
- Medium (23" High)** SP3-30514 n/c
- Tall (26" High)** SP3-30516 n/c
- Push Handle Extension-Quick Release** SP3-30517 **\$100**
- One Piece Push Handle Extension - Quick Release** SP3-30526 **\$125**
- One Piece Push Button Angle Adjustable Stroller Handle** SP3-30527 **\$175**

Available on chairs 14"W and up.
Note: Chair will not fold with this stroller handle.

Back canes can achieve an additional height of 6.5" when quick release push handle extensions are ordered.

Swing Away Hangers

- Omit Hangers & Footplates** SP3-33000 **(\$30)**
- Spring Loaded Pop Up** SP3-30272 n/c
- Latch Style Swing Away** SP3-30278 n/c
Not available in reverse configuration.

Specify:

- 90 Degree Hangers*** SP3-33334 n/c
- 75 Degree Hangers*** SP3-33335 n/c

*See Configuration Chart



- Elevating Legrest (E0990)** SP3-30279 **\$235**
Not available with latch style.

Footrest Options



- Omit Footplates Only** SP3-36404 **(\$15)**
- Composite Footplate - Flip Up** SP3-17070 n/c
- Kids Footplate (K0040) Flip Up Footplate—Adjustable Angle** SP3-36020 **\$150**
For chairs 14" wide and under
- Adult Footplate 15"W & Over (K0040) Flip Up Footplate—Adjustable Angle w/ Heel Loop** SP3-36021 **\$150**
For chairs 15" wide and over.

Footrest Options (Continued)

- Omit Heel Loop** SP3-36023 n/c
- Heel Loop (E0951)** SP3-36022 **\$ 50**
- Footboard - Removable, Angle Adj.** SP3-36029 **\$150**
- Hanger Extensions - .5—4.5" Long** SP3-36402 **\$ 55**
- Hanger Extensions - 1.5—6.5" Long** SP3-36403 **\$ 65**
1.5" - 6.5" long only available with 75 degree hangers.

Configuration	Knee to Heel Range (from seat rail)
Hanger with High Mount Clamps	4" - 9"
Hanger Only	6" - 9"
Hanger with .5" - 4.5" Extension	to 13.5"
Hanger with 1.5" - 6.5" Extension	to 15.5"

- High Mount Clamp (K0037)** SP3-36407 **\$100**
Not available on 75 degree hangers.
- Knee to Heel Measurement** _____ **Inches**

Casters

- Omit Caster Wheels** SP3-34333 n/c
- 3" Poly Casters** SP3-34888 **\$ 55**
- 3" Light Up Casters** SP3-34111 **\$ 55**
3" casters in combination with 20" wheels are not available.
3" casters recommended for low frame configurations.
- 4" Poly Casters** SP3-40000 **\$ 65**
- 4" Light Up Casters** SP3-41000 **\$ 99**
- 4" x 1.4" Poly Aluminum Mag** SP3-43000 **\$ 90**
- 5" Poly Casters** SP3-34555 n/c
- 5" x 1.5" Poly Casters** SP3-34222 **\$ 85**
- 5" x 1.4" Poly Aluminum Mag** SP3-42000 **\$ 90**
- 6" Poly Casters** SP3-34666 n/c
- 6" x 1.5" Poly Casters** SP3-34999 **\$ 85**
- 6" x 2" Poly Casters with Tread** SP3-34444 **\$ 90**
- 6" Pneumatic Casters** SP3-34667 **\$ 75**
- 7" Poly Casters** SP3-34777 **\$ 75**
- 7" x 1.5" Poly Casters** SP3-34000 **\$116**
- 7" x 2" Poly Casters** SP3-34200 **\$116**
- Frog Legs (E1015)** SP3-35000 **\$325**
Not available with 6" x 2", 7", 7" x 1.5", & 7" x 2" casters.
- Caster Stem Bolt—.75" Longer** SP3-34778 n/c
- Caster Stem Bolt—1.5" Longer** SP3-34779 n/c

Date: _____ PO#: _____

Marked For _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Anti-Tip Tubes

- Anti-Tip Tubes—Fixed** (E0971) SP3-31600 **\$ 95**
- Anti-Tip Tubes—Swing Away** (E0971) SP3-31602 **\$ 95**
Can not be used with reverse configuration.

Wheels & Wheel Accessories

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** SP3-31524 **(\$55)**

Mag Wheel



Spoke Wheel



Mag Wheels

- 12" Pneumatic Mag** SP3-31525 **n/c**
- 12" Polyurethane Mag** SP3-31526 **\$ 35**
- 12" Semi Pneumatic Mag** SP3-31523 **\$165**
Handrims are not available with 12" wheels.

- 16" Pneumatic Mag** SP3-31530 **\$ 40**
- 18" Pneumatic Mag** SP3-31533 **\$ 40**

- 20" Pneumatic Mag** SP3-31527 **n/c**
- 22" Pneumatic Mag** SP3-31528 **n/c**
- 24" Pneumatic Mag** SP3-31529 **n/c**

Aluminum handrims are included on 20", 22" and 24" wheels.

To omit handrims, please check below.

24" wheels are not available with small frame lengths.

Spoke Wheels

- 20" Pneumatic Spoke** SP3-31540 **\$200**
- 22" Pneumatic Spoke** SP3-31541 **\$200**
- 24" Pneumatic Spoke** SP3-31542 **\$200**

Aluminum handrims are included on 20", 22" and 24" wheels.

To omit handrims, please check below.

24" wheels are not available with small frame lengths.

Plastic coated handrims are not available with spoke wheels.

Wheel Accessories

- Omit Handrims** SP3-31006 **(\$20)**
On 20", 22" and 24" mag or spoke wheels.
- Aluminum Handrims 16" & 18" (pr.)** SP3-31007 **\$166**
- Plastic Coated Handrims (pr.)** SP3-31005 **\$117**
Plastic coated handrims are not available with spoke wheels.

Natural Fit Handrims

Please contact your customer service representative for details.

- Airless Inserts (pr.)** (E2213) SP3-31000 **\$ 75**
For pneumatic tires only.

Wheel Locks



Specify lever type:

- Push to Lock** SP3-31501 **n/c**
- Pull to Lock** SP3-31506 **n/c**
- Rear Mount Push to Lock** SP3-31507 **n/c**
- Rear Mount Pull to Lock** SP3-31508 **n/c**
- Wheel Lock Handle Ext.** (E0961) SP3-31801 **\$60**

Armrests

- Omit Flip Back Armrests Back Cane Mount** SP3-38209 **(\$20)**
- Flip Back Armrests—13" Long** SP3-38204 **n/c**
- Flip Back Armrests—16" Long** SP3-38202 **n/c**
- Locking Flip Back Armrests—13" Long** (E0973) SP3-38201 **\$235**

- Locking Flip Back Armrests—16" Long** (E0973) SP3-38200 **\$235**

- Flip Back Armrests—Telescoping** SP3-38208 **\$100**
- Locking Flip Back Armrests—Telescoping** SP3-38207 **\$260**

Frame Mount

- Height Adjustable Padded Armrests** (E0973) SP3-38203 **\$235**



Flip Back Armrest



Height Adjustable Padded Armrest



Lock-


- Desk Length Flat Armrest Pad** SP3-38213 **n/c**
- Full Length Flat Armrest Pad** SP3-38214 **n/c**
Only available on 16" armrests.
- Seat Rail to Top of Armrest** _____ **Inches**

Date: _____ PO#: _____

Marked For _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Seating Options



- Rigid Seat Base** SP3-39006 **\$100**
- Removable Seat Base** SP3-39012 **\$150**
- Seat Rail to Floor Height** _____ **Inches**

Accessories

- Clothing Guards (Fabric)** SP3-CGF **\$ 70**
- Clothing Guards (Quick Release - Rigid)** SP3-CGR **\$120**
- Medical Necessities Pouch** SP3-23207 **\$ 67**
Note: Black nylon.
- Footrest Hanger Pads - Pair** SFDI-9HP **\$ 90**
Note: Wrap around pads made from black neoprene.
- Frame Pads - Pair** 10911-2 **\$ 80**
Note: Wrap around pads made from black neoprene.
Note: Typically used to pad exposed back post tubes.
- Freedom Backpack** SP3-10900 **\$ 55**
Note: Back nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom Backpack.
- Freedom Canopy - Flip Back** SP3-10905 **\$125**
Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.
Note: See offer for FREE Backpack and Canopy on last page of order form.

Notes

Wheel & Caster Matrix

Floor to Front of Seat Rail	Caster Size	Wheel Size	Arc Over Rail
11.5"	3"	16"	4.5"
12"	3"	16"	4"
12.5"	3"	16"	3.5"
12.5"	3"	18"	5.5"
13"	3"	16"	3"
13"	3"	18"	5"
13.5"	3"	18"	4.5"
14"	3"	18"	4"
14.5"	3"	20"	5.5"
15"	3"	20"	5"
14"	5"	20"	6 "
14.5"	5"	18"	3.5"
14"	5"	20"	6"
14.5"	5"	18"	3.5"
14.5"	5"	20"	5.5"
15"	5"	20"	5"
15"	5"	22"	7"
15.5"	5"	20"	4.5"
15.5"	5"	22"	4.5"
16"	5"	20"	4"
16"	5"	22"	6"
15.5"	6"	20"	4.5"
15.5"	6"	22"	6.5"
16"	6"	20"	4"
16"	6"	22"	6"
16"	6"	24"	8"
16.5"	6"	22"	5.5"
16.5"	6"	24"	7.5"
17	6	22	5
17"	6"	24"	7"
16.5"	7"	22"	5.5"
16.5"	7"	24"	7.5"
17"	7"	22"	5"
17"	7"	24"	7"

Maximum Rider Weight—250 lbs.
Frame weight is approximately 23 lbs.

Date: _____ PO#: _____

Marked For: _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***

Freedom Designs & Pindot

Are Pleased To Announce
"A Sizzling Summer Sensation"
And For a Limited Time It's...

FREE * FREE * FREE

When ordering a Freedom
Wheelchair and a
complete Freedom or Pindot
seating system including seat,
back, and attachment hardware,
simply provide the seating Quote
or seating PO number
below and receive a
Freedom Canopy and
Freedom Backpack

*** FREE ***



Freedom Backpack with Complete Freedom Seating System SP3-SYS-10900 n/c

Note: Black Nylon only with Freedom Designs logo.

Note: Monogramming is not available on Freedom backpack.

Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

Provide Freedom/PinDot Seating Quote/PO # _____



Freedom Canopy - Flip-Back with Complete Freedom Seating System SP3-SYS-10905 n/c

Note: Canopy assembly includes removable mounting hardware.

Note: Canopy is available in black only with Freedom Designs logo.

Note: Canopy includes rear flap with zippered storage pocket.

Note: Canopy includes clear-view window.

Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

Provide Freedom/PinDot Seating Quote/PO # _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***