



Our kids are our passion.™

# Freedom SP3—Mini Series

2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

**Bill to:**  
 Date: \_\_\_\_\_ PO#: \_\_\_\_\_  
 Company: \_\_\_\_\_  
 Buyer: \_\_\_\_\_  
 Marked For: \_\_\_\_\_

**Ship to:**  
 Company: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type		
<input type="radio"/> <b>Freedom SP3 Mini Folding</b>	SPM-30270	<b>\$2315</b>
(E1238) or (K0005) HCPCS codes for frames 10"W—18"W		
<input type="radio"/> With double seat rail for upholstery & seating. <input type="radio"/> With folding single seat rail for seating only. *		
Folding Frame With Seat Sling 	Folding Frame For Seating Only 	
<input type="radio"/> <b>Freedom SP3 Mini Rigid</b>	SPM-30270R	<b>\$2315</b>
(E1237) or (K0005) HCPCS codes for frames 10"W—18"W		
<input type="radio"/> With double seat rail for upholstery & seating. <input type="radio"/> With folding single seat rail for seating only. *		
Rigid Frame With Seat Sling 	Rigid Frame For Seating Only 	
*Seat to floor height - 1 5/8" less *Lighter weight		
Transport Option		
<input type="radio"/> <b>Transport Package</b>	SPM-99998	<b>\$150</b>
Available with FDI WC19 approved seating or on frame only.		
<input type="radio"/> <b>No Transport Package</b>	SPM-99999-1	n/c
Configuration Set-up		
<input type="radio"/> <b>Standard—Rear Wheel Drive</b>	SPM-STD	n/c
<input type="radio"/> <b>Reverse Configuration—Front Wheel Drive</b>	SPM-REV	n/c
Standard Configuration 	Reverse Configuration 	

Frame Colors	
<input type="radio"/> Sugar Cookie	<input type="radio"/> Pacific Blue
<input type="radio"/> Neon Yellow	<input type="radio"/> Blue Blast
<input type="radio"/> Yellow	<input type="radio"/> Electric Blue
<input type="radio"/> Bright Pink	<input type="radio"/> Galaxy Blue
<input type="radio"/> Pink Parfait	<input type="radio"/> Passion Purple
<input type="radio"/> Sugar Plum	<input type="radio"/> Octane Orange
<input type="radio"/> Vigorous Violet	<input type="radio"/> Red Sunset
<input type="radio"/> Candy Violet	<input type="radio"/> Really Red
<input type="radio"/> Candy Raspberry	<input type="radio"/> Fiery Red
<input type="radio"/> Extreme Green	<input type="radio"/> Twilight
<input type="radio"/> Sparkle Green	<input type="radio"/> Stormy Gray
<input type="radio"/> Mystic Teal	<input type="radio"/> Black
<input type="radio"/> Blizzard Blue	

Frame Width		
Specify Frame Width:		
<input type="radio"/> <b>10" Wide</b> —Available without sling only.	SPM-30250	n/c
<input type="radio"/> <b>11" Wide</b>	SPM-30210	n/c
<input type="radio"/> <b>12" Wide</b>	SPM-30220	n/c
<input type="radio"/> <b>13" Wide</b>	SPM-30230	n/c
<input type="radio"/> <b>14" Wide</b>	SPM-30240	n/c
<input type="radio"/> <b>15" Wide</b>	SPM-30290	n/c
<input type="radio"/> <b>16" Wide</b>	SPM-30260	n/c
<input type="radio"/> <b>17" Wide</b>	SPM-30330	<b>\$ 75</b>
<input type="radio"/> <b>18" Wide</b>	SPM-30280	<b>\$ 75</b>

Frame Length		
Specify Frame Length:		
<input type="radio"/> <b>Small (9" - 12")</b>	SPM-30509	n/c
<input type="radio"/> <b>9" Deep</b>	SPM- 9	n/c
<input type="radio"/> <b>10" Deep</b>	SPM-10	n/c
<input type="radio"/> <b>11" Deep</b>	SPM-11	n/c
<input type="radio"/> <b>12" Deep</b>	SPM-12	n/c
<input type="radio"/> <b>Long (11" - 14")</b>	SPM-30510	<b>\$ 75</b>
<input type="radio"/> <b>11" Deep</b>	SPM-11	n/c
<input type="radio"/> <b>12" Deep</b>	SPM-12	n/c
<input type="radio"/> <b>13" Deep</b>	SPM-13	n/c
<input type="radio"/> <b>14" Deep</b>	SPM-14	n/c

Sling Seat		
<input type="radio"/> <b>Omit Seat Sling</b>	SPM-39004	<b>(\$15)</b>
<input type="radio"/> <b>Seat Sling</b>	SPM-39011	n/c
Specify Sling Length: Front of seat rail to back upholstery.		

<input type="radio"/> <b>Seat Sling Length:</b>		
<input type="radio"/> <b>9" Deep</b>	SPM-30409	n/c
<input type="radio"/> <b>10" Deep</b>	SPM-30500	n/c
<input type="radio"/> <b>11" Deep</b>	SPM-30501	n/c
<input type="radio"/> <b>12" Deep</b>	SPM-30502	n/c
<input type="radio"/> <b>13" Deep</b>	SPM-30503	n/c
<input type="radio"/> <b>14" Deep</b>	SPM-30504	n/c

Cushions		
<input type="radio"/> <b>1" Cushion</b>	SPM-39003	<b>\$ 85</b>
<input type="radio"/> <b>2" Cushion</b>	SPM-39007	<b>\$ 95</b>
Includes 4# foam & black knit cover.		

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles  
Hand Grip Bend Options**

- Traditional Bend Push Handles** SPM-30518 n/c
- Ergo Bend Push Handles** SPM-30519 n/c



Traditional Push Handle

Ergo Push Handle

Specify height from rear seat rail to top of bend on cane:

- Short (19" High)** SPM-30515 n/c
- Medium (23" High)** SPM-30514 n/c
- Tall (26" High)** SPM-30516 n/c

- Push Handle Extension - Quick Release** SPM-30517 \$100
- One Piece Push Handle Extension - Quick Release** SPM-30526 \$125
- One Piece Push Button Angle Adjustable Stroller Handle** SPM-30527 \$175

Available on chairs 14"W and up.  
Note: Chair will not fold with this stroller handle.

Back canes can achieve an additional height of 6.5" when quick release push handle extensions are ordered.

**Sling Backs**

- Omit Sling Back** SPM-39005 (\$15)
- Sling Back\*** SPM-39015 n/c
- Tension Adjustable Back\*** SPM-39008 \$231

Specify Sling Back Height

(Top of seat to top of back upholstery)

- 10" High** SPM-30910 n/c
- 11" High** SPM-30911 n/c
- 12" High** SPM-30912 n/c
- 13" High** SPM-30913 n/c
- 14" High** SPM-30914 n/c
- 15" High** SPM-30915 n/c
- 16" High** SPM-30916 n/c
- 17" High** SPM-30917 n/c
- 18" High** SPM-30918 n/c

\*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery

**Swing Away Hangers**

- Omit Hangers & Footplates** SPM-33000 (\$30)
- Spring Loaded Pop-Up** (pictured below) SPM-30272 n/c
- Latch Style Swing Away** SPM-30278 n/c  
*Not available in reverse configuration.*

Specify:

- 90 Degree Hangers\*** SPM-33334 n/c
  - 75 Degree Hangers\*** SPM-33335 n/c
- \*See Configuration Chart



90° Latch Style Hanger

75° Pop-Up Style Hanger

**Footrest Options**



Angle Adjustable Flip Up Footplate

Removable Footboard

- Omit Footplates Only** SPM-36404 (\$15)
- Composite Footplate - Flip Up** SPM-17070 n/c
- Kids Footplate (K0040) Flip Up Footplate—Adjustable Angle** SPM-36020 \$150  
*For chairs 14" wide and under.*
- Adult Footplate 15"W & Over (K0040) Flip Up Footplate—Adjustable Angle w/ Heel Loop** SPM-36021 \$150  
*For chairs 15" wide and over.*
- Omit Heel Loop** SPM-36023 n/c
- Heel Loop (E0951)** SPM-36022 \$ 50
- Footboard—Removable, Angle Adj.** SPM-36029 \$150
- Hanger Extensions - .5—4.5" Long** SPM-36402 \$ 55
- Hanger Extensions - 1.5—6.5" Long** SPM-36403 \$ 65  
*1.5" - 6.5" long only available with 75 degree hangers.*

Configuration	Knee to Heel Range (from seat rail)
Hanger with High Mount Clamps	4" - 9"
Hanger Only	6" - 9"
Hanger with .5" - 4.5" Extension	to 13.5"
Hanger with 1.5" - 6.5" Extension	to 15.5"

- High Mount Clamp (K0037)** SPM-36407 \$100  
*Not available on 75 degree hangers.*
- Knee to Heel Measurement** \_\_\_\_\_ Inches

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Casters**

- Omit Caster Wheels** SPM-34333 n/c
- 3" Poly Casters** SPM-34888 \$ 55
- 3" Light Up Casters** SPM-34111 \$ 55  
3" casters in combination with 20" wheels are not available.  
3" casters recommended for low frame configurations.
- 4" Poly Casters** SPM-40000 \$ 65
- 4" Light Up Casters** SPM-41000 \$ 99
- 4" x 1.4" Poly Aluminum Mag** SPM-43000 \$ 90
- 5" Poly Casters** SPM-34555 n/c
- 5" x 1.5" Poly Casters** SPM-34222 \$ 85
- 5" x 1.4" Poly Aluminum Mag** SPM-42000 \$ 90
- 6" Poly Casters** SPM-34666 n/c
- 6" x 1.5" Poly Casters** SPM-34999 \$ 85
- 6" x 2" Poly Casters with Tread** SPM-34444 \$ 90
- 6" Pneumatic Casters** SPM-34667 n/c
- Frog Legs (E1015)** SPM-35000 \$325  
Not available with 6" x 2" casters.
- Caster Stem Bolt—.75" Longer** SPM-34778 n/c
- Caster Stem Bolt—1.5" Longer** SPM-34779 n/c



**Anti-Tip Tubes**

- Anti-Tip Tubes—Fixed (E0971)** SPM-31600 \$ 95
- Anti-Tip Tubes—Swing Away (E0971)** SPM-31602 \$ 95  
Can not be used with reverse configuration.

**Wheels & Wheel Accessories**

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** SPM-31524 (\$55)  
Mag Wheel  Spoke Wheel 

**Mag Wheels**

- 12" Pneumatic Mag** SPM-31525 n/c
- 12" Polyurethane Mag** SPM-31526 \$ 35
- 12" Semi Pneumatic Mag** SPM-31523 \$165  
Handrims not available with 12" wheels.
- 16" Pneumatic Mag** SPM-31530 \$ 40
- 18" Pneumatic Mag** SPM-31533 \$ 40
- 20" Pneumatic Mag** SPM-31527 n/c  
Available on long frame only.  
Aluminum handrims are included on 20" mag wheels.  
To omit handrims, please check below

**Spoke Wheels**

- 20" Pneumatic Spoke** SPM-31540 \$200  
Available on long frame only.  
Aluminum handrims are included on 20" spoke wheels.  
Plastic coated handrims are not available with spoke wheels.

**Wheel Accessories**

- Omit Handrims** SPM-31006 (\$20)  
On 20" mag or spoke wheels.
- Aluminum Handrims 16" & 18" (pr.)** SPM-31007 \$166
- Plastic Coated Handrims (pr.)** SPM-31005 \$117  
Plastic coated handrims are not available with spoke wheels.
- Natural Fit Handrims**  
Please contact your customer service representative for details.
- Airless Inserts (pr.) (E2213)** SPM-31000 \$ 75  
For pneumatic tires only.

**Wheel Locks**



Specify lever type:

- Push to Lock** SPM-31501 n/c
- Pull to Lock** SPM-31506 n/c
- Rear Mount Push to Lock** SPM-31507 n/c
- Rear Mount Pull to Lock** SPM-31508 n/c
- Wheel Lock Handle Ext. (E0961)** SPM-31801 \$ 60

**Armrests**

- Omit Flip Back Armrest** SPM-38209 (\$20)
- Back Cane Mount**
- Flip Back Armrests—13" Long** SPM-38204 n/c
- Flip Back Armrests—16" Long** SPM-38202 n/c
- Locking Flip Back Armrests—13" Long (E0973)** SPM-38201 \$235
- Locking Flip Back Armrests—16" Long (E0973)** SPM-38200 \$235
- Flip Back Armrests—Telescoping** SPM-38208 \$100
- Locking Flip Back Armrests—Telescoping** SPM-38207 \$260

**Frame Mount**

- Height Adjustable Padded Armrest (E0973)** SPM-38203 \$235



Flip Back Armrest



Height Adjustable Padded Armrest




Locking Flip Back Armrest

- Desk Length Flat Armrest Pad** SPM-38213 n/c
- Full Length Flat Armrest Pad** SPM-38214 n/c  
Only available on 16" armrests.
- Seat Rail to Top of Armrest** \_\_\_\_\_ Inches

Date: \_\_\_\_\_ PO#: \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Seating Options**



- Rigid Seat Base** SPM-39006 **\$100**
- Removable Seat Base** SPM-39012 **\$150**
- Seat Rail to Floor Height** \_\_\_\_\_ **Inches**

**Accessories**

- Clothing Guards (Fabric)** SPM-CGF **\$ 70**
- Clothing Guards (Quick Release - Rigid)** SPM-CGR **\$120**
- Medical Necessities Pouch** SPM-23207 **\$ 67**  
*Note: Black nylon.*
- Footrest Hanger Pads - Pair** SFDI-9HP **\$ 90**  
*Note: Wrap around pads made from black neoprene.*
- Frame Pads - Pair** 10911-2 **\$ 80**  
*Note: Wrap around pads made from black neoprene.*  
*Note: Typically used to pad exposed back post tubes.*
- Freedom Backpack** SPM-10900 **\$ 55**  
*Note: Back nylon only with Freedom Designs logo.*  
*Note: Monogramming is not available on Freedom Backpack.*
- Freedom Canopy - Flip Back** SPM-10905 **\$125**  
*Note: Canopy assembly includes removable mounting hardware.*  
*Note: Canopy is available in black only with Freedom Designs logo.*  
*Note: Canopy includes rear flap with zippered storage pocket.*  
*Note: Canopy includes clear-view window.*  
  
*Note: See offer for ~~FREE~~ Backpack and Canopy on last page of order form.*

**Wheel & Caster Matrix**

Floor to Front of Seat Rail	Caster Size	Wheel Size	Arc Over Rail	Fork Size
10.5"	3"	16"	4.5"	Short
11.5"	3"	18"	6.5"	Short
12.5"	5"	18"	5.5"	Long
12.5"	5"	20"	7.5"	Long
13"	5"	18"	5"	Long
13"	5"	20"	7"	Long
13.5"	5"	18"	4.5"	Long
13.5"	5"	20"	6.5"	Long

**Maximum Rider Weight—200 lbs.**  
**Frame weight is approximately 23 lbs.**

**Notes:**

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

**\* FREE \* FREE \* FREE \* FREE \* FREE \* FREE \***

# Freedom Designs & Pindot

Are Pleased To Announce  
*"A Sizzling Summer Sensation"*  
And For a Limited Time It's...

**FREE \* FREE \* FREE**

When ordering a Freedom  
Wheelchair and a  
complete Freedom or Pindot  
seating system including seat,  
back, and attachment hardware,  
simply provide the seating Quote  
or seating PO number  
below and receive a  
Freedom Canopy and  
Freedom Backpack

**\* FREE \***



**Freedom Backpack with Complete Freedom Seating System** SPM-SYS-10900 n/c

**Note:** Black Nylon only with Freedom Designs logo.  
**Note:** Monogramming is not available on Freedom backpack.  
**Note:** Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

Provide Freedom/PinDot Seating Quote/PO # \_\_\_\_\_



**Freedom Canopy - Flip-Back with Complete Freedom Seating System** SPM-SYS-10905 n/c

**Note:** Canopy assembly includes removable mounting hardware.  
**Note:** Canopy is available in black only with Freedom Designs logo.  
**Note:** Canopy includes rear flap with zippered storage pocket.  
**Note:** Canopy includes clear-view window.  
**Note:** Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

Provide Freedom/PinDot Seating Quote/PO # \_\_\_\_\_

**\* FREE \* FREE \* FREE \* FREE \* FREE \* FREE \***