



**Like no other
tilt-in-space
wheelchairs
in the world...**

Freedom[®] NXT[™]

Generation NEXT

Tilt-in-space Mobility:

- **Rigid**
- **Folding**
- **Specialized**



Why?

Since its debut in March 2011, the NXT has proven that superior strength, durability and reliability don't have to come at the cost of weight, adjustability and portability.

The NXT is easy to fold, lift and stow—and in its lightest configuration,

weighs in at only 22 lbs. The robust design has a rigid-like feel, a maximum weight capacity of 225 lbs., adjusts 2" in width and 4" in depth for growth and still fits into the trunk of a small car when folded.

The NXT is unique in that it is the first product ever to have a true "specialized" positioning seat frame. With the LEAF & MAPS options, it can fit a seated posture at the hips and at the knees and address foot and ankle positioning challenges.

Even with all of these special features, there is no impact on the overall size of the NXT.

The NXT is truly the first fully functional positioning frame that remains lightweight, compact and portable and wherever you take it, you can rest assured it won't let you down!



NXT "Generation NEXT" Folding



NXT "Generation NEXT" Rigid

- Streamlined and aesthetically pleasing
- Very light 22 lb. transfer weight
- Compact to fit into small spaces
- Folding & Rigid frame choices
- Folding frame has "rigid-like" feel
- Short, medium & long frame choices
- Rider weight capacity of 225 lbs.
- 5° anterior to 45° posterior tilt
- Easy to adjust for growing kids
- Pediatric and adult sizes
- 21 colors to choose from
- Eleven sizes in 1" increments
- Most sizes grow 2" in frame width
- All sizes grow 4" in seat depth
- Back angle adjusts from 80° to 120°
- 2", 3" and 4" seat frame extensions
- Anti-rotation seat frame
- Anti-rotation cross brace
- Welded caster journals
- Flared front base frame
- Seat tube leveling guides
- Lightweight front fork
- Aluminum front casters
- Universal axle plate
- Patented dual locking safety axle
- NEW** patented **Lower Extremity Accommodation Frame LEAF** accommodates abduction, adduction or windswept legs
- NEW** patented **Multi Axis Positioning System footrest MAPS** rotates and accommodates inversion, eversion, plantar and dorsiflexion
- NEW** **Seat Rail Only seat frame option SRO** used when footrest hanger attachments are not necessary
- NEW** asymmetric seat frame adjustment
- NEW** 5° anterior tilt frame setting
- NEW** Seating Dynamics Rocker Back option
- Single or dual hand lever activated tilt
- Foot activated tilt
- Foot activated wheel lock
- NEW** attendant control hand lever activated tilt, brake and wheel lock
- Pop Up or Swing Away footrest hangers
- NEW** length adjustable footrest extension tubes
- Combination anti-tipper & curb assist leverage bar
- Crash tested with Freedom Seating
- Transit tie downs mounted to base frame
- Vent and battery tray capability
- Oxygen tank holder
- Safety light reflecting coating
- Light up front casters
- Freedom's "Future Fit" program



*Like no other
tilt-in-space
wheelchairs in
the world...*

Freedom: for Kids (and other people too)



NXT "Generation NEXT" Specialized with LEAF & MAPS



Our kids are our passion.™

Freedom Designs Incorporated

2241 N. Madera Road, Simi Valley, CA 93065

phone: 800-331-8551

fax: 888-582-1509

e: CustomerService@freedomdesigns.com

www.freedomdesigns.com

Freedom Designs Incorporated (FDI) maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without notification. Check with FDI for the latest information.



Copyright © 2012
Freedom Designs, Inc.
All rights reserved.
Printed in the USA.
WBNT45R 92012



Manufactured in the U.S.A.
Patent pending.