



**NXT<sup>TM</sup>Mini**  
Generation NEXT

**Pediatric**  
**Manual Wheelchair**



2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: \_\_\_\_\_ Buyer: \_\_\_\_\_  
Company: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Email: \_\_\_\_\_

Quote Order PO# \_\_\_\_\_  
Location: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail																																																																																					
<p align="center"><b>Transport Optional Base</b></p> <p><b>NXT-Mini Pediatric Folding</b> (E1234) <b>NXT-30330</b> <b>\$2800</b> <i>Transport Optional</i> <b>Note:</b> Transportation Option selected separately below.</p> <p><b>NXT-Mini Pediatric Rigid</b> (E1233) <b>NXT-30340</b> <b>\$2700</b> <i>Transport Optional</i> <b>Note:</b> Transportation Option selected separately below. <b>Note:</b> Seat frame is offered as a Seat Rail Only for use with a seating system mounted footrest, or with a fixed 70°, 80°, or 90° bend legrest insert that is manufactured with holes located in 1/2" increments for footrest hardware attachment and footrest plate adjustment. <b>Note:</b> Not available with any Pop-Up or Latch swing away or elevating footrest hangers. <b>Note:</b> Freedom Future Fit Program includes folding cross brace &amp; footplates for growing frame WIDTH only. <b>Note:</b> ANY DEPTH growth of NXT-MINI requires conversion to a SMALL frame NXT. Contact Customer Service for pricing and availability.</p> <p align="center"><b>Transportation Option</b></p> <p><b>Note:</b> Used when Transport Optional frame is selected above. <b>Transportation Package</b> <b>NXT-99998</b> <b>\$ 250</b> <b>Note:</b> Must be installed by Freedom Designs manufacturing. <b>Note:</b> Not available with all vent and battery trays. <b>Note:</b> Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved. <b>Note:</b> Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.</p>			<p align="center"><b>Frame Colors</b></p> <table border="0"> <tr><td>Sugar Cookie</td><td>n/c</td><td>Pacific Blue</td><td>n/c</td></tr> <tr><td>Neon Yellow</td><td>n/c</td><td>Blue Blast</td><td>n/c</td></tr> <tr><td>Yellow</td><td>n/c</td><td>Electric Blue</td><td>n/c</td></tr> <tr><td>Bright Pink</td><td>n/c</td><td>Galaxy Blue</td><td>n/c</td></tr> <tr><td>Pink Parfait</td><td>n/c</td><td>Passion Purple</td><td>n/c</td></tr> <tr><td>Sugar Plum</td><td>n/c</td><td>Octane Orange</td><td>n/c</td></tr> <tr><td>Candy Violet</td><td>n/c</td><td>Red Sunset</td><td>n/c</td></tr> <tr><td>Candy Raspberry</td><td>n/c</td><td>Really Red</td><td>n/c</td></tr> <tr><td>Extreme Green</td><td>n/c</td><td>Fiery Red</td><td>n/c</td></tr> <tr><td>Sparkle Green</td><td>n/c</td><td>Twilight</td><td>n/c</td></tr> <tr><td>Mystic Teal</td><td>n/c</td><td>Stormy Gray</td><td>n/c</td></tr> <tr><td></td><td></td><td>Black</td><td>n/c</td></tr> </table> <p align="center"><b>Frame Dimensions</b></p> <p align="center"><b>Specify Frame Width Setting</b></p> <table border="0"> <tr><td>10"</td><td>11"</td><td>12"</td><td>NXT-30250</td><td>n/c</td></tr> <tr><td>12"</td><td>13"</td><td>14"</td><td>NXT-30220</td><td>n/c</td></tr> </table> <p><b>Note:</b> Frame widths measured across and from outside of seat frame tubes. <b>Note:</b> NXT-MINI is offered as a Seat Rail Only for use with a seating system mounted footrest, or with a fixed 75° or 85° bend legrest insert that is manufactured with holes located in 1/2" increments for footrest hardware attachment and footrest plate adjustment. <b>Note:</b> Freedom Future Fit Program applies to NXT-Mini size NXT-30250 ONLY. No Freedom Future Fit Program is available for NXT-Mini size NXT-30220. ANY growth greater than 14" in seat width and 14" in seat depth requires a conversion to a small frame NXT. Contact Freedom Customer Service for assistance in pricing and ordering conversion components.</p> <p align="center"><b>Back Post Locations</b></p> <p align="center"><b>Specify Frame Type and Frame Length</b></p> <p><b>Note:</b> Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.</p> <table border="0"> <tr><td><b>Mini Frame</b></td><td><b>NXT-30510</b></td><td><b>n/c</b></td></tr> <tr><td>8"</td><td>NXT-30500XTP</td><td>n/c</td></tr> <tr><td>9"</td><td>NXT-30501XTP</td><td>n/c</td></tr> <tr><td>10"</td><td>NXT-30502XTP</td><td>n/c</td></tr> <tr><td>11"</td><td>NXT-30503XTP</td><td>n/c</td></tr> <tr><td>12"</td><td>NXT-30504XTP</td><td>n/c</td></tr> <tr><td>* 13" Must select 2" extension</td><td>NXT-30505XTP</td><td>n/c</td></tr> <tr><td>* 14" Must select 2" extension</td><td>NXT-30506XTP</td><td>n/c</td></tr> <tr><td>2" Extension - Pair (K0108)</td><td>NXT-30525</td><td>\$ 26</td></tr> </table> <p><b>Specify Back Post Location</b> (if other than specified above) _____" <b>Note:</b> Floor to seat frame height range is 15" to 19". <b>Note:</b> Must choose an extension for 13" and 14" frame lengths.</p> <p><b>Note:</b> NXT Mini has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 6" + back / hardware thickness of 2" = back post location set at 8".</p>			Sugar Cookie	n/c	Pacific Blue	n/c	Neon Yellow	n/c	Blue Blast	n/c	Yellow	n/c	Electric Blue	n/c	Bright Pink	n/c	Galaxy Blue	n/c	Pink Parfait	n/c	Passion Purple	n/c	Sugar Plum	n/c	Octane Orange	n/c	Candy Violet	n/c	Red Sunset	n/c	Candy Raspberry	n/c	Really Red	n/c	Extreme Green	n/c	Fiery Red	n/c	Sparkle Green	n/c	Twilight	n/c	Mystic Teal	n/c	Stormy Gray	n/c			Black	n/c	10"	11"	12"	NXT-30250	n/c	12"	13"	14"	NXT-30220	n/c	<b>Mini Frame</b>	<b>NXT-30510</b>	<b>n/c</b>	8"	NXT-30500XTP	n/c	9"	NXT-30501XTP	n/c	10"	NXT-30502XTP	n/c	11"	NXT-30503XTP	n/c	12"	NXT-30504XTP	n/c	* 13" Must select 2" extension	NXT-30505XTP	n/c	* 14" Must select 2" extension	NXT-30506XTP	n/c	2" Extension - Pair (K0108)	NXT-30525	\$ 26
Sugar Cookie	n/c	Pacific Blue	n/c																																																																																							
Neon Yellow	n/c	Blue Blast	n/c																																																																																							
Yellow	n/c	Electric Blue	n/c																																																																																							
Bright Pink	n/c	Galaxy Blue	n/c																																																																																							
Pink Parfait	n/c	Passion Purple	n/c																																																																																							
Sugar Plum	n/c	Octane Orange	n/c																																																																																							
Candy Violet	n/c	Red Sunset	n/c																																																																																							
Candy Raspberry	n/c	Really Red	n/c																																																																																							
Extreme Green	n/c	Fiery Red	n/c																																																																																							
Sparkle Green	n/c	Twilight	n/c																																																																																							
Mystic Teal	n/c	Stormy Gray	n/c																																																																																							
		Black	n/c																																																																																							
10"	11"	12"	NXT-30250	n/c																																																																																						
12"	13"	14"	NXT-30220	n/c																																																																																						
<b>Mini Frame</b>	<b>NXT-30510</b>	<b>n/c</b>																																																																																								
8"	NXT-30500XTP	n/c																																																																																								
9"	NXT-30501XTP	n/c																																																																																								
10"	NXT-30502XTP	n/c																																																																																								
11"	NXT-30503XTP	n/c																																																																																								
12"	NXT-30504XTP	n/c																																																																																								
* 13" Must select 2" extension	NXT-30505XTP	n/c																																																																																								
* 14" Must select 2" extension	NXT-30506XTP	n/c																																																																																								
2" Extension - Pair (K0108)	NXT-30525	\$ 26																																																																																								
<p align="center"><b>Transport Included Base</b></p> <p><b>NXT-Mini Pediatric Folding</b> (E1234) <b>NXT-30331</b> <b>\$2950</b> <i>Transport Included</i> <b>NXT-Mini Pediatric Rigid</b> (E1233) <b>NXT-30341</b> <b>\$2850</b> <i>Transport Included</i> <b>Note:</b> Seat frame is offered as a Seat Rail Only for use with a seating system mounted footrest, or with a fixed 70°, 80°, or 90° bend legrest insert that is manufactured with holes located in 1/2" increments for footrest hardware attachment and footrest plate adjustment. <b>Note:</b> Not available with any Pop-Up or Latch swing away or elevating footrest hangers. <b>Note:</b> Freedom Future Fit Program includes folding cross brace &amp; footplates for growing frame WIDTH only. <b>Note:</b> ANY DEPTH growth of NXT-MINI requires conversion to a SMALL frame NXT. Contact Customer Service for pricing and availability.</p>			<p align="center"><b>Tilt Activations</b></p> <table border="0"> <tr><td><b>Dual Lever Hand Activated Tilt</b></td><td><b>NXT-30272</b></td><td><b>n/c</b></td></tr> <tr><td><b>Foot Activated Tilt</b></td><td><b>NXT-30273</b></td><td><b>n/c</b></td></tr> <tr><td><b>Note:</b> Foot pedal installed on right rear side frame.</td><td></td><td></td></tr> <tr><td><b>Combination Hand Activated Tilt, Brake and Wheel Lock System</b></td><td><b>NXT-30274</b></td><td><b>n/c</b></td></tr> <tr><td><b>Note:</b> Includes manual wheel locks to be used with public transportation.</td><td></td><td></td></tr> </table> <p align="center"><b>Vent &amp; Battery Trays</b></p> <p><b>Note:</b> All vent &amp; battery trays not available with NXT-MINI. <b>Note:</b> All vent and battery trays not with transport option.</p>			<b>Dual Lever Hand Activated Tilt</b>	<b>NXT-30272</b>	<b>n/c</b>	<b>Foot Activated Tilt</b>	<b>NXT-30273</b>	<b>n/c</b>	<b>Note:</b> Foot pedal installed on right rear side frame.			<b>Combination Hand Activated Tilt, Brake and Wheel Lock System</b>	<b>NXT-30274</b>	<b>n/c</b>	<b>Note:</b> Includes manual wheel locks to be used with public transportation.																																																																								
<b>Dual Lever Hand Activated Tilt</b>	<b>NXT-30272</b>	<b>n/c</b>																																																																																								
<b>Foot Activated Tilt</b>	<b>NXT-30273</b>	<b>n/c</b>																																																																																								
<b>Note:</b> Foot pedal installed on right rear side frame.																																																																																										
<b>Combination Hand Activated Tilt, Brake and Wheel Lock System</b>	<b>NXT-30274</b>	<b>n/c</b>																																																																																								
<b>Note:</b> Includes manual wheel locks to be used with public transportation.																																																																																										
<p align="center"><b>NXT-MINI Rider Weight Capacity is 75 lbs.</b></p>			<p align="center">Page 1 of 4 <b>NXT-MINI Pediatric</b> €€</p>																																																																																							

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

**Description** **Part No.** **Retail**

**Description** **Part No.** **Retail**

**Back Post Tubes**

**Specify Push Handle, Back Post and Extensions**

**Standard 90° Bend Push Handles and Back Post Tubes** NXT-30518 n/c

**Ergo Bend Push Handles and Back Post Tubes** NXT-30519 n/c

**Note:** All push/stroller handle extensions not recommended with Ergo Bend.

**Back Post Height**

**Note:** Back Posts are measured from top of seat rail to top of handle bend.

**Short Push Handle - 20" high** NXT-30515 n/c

**Medium Push Handle - 23" high** NXT-30514 n/c

**Tall Push Handle - 26" high** NXT-30516 n/c

**Stroller Handle Extensions**

**Note:** All stroller handle ext. add approximately 6.5" to push handle height.

**Note:** All stroller/push handle extensions except Individual include a black neoprene wrap around pad.



-30560



-30570



-30580



-30517

**Push Handle Extension - Width Adjustable** NXT-30560 n/c

**Push Button, Angle-Adjustable Folding Stroller Handle Extension** NXT-30570 \$ 112

**Note:** Available 10W-15W wheelchair frames only.

**Note:** Not removable, but allows frame to fold.

**Note:** Mounts direct to back post push handle

**Push Button, Angle-Adjustable, Folding Push Handle Extension** NXT-30580 \$ 112

**Note:** Available 10W-15W wheelchair frames only.

**Note:** Not removable, but allows frame to fold.

**Note:** Mounts direct to straight back post tubes.

**Individual Stroller Handle Extensions Removable - Pair** NXT-30517 \$ 102

**Specify Knee to Heel**

3"	3.5"	4"	4.5"	5"	5.5"	6"
6.5"	7"	7.5"	Requires 4" Legrest Extension			
8"	8.5"	9"	Requires 4" Legrest Extension			
9.5"	10"	Requires 4" Legrest Extension				

**Legrests, Footrests and Extensions**

**Legrests**

NXT-Mini offers a series of fixed legrest tube options which insert into the seat rail. Footplate mounting holes are placed at 1/2" increments for height adjustment. Footplate range is measured from the top of the seat rail to the top of the footplate. Range of adjustment is 3" to 10" with attachment hardware and optional 4" legrest extension.

**Seat Rail Only - Mini** NXT-30326 n/c

**Note:** For use with seat mounted footrest.

**75° Fixed Legrest Insert - Pair** NXT-30370 n/c

**Note:** Bottom mounting hole usable after installation of 4" extension tube only.

**Note:** 75° & 85° fixed legrest use an interchangeable insert tube.

**85° Fixed Legrest Insert - Pair** NXT-30371 n/c

**Note:** Bottom mounting hole usable after installation of 4" extension tube only.

**Note:** 75° & 85° fixed legrest use an interchangeable insert tube.

**Note:** Seat rail must be moved 1.5" forward for base frame clearance when installing 4" extensions.

**15° Abduction Seat Rail Insert with 75° Fixed Legrest Extensions** NXT-30327 n/c

**Note:** Features a fixed 15 degree lateral bend in the seat rail insert to allow for abducted lower extremity positioning.

**Note:** When used with 15° Fixed Abduction Seat Rail and the 75° Fixed Legrest Extension, the footplate is mounted in a 15° abducted position.

**Note:** Includes a 75° fixed legrest extensions only.

**Note:** Not available with Bolteron Flip-Up Footboard.

**Note:** Not available with 4" length adjustable extensions.

**Footplates & Footboards**

**Omit Footrest Plates** NXT-36404

**Composite Flip-Up** NXT-17070 n/c

**MINI Aluminum Angle Adjustable, Flip-up** (K0040) NXT-36019 \$ 160

**Mid-line Mounted Footplates** NXT-36025 n/c

**Note:** When used with 15° Fixed Abduction Seat Rail and the 75° Fixed Legrest Extension, the footplate is mounted in a mid-line position.

**Abduction Mounted Footplates** NXT-36026 n/c

**Note:** When used with 15° Fixed Abduction Seat Rail and the 75° Fixed Legrest Extension, the footplate is mounted in a 15° abducted position.

**Note:** Only available with the 15° Fixed Abduction Seat Rail and the 75° Fixed Legrest Extension.

**Heel Loops** (E0951) NXT-36022 \$ 41

**Calf Strap - Webbing** (K0038) NXT-36024 \$ 41

**Bolteron Angle Adjustable Flip-Up Footboard** (K0040) NXT-36030 \$ 160

**Note:** Not available with 15° Fixed Abduction Seat Rail.

**Legrest Extensions**

**4" Mini Legrest Extensions** NXT-30289 n/c

**Note:** Seat rail must be moved 1.5" forward for base frame clearance when installing 4" extensions.

**Note:** Not available with 15° Fixed Abduction Seat Rail.

**Specify Knee to Heel Measurement \_\_\_\_\_"**

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Armrests and Armrest Pads</b>			<b>Rear Wheels - Tires - Accessories (continued)</b>		
<b>8.5" Long Flip-Back</b>	NXT-38218	n/c	<b>Spoke Wheels</b>		
<b>10.5" Long Flip-Back</b>	NXT-38204	n/c	<b>16" Gray Pneumatic Tire &amp; Silver Spoke Aluminum Wheel</b>	NXT-31534	n/c
<b>10.5" to 13.5" Long Flip-Back, Telescoping</b> (E0973)	NXT-38208	\$ 240	Note: Rim, hub & spokes are anodized silver.		
Note: All Flip-Back armrests are mounted to backpost tubes.			<b>16" Airless Inserts</b> (E2213)	NXT-31020	\$ 77
Note: 8.5" armrests not available with desk and full length armrest pads.			<b>16" Black Pneumatic Tire &amp; Black Spoke Aluminum Wheel</b>	NXT-31535	n/c
Note: Flip-Back armrests length are measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.			Note: Rim, hub & spokes are anodized black.		
Note: Armrests measure 6" to 11" from seat rail with foam wrap arm pads installed.			<b>16" Airless Inserts</b> (E2213)	NXT-31020	\$ 77
Note: Armrests measure 7" to 12" from seat rail with desk length arm pads installed.			<b>18" Gray Pneumatic Tire &amp; Silver Spoke Aluminum Wheel</b>	NXT-31537	n/c
			Note: Rim, hub & spokes are anodized silver.		
<b>Foam Wraps - Standard</b>	NXT-38216	n/c	<b>18" Airless Inserts</b> (E2213)	NXT-31021	\$ 77
<b>Desk Length Pads</b>	NXT-38213	n/c	<b>18" Black Pneumatic Tire &amp; Black Spoke Aluminum Wheel</b>	NXT-31538	n/c
Note: 8.5" armrests not available with desk and full length armrest pads.			Note: Rim, hub & spokes are anodized black.		
<b>Specify Seat Rail to Top of Armrest Measurement</b> _____ "			<b>18" Airless Inserts</b> (E2213)	NXT-31021	\$ 77
<b>Front Wheels - Casters</b>			<b>Wheel Locks</b>		
<b>3" Polyurethane</b>	NXT-34888	\$ 56	<b>Push to Lock - Standard Extension Handles</b>	NXT-31501 (E0961) NXT-31801	n/c \$ 61
<b>3" Polyurethane Light-Up</b>	NXT-34111	\$ 56	<b>Foot Activated Wheel Lock</b> (E2228)	NXT-31508	n/c
Note: 3" casters recommended for low seat to floor heights.			Note: Includes manual wheel locks to be used with public transportation.		
<b>4" Polyurethane</b>	NXT-40000	\$ 66	Note: Not available with 12" poly plastic mag wheels.		
<b>4" Polyurethane Light-Up</b>	NXT-41000	\$ 101	<b>Hand Activated Brake and Wheel Lock System</b> (E2228)	NXT-30275	n/c
<b>4" x 1.4" Poly, Aluminum Mag</b>	NXT-43000	\$ 92	Note: Includes manual wheel locks to be used with public transportation.		
<b>5" x 1" Poly, Plastic Mag</b>	NXT-34555	n/c	Note: Not available with 12" mag wheels with low pro poly tires.		
<b>5" Polyurethane Light-Up</b>	NXT-35050	\$ 101	Note: Used to control direction, speed, brake and lock wheels from strolling position.		
<b>5" x 1.5" Poly, Plastic Mag</b>	NXT-34222	\$ 87	<b>Combination Hand Activated Tilt, Brake and Wheel Lock System</b> (E2228)	NXT-30277	n/c
<b>5" x 1.4" Poly, Aluminum Mag</b>	NXT-42000	\$ 92	Note: Includes manual wheel locks to be used with public transportation.		
<b>6" x 1" Poly, Plastic Mag</b>	NXT-34666	n/c	Note: Not available with 12" poly plastic mag wheels.		
<b>6" x 1.5" Poly, Plastic Mag</b>	NXT-34999	\$ 87	Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.		
<b>6" x 1.4" Poly, Aluminum Mag</b>	NXT-44000	\$ 92			
<b>6" x 2" Poly with Tread</b>	NXT-34444	\$ 92			
<b>6" x 1" Pneumatic Plastic Mag</b> (E2214/E2215)	NXT-34667	\$ 77			
<b>7" x 1" Poly, Plastic Mag</b>	NXT-34777	\$ 77			
<b>7" x 2" Polyurethane</b>	NXT-34200	\$ 118			
Note: All 7" casters not available with all 12" rear wheels.					
<b>Rear Wheels - Tires - Accessories</b>			<b>Rear Anti-Tipper Options</b>		
Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.			<b>Swing-Away Anti-Tip Assembly Foot Activated</b> (E0971)	NXT-31605	\$ 153
<b>Mag Wheels</b>			<b>Fixed Anti-Tip Assembly</b> (E0971)	NXT-31606	\$ 97
<b>12" Low Profile Poly Mag</b> (E2218)	NXT-31526	\$ 36			
Note: Not available with any foot/hand activated brake/wheel lock system.					
Note: Not available with airless inserts.					
<b>12" Semi Pneumatic Plastic Mag</b> (K0069)	NXT-31523	\$ 168			
Note: Not available with airless inserts.					
<b>12" Pneumatic Plastic Mag</b>	NXT-31525	n/c			
Note: All 12" wheels not available with all handrims.					
Note: All 12" wheels not available with all 7" casters.					
<b>12" Airless Inserts</b> (E2213)	NXT-31019	\$ 77			
<b>16" Pneumatic Plastic Mag</b>	NXT-31530	\$ 41			
Note: 16" mag wheels not available with plastic coated handrims.					
<b>16" Airless Inserts</b> (E2213)	NXT-31020	\$ 77			
<b>18" Pneumatic Plastic Mag</b>	NXT-31533	\$ 41			
Note: 18" mag wheels not available with plastic coated handrims.					
<b>18" Airless Inserts</b> (E2213)	NXT-31021	\$ 77			

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

**Description Part No. Retail**

**Description Part No. Retail**

**Seat Pan Options**

- Fixed Solid Seat Pan - Bolteron** (E2231) NXT-39010 \$ 150  
**Note:** Not available with folding frame.
- Fixed Solid Seat Pan - Aluminum** (E2231) NXT-39030 \$ 170  
**Note:** Fixed Seat Pan raises the effective seat to floor height by 1/2".  
**Note:** Not available with folding frame.
- Removable Solid Seat Pan - Bolteron** (E2231) NXT-39012 \$ 170
- Removable Solid Seat Pan - Aluminum** (E2231) NXT-39032 \$ 255  
**Note:** Removable Seat Pan raises the effective seat to floor height by 5/8".  
**Note:** Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.  
**Note:** Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.  
**Note:** All seat pans include 2" x 2" notches for growth adjustment.  
**Note:** All seat pans made to fit frame dimensions unless specified below.

Seat Width \_\_\_\_\_" X Seat Depth \_\_\_\_\_"

**Accessories**

- Medical Necessities Pouch** NXT-23207 \$ 68  
**Note:** Black nylon.
- Footrest Hanger Pads - Pair** SFDI-9HP \$ 92  
**Note:** Wrap-around pads made from black neoprene.
- Frame Pads - Pair** 10911-2 \$ 82  
**Note:** Wrap-around pads made from black neoprene.  
**Note:** Typically used to pad exposed back post tubes.
- Freedom Backpack** NXT-10900 \$ 56  
**Note:** Black Nylon only with Freedom Designs logo.  
**Note:** Monogramming is not available on Freedom Backpack.
- Freedom Backpack with Complete Freedom Seating System** NXT-SYS-10900 n/c  
**Note:** Black Nylon only with Freedom Designs logo.  
**Note:** Monogramming is not available on Freedom Backpack.  
**Note:** Complete Freedom Designs seating system including solid seat, solid back and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.
- Provide Freedom/PinDot seating quote/PO #** \_\_\_\_\_
- Freedom Canopy - Flip-Back** NXT-10905 \$ 128  
**Note:** Canopy assembly includes removable mounting hardware.  
**Note:** Canopy is available in black only with Freedom Designs logo.  
**Note:** Canopy includes rear flap with zippered storage pocket.  
**Note:** Canopy includes clear-view window.
- Freedom Canopy - Flip-Back with Complete Freedom Seating System** NXT-SYS-10905 n/c  
**Note:** Canopy assembly includes removable mounting hardware.  
**Note:** Canopy is available in black only with Freedom Designs logo.  
**Note:** Canopy includes rear flap with zippered storage pocket.  
**Note:** Canopy includes clear-view window.  
**Note:** Complete Freedom Designs seating system including solid seat, solid back and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.
- Provide Freedom/PinDot seating quote/PO #** \_\_\_\_\_

**Anterior Support Surfaces – (Trays)**

- 3/8" Lexan Tray - Channel Bracket Hardware** (E0950) 08-2000 \$ 271
- 3/8" Lexan Tray - Adjustable Knob Hardware & Ribs** (E0950) 08-6000 \$ 219
- 1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs** (E0950) 08-13000 \$ 183
- 1/4" Lexan Tray - Hook & Loop Straps** (E0950) 08-16000 \$ 106

**HCPCS CODES** are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Choose NXT Seat Rail to Floor Height \_\_\_\_\_"**

- Note:** Use chart to determine seat rail to floor height. Seat rail to floor height is determined by frame type, rear wheels, front caster size and tilt pivot height.
- Note:** Caster sizes vary slightly in height by type.
- Note:** Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".
- Note:** Seat pan choice increases seat rail to floor height 3/8" to 3/4".

**NXT Wheel / Caster Seat Rail to Floor Heights**

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Mini	15-15.5"	16"	6"	Mini	18"
12"	4"	Mini	16"	16"	7"	Mini	19"
12"	5"	Mini	18"	18"	3"	Mini	15-15.5"
12"	6"	Mini	18.5"	18"	4"	Mini	17.5-19"
16"	3"	Mini	15.5-16"	18"	5"	Mini	18-19"
16"	4"	Mini	16.5-17"	18"	6"	Mini	18"
16"	5"	Mini	17"	18"	7"	Mini	N/A