



Our kids are our passion.™

2241 N Madera Rd

NXT Mini
Generation NEXT

Simi Valley, CA 93065

Pediatric Manual Wheelchair Transport Package Included

Phone: (800) 331-8551 Fax: (888) 582-1509

www.freedomdesigns.com



Date: _____ Buyer: _____
 Company: _____
 Marked For: _____
 Email: _____

Quote Order PO# _____
 Location: _____
 Phone: _____
 Fax: _____

Description	Part No.	Retail
-------------	----------	--------

NXT Pediatric Frame Options

- NXT-MINI Pediatric Folding** (E1234) NXT-30331 **\$2950**
Transport Included
Note: Seat frame offered as seat rail only for use with seating system mounted footrest, or for use with 75° or 85° fixed standard legrests or with 15° abduction frame insert with 75° fixed legrest.
Note: Legrests manufactured with holes located in 1/2" increments for footrest hardware attachment and footrest plate adjustment.
Note: Not available with any Pop-Up or Latch swing away or elevating footrest hangers.
Note: Freedom Future Fit Program includes folding cross brace and footplates for growing frame WIDTH only.
Note: ANY DEPTH growth of NXT-MINI requires conversion to a SMALL frame NXT. Contact Customer Service for pricing and availability.
Note: All vent and battery trays not available with NXT-MINI.

- NXT-MINI Pediatric Rigid** (E1233) NXT-30341 **\$2850**
Transport Included
Note: Seat frame offered as seat rail only for use with seating system mounted footrest, or for use with 75° or 85° fixed standard legrests or with 15° abduction frame insert with 75° fixed legrest.
Note: Legrests manufactured with holes located in 1/2" increments for footrest hardware attachment and footrest plate adjustment.
Note: Not available with any Pop-Up or Latch swing away or elevating footrest hangers.
Note: Freedom Future Fit Program includes rigid cross tubes and footplates for growing frame WIDTH only.
Note: ANY DEPTH growth of NXT-MINI requires conversion to a SMALL frame NXT. Contact Customer Service for pricing and availability.
Note: All vent and battery trays not available with NXT-MINI.

Transportation Option

- Transportation Package Included** NXT-99999-1 **n/c**
Note: Must be installed by Freedom Designs manufacturing.
Note: Not available as parts.
Note: Not available with all vent and battery trays.
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.
Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit

Tilt Activations

- Dual Lever Hand Activated Tilt** NXT-30272 **n/c**
- Foot Activated Tilt** NXT-30273 **n/c**
Note: Foot pedal installed on right rear side frame.
- Combination Dual Lever Hand Activated Tilt, Brake and Wheel Lock System** NXT-30277 **n/c**
Note: Includes manual wheel locks to be used with public transportation.

Vent & Battery Trays

Note: All vent & battery trays not available with NXT-MINI.
Note: All vent and battery trays not with transport option.

NXT-MINI Rider Weight Capacity is 75 lbs.

Description	Part No.	Retail
-------------	----------	--------

Frame Colors

- Sugar Cookie **n/c**
- Neon Yellow **n/c**
- Yellow **n/c**
- Bright Pink **n/c**
- Pink Parfait **n/c**
- Sugar Plum **n/c**
- Vigorous Violet **n/c**
- Candy Violet **n/c**
- Candy Raspberry **n/c**
- Extreme Green **n/c**
- Sparkle Green **n/c**
- Mystic Teal **n/c**
- Blizzard Blue **n/c**
- Pacific Blue **n/c**
- Blue Blast **n/c**
- Electric Blue **n/c**
- Galaxy Blue **n/c**
- Passion Purple **n/c**
- Octane Orange **n/c**
- Red Sunset **n/c**
- Really Red **n/c**
- Fiery Red **n/c**
- Twilight **n/c**
- Stormy Gray **n/c**
- Black **n/c**

Frame Dimensions

Specify Frame Width Setting

- 10" 11" 12" **NXT-30250** **n/c**
Note: NXT-MINI frame is width adjustable from 10" to 12" wide.
Note: Frame widths measured across and from outside of seat frame tubes.
Note: Seat frame offered as seat rail only for use with a seating system mounted footrest, or for use with 75° or 85° fixed standard legrests or with 15° abduction frame insert with 75° fixed legrest.
Note: With seating, positioning accessories, and attachment hardware MINI can accommodate as little as 5" seat width and 6" seat depth.

Back Post Locations

Specify Frame Type and Frame Length

- Note:** Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.
- Mini Frame** **NXT-30510** **n/c**
 - 8" **NXT-30500XTP** **n/c**
 - 9" **NXT-30501XTP** **n/c**
 - 10" **NXT-30502XTP** **n/c**
 - 11" **NXT-30503XTP** **n/c**
 - 12" **NXT-30504XTP** **n/c**
 - * 13" Must select 2" extension **NXT-30505XTP** **n/c**
 - * 14" Must select 2" extension **NXT-30506XTP** **n/c**
 - 2" Extension - Pair (K0108) **NXT-30525** **\$ 25**

- Specify Back Post Location** (if other than specified above) _____"
Note: Floor to seat frame height range is 15" to 19".
Note: Must choose 2" extension for 13" and 14" frame lengths.
Note: With seating, positioning accessories, and attachment hardware MINI can accommodate as little as 6" seat width and 6" seat depth.

Note: NXT Mini has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 6" + back / hardware thickness of 2" = back post location set at 8".

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Back Post Tubes

Specify Push Handle, Back Post and Extensions

Note: Back Posts are measured from top of seat rail to top of handle bend.

- Standard 90° Bend Push Handles and Back Post Tubes** NXT-30518 n/c

Standard Fixed Height

- Short Push Handle - 20" high** NXT-30515 n/c
- Medium Push Handle - 23" high** NXT-30514 n/c
- Tall Push Handle - 26" high** NXT-30516 n/c

OR

Back Post Adjustable Height

- Adjustable Height 90° Bend Push Handle - Short 18" to 24" High**  NXT-30533 \$ 150
- Adjustable Height 90° Bend Push Handle - Tall 20" to 26" High**  NXT-30534 \$ 150

Note: Adjustable height back posts includes Standard 90° bend.

Note: Handles quickly adjust in height with a range of 6".

Note: ALL stroller / push handle extensions not available with adjustable height back posts.

Stroller Handle Extensions

Note: All stroller handle ext. add approximately 6.5" to push handle height.

Note: All stroller / push handle extensions not available with adjustable height back posts.

- Push Handle Extension - Width Adjustable Removable** NXT-30513 n/c
- Stroller Handle Extension - Width Adjustable Removable** NXT-30526 \$ 125
- Push Button Angle Adjustable Stroller Handle Extension (PB)** NXT-30536 \$ 175
 Note: PB stroller handles not available with frames less than 12" wide.
 Note: PB stroller handles are not removable and not with folding frame.
 Note: All stroller/push handle extensions include black neoprene wrap around pad.
- Individual Stroller Handle Extensions Removable - Pair** NXT-30517 \$ 100

Specify Knee to Heel

- 3" 3.5" 4" 4.5" 5" 5.5" 6"
- 6.5" 7" 7.5" Requires 4" Legrest Extension
- 8" 8.5" 9" Requires 4" Legrest Extension
- 9.5" 10" Requires 4" Legrest Extension

Legrests, Footrests and Extensions

Legrests

NXT-Mini offers a series of fixed legrest options which insert into the seat rail. Footplate mounting holes are placed at 1/2" increments for height adjustment. Footplate range is measured from the top of the seat rail to the top of the footplate. Range of adjustment is 3" to 10" with attachment hardware and optional 4" legrest extension.

- 75° Fixed Legrest - Pair**  NXT-30370 n/c
 Note: Bottom 2 mounting holes usable only with 4" extension.
 Note: 75° & 85° standard legrests are interchangeable.
 Note: Always with Mid-Line mounted footrest plates.
- 85° Fixed Legrest - Pair**  NXT-30371 n/c
 Note: Bottom 2 mounting holes usable only with 4" extension.
 Note: 75° & 85° standard legrests are interchangeable.
 Note: Always with mid-line mounted footrest plates.
 Note: Seat rail must be moved 1.5" forward for base frame clearance when installing 4" extensions.
- 15° Abduction Seat Rail Insert with 75° Fixed Legrest** NXT-30327 \$ 150
 Note: Features a fixed 15° lateral bend in the seat rail insert to allow for abducted lower extremity positioning.
 Note: Footplates adjust for 15° abduction seat rail or for mid-line positioning.
 Note: Not available with swing-away or pop-up legrest hangers.
 Note: Not available with Bolteron flip-up footboard.
 Note: Not available with 4" length adjustable extensions.
- Seat Rail Only - Mini** NXT-30326 n/c
 Note: For use with seat mounted footrest.

Footplates & Footboards

- Composite Flip-Up** NXT-17070 n/c
- MINI Aluminum Angle Adjustable, Flip-up**  (K0040) NXT-36019 \$ 150
- Mid-line Mounted Footplates**  NXT-36025 n/c
 Note: Always with fixed standard legrests.
- Abduction Mounted Footplates**  NXT-36026 n/c
 Note: Only with abduction seat rail with 75° legrest.
 Note: Set externally rotated 15°.
 Note: Can be adjusted for mid-line foot positioning.
- Heel Loops** (E0951) NXT-36022 \$ 50
- Calf Strap - Webbing** (K0038) NXT-36024 \$ 40
- Bolteron Angle Adjustable Flip-Up Footboard**  (K0040) NXT-36030 \$ 150
 Note: Not available with 15° Abduction Seat Rail Insert.



Legrest Extensions

- 4" Mini Legrest Extensions**  NXT-30289 n/c
 Note: Seat rail must be moved 1.5" forward for base frame clearance when installing 4" extensions.
- Specify Knee to Heel Measurement _____"**

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Armrests and Armrest Pads

- 8.5" Long Flip-Back**  NXT-38218 n/c
 - 10.5" Long Flip-Back**  NXT-38204 n/c
- Note:** All Flip-Back armrests are mounted to backpost tubes.
Note: 8.5" armrests not available with desk and full length armrest pads.
Note: Flip-Back armrests length are measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.
Note: Armrests measure 6" to 11" from seat rail with foam wrap arm pads installed.
Note: Armrests measure 7" to 12" from seat rail with desk length arm pads installed.
- Foam Wraps - Standard** NXT-30816 n/c
 - Desk Length Pads** NXT-38213 n/c
- Note:** 8.5" armrests not available with desk and full length armrest pads.
- Specify Seat Rail to Top of Armrest Measurement** _____ "

Front Wheels - Casters

- 3" Polyurethane** NXT-34888 \$ 55
 - 3" Polyurethane Light-Up** NXT-34111 \$ 55
- Note:** 3" casters recommended for low seat to floor heights.
- 4" Polyurethane** NXT-40000 \$ 65
 - 4" Polyurethane Light-Up** NXT-41000 \$ 99
 - 4" x 1.4" Poly, Aluminum Mag** NXT-43000 \$ 90
- 5" x 1" Poly, Plastic Mag** NXT-34555 n/c
 - 5" Polyurethane Light-Up** NXT-35050 \$ 99
 - 5" x 1.5" Poly, Plastic Mag** NXT-34222 \$ 85
 - 5" x 1.4" Poly, Aluminum Mag** NXT-42000 \$ 90
- 6" x 1" Poly, Plastic Mag** NXT-34666 n/c
 - 6" x 1.5" Poly, Plastic Mag** NXT-34999 \$ 85
 - 6" x 1.4" Poly, Aluminum Mag** NXT-44000 \$ 90
 - 6" x 2" Poly with Tread** NXT-34444 \$ 90
 - 6" x 1" Pneumatic Plastic Mag** (E2214/ E2215) NXT-34667 \$ 75
 - 7" x 1" Poly, Plastic Mag** NXT-34777 \$ 75
 - 7" x 1.5" Poly, Plastic Mag** NXT-34000 \$ 116
 - 7" x 2" Polyurethane** NXT-34200 \$ 116
- Note:** All 7" casters not available with all 12" rear wheels.

Rear Wheels - Tires - Accessories

Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.

Mag Wheels


- 12" Low Profile Poly Mag** (E2218) NXT-31526 \$ 35
- Note:** Not available with any foot/hand activated brake/wheel lock system.
Note: Not available with airless inserts.
- 12" Semi Pneumatic Plastic Mag** NXT-31523 \$ 165
- Note:** Not available with airless inserts.
- 12" Pneumatic Plastic Mag** NXT-31525 n/c
- Note:** All 12" wheels not available with all handrims.
Note: All 12" wheels not available with all 7" casters.
- 12" Airless Inserts** (E2213) NXT-31019 \$ 75
- 16" Pneumatic Plastic Mag** NXT-31530 \$ 40
- Note:** 16" mag wheels not available with plastic coated handrims.
- 16" Airless Inserts** (E2213) NXT-31020 \$ 75
 - 16" Aluminum Handrims** NXT-31014 \$ 166
- 18" Pneumatic Plastic Mag** NXT-31533 \$ 40
- Note:** 18" mag wheels not available with plastic coated handrims.
- 18" Airless Inserts** (E2213) NXT-31021 \$ 75
 - 18" Aluminum Handrims** NXT-31015 \$ 166

Rear Wheels - Tires - Accessories (continued)
Spoke Wheels



Spoke Wheels Compatible with Drum Braking Systems
AVAILABLE SUMMER 2014

- 16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel**  NXT-31534 n/c
- Note:** Rim, hub, & spokes are anodized silver.
- 16" Silver Aluminum Handrim**  NXT-31030 \$ 75
 - 16" Black Plastic Coated Handrim**  NXT-31032 \$ 117
 - 16" Airless Inserts** (E2213) NXT-31020 \$ 75
- 16" Black Pneumatic Tire & Black Spoke Aluminum Wheel**  NXT-31535 n/c
- Note:** Rim, hub, & spokes are anodized black.
- 16" Black Aluminum Handrim**  NXT-31031 \$ 75
 - 16" Black Plastic Coated Handrim**  NXT-31032 \$ 117
 - 16" Airless Inserts** (E2213) NXT-31020 \$ 75
- 18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel**  NXT-31537 n/c
- Note:** Rim, hub, & spokes are anodized silver.
- 18" Silver Aluminum Handrim**  NXT-31034 \$ 75
 - 18" Black Plastic Coated Handrim**  NXT-31036 \$ 117
 - 18" Airless Inserts** (E2213) NXT-31021 \$ 75
- 18" Black Pneumatic Tire & Black Spoke Aluminum Wheel**  NXT-31538 n/c
- Note:** Rim, hub, & spokes are anodized black.
- 18" Black Aluminum Handrim**  NXT-31035 \$ 75
 - 18" Black Plastic Coated Handrim**  NXT-31036 \$ 117
 - 18" Airless Inserts** (E2213) NXT-31021 \$ 75
- Note:** 16" & 18" spoke wheels do not include handrims.

Wheel Locks

- Push to Lock - Standard** NXT-31501 n/c
 - Extension Handles** (E0961) NXT-31801 \$ 60
- Foot Activated Wheel Lock** (E2228) NXT-31508 n/c
- Note:** Includes manual wheel locks to be used with public transportation.
Note: Not available with 12" poly plastic mag wheels.
- Hand Activated Brake and Wheel Lock System** (E2228) NXT-30275 n/c
- Note:** Includes manual wheel locks to be used with public transportation.
Note: Not available with 12" mag wheels with low pro poly tires.
Note: Used to control direction, speed, brake and lock wheels from strolling position.
- Combination Dual Lever Hand Activated Tilt, and Wheel Lock System**  (E2228) NXT-30277 n/c
- Note:** Includes manual wheel locks to be used with public transportation.
Note: Not available with 12" poly plastic mag wheels.
Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.

Rear Anti-Tipper Options

- Swing-Away Anti-Tip Assembly Foot Activated**  (E0971) NXT-31605 \$ 150
- Fixed Anti-Tip Assembly**  (E0971) NXT-31606 \$ 95

NXT-MINI Rider Weight Capacity is 75 lbs.

Page 3 of 4 **NXT-MINI Pediatric Transport Included**

CE

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
Solid Seat Options		
<input type="radio"/> Fixed Solid Seat Pan - Bolteron (E2231) NXT-39010 \$ 100 Note: Not available with folding frame. Note: All seat pans include 2" x 2" notches for growth adjustment. Note: All Bolteron seat pans made to fit frame dimensions unless specified below.		
<input type="radio"/> Removable Solid Seat Pan - Bolteron (E2231) NXT-39012 \$ 150 Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware. Note: All seat pans include 2" x 2" notches for growth adjustment. Note: All Bolteron seat pans made to fit frame dimensions unless specified below.		
<input type="radio"/> Width _____ " X Depth _____ "		

Choose NXT Seat Rail to Floor Height _____"
Note: Use chart to determine seat rail to floor height. Seat rail to floor height is determined by frame type, rear wheels, front caster size and tilt pivot height.
Note: Caster sizes vary slightly in height by type.
Note: Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".
Note: Seat pan choice increases seat rail to floor height 3/8" to 3/4".

NXT Wheel / Caster Seat Rail to Floor Heights							
Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Mini	15-15.5"	16"	6"	Mini	18"
12"	4"	Mini	16"	16"	7"	Mini	19"
12"	5"	Mini	18"	18"	3"	Mini	15-15.5"
12"	6"	Mini	18.5"	18"	4"	Mini	17.5-19"
16"	3"	Mini	15.5-16"	18"	5"	Mini	18-19"
16"	4"	Mini	16.5-17"	18"	6"	Mini	18"
16"	5"	Mini	17"	18"	7"	Mini	N/A

AVAILABLE SUMMER 2014		
<input type="radio"/> Adjustable, Fixed Solid Seat Pan - Aluminum (E2231) NXT-39030 \$ 150 Note: Not available with folding frame. Note: Adjustable Solid Seat Pan adjusts 2" in width and depth to grow with NXT frame. Note: All seat pans include 2" x 2" notches for growth adjustment. Note: All aluminum seat pans available in standard sizes to fit frame dimensions.		
<input type="radio"/> Adjustable, Removable Solid Seat Pan - Aluminum (E2231) NXT-39032 \$ 250 Note: Adjustable Solid Seat Pan adjusts 2" in width and depth to grow with NXT frame. Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware. Note: All seat pans include 2" x 2" notches for growth adjustment. Note: All aluminum seat pans available in standard sizes to fit frame dimensions.		

Freedom Designs & Pindot
 Are Pleased To Announce
"A Sizzling Summer Sensation"
 And For a Limited Time It's...

FREE * FREE * FREE

When ordering a Freedom
 ...Wheelchair and a
 complete Freedom or Pindot
 seating system including seat,
 back, and attachment hardware,
 simply provide the seating Quote or
 seating PO number
 below and receive a
**Freedom Canopy and
 Freedom Backpack**

*** FREE ***

Accessories		
<input type="radio"/> Medical Necessities Pouch NXT-23207 \$ 67 Note: Black nylon.		
<input type="radio"/> Footrest Hanger Pads - Pair SFDI-9HP \$ 90 Note: Wrap-around pads made from black neoprene.		
<input type="radio"/> Frame Pads - Pair 10911-2 \$ 80 Note: Wrap-around pads made from black neoprene. Note: Typically used to pad exposed back post tubes.		
<input type="radio"/> Freedom Backpack (New) NXT-10900 \$ 55 Note: Black Nylon only with Freedom Designs logo. Note: Monogramming is not available on Freedom Backpack.		
<input type="radio"/> Freedom Canopy - Flip-Back (New) NXT-10905 \$ 125 Note: Canopy assembly includes removable mounting hardware. Note: Canopy is available in black only with Freedom Designs logo. Note: Canopy includes rear flap with zippered storage pocket. Note: Canopy includes clear-view window.		

Freedom Backpack with Complete Freedom Seating System (New) NXT-SYS-10900 **n/c**
Note: Black Nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom backpack.
Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

Provide Freedom/PinDot Seating Quote/PO # _____

Freedom Canopy - Flip-Back with Complete Freedom Seating System (New) NXT-SYS-10905 **n/c**
Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.
Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

Provide Freedom/PinDot Seating Quote/PO # _____

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.