



Our kids are our passion.™

# NXT-SC

Generation NEXT

## Adult Manual Wheelchair



Yes, you can.

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: \_\_\_\_\_ Buyer: \_\_\_\_\_  
Company: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Email: \_\_\_\_\_

Quote Order PO# \_\_\_\_\_  
Location: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Transport Optional Base</b>			<b>Vent &amp; Battery Trays</b>		
<b>NXT-SC Adult Folding</b> (E1161) NXT-30332		<b>\$3465</b>	<b>Note:</b> All vent and battery trays not with transport.		
<b>Transport Optional</b>			<b>Note:</b> All vent trays on NXT only with Dual Lever Hand Activated Tilt.		
<b>Note:</b> Transportation Option selected separately below.			<b>Note:</b> All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.		
<b>Note:</b> Includes Freedom Future Fit Program.			<b>Vent Tray - Fixed</b> (E1029) 11-7000		<b>\$ 439</b>
<b>NXT-SC Adult Rigid</b> (E1161) NXT-30342		<b>\$3465</b>	<b>Battery Tray - Fixed</b> (K0108) 11-8000		<b>\$ 429</b>
<b>Transport Optional</b>			<b>Quick Release Tray Hardware</b> (E1028) 563		<b>\$ 272</b>
<b>Note:</b> Transportation Option selected separately below.			<b>Note:</b> For LP10 style vent and battery tray on frames 16" wide and greater.		
<b>Note:</b> Includes Freedom Future Fit Program.			<b>Note:</b> Vent tray is 4" narrower than frame width and 16" deep.		
<b>Transportation Option</b>			<b>Note:</b> Battery tray is 4" narrower than frame width and 10" deep.		
<b>Note:</b> Used when Transport Optional frame is selected above.			<b>Frame Colors</b>		
<b>Transportation Package</b> NXT-99998		<b>\$ 250</b>	Sugar Cookie	n/c	Pacific Blue n/c
<b>Note:</b> Must be installed by Freedom Designs manufacturing.			Neon Yellow	n/c	Blue Blast n/c
<b>Note:</b> Not available with all vent and battery trays.			Yellow	n/c	Electric Blue n/c
<b>Note:</b> Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.			Bright Pink	n/c	Galaxy Blue n/c
<b>Note:</b> Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.			Pink Parfait	n/c	Passion Purple n/c
<b>Transport Included Base</b>			Sugar Plum	n/c	Octane Orange n/c
<b>NXT-SC Adult Folding</b> (E1161) NXT-30333		<b>\$3715</b>	Candy Violet	n/c	Red Sunset n/c
<b>Transport Included</b>			Candy Raspberry	n/c	Really Red n/c
<b>Note:</b> Includes Freedom Future Fit Program.			Extreme Green	n/c	Fiery Red n/c
<b>Note:</b> Not available with all vent and battery trays.			Sparkle Green	n/c	Twilight n/c
<b>NXT-SC Adult Rigid</b> (E1161) NXT-30343		<b>\$3715</b>	Mystic Teal	n/c	Stormy Grey n/c
<b>Transport Included</b>				n/c	Black n/c
<b>Note:</b> Includes Freedom Future Fit Program.			<b>Frame Dimensions</b>		
<b>Note:</b> Not available with all vent and battery trays.			<b>Specify Frame Width Setting</b>		
<b>Seat Rail Options</b>			14"	15"	16" Adjusts 14" to 16" NXT-30240 n/c
<b>Note:</b> Not available with ALL Contracture, 2D & 3D Dynamic footrest systems.			15"	16"	17" Adjusts 15" to 17" NXT-30290 n/c
<b>NXT LEAF Seat Frame- Pair</b> (K0108) NXT-30324		<b>\$ 714</b>	16"	17"	18" Adjusts 16" to 18" NXT-30260 n/c
<b>NXT LEAF Seat Frame- Left</b> (K0108) NXT-30324-L		<b>\$ 357</b>	17"	18"	Adjusts 17" to 18" NXT-30261 n/c
<b>NXT LEAF Seat Frame- Right</b> (K0108) NXT-30324-R		<b>\$ 357</b>	18"		Not width adjustable NXT-30262 n/c
<b>Note:</b> LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.			19"		NXT-30263 n/c
<b>Note:</b> Freedom customized seats can be ordered to fit pre-configured LEAF frame models.			20"	(E2201) NXT-30264	<b>\$ 400</b>
<b>Note:</b> No seating system is provided with NXT frames. A seating system must be purchased separately.			<b>Note:</b> Frame widths measured across and from outside of seat frame tubes.		
<b>NXT SRO Seat Frame- Pair</b> (K0108) NXT-30325		<b>\$ 306</b>	<b>Note:</b> Vent and battery trays only with 14" and wider frames.		
<b>NXT SRO Seat Frame- Left</b> (K0108) NXT-30325-L		<b>\$ 153</b>			
<b>NXT SRO Seat Frame- Right</b> (K0108) NXT-30325-R		<b>\$ 153</b>			
<b>Note:</b> SRO is manufactured without standard footrest hanger attachments.					
<b>Note:</b> SRO is used when customer attaches leg/footrest(s) directly to the underside of the solid seat.					
<b>Tilt Activations</b>					
<b>Dual Lever Hand Activated Tilt</b>	NXT-30272	n/c			
<b>Foot Activated Tilt</b>	NXT-30273	n/c			
<b>Note:</b> Foot pedal installed on right rear side frame.					
<b>Note:</b> Not with ALL vent and battery trays.					

NXT Rider Weight Capacity is 225 lbs.

Page 1 of 7

NXT-SC Adult

CE

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

**Back Post Locations**  
Specify Frame Type and Frame Length

**Notes:**

**Note:** Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

<b>Medium Frame &amp; Short Tilt Bracket</b>	NXT-30512	n/c
14"	NXT-30506	n/c
15"	NXT-30521	n/c
16"	NXT-30522	n/c
17"	NXT-30523	n/c
18"	NXT-30524	n/c
19" Must select 2" extension	NXT-30527	n/c
20" Must select 2" extension (E2203)	NXT-30528	\$ 500
21" Must select 3" extension (E2203)	NXT-30529	\$ 500
2" Extension - Pair (K1008)	NXT-SC-30525	std
3" Extension - Pair (K1008)	NXT-SC-30530	std

**Specify Back Post Location** (if other than specified above) \_\_\_\_\_ "

**Note:** Floor to seat frame height range is 12.5" to 16.5".

**Note:** Must choose an extension for 19", 20", and 21" frame lengths.

**Note:** 22" & 24" wheels not available with medium frame and short tilt bracket.

<b>Medium Frame &amp; Tall Tilt Bracket</b>	NXT-30540	n/c
14"	NXT-30506TP	n/c
15"	NXT-30521TP	n/c
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
19" Must select 2" extension	NXT-30527TP	n/c
20" Must select 2" extension (E2203)	NXT-30528TP	\$ 500
21" Must select 3" extension (E2203)	NXT-30529TP	\$ 500
22" Must select 4" extension (E2204)	NXT-30535TP	\$ 900
2" Extension - Pair (K1008)	NXT-SC-30525	std
3" Extension - Pair (K1008)	NXT-SC-30530	std
4" Extension - Pair (K1008)	NXT-SC-30531	std

**Specify Back Post Location** (if other than specified above) \_\_\_\_\_ "

**Note:** Tall tilt bracket increases floor to seat rail height 1.5".

**Note:** Floor to seat frame height range is 14" to 18".

**Note:** Must choose an extension for 19", 20", 21", and 22" frame lengths.

<b>Long Frame &amp; Tall Tilt Bracket</b>	NXT-30520	n/c
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
19"	NXT-30527TP	n/c
20" (E2203)	NXT-30528TP	\$ 500
21" (E2203)	NXT-30529TP	\$ 500
22" (E2204)	NXT-30535TP	\$ 900
23" Must select 2" extension (E2204)	NXT-30537TP	\$ 900
24" Must select 2" extension (E2204)	NXT-30538TP	\$ 900
25" Must select 3" extension (E2204)	NXT-30542TP	\$ 900
26" Must select 4" extension (E2204)	NXT-30543TP	\$ 900
2" Extension - Pair (K1008)	NXT-SC-30525	std
3" Extension - Pair (K1008)	NXT-SC-30530	std
4" Extension - Pair (K1008)	NXT-SC-30531	std

**Note:** Tall tilt bracket increases floor to seat rail height 1.5".

**Note:** Floor to seat frame height range is 14" to 18".

**Note:** Must choose an extension for 23", 24", 25", and 26" frame lengths.


**Specify Back Post Location** (if other than specified above) \_\_\_\_\_ "

**5° Anterior Tilt Seat Frame** NXT-30276 n/c

**Note:** Limits posterior tilt to 40°, includes split collars for adjustment.

**Note:** Medium with Tall Tilt and Long frames only.

**Note:** NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14".



**Back Post Tubes**

**Specify Push Handle, Back Post and Extensions**

**Standard 90° Bend Push Handles and Back Post Tubes** NXT-30518 n/c

**Ergo Bend Push Handles and Back Post Tubes** NXT-30519 n/c

**Note:** All push/stroller handle extensions not recommended with Ergo Bend.

**Seating Dynamics Rocker Back** (K0108) NXT-30532 \$1734

**Note:** Angle adjustable from 72° to 108° in 6° increments.

**Note:** Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.

**Note:** Locks into a non-dynamic position when necessary.

**Note:** Does not allow the back post tubes to fold forward.

**Back Post Height**

**Note:** Back Posts are measured from top of seat rail to top of handle bend.


**Short Push Handle - 20" high** NXT-30515 n/c

**Medium Push Handle - 23" high** NXT-30514 n/c

**Tall Push Handle - 26" high** NXT-30516 n/c

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Back Post Tubes (continued)</b>			<b>Legrests, Footrests and Extensions(continued)</b>		
<b>Stroller Handle Extensions</b>			<b>Footplates &amp; Footboards</b>		
<p><b>Note:</b> All stroller handle ext. add approximately 6.5" to push handle height.  <b>Note:</b> All stroller/push handle extensions except Individual include a black neoprene wrap around pad.</p>			<p><b>Omit Footrest Plates</b> NXT-36404</p> <p><b>Composite Flip-Up</b> NXT-17070 <b>n/c</b></p> <p><b>Kids Aluminum Angle Adjustable, Flip-up</b> (K0040) NXT-36020 <b>\$ 160</b>  <b>Note:</b> For frames 14" wide and narrower.</p> <p><b>Heel Loops</b> (E0951) NXT-36022 <b>\$ 41</b></p> <p><b>Calf Strap - Webbing</b> (K0038) NXT-36024 <b>\$ 41</b></p> <p><b>Adult Aluminum Angle Adj. Flip-up w/ Heel Loops</b> (K0040/E0951) NXT-36021 <b>\$ 180</b>  <b>Note:</b> For frames 15" wide and greater.</p> <p><b>Omit Heel Loops</b> NXT-36023 <b>n/c</b></p> <p><b>Aluminum Angle Adjustable Footboard Removable</b> (K0040) NXT-36029 <b>\$ 160</b></p>		
			<p style="text-align: center;"><b>MAPS Footplate Mounting Hardware</b></p> <p><b>MAPS</b> (Multi Axis Positioning System) - Pair NXT-36040 <b>\$ 306</b></p> <p><b>MAPS</b> (Multi Axis Positioning System) - Left NXT-36040-L <b>\$ 153</b></p> <p><b>MAPS</b> (Multi Axis Positioning System) - Right NXT-36040-R <b>\$ 153</b>  <b>Note:</b> MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.  <b>Note:</b> Must choose aluminum footplates.  <b>Note:</b> Not available on 10" &amp; 11" wide frames.  <b>Note:</b> MAPS not available with composite footplates or aluminum footboard.</p>		
<p><b>Push Handle Extension - Width Adjustable</b> NXT-30560 <b>n/c</b>  <b>Note:</b> Not removable, but allows frame to fold.  <b>Note:</b> Mounts direct to back post push handle.</p> <p><b>Push Button, Angle-Adjustable Folding Stroller Handle Extension</b> NXT-SC-30570 <b>std</b>  <b>Note:</b> Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.  <b>Note:</b> Not removable, but allows frame to fold.  <b>Note:</b> Mounts direct to back post push handle.</p> <p><b>Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty</b> NXT-30575 <b>\$ 179</b>  <b>Note:</b> Available 10W-22W wheelchair frames.  <b>Note:</b> Not removable, but allows frame to fold.  <b>Note:</b> Same as above but with metal hinge mechanism.  <b>Note:</b> Mounts direct to back post push handle.</p> <p><b>Push Button, Angle-Adjustable, Folding Push Handle Extension</b> NXT-SC-30580 <b>std</b>  <b>Note:</b> Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.  <b>Note:</b> Not removable, but allows frame to fold.  <b>Note:</b> Mounts direct to straight back post tubes.</p> <p><b>Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty</b> NXT-30585 <b>\$ 179</b>  <b>Note:</b> Available 10W-22W wheelchair frames.  <b>Note:</b> Not removable, but allows frame to fold.  <b>Note:</b> Same as above but with metal hinge mechanism.  <b>Note:</b> Mounts direct to straight back post tubes.</p> <p><b>Individual Stroller Handle Extensions Removable - Pair</b> NXT-30517 <b>\$ 102</b></p>			<p style="text-align: center;"><b>Legrest Extensions</b></p> <p><b>4" Short Length Adjustable Extensions</b> NXT-30285 <b>n/c</b></p> <p><b>6" Medium Length Adjustable Extensions</b> NXT-30286 <b>n/c</b></p> <p><b>8" Long Length Adjustable Extensions</b> NXT-30287 <b>n/c</b>  <b>Note:</b> Not available with all 90° hangers.</p> <p><b>4.5" Heavy Duty ELR Extensions</b> NXT-SC-36401 <b>std</b></p> <p><b>4.5" Heavy Duty Extensions</b> NXT-SC-36402 <b>std</b>  <b>Note:</b> Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p> <p><b>6.5" Heavy Duty Extensions</b> NXT-SC-36403 <b>std</b>  <b>Note:</b> Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p> <p><b>High Mount Brackets</b> NXT-SC-36407 <b>std</b>  <b>Note:</b> Includes flip-up mounting hardware for future growth adjustment.  <b>Note:</b> Not available with all 75° hangers.</p> <p><b>Specify Knee to Heel Measurement _____"</b> (see matrix page)</p>		
<b>Legrests, Footrests and Extensions</b>			<b>MAPS Footplate Mounting Hardware</b>		
<b>Legrests</b>			<b>Legrest Extensions</b>		
<p><b>Omit Legrest Hangers &amp; Footrests</b> NXT-33000</p> <p><b>75° Pop-Up Swing Away Hangers</b> NXT-30280 <b>n/c</b>  <b>Note:</b> Not available with high mount bracket.</p> <p><b>90° Pop-Up Swing Away Hangers</b> NXT-30281 <b>n/c</b>  <b>Note:</b> Not available with 8" adjustable extensions.</p> <p><b>90° Pop-Up Swing Away Elevating Legrests (ELRs)</b> NXT-30279 <b>\$ 255</b>  <b>Note:</b> Not available with all length adjustable extensions.  <b>Note:</b> Not available with 6.5" heavy duty extensions.</p> <p><b>75° Latch Style Swing Away Hangers</b> NXT-30282 <b>n/c</b>  <b>Note:</b> Not available with high mount bracket.</p> <p><b>90° Latch Style Swing Away Hangers</b> NXT-30283 <b>n/c</b>  <b>Note:</b> Not available with 8" adjustable extensions.</p>			<p style="text-align: center;"><b>Continue To Next Page</b></p>		
<b>NXT Rider Weight Capacity is 225 lbs.</b>			<b>NXT-SC Adult</b>		

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<h3>Legrests, Footrests and Extensions</h3> <h4>Specialized Footrest Systems</h4> <p><b>Note:</b> ALL legrest measurements are from the seat rail to footplate.  <b>Note:</b> Not available with LEAF or SRO seat rail options.</p>			<h3>Armrests and Armrest Pads (continued)</h3>		
<h4>Super Contracture Footrest Hanger (K0108)</h4>			<h4>10.5" to 13.5" Long Flip-Back, Telescoping (E0973) NXT-38208 \$ 240</h4> <p><b>Note:</b> Full length pad on 13.5" or longer armrests only.</p>		
Short Pair - 11" to 15"	NXT-30291-2	\$3570	10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973) NXT-38210		\$ 265
Long Pair - 13" to 18"	NXT-30292-2	\$3570	17" to 20" Long Flip-Back Locking, Telescoping (E0973) NXT-38215		\$ 265
<p><b>Note:</b> Super Contracture for frames 15" wide and greater.  <b>Note:</b> Super Contracture is both swing-away and removable.  <b>Note:</b> Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.  <b>Note:</b> Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar / dorsiflexion positioning.  <b>Note:</b> Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.  <b>Note:</b> Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.  <b>Note:</b> Super Contracture position may limit ankle joint/ footplate adjustment range.</p>			<p>Foam Wraps - Standard NXT-38216 n/c                  Desk Length Pads NXT-38213 n/c                  Full Length Pads NXT-38214 n/c  <b>Note:</b> 8.5" armrests not available with desk and full length armrest pads.  <b>Note:</b> Full length pad on 13.5" or longer armrests only.</p>		
<h4>2D Dynamic Footrest Hanger (K0108)</h4>			<h4>Frame Mounted Armrests</h4>		
Short Pair - 14" to 16"	NXT-30293-2	\$4080	Single Post Fixed Height, Removable - Short (K0015) NXT-38222		\$ 225
Long Pair - 16" to 18"	NXT-30294-2	\$4080	Single Post Fixed Height, Removable - Tall (K0015) NXT-38224		\$ 225
<p><b>Note:</b> 2D Dynamic for frames 15" wide and greater.  <b>Note:</b> 2D Dynamic is both swing-away and removable.  <b>Note:</b> 2D Dynamic has adjustable knee joints that have 48 positions that adjust in 7.5° increments.  <b>Note:</b> 2D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.  <b>Note:</b> 2D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.  <b>Note:</b> 2D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.  <b>Note:</b> Frame and/or front caster choice may limit 2D Dynamic adjustment range and/or cause footplate/front caster interference.</p>			<p><b>Note:</b> Single Post Fixed Height arms include a removable side guard.  <b>Note:</b> Adjustable in height with tools only. See matrix for height ranges.</p>		
<h4>3D Dynamic Footrest Hanger (K0108)</h4>			<h4>Single Post Height Adjustable, Removable - Short (E0973) NXT-38217 \$ 265</h4>		
Short Pair - 13" to 15"	NXT-30295-2	\$4590	Single Post Height Adjustable, Removable - Tall (E0973) NXT-38203		\$ 265
Long Pair - 15" to 17"	NXT-30296-2	\$4590	Dual Post Fixed Height Removable (E0973) NXT-38205		n/c
<p><b>Note:</b> 3D Dynamic for frames 15" wide and greater.  <b>Note:</b> 3D Dynamic is both swing-away and removable.  <b>Note:</b> 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.  <b>Note:</b> 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.  <b>Note:</b> 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.  <b>Note:</b> 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.  <b>Note:</b> Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.</p>			<p><b>Note:</b> All dual post armrests not available with 22" or 24" wheels.                  Dual Post Flip-Back Height Adjustable, Removable (E0973) NXT-38206 \$ 275  <b>Note:</b> All dual post armrests not available with 22" or 24" wheels.</p>		
<h4>Armrests and Armrest Pads</h4>			<h4>Desk Length Pads (NXT-38213) n/c</h4>		
Omit Armrests	NXT-38209		Full Length Pads (NXT-38214) n/c		
<h4>Back Post Mounted Armrests</h4>			<h4>Specify Seat Rail to Top of Armrest Measurement _____ "</h4>		
<p><b>Note:</b> Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.</p>			<h4>Front Wheels - Casters</h4>		
10.5" Long Flip-Back	NXT-38204	n/c	Omit Caster Wheels	NXT-34333	
13.5" Long Flip-Back	NXT-38202	n/c	3" Polyurethane	NXT-34888	\$ 56
<p><b>Note:</b> Full length pad on 13.5" or longer armrests only.</p>			3" Polyurethane Light-Up	NXT-34111	\$ 56
10.5" Long Flip-Back Locking (E0973) NXT-38201		\$ 245	<p><b>Note:</b> 3" casters recommended for low seat to floor heights.</p>		
13.5" Long Flip-Back Locking (E0973) NXT-38200		\$ 245	4" Polyurethane	NXT-40000	\$ 66
<p><b>Note:</b> Full length pad on 13.5" or longer armrests only.</p>			4" Polyurethane Light-Up	NXT-41000	\$ 101
			4" x 1.4" Poly, Aluminum Mag	NXT-43000	\$ 92
			5" x 1" Poly, Plastic Mag	NXT-34555	n/c
			5" Polyurethane Light-Up	NXT-35050	\$ 101
			5" x 1.5" Poly, Plastic Mag	NXT-SC-34222	std
			5" x 1.4" Poly, Aluminum Mag	NXT-42000	\$ 92
			6" x 1" Poly, Plastic Mag	NXT-34666	n/c
			6" x 1.5" Poly, Plastic Mag	NXT-SC-34999	std
			6" x 1.4" Poly, Aluminum Mag	NXT-44000	\$ 92
			6" x 2" Poly with Tread	NXT-34444	\$ 92
			6" x 1" Pneumatic Plastic Mag (E2214/E2215)	NXT-34667	\$ 77
			7" x 1" Poly, Plastic Mag	NXT-SC-34777	std
			7" x 2" Polyurethane	NXT-34200	\$ 118
			<p><b>Note:</b> All 7" casters not available with all 12" rear wheels.</p>		

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail
<b>Rear Wheels - Tires - Accessories</b>		
<b>Note:</b> All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.		
<b>Omit Rear Wheels and Axles</b>	NXT-31524	
<b>Mag Wheels</b>		
<b>12" Low Profile Poly Mag</b> (E2218)	NXT-31526	\$ 36
<b>Note:</b> Not available with any foot/hand activated brake/wheel lock system.		
<b>Note:</b> Not available with airless inserts.		
<b>Note:</b> All 12" wheels not available with all 7" casters.		
<b>12" Semi Pneumatic Plastic Mag</b> (K0069)	NXT-31523	\$ 168
<b>Note:</b> Not available with airless inserts.		
<b>Note:</b> All 12" wheels not available with all 7" casters.		
<b>12" Pneumatic Plastic Mag</b>	NXT-31525	n/c
<b>12" Airless Inserts</b> (E2213)	NXT-31019	\$ 77
<b>Note:</b> All 12" wheels not available with all 7" casters.		
<b>16" Pneumatic Plastic Mag</b>	NXT-SC-31530	std
<b>16" Airless Inserts</b> (E2213)	NXT-31020	\$ 77
<b>16" Aluminum Handrims</b>	NXT-31014	\$ 169
<b>Note:</b> 16" mag wheels not available with plastic coated handrims.		
<b>18" Pneumatic Plastic Mag</b>	NXT-SC-31533	std
<b>18" Airless Inserts</b> (E2213)	NXT-31021	\$ 77
<b>18" Aluminum Handrims</b>	NXT-31015	\$ 169
<b>Note:</b> 18" mag wheels not available with plastic coated handrims.		
<b>20" Pneumatic Plastic Mag with Aluminum Handrim</b>	NXT-31527	n/c
<b>20" Airless Inserts</b> (E2213)	NXT-31022	\$ 77
<b>20" Plastic Coated Handrims</b>	NXT-31016	\$ 119
<b>22" Pneumatic Plastic Mag with Aluminum Handrim</b>	NXT-31528	n/c
<b>22" Airless Inserts</b> (E2213)	NXT-31023	\$ 77
<b>22" Plastic Coated Handrims</b>	NXT-31017	\$ 119
<b>Note:</b> All 22" & 24" wheels not available with short tilt bracket.		
<b>Note:</b> All 22" & 24" wheels not available with dual post armrests.		
<b>24" Pneumatic Plastic Mag with Aluminum Handrim</b>	NXT-31529	n/c
<b>24" Airless Inserts</b> (E2213)	NXT-31024	\$ 77
<b>24" Plastic Coated Handrims</b>	NXT-31018	\$ 119
<b>Note:</b> All 22" & 24" wheels not available with short tilt bracket.		
<b>Note:</b> All 22" & 24" wheels not available with dual post armrests.		
<b>Omit Handrims from 20", 22" or 24" Wheels</b>	NXT-31006	
<b>Spoke Wheels</b>		
<b>16" Grey Pneumatic Tire &amp; Silver Spoke Aluminum Wheel</b>	NXT-31534	n/c
<b>Note:</b> Rim, hub & spokes are anodized silver.		
<b>16" Airless Inserts</b> (E2213)	NXT-31020	\$ 77
<b>16" Silver Aluminum Handrim</b>	NXT-31030	\$ 77
<b>16" Black Plastic Coated Handrim</b>	NXT-31032	\$ 77
<b>Note:</b> 16" & 18" spoke wheels do not include handrims.		
<b>16" Black Pneumatic Tire &amp; Black Spoke Aluminum Wheel</b>	NXT-31535	n/c
<b>Note:</b> Rim, hub & spokes are anodized black.		
<b>16" Airless Inserts</b> (E2213)	NXT-31020	\$ 77
<b>Note:</b> 16" & 18" spoke wheels do not include handrims.		
<b>16" Black Aluminum Handrim</b>	NXT-31031	\$ 77
<b>16" Black Plastic Coated Handrim</b>	NXT-31032	\$ 77

Description	Part No.	Retail
<b>Rear Wheels - Tires - Accessories (continued)</b>		
<b>Note:</b> 16" & 18" spoke wheels do not include handrims.		
<b>18" Grey Pneumatic Tire &amp; Silver Spoke Aluminum Wheel</b>	NXT-31537	n/c
<b>Note:</b> Rim, hub & spokes are anodized silver.		
<b>18" Airless Inserts</b> (E2213)	NXT-31021	\$ 77
<b>18" Silver Aluminum Handrim</b>	NXT-31034	\$ 77
<b>18" Black Plastic Coated Handrim</b>	NXT-31036	\$ 77
<b>Note:</b> 16" & 18" spoke wheels do not include handrims.		
<b>18" Black Pneumatic Tire &amp; Black Spoke Aluminum Wheel</b>	NXT-31538	n/c
<b>Note:</b> Rim, hub & spokes are anodized black.		
<b>18" Airless Inserts</b> (E2213)	NXT-31021	\$ 77
<b>18" Black Aluminum Handrim</b>	NXT-31035	\$ 77
<b>18" Black Plastic Coated Handrim</b>	NXT-31036	\$ 77
<b>Note:</b> 16" & 18" spoke wheels do not include handrims.		
<b>20" Pneumatic Spoke Wheel with Aluminum Handrim</b>	NXT-31540	\$ 204
<b>20" Airless Inserts</b> (E2213)	NXT-31022	\$ 77
<b>Note:</b> Not available with foot/hand activated wheel locks.		
<b>Note:</b> 20", 22" & 24" spoke wheels not available plastic coated handrims.		
<b>22" Pneumatic Spoke Wheel with Aluminum Handrim</b>	NXT-31541	\$ 204
<b>22" Airless Inserts</b> (E2213)	NXT-31023	\$ 77
<b>Note:</b> Not available with foot/hand activated wheel locks.		
<b>Note:</b> All 22" & 24" wheels not available with short tilt bracket.		
<b>Note:</b> All 22" & 24" wheels not available with dual post armrests.		
<b>Note:</b> 20", 22" & 24" spoke wheels not available plastic coated handrims.		
<b>24" Pneumatic Spoke Wheel with Aluminum Handrim</b>	NXT-31542	\$ 204
<b>24" Airless Inserts</b> (E2213)	NXT-31024	\$ 77
<b>Note:</b> Not available with foot/hand activated wheel locks.		
<b>Note:</b> All 22" & 24" wheels not available with short tilt bracket.		
<b>Note:</b> All 22" & 24" wheels not available with dual post armrests.		
<b>Note:</b> 20", 22" & 24" spoke wheels not available plastic coated handrims.		
<b>Omit Handrims from 20", 22" or 24" Wheels</b>	NXT-31006	

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Wheel Locks</b>			<b>Accessories</b>		
<b>Push to Lock - Standard</b>	NXT-31501	n/c	<b>Medical Necessities Pouch</b>	NXT-23207	\$ 68
<b>Pull to Lock</b>	NXT-31506	n/c	<b>Note:</b> Black nylon.		
<b>Extension Handles</b> (E0961)	NXT-31801	\$ 61	<b>Footrest Hanger Pads - Pair</b>	SFDI-9HP	\$ 92
			<b>Note:</b> Wrap-around pads made from black neoprene.		
<b>Foot Activated Drum Brake Wheel Lock</b> (E2228)	NXT-SC-31508	\$ 900	<b>Frame Pads - Pair</b>	10911-2	\$ 82
<b>Note:</b> Includes manual wheel locks to be used with public transportation.			<b>Note:</b> Wrap-around pads made from black neoprene.		
<b>Note:</b> Not available with 12" mag wheels with low pro poly tires.			<b>Note:</b> Typically used to pad exposed back post tubes.		
<b>Note:</b> Not available with ALL vent and battery trays.			<b>Flip-Down Rigidizing Bar</b>	NXT-23210	\$ 214
<b>Hand Activated Brake and Wheel Lock System</b> (E2228)	NXT-SC-30275	\$ 900	<b>Note:</b> Mounts to back post tubes to add rigidity to frame.		
<b>Note:</b> Includes manual wheel locks to be used with public transportation.			<b>Note:</b> Does NOT adjust in width to grow with frame.		
<b>Note:</b> Not available with 12" mag wheels with low pro poly tires.			<b>Freedom Backpack</b>	NXT-10900	\$ 56
<b>Note:</b> Used to control direction, speed, brake and lock wheels from strolling position.			<b>Note:</b> Black Nylon only with Freedom Designs logo.		
			<b>Note:</b> Monogramming is not available on Freedom Backpack.		
<b>Rear Anti-Tippers</b>			<b>Freedom Canopy - Flip-Back</b>	NXT-10905	\$ 128
<b>Curb Assist Anti-Tip Assembly Flip-Up, Adjustable</b> (E0971)	NXT-31600	\$ 97	<b>Note:</b> Canopy assembly includes removable mounting hardware.		
<b>Note:</b> Pair includes one anti-tip tube with curb assist feature.			<b>Note:</b> Canopy is available in black only with Freedom Designs logo.		
<b>Standard Anti-Tip Assembly Flip-Up, Adjustable</b> (E0971)	NXT-31603	\$ 97	<b>Note:</b> Canopy includes rear flap with zippered storage pocket.		
			<b>Note:</b> Canopy includes clear-view window.		
<b>Seat Pan Options</b>			<b>Anterior Support Surfaces – (Trays)</b>		
<b>Fixed Solid Seat Pan - Bolteron</b> (E2231)	NXT-39010	\$ 150	<b>3/8" Lexan Tray - Channel Bracket Hardware</b> (E0950)	08-2000	\$ 271
<b>Note:</b> Not available with folding frame.			<b>3/8" Lexan Tray - Adjustable Knob Hardware &amp; Ribs</b> (E0950)	08-6000	\$ 219
<b>Fixed Solid Seat Pan - Aluminum</b> (E2231)	NXT-39030	\$ 170	<b>1/4" Lexan Tray - EZ Lite Channel Hardware &amp; Ribs</b> (E0950)	08-13000	\$ 183
<b>Note:</b> Fixed Seat Pan raises the effective seat to floor height by 1/2".			<b>1/4" Lexan Tray - Hook &amp; Loop Straps</b> (E0950)	08-16000	\$ 106
<b>Note:</b> Not available with folding frame.					
<b>Removable Solid Seat Pan - Bolteron</b> (E2231)	NXT-39012	\$ 170			
<b>Removable Solid Seat Pan - Aluminum</b> (E2231)	NXT-39032	\$ 255			
<b>Note:</b> Removable Seat Pan raises the effective seat to floor height by 5/8".					
<b>Note:</b> Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.					
<b>Note:</b> Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.					
<b>Note:</b> All seat pans include 2" x 2" notches for growth adjustment.					
<b>Note:</b> All seat pans made to fit frame dimensions unless specified below.					
Seat Width _____" X Seat Depth _____"					

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

**NXT Footrest / Extension Lengths**

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Choose NXT Seat to Floor Height \_\_\_\_\_"**

**Note:** Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height.

**Note:** Caster sizes vary slightly in height by type.

**Note:** Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".

**NXT MAPS Footrest / Extension Lengths**

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

**Note:** Length from top of footrest plate to top of seat rail.  
**Note:** Without extensions bottom hanger hole not accessible for hardware.  
**Note:** Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

**Freedom Designs Armrest Heights**

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

**Note:** All heights from top of seat frame to top of armrest pad.

**NXT Wheel / Caster Seat to Floor Heights**

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5-13"	20"	3"	Tall	N/A
12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"
12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"
12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"
12"	5"	Short	15.5"	20"	5"	Tall	16-17"
12"	5"	Tall	17"	20"	6"	Short	15.5-16"
12"	6"	Short	16"	20"	6"	Tall	17-17.5"
12"	6"	Tall	17.5"	20"	7"	Short	16.5"
16"	3"	Short	12.5-13"	20"	7"	Tall	18"
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A
16"	5"	Short	14.5"	22"	4"	Tall	16.5"
16"	5"	Tall	16"	22"	5"	Short	N/A
16"	6"	Short	16"	22"	5"	Tall	16.5-17"
16"	6"	Tall	17.5"	22"	6"	Short	N/A
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"
16"	7"	Tall	18"	22"	7"	Short	N/A
18"	3"	Short	12.5-13"	22"	7"	Tall	18"
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A
18"	5"	Tall	16-17"	24"	5"	Short	N/A
18"	6"	Short	15.5"	24"	5"	Tall	N/A
18"	6"	Tall	17"	24"	6"	Short	N/A
18"	7"	Short	N/A	24"	6"	Tall	17.5"
18"	7"	Tall	N/A	24"	7"	Short	N/A
20"	3"	Short	N/A	24"	7"	Tall	18"