



Our kids are our passion.™

NXT-SE

Generation NEXT

Pediatric Manual Wheelchair



Yes, you can.®

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: _____ Buyer: _____
Company: _____
Marked For: _____
Email: _____

Quote Order PO# _____
Location: _____
Phone: _____
Fax: _____

Description	Part No.	Retail	Description	Part No.	Retail
NXT Pediatric Frame Options			Vent & Battery Trays		
NXT-SE Pediatric Folding Transport Optional (E1234) NXT-30304		\$2650	Note: All vent and battery trays not with transport.		
Note: Transportation Option selected separately below.			Note: All vent and battery trays not with short frames.		
Note: SE frames <u>Do Not</u> include Freedom Future Fit Program.			Note: All vent trays on NXT only with Dual Lever Hand Activated Tilt.		
NXT-SE Pediatric Rigid Transport Optional (E1233) NXT-30314		\$2525	Note: All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.		
Note: Transportation Option selected separately below.			Vent Tray - Fixed (E1029) 11-7000		\$ 439
Note: SE frames <u>Do Not</u> include Freedom Future Fit Program.			Battery Tray - Fixed (K0108) 11-8000		\$ 429
Transportation Option			Quick Release Tray Hardware (E1028) 563		\$ 272
Note: Used when Transport Optional frame is selected above.			Note: For LP10 style vent and battery tray on frames 16" wide and greater.		
Transportation Package - SE NXT-99997		\$ 300	Note: Vent tray is 4" narrower than frame width and 16" deep.		
Note: Must be installed by Freedom Designs manufacturing.			Note: Battery tray is 4" narrower than frame width and 10" deep.		
Note: Not available with all vent and battery trays.			Frame Colors		
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.			Sugar Cookie	n/c	Pacific Blue n/c
Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.			Neon Yellow	n/c	Blue Blast n/c
Seat Rail Options			Yellow	n/c	Electric Blue n/c
Note: Not available with ALL Contracture, 2D & 3D Dynamic footrest systems.			Bright Pink	n/c	Galaxy Blue n/c
NXT LEAF Seat Frame- Pair (K0108) NXT-30324		\$ 714	Pink Parfait	n/c	Passion Purple n/c
NXT LEAF Seat Frame- Left (K0108) NXT-30324-L		\$ 357	Sugar Plum	n/c	Octane Orange n/c
NXT LEAF Seat Frame- Right (K0108) NXT-30324-R		\$ 357	Candy Violet	n/c	Red Sunset n/c
Note: LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.			Candy Raspberry	n/c	Really Red n/c
Note: Freedom customized seats can be ordered to fit pre-configured LEAF frame models.			Extreme Green	n/c	Fiery Red n/c
Note: No seating system is provided with NXT frames. A seating system must be purchased separately.			Sparkle Green	n/c	Twilight n/c
NXT SRO Seat Frame- Pair (K0108) NXT-30325		\$ 306	Mystic Teal	n/c	Stormy Gray n/c
NXT SRO Seat Frame- Left (K0108) NXT-30325-L		\$ 153			Black n/c
NXT SRO Seat Frame- Right (K0108) NXT-30325-R		\$ 153	Frame Dimensions		
Note: SRO is manufactured without standard footrest hanger attachments.			Specify Frame Width Setting		
Note: SRO is used when customer attaches leg/footrest(s) directly to the underside of the solid seat.			10"	11"	12" Adjusts 10" to 12" NXT-30250 n/c
Tilt Activations			11"	12"	13" Adjusts 11" to 13" NXT-30210 n/c
Dual Lever Hand Activated Tilt NXT-30272		n/c	12"	13"	14" Adjusts 12" to 14" NXT-30220 n/c
Note: Must be selected for ALL vent and battery trays.			13"	14"	15" Adjusts 13" to 15" NXT-30230 n/c
Foot Activated Tilt NXT-30273		n/c	14"	15"	16" Adjusts 14" to 16" NXT-30240 n/c
Note: Foot pedal installed on right rear side frame.			15"	16"	17" Adjusts 15" to 17" NXT-30290 n/c
Note: Not with ALL vent and battery trays.			16"	17"	18" Adjusts 16" to 18" NXT-30260 n/c
Combination Dual Hand Activated Tilt, Brake and Wheel Lock System NXT-30277		n/c	17"	18"	Adjusts 17" to 18" NXT-30261 n/c
Note: Includes manual wheel locks to be used with public transportation.			18"		Not width adjustable NXT-30262 n/c
Note: Not available with 12" mag wheels with low pro poly tires.			19"		Not adjustable, not w/ short frame NXT-30263 n/c
			20"		Not adjustable, not w/ short frame (E2201) NXT-30264 \$ 400
			Note: Frame widths measured across and from outside of seat frame tubes.		
			Note: Seat frame widths of 19" and 20" not available with short frame.		
			Note: Vent and battery trays only with 14" and wider frames.		

NXT Rider Weight Capacity is 225 lbs.

Page 1 of 7

NXT-SE Pediatric

CE

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

Back Post Locations
Specify Frame Type and Frame Length

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

Short Frame & Short Tilt Bracket	NXT-30511	n/c
12"	NXT-30504	n/c
13"	NXT-30505	n/c
14"	NXT-30506	n/c
15"	NXT-30521	n/c
16"	NXT-30522	n/c
17" Must select 2" extension	NXT-30523	n/c
18" Must select 2" extension	NXT-30524	n/c
2" Extension - Pair (K1008)	NXT-30525	\$ 26

Specify Back Post Location (if other than specified above) _____ "

- Note:** Floor to seat frame height range is 12.5 to 16.5".
- Note:** Seat frame widths of 19" and 20" not available with short frame.
- Note:** Must choose an extension for 17" and 18" frame lengths.
- Note:** 22" & 24" wheels not available with short frame and short tilt bracket.
- Note:** All vent and battery trays not available with short frame.

Medium Frame & Short Tilt Bracket	NXT-30512	n/c
14"	NXT-30506	n/c
15"	NXT-30521	n/c
16"	NXT-30522	n/c
17"	NXT-30523	n/c
18"	NXT-30524	n/c
19" Must select 2" extension	NXT-30527	n/c
20" Must select 2" extension (E2203)	NXT-30528	\$ 500
21" Must select 3" extension (E2203)	NXT-30529	\$ 500
2" Extension - Pair (K1008)	NXT-30525	\$ 26
3" Extension - Pair (K1008)	NXT-30530	\$ 26

Specify Back Post Location (if other than specified above) _____ "

- Note:** Floor to seat frame height range is 12.5" to 16.5".
- Note:** Must choose an extension for 19", 20", and 21" frame lengths.
- Note:** 22" & 24" wheels not available with medium frame and short tilt bracket.

Medium Frame & Tall Tilt Bracket	NXT-30540	n/c
14"	NXT-30506TP	n/c
15"	NXT-30521TP	n/c
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
19" Must select 2" extension	NXT-30527TP	n/c
20" Must select 2" extension (E2203)	NXT-30528TP	\$ 500
21" Must select 3" extension (E2203)	NXT-30529TP	\$ 500
22" Must select 4" extension (E2204)	NXT-30535TP	\$ 900
2" Extension - Pair (K1008)	NXT-30525	\$ 26
3" Extension - Pair (K1008)	NXT-30530	\$ 26
4" Extension - Pair (K1008)	NXT-30531	\$ 26

Specify Back Post Location (if other than specified above) _____ "

- Note:** Tall tilt bracket increases floor to seat rail height 1.5".
- Note:** Floor to seat frame height range is 14" to 18".
- Note:** Must choose an extension for 19", 20", 21", and 22" frame lengths.

Back Post Locations (continued)
Specify Frame Type and Frame Length

Pediatric Long Frame & Tall Tilt	(E1298) NXT-30520P	\$ 153
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
19"	NXT-30527TP	n/c
20"	(E2203) NXT-30528TP	\$ 500
21"	(E2203) NXT-30529TP	\$ 500
22"	(E2204) NXT-30535TP	\$ 900
23" Must select 2" extension	(E2204) NXT-30537TP	\$ 900
24" Must select 2" extension	(E2204) NXT-30538TP	\$ 900
25" Must select 3" extension	(E2204) NXT-30542TP	\$ 900
26" Must select 4" extension	(E2204) NXT-30543TP	\$ 900
2" Extension - Pair (K1008)	NXT-30525	\$ 26
3" Extension - Pair (K1008)	NXT-30530	\$ 26
4" Extension - Pair (K1008)	NXT-30531	\$ 26

- Note:** Tall tilt bracket increases floor to seat rail height 1.5".
- Note:** Floor to seat frame height range is 14" to 18".
- Note:** Must choose extensions for 23", 24", 25" and 26" frame lengths.

Specify Back Post Location (if other than specified above) _____ "

5° Anterior Tilt Seat Frame	NXT-30276	n/c
------------------------------------	-----------	-----

- Note:** Limits posterior tilt to 40°, includes split collars for adjustment.
- Note:** Medium with Tall Tilt and Long frames only.
- Note:** NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14".

Back Post Tubes

Specify Push Handle, Back Post and Extensions		
Standard 90° Bend Push Handles and Back Post Tubes	NXT-30518	n/c
Ergo Bend Push Handles and Back Post Tubes	NXT-30519	n/c
Seating Dynamics Rocker Back	(K0108) NXT-30532	\$1734

- Note:** All push/stroller handle extensions not recommended with Ergo Bend.
- Note:** Angle adjustable from 72° to 108° in 6° increments.
- Note:** Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.
- Note:** Locks into a non-dynamic position when necessary.
- Note:** Does not allow the back post tubes to fold forward.

Back Post Height

Note: Back Posts are measured from top of seat rail to top of handle bend.

Short Push Handle - 20" high	NXT-30515	n/c
Medium Push Handle - 23" high	NXT-30514	n/c
Tall Push Handle - 26" high	NXT-30516	n/c

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

Stroller Handle Extensions

Note: All stroller handle ext. add approximately 6.5" to push handle height.
Note: All stroller/push handle extensions except individual include a black neoprene wrap around pad.



Push Handle Extension - Width Adjustable NXT-30560 n/c

Note: Not removable, but allows frame to fold.
Note: Mounts direct to back post push handle.

Push Button, Angle-Adjustable Folding Stroller Handle Extension NXT-30570 \$ 112

Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.
Note: Not removable, but allows frame to fold.
Note: Mounts direct to back post push handle.

Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty NXT-30575 \$ 179

Note: Available 10W-22W wheelchair frames.
Note: Not removable, but allows frame to fold.
Note: Same as above but with metal hinge mechanism.
Note: Mounts direct to back post push handle.

Push Button, Angle-Adjustable, Folding Push Handle Extension NXT-30580 \$ 112

Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.
Note: Not removable, but allows frame to fold.
Note: Mounts direct to straight back post tubes.

Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty NXT-30585 \$ 179

Note: Available 10W-22W wheelchair frames.
Note: Not removable, but allows frame to fold.
Note: Same as above but with metal hinge mechanism.
Note: Mounts direct to straight back post tubes.

Individual Stroller Handle Extensions Removable - Pair NXT-30517 \$ 102

Legrests, Footrests and Extensions Legrests

Omit Legrest Hangers & Footrests NXT-33000

75° Pop-Up Swing Away Hangers NXT-30280 n/c
Note: Not available with high mount bracket.

90° Pop-Up Swing Away Hangers NXT-30281 n/c
Note: Not available with 8" adjustable extensions.

90° Pop-Up Swing Away Elevating Legrests (ELRs) (E0990) NXT-30279 \$ 255
Note: Not available with all length adjustable extensions.
Note: Not available with 6.5" heavy duty extensions.

75° Latch Style Swing Away Hangers NXT-30282 n/c
Note: Not available with high mount bracket.

90° Latch Style Swing Away Hangers NXT-30283 n/c
Note: Not available with 8" adjustable extensions.

Legrests, Footrests and Extensions(continued) Footplates & Footboards

Omit Footrest Plates NXT-36404

Composite Flip-Up NXT-17070 n/c

Kids Aluminum Angle Adjustable, Flip-up (K0040) NXT-36020 \$ 160
Note: For frames 14" wide and narrower.

Heel Loops (E0951) NXT-36022 \$ 41

Calf Strap - Webbing (K0038) NXT-36024 \$ 41

Adult Aluminum Angle Adj. Flip-up w/ Heel Loops (K0040/ E0951) NXT-36021 \$ 180
Note: For frames 15" wide and greater.

Omit Heel Loops NXT-36023 n/c

Aluminum Angle Adjustable Footboard Removable (K0040) NXT-36029 \$ 160

MAPS Footplate Mounting Hardware

MAPS (Multi Axis Positioning System) - Pair NXT-36040 \$ 306

MAPS (Multi Axis Positioning System) - Left NXT-36040-L \$ 153

MAPS (Multi Axis Positioning System) - Right NXT-36040-R \$ 153

Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.

Note: Must choose aluminum footplates.

Note: Not available on 10" & 11" wide frames.

Note: MAPS not available with composite footplates or aluminum footboard.

Legrest Extensions

4" Short Length Adjustable Extensions NXT-30285 n/c

6" Medium Length Adjustable Extensions NXT-30286 n/c

8" Long Length Adjustable Extensions NXT-30287 n/c

Note: Not available with all 90° hangers.

4.5" Heavy Duty ELR Extensions NXT-36401 \$ 56

4.5" Heavy Duty Extensions NXT-36402 \$ 56
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

6.5" Heavy Duty Extensions NXT-36403 \$ 66
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

High Mount Brackets NXT-36407 \$ 102
Note: Includes flip-up mounting hardware for future growth adjustment.
Note: Not available with all 75° hangers.

Specify Knee to Heel Measurement _____" (see matrix page)

Continue To Next Page

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
<h3>Legrests, Footrests and Extensions</h3> <h4>Specialized Footrest Systems</h4> <p>Note: ALL legrest measurements are from the seat rail to footplate. Note: Not available with LEAF or SRO seat rail options.</p>			<h3>Armrests and Armrest Pads (continued)</h3>		
<h4>Super Contracture Footrest Hanger (K0108)</h4> <p>Short Pair - 11" to 15" NXT-30291-2 \$3570 Long Pair - 13" to 18" NXT-30292-2 \$3570</p> <p>Note: Super Contracture for frames 15" wide and greater. Note: Super Contracture is both swing-away and removable. Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments. Note: Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar / dorsiflexion positioning. Note: Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment. Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference. Note: Super Contracture position may limit ankle joint/ footplate adjustment range.</p>			<p>10.5" to 13.5" Long Flip-Back, Telescoping (E0973) NXT-38208 \$ 240 Note: Full length pad on 13.5" or longer armrests only. 10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973) NXT-38210 \$ 265 Note: Full length pad on 13.5" or longer armrests only. 17" to 20" Long Flip-Back Locking, Telescoping (E0973) NXT-38215 \$ 265</p>		
<h4>2D Dynamic Footrest Hanger (K0108)</h4> <p>Short Pair - 14" to 16" NXT-30293-2 \$4080 Long Pair - 16" to 18" NXT-30294-2 \$4080</p> <p>Note: 2D Dynamic for frames 15" wide and greater. Note: 2D Dynamic is both swing-away and removable. Note: 2D Dynamic has adjustable knee joints that have 48 positions that adjust in 7.5° increments. Note: 2D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel. Note: 2D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning. Note: 2D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment. Note: Frame and/or front caster choice may limit 2D Dynamic adjustment range and/or cause footplate/front caster interference.</p>			<p>Foam Wraps - Standard NXT-38216 n/c Desk Length Pads NXT-38213 n/c Full Length Pads NXT-38214 n/c Note: 8.5" armrests not available with desk and full length armrest pads. Note: Full length pad on 13.5" or longer armrests only.</p>		
<h4>3D Dynamic Footrest Hanger (K0108)</h4> <p>Short Pair - 13" to 15" NXT-30295-2 \$4590 Long Pair - 15" to 17" NXT-30296-2 \$4590</p> <p>Note: 3D Dynamic for frames 15" wide and greater. Note: 3D Dynamic is both swing-away and removable. Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel. Note: 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel. Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning. Note: 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment. Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.</p>			<h4>Frame Mounted Armrests</h4> <p>Single Post Fixed Height, Removable - Short (K0015) NXT-38222 \$ 225 Single Post Fixed Height, Removable - Tall (K0015) NXT-38224 \$ 225</p> <p>Note: Single Post Fixed Height arms include a removable side guard. Note: Adjustable in height with tools only. See matrix for height ranges.</p> <p>Single Post Height Adjustable, Removable - Short (E0973) NXT-38217 \$ 265 Single Post Height Adjustable, Removable - Tall (E0973) NXT-38203 \$ 265</p> <p>Dual Post Fixed Height Removable NXT-38205 n/c Note: All dual post armrests not available with short frame. Note: All dual post armrests not available with 22" or 24" wheels. Dual Post Flip-Back Height Adjustable, Removable (E0973) NXT-38206 \$ 275 Note: All dual post armrests not available with short frame. Note: All dual post armrests not available with 22" or 24" wheels.</p>		
<h4>Armrests and Armrest Pads</h4> <p>Omit Armrests NXT-38209</p> <p>Back Post Mounted Armrests</p> <p>Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.</p> <p>10.5" Long Flip-Back NXT-38204 n/c Note: Full length pad on 13.5" or longer armrests only. 13.5" Long Flip-Back NXT-38202 n/c Note: Full length pad on 13.5" or longer armrests only.</p> <p>10.5" Long Flip-Back Locking (E0973) NXT-38201 \$ 245 Note: Full length pad on 13.5" or longer armrests only. 13.5" Long Flip-Back Locking (E0973) NXT-38200 \$ 245 Note: Full length pad on 13.5" or longer armrests only.</p>			<p>Desk Length Pads NXT-38213 n/c Full Length Pads NXT-38214 n/c</p> <p>Specify Seat Rail to Top of Armrest Measurement _____ "</p> <h4>Front Wheels - Casters</h4> <p>Omit Caster Wheels NXT-34333</p> <p>3" Polyurethane NXT-34888 \$ 56 3" Polyurethane Light-Up NXT-34111 \$ 56 Note: 3" casters recommended for low seat to floor heights.</p> <p>4" Polyurethane NXT-40000 \$ 66 4" Polyurethane Light-Up NXT-41000 \$ 101 4" x 1.4" Poly, Aluminum Mag NXT-43000 \$ 92</p> <p>5" x 1" Poly, Plastic Mag NXT-34555 n/c 5" Polyurethane Light-Up NXT-35050 \$ 101 5" x 1.5" Poly, Plastic Mag NXT-34222 \$ 87 5" x 1.4" Poly, Aluminum Mag NXT-42000 \$ 92</p> <p>6" x 1" Poly, Plastic Mag NXT-34666 n/c 6" x 1.5" Poly, Plastic Mag NXT-34999 \$ 87 6" x 1.4" Poly, Aluminum Mag NXT-44000 \$ 92 6" x 2" Poly with Tread NXT-34444 \$ 92 6" x 1" Pneumatic Plastic Mag (E2214/E2215) NXT-34667 \$ 77</p> <p>7" x 1" Poly, Plastic Mag NXT-34777 \$ 77 7" x 2" Polyurethane NXT-34200 \$ 118 Note: All 7" casters not available with all 12" rear wheels.</p>		

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
Rear Wheels - Tires - Accessories		
Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.		
Omit Rear Wheels and Axles	NXT-31524	
Mag Wheels		
12" Low Profile Poly Mag (E2218)	NXT-31526	\$ 36
Note: Not available with any foot/hand activated brake/wheel lock system.		
Note: Not available with airless inserts.		
Note: All 12" wheels not available with all 7" casters.		
12" Semi Pneumatic Plastic Mag (K0069)	NXT-31523	\$ 168
Note: Not available with airless inserts.		
Note: All 12" wheels not available with all 7" casters.		
12" Pneumatic Plastic Mag	NXT-31525	n/c
12" Airless Inserts (E2213)	NXT-31019	\$ 77
Note: All 12" wheels not available with all handrims.		
Note: All 12" wheels not available with all 7" casters.		
16" Pneumatic Plastic Mag	NXT-31530	\$ 41
16" Airless Inserts (E2213)	NXT-31020	\$ 77
16" Aluminum Handrims	NXT-31014	\$ 169
Note: 16" mag wheels not available with plastic coated handrims.		
18" Pneumatic Plastic Mag	NXT-31533	\$ 41
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Aluminum Handrims	NXT-31015	\$ 169
Note: 18" mag wheels not available with plastic coated handrims.		
20" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31527	n/c
20" Airless Inserts (E2213)	NXT-31022	\$ 77
20" Plastic Coated Handrims	NXT-31016	\$ 119
22" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31528	n/c
22" Airless Inserts (E2213)	NXT-31023	\$ 77
22" Plastic Coated Handrims	NXT-31017	\$ 119
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
24" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31529	n/c
24" Airless Inserts (E2213)	NXT-31024	\$ 77
24" Plastic Coated Handrims	NXT-31018	\$ 119
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	
Spoke Wheels		
16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	NXT-31534	n/c
Note: Rim, hub & spokes are anodized silver.		
16" Airless Inserts (E2213)	NXT-31020	\$ 77
16" Silver Aluminum Handrim	NXT-31030	\$ 77
16" Black Plastic Coated Handrim	NXT-31032	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
16" Black Pneumatic Tire & Black Spoke Aluminum Wheel	NXT-31535	n/c
Note: Rim, hub & spokes are anodized black.		
16" Airless Inserts (E2213)	NXT-31020	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
16" Black Aluminum Handrim	NXT-31031	\$ 77
16" Black Plastic Coated Handrim	NXT-31032	\$ 77

Description	Part No.	Retail
Rear Wheels - Tires - Accessories (continued)		
Note: 16" & 18" spoke wheels do not include handrims.		
18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	NXT-31537	n/c
Note: Rim, hub & spokes are anodized silver.		
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Silver Aluminum Handrim	NXT-31034	\$ 77
18" Black Plastic Coated Handrim	NXT-31036	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
18" Black Pneumatic Tire & Black Spoke Aluminum Wheel	NXT-31538	n/c
Note: Rim, hub & spokes are anodized black.		
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Black Aluminum Handrim	NXT-31035	\$ 77
18" Black Plastic Coated Handrim	NXT-31036	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
20" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31540	\$ 204
20" Airless Inserts (E2213)	NXT-31022	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
22" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31541	\$ 204
22" Airless Inserts (E2213)	NXT-31023	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
24" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31542	\$ 204
24" Airless Inserts (E2213)	NXT-31024	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
Wheel Locks			Accessories		
Push to Lock - Standard	NXT-31501	n/c	Medical Necessities Pouch	NXT-23207	\$ 68
Pull to Lock	NXT-31506	n/c	Note: Black nylon.		
Extension Handles (E0961)	NXT-31801	\$ 61	Footrest Hanger Pads - Pair	SFDI-9HP	\$ 92
			Note: Wrap-around pads made from black neoprene.		
Foot Activated Wheel Lock (E2228)	NXT-31508	n/c	Frame Pads - Pair	10911-2	\$ 82
Note: Includes manual wheel locks to be used with public transportation.			Note: Wrap-around pads made from black neoprene.		
Note: Not available with 12" mag wheels with low pro poly tires.			Note: Typically used to pad exposed back post tubes.		
Note: Not available with ALL vent and battery trays.			Flip-Down Rigidizing Bar	NXT-23210	\$ 214
Hand Activated Brake and Wheel Lock System (E2228)	NXT-30275	n/c	Note: Mounts to back post tubes to add rigidity to frame.		
Note: Includes manual wheel locks to be used with public transportation.			Note: Does NOT adjust in width to grow with frame.		
Note: Not available with 12" mag wheels with low pro poly tires.			Freedom Backpack	NXT-10900	\$ 56
Note: Used to control direction, speed, brake and lock wheels from strolling position.			Note: Black Nylon only with Freedom Designs logo.		
Note: Not available with all vent and battery trays.			Note: Monogramming is not available on Freedom Backpack.		
Combination Dual Hand Activated Tilt, Brake and Wheel Lock System (E2228)	NXT-30277	n/c	Freedom Canopy - Flip-Back	NXT-10905	\$ 128
Note: Includes manual wheel locks to be used with public transportation.			Note: Canopy assembly includes removable mounting hardware.		
Note: Not available with 12" mag wheels with low pro poly tires.			Note: Canopy is available in black only with Freedom Designs logo.		
Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.			Note: Canopy includes rear flap with zippered storage pocket.		
Note: Not available with all vent and battery trays.			Note: Canopy includes clear-view window.		
Rear Anti-Tippers			Anterior Support Surfaces – (Trays)		
Curb Assist Anti-Tip Assembly Flip-Up, Adjustable (E0971)	NXT-31600	\$ 97	3/8" Lexan Tray - Channel Bracket Hardware (E0950)	08-2000	\$ 271
Note: Pair includes one anti-tip tube with curb assist feature.			3/8" Lexan Tray - Adjustable Knob Hardware & Ribs (E0950)	08-6000	\$ 219
Standard Anti-Tip Assembly Flip-Up, Adjustable (E0971)	NXT-31603	\$ 97	1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs (E0950)	08-13000	\$ 183
			1/4" Lexan Tray - Hook & Loop Straps (E0950)	08-16000	\$ 106
Seat Pan Options					
Fixed Solid Seat Pan - Bolteron (E2231)	NXT-39010	\$ 150			
Note: Not available with folding frame.					
Fixed Solid Seat Pan - Aluminum (E2231)	NXT-39030	\$ 170			
Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".					
Note: Not available with folding frame.					
Removable Solid Seat Pan - Bolteron (E2231)	NXT-39012	\$ 170			
Removable Solid Seat Pan - Aluminum (E2231)	NXT-39032	\$ 255			
Note: Removable Seat Pan raises the effective seat to floor height by 5/8".					
Note: Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.					
Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.					
Note: All seat pans include 2" x 2" notches for growth adjustment.					
Note: All seat pans made to fit frame dimensions unless specified below.					
Seat Width _____ " X Seat Depth _____ "					

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

NXT Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Choose NXT Seat to Floor Height _____"

Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height.

Note: Caster sizes vary slightly in height by type.

Note: Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".

NXT MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.
Note: Without extensions bottom hanger hole not accessible for hardware.
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

NXT Wheel / Caster Seat to Floor Heights

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5-13"	20"	3"	Tall	N/A
12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"
12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"
12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"
12"	5"	Short	15.5"	20"	5"	Tall	16-17"
12"	5"	Tall	17"	20"	6"	Short	15.5-16"
12"	6"	Short	16"	20"	6"	Tall	17-17.5"
12"	6"	Tall	17.5"	20"	7"	Short	16.5"
16"	3"	Short	12.5-13"	20"	7"	Tall	18"
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A
16"	5"	Short	14.5"	22"	4"	Tall	16.5"
16"	5"	Tall	16"	22"	5"	Short	N/A
16"	6"	Short	16"	22"	5"	Tall	16.5-17"
16"	6"	Tall	17.5"	22"	6"	Short	N/A
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"
16"	7"	Tall	18"	22"	7"	Short	N/A
18"	3"	Short	12.5-13"	22"	7"	Tall	18"
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A
18"	5"	Tall	16-17"	24"	5"	Short	N/A
18"	6"	Short	15.5"	24"	5"	Tall	N/A
18"	6"	Tall	17"	24"	6"	Short	N/A
18"	7"	Short	N/A	24"	6"	Tall	17.5"
18"	7"	Tall	N/A	24"	7"	Short	N/A
20"	3"	Short	N/A	24"	7"	Tall	18"