



Our kids are our passion.™

NXT™

Generation NEXT

Adult Manual Wheelchair



Yes, you can.™

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: _____ Buyer: _____
Company: _____
Marked For: _____
Email: _____

Quote Order PO# _____
Location: _____
Phone: _____
Fax: _____

Description	Part No.	Retail
Transport Optional Base		
NXT Adult Folding <i>Transport Optional</i>	(E1161) NXT-30302	\$3025
Note: Transportation Option selected separately below. Note: Includes Freedom Future Fit Program.		
NXT Adult Rigid <i>Transport Optional</i>	(E1161) NXT-30312	\$3025
Note: Transportation Option selected separately below. Note: Includes Freedom Future Fit Program.		
Transportation Option		
Note: Used when Transport Optional frame is selected above.		
Transportation Package	NXT-99998	\$ 250
Note: Must be installed by Freedom Designs manufacturing. Note: Not available with all vent and battery trays. Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved. Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.		

Transport Included Base		
NXT Adult Folding <i>Transport Included</i>	(E1161) NXT-30303	\$3175
Note: Includes Freedom Future Fit Program. Note: Not available with all vent and battery trays.		
NXT Adult Rigid <i>Transport Included</i>	(E1161) NXT-30313	\$3175
Note: Includes Freedom Future Fit Program. Note: Not available with all vent and battery trays.		

Seat Rail Options			
Note: Not available with ALL Contracture, 2D & 3D Dynamic footrest systems.			
NXT LEAF Seat Frame- Pair	(K0108) NXT-30324	\$ 714	
NXT LEAF Seat Frame- Left	(K0108) NXT-30324-L	\$ 357	
NXT LEAF Seat Frame- Right	(K0108) NXT-30324-R	\$ 357	
Note: LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.			
Note: Freedom customized seats can be ordered to fit pre-configured LEAF frame models.			
Note: No seating system is provided with NXT frames. A seating system must be purchased separately.			
NXT SRO Seat Frame- Pair	(K0108) NXT-30325	\$ 306	
NXT SRO Seat Frame- Left	(K0108) NXT-30325-L	\$ 153	
NXT SRO Seat Frame- Right	(K0108) NXT-30325-R	\$ 153	
Note: SRO is manufactured without standard footrest hanger attachments.			
Note: SRO is used when customer attaches leg/footrest(s) directly to the underside of the solid seat.			

Tilt Activations			
Dual Lever Hand Activated Tilt	NXT-30272	n/c	
Foot Activated Tilt	NXT-30273	n/c	
Note: Foot pedal installed on right rear side frame. Note: Not with ALL vent and battery trays.			
Combination Dual Hand Activated Tilt, Brake and Wheel Lock System	NXT-30277	n/c	
Note: Includes manual wheel locks to be used with public transportation. Note: Not available with 12" mag wheels with low pro poly tires.			

Description	Part No.	Retail
Vent & Battery Trays		
Note: All vent and battery trays not with transport.		
Note: All vent trays on NXT only with Dual Lever Activated Tilt.		
Note: All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.		
Vent Tray - Fixed	(E1029) 11-7000	\$ 439
Battery Tray - Fixed	(K0108) 11-8000	\$ 429
Quick Release Tray Hardware	(E1028) 563	\$ 272
Note: For LP10 style vent and battery tray on frames 16" wide and greater.		
Note: Vent tray is 4" narrower than frame width and 16" deep.		
Note: Battery tray is 4" narrower than frame width and 10" deep.		

Frame Colors			
Sugar Cookie	n/c	Pacific Blue	n/c
Neon Yellow	n/c	Blue Blast	n/c
Yellow	n/c	Electric Blue	n/c
Bright Pink	n/c	Galaxy Blue	n/c
Pink Parfait	n/c	Passion Purple	n/c
Sugar Plum	n/c	Octane Orange	n/c
Candy Violet	n/c	Red Sunset	n/c
Candy Raspberry	n/c	Really Red	n/c
Extreme Green	n/c	Fiery Red	n/c
Sparkle Green	n/c	Twilight	n/c
Mystic Teal	n/c	Stormy Gray	n/c
		Black	n/c

Frame Dimensions					
Specify Frame Width Setting					
14"	15"	16"	Adjusts 14" to 16"	NXT-30240	n/c
15"	16"	17"	Adjusts 15" to 17"	NXT-30290	n/c
16"	17"	18"	Adjusts 16" to 18"	NXT-30260	n/c
17"	18"		Adjusts 17" to 18"	NXT-30261	n/c
18"			Not width adjustable	NXT-30262	n/c
19"				NXT-30263	n/c
20"				NXT-30264	\$ 400
Note: Frame widths measured across and from outside of seat frame tubes.					
Note: Vent and battery trays only with 14" and wider frames.					

NXT Rider Weight Capacity is 225 lbs.

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Back Post Locations
Specify Frame Type and Frame Length

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

Medium Frame & Short Tilt Bracket	NXT-30512	n/c
14"	NXT-30506	n/c
15"	NXT-30521	n/c
16"	NXT-30522	n/c
17"	NXT-30523	n/c
18"	NXT-30524	n/c
* 19" Must select 2" extension	NXT-30527	n/c
* 20" Must select 2" extension (E2203)	NXT-30528	\$ 500
* 21" Must select 3" extension (E2203)	NXT-30529	\$ 500
2" Extension - Pair (K1008)	NXT-30525	\$ 26
3" Extension - Pair (K1008)	NXT-30530	\$ 26

Specify Back Post Location (if other than specified above) _____ "

Note: Floor to seat frame height range is 12.5" to 16.5".

* **Note:** Must choose an extension for 19", 20", and 21" frame lengths.

Note: 22" & 24" wheels not available with medium frame and short tilt bracket.

Medium Frame & Tall Tilt Bracket	NXT-30540	n/c
14"	NXT-30506TP	n/c
15"	NXT-30521TP	n/c
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
* 19" Must select 2" extension	NXT-30527TP	n/c
* 20" Must select 2" extension (E2203)	NXT-30528TP	\$ 500
* 21" Must select 3" extension (E2203)	NXT-30529TP	\$ 500
* 22" Must select 4" extension (E2204)	NXT-30535TP	\$ 900
2" Extension - Pair (K1008)	NXT-30525	\$ 26
3" Extension - Pair (K1008)	NXT-30530	\$ 26
4" Extension - Pair (K1008)	NXT-30531	\$ 26

Specify Back Post Location (if other than specified above) _____ "

Note: Tall tilt bracket increases floor to seat rail height 1.5".

Note: Floor to seat frame height range is 14" to 18".

* **Note:** Must choose an extension for 19", 20", 21", and 22" frame lengths.

Long Frame & Tall Tilt Bracket	NXT-30520	n/c
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
19"	NXT-30527TP	n/c
20" (E2203)	NXT-30528TP	\$ 500
21" (E2203)	NXT-30529TP	\$ 500
22" (E2204)	NXT-30535TP	\$ 900
* 23" Must select 2" extension (E2204)	NXT-30537TP	\$ 900
* 24" Must select 2" extension (E2204)	NXT-30538TP	\$ 900
* 25" Must select 3" extension (E2204)	NXT-30542TP	\$ 900
* 26" Must select 4" extension (E2204)	NXT-30543TP	\$ 900
2" Extension - Pair (K1008)	NXT-30525	\$ 26
3" Extension - Pair (K1008)	NXT-30530	\$ 26
4" Extension - Pair (K1008)	NXT-30531	\$ 26

Note: Tall tilt bracket increases floor to seat rail height 1.5"

Note: Floor to seat frame height range is 14" to 18".

* **Note:** Must choose an extension for 23", 24", 25", and 26" frame lengths.

Specify Back Post Location (if other than specified above) _____ "

5° Anterior Tilt Seat Frame NXT-30276 n/c

Note: Limits posterior tilt to 40°, includes split collars for adjustment.

Note: Medium with Tall Tilt and Long frames only.

Note: NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14".

Notes:

Back Post Tubes

Specify Push Handle, Back Post and Extensions

Standard 90° Bend Push Handles and Back Post Tubes NXT-30518 n/c

Ergo Bend Push Handles and Back Post Tubes NXT-30519 n/c

Note: All push/stroller handle extensions not recommended with Ergo Bend.

Seating Dynamics Rocker Back (K0108) NXT-30532 \$1734

Note: Angle adjustable from 72° to 108° in 6° increments.

Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.

Note: Locks into a non-dynamic position when necessary.

Note: Does not allow the back post tubes to fold forward.

Back Post Height

Note: Back Posts are measured from top of seat rail to top of handle bend.

Short Push Handle - 20" high NXT-30515 n/c





Medium Push Handle - 23" high NXT-30514 n/c

Tall Push Handle - 26" high NXT-30516 n/c

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
Back Post Tubes (continued)		
Stroller Handle Extensions		
<p>Note: All stroller handle ext. add approximately 6.5" to push handle height. Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.</p>		
		
-30560	-30570/5	-30580/5
		-30517
Push Handle Extension - Width Adjustable	NXT-30560	n/c
<p>Note: Not removable, but allows frame to fold. Note: Mounts direct to back post push handle.</p>		
Push Button, Angle-Adjustable Folding Stroller Handle Extension	NXT-30570	\$ 112
<p>Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below. Note: Not removable, but allows frame to fold. Note: Mounts direct to back post push handle.</p>		
Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty	NXT-30575	\$ 179
<p>Note: Available 10W-22W wheelchair frames. Note: Not removable, but allows frame to fold. Note: Same as above but with metal hinge mechanism. Note: Mounts direct to back post push handle.</p>		
Push Button, Angle-Adjustable, Folding Push Handle Extension	NXT-30580	\$ 112
<p>Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below. Note: Not removable, but allows frame to fold. Note: Mounts direct to straight back post tubes.</p>		
Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty	NXT-30585	\$ 179
<p>Note: Available 10W-22W wheelchair frames. Note: Not removable, but allows frame to fold. Note: Same as above but with metal hinge mechanism. Note: Mounts direct to straight back post tubes.</p>		
Individual Stroller Handle Extensions Removable - Pair	NXT-30517	\$ 102
Legrests, Footrests and Extensions		
Legrests		
Omit Legrest Hangers & Footrests	NXT-33000	
75° Pop-Up Swing Away Hangers	NXT-30280	n/c
Note: Not available with high mount bracket.		
90° Pop-Up Swing Away Hangers	NXT-30281	n/c
Note: Not available with 8" adjustable extensions.		
90° Pop-Up Swing Away Elevating Legrests (ELRs)	(E0990) NXT-30279	\$ 255
<p>Note: Not available with all length adjustable extensions. Note: Not available with 6.5" heavy duty extensions.</p>		
75° Latch Style Swing Away Hangers	NXT-30282	n/c
Note: Not available with high mount bracket.		
90° Latch Style Swing Away Hangers	NXT-30283	n/c
Note: Not available with 8" adjustable extensions.		

Description	Part No.	Retail
Legrests, Footrests and Extensions(continued)		
Footplates & Footboards		
Omit Footrest Plates	NXT-36404	
Composite Flip-Up	NXT-17070	n/c
Kids Aluminum Angle Adjustable, Flip-up	(K0040) NXT-36020	\$ 160
Note: For frames 14" wide and narrower.		
Heel Loops	(E0951) NXT-36022	\$ 41
Calf Strap - Webbing	(K0038) NXT-36024	\$ 41
Adult Aluminum Angle Adj. Flip-up w/ Heel Loops	(K0040/ E0951) NXT-36021	\$ 180
Note: For frames 15" wide and greater.		
Omit Heel Loops	NXT-36023	n/c
Aluminum Angle Adjustable Footboard Removable	(K0040) NXT-36029	\$ 160
MAPS Footplate Mounting Hardware		
MAPS (Multi Axis Positioning System) - Pair	NXT-36040	\$ 306
MAPS (Multi Axis Positioning System) - Left	NXT-36040-L	\$ 153
MAPS (Multi Axis Positioning System) - Right	NXT-36040-R	\$ 153
<p>Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning. Note: Must choose aluminum footplates. Note: Not available on 10" & 11" wide frames. Note: MAPS not available with composite footplates or aluminum footboard.</p>		
Legrest Extensions		
4" Short Length Adjustable Extensions	NXT-30285	n/c
6" Medium Length Adjustable Extensions	NXT-30286	n/c
8" Long Length Adjustable Extensions	NXT-30287	n/c
Note: Not available with all 90° hangers.		
4.5" Heavy Duty ELR Extensions	NXT-36401	\$ 56
4.5" Heavy Duty Extensions	NXT-36402	\$ 56
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.		
6.5" Heavy Duty Extensions	NXT-36403	\$ 66
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.		
High Mount Brackets	NXT-36407	\$ 102
<p>Note: Includes flip-up mounting hardware for future growth adjustment. Note: Not available with all 75° hangers.</p>		
<p>Specify Knee to Heel Measurement _____" (see matrix page)</p>		
Continue To Next Page		

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
<h3>Legrests, Footrests and Extensions</h3> <h4>Specialized Footrest Systems</h4> <p>Note: ALL legrest measurements are from the seat rail to footplate. Note: Not available with LEAF or SRO seat rail options.</p>			<h3>Armrests and Armrest Pads (continued)</h3>		
<h4>Super Contracture Footrest Hanger (K0108)</h4> <p>Short Pair - 11" to 15" NXT-30291-2 \$3570 Long Pair - 13" to 18" NXT-30292-2 \$3570</p> <p>Note: Super Contracture for frames 15" wide and greater. Note: Super Contracture is both swing-away and removable. Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments. Note: Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar / dorsiflexion positioning. Note: Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment. Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference. Note: Super Contracture position may limit ankle joint/ footplate adjustment range.</p>			<p>10.5" to 13.5" Long Flip-Back, Telescoping (E0973) NXT-38208 \$ 240 Note: Full length pad on 13.5" or longer armrests only. 10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973) NXT-38210 \$ 265 Note: Full length pad on 13.5" or longer armrests only. 17" to 20" Long Flip-Back Locking, Telescoping (E0973) NXT-38215 \$ 265</p> <p>Foam Wraps - Standard NXT-38216 n/c Desk Length Pads NXT-38213 n/c Full Length Pads NXT-38214 n/c Note: 8.5" armrests not available with desk and full length armrest pads. Note: Full length pad on 13.5" or longer armrests only.</p>		
<h4>2D Dynamic Footrest Hanger (K0108)</h4> <p>Short Pair - 14" to 16" NXT-30293-2 \$4080 Long Pair - 16" to 18" NXT-30294-2 \$4080</p> <p>Note: 2D Dynamic for frames 15" wide and greater. Note: 2D Dynamic is both swing-away and removable. Note: 2D Dynamic has adjustable knee joints that have 48 positions that adjust in 7.5° increments. Note: 2D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel. Note: 2D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning. Note: 2D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment. Note: Frame and/or front caster choice may limit 2D Dynamic adjustment range and/or cause footplate/front caster interference.</p>			<h4>Frame Mounted Armrests</h4> <p>Single Post Fixed Height, Removable - Short (K0015) NXT-38222 \$ 225 Single Post Fixed Height, Removable - Tall (K0015) NXT-38224 \$ 225</p> <p>Note: Single Post Fixed Height arms include a removable side guard. Note: Adjustable in height with tools only. See matrix for height ranges.</p> <p>Single Post Height Adjustable, Removable - Short (E0973) NXT-38217 \$ 265 Single Post Height Adjustable, Removable - Tall (E0973) NXT-38203 \$ 265</p> <p>Dual Post Fixed Height Removable NXT-38205 n/c Note: All dual post armrests not available with 22" or 24" wheels. Dual Post Flip-Back Height Adjustable, Removable (E0973) NXT-38206 \$ 275 Note: All dual post armrests not available with 22" or 24" wheels.</p>		
<h4>3D Dynamic Footrest Hanger (K0108)</h4> <p>Short Pair - 13" to 15" NXT-30295-2 \$4590 Long Pair - 15" to 17" NXT-30296-2 \$4590</p> <p>Note: 3D Dynamic for frames 15" wide and greater. Note: 3D Dynamic is both swing-away and removable. Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel. Note: 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel. Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning. Note: 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment. Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.</p>			<p>Desk Length Pads NXT-38213 n/c Full Length Pads NXT-38214 n/c</p> <p>Specify Seat Rail to Top of Armrest Measurement _____ "</p>		
<h3>Armrests and Armrest Pads</h3> <p>Omit Armrests NXT-38209</p> <h4>Back Post Mounted Armrests</h4> <p>Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.</p> <p>10.5" Long Flip-Back NXT-38204 n/c Note: Full length pad on 13.5" or longer armrests only. 13.5" Long Flip-Back NXT-38202 n/c Note: Full length pad on 13.5" or longer armrests only.</p> <p>10.5" Long Flip-Back Locking (E0973) NXT-38201 \$ 245 Note: Full length pad on 13.5" or longer armrests only. 13.5" Long Flip-Back Locking (E0973) NXT-38200 \$ 245 Note: Full length pad on 13.5" or longer armrests only.</p>			<h4>Front Wheels - Casters</h4> <p>Omit Caster Wheels NXT-34333</p> <p>3" Polyurethane NXT-34888 \$ 56 3" Polyurethane Light-Up NXT-34111 \$ 56 Note: 3" casters recommended for low seat to floor heights.</p> <p>4" Polyurethane NXT-40000 \$ 66 4" Polyurethane Light-Up NXT-41000 \$ 101 4" x 1.4" Poly, Aluminum Mag NXT-43000 \$ 92</p> <p>5" x 1" Poly, Plastic Mag NXT-34555 n/c 5" Polyurethane Light-Up NXT-35050 \$ 101 5" x 1.5" Poly, Plastic Mag NXT-34222 \$ 87 5" x 1.4" Poly, Aluminum Mag NXT-42000 \$ 92</p> <p>6" x 1" Poly, Plastic Mag NXT-34666 n/c 6" x 1.5" Poly, Plastic Mag NXT-34999 \$ 87 6" x 1.4" Poly, Aluminum Mag NXT-44000 \$ 92 6" x 2" Poly with Tread NXT-34444 \$ 92 6" x 1" Pneumatic Plastic Mag (E2214/E2215) NXT-34667 \$ 77</p> <p>7" x 1" Poly, Plastic Mag NXT-34777 \$ 77 7" x 2" Polyurethane NXT-34200 \$ 118 Note: All 7" casters not available with all 12" rear wheels.</p>		

NXT Rider Weight Capacity is 225 lbs.

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
Rear Wheels - Tires - Accessories		
Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.		
Omit Rear Wheels and Axles	NXT-31524	
Mag Wheels		
12" Low Profile Poly Mag (E2218)	NXT-31526	\$ 36
Note: Not available with any foot/hand activated brake/wheel lock system.		
Note: Not available with airless inserts.		
Note: All 12" wheels not available with all 7" casters.		
12" Semi Pneumatic Plastic Mag (K0069)	NXT-31523	\$ 168
Note: Not available with airless inserts.		
Note: All 12" wheels not available with all 7" casters.		
12" Pneumatic Plastic Mag	NXT-31525	n/c
12" Airless Inserts (E2213)	NXT-31019	\$ 77
Note: All 12" wheels not available with all handrims.		
Note: All 12" wheels not available with all 7" casters.		
16" Pneumatic Plastic Mag	NXT-31530	\$ 41
16" Airless Inserts (E2213)	NXT-31020	\$ 77
16" Aluminum Handrims	NXT-31014	\$ 169
Note: 16" mag wheels not available with plastic coated handrims.		
18" Pneumatic Plastic Mag	NXT-31533	\$ 41
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Aluminum Handrims	NXT-31015	\$ 169
Note: 18" mag wheels not available with plastic coated handrims.		
20" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31527	n/c
20" Airless Inserts (E2213)	NXT-31022	\$ 77
20" Plastic Coated Handrims	NXT-31016	\$ 119
22" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31528	n/c
22" Airless Inserts (E2213)	NXT-31023	\$ 77
22" Plastic Coated Handrims	NXT-31017	\$ 119
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
24" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31529	n/c
24" Airless Inserts (E2213)	NXT-31024	\$ 77
24" Plastic Coated Handrims	NXT-31018	\$ 119
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	
Spoke Wheels		
16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	NXT-31534	n/c
Note: Rim, hub & spokes are anodized silver.		
16" Airless Inserts (E2213)	NXT-31020	\$ 77
16" Silver Aluminum Handrim	NXT-31030	\$ 77
16" Black Plastic Coated Handrim	NXT-31032	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
16" Black Pneumatic Tire & Black Spoke Aluminum Wheel	NXT-31535	n/c
Note: Rim, hub & spokes are anodized black.		
16" Airless Inserts (E2213)	NXT-31020	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
16" Black Aluminum Handrim	NXT-31031	\$ 77
16" Black Plastic Coated Handrim	NXT-31032	\$ 77

Description	Part No.	Retail
Rear Wheels - Tires - Accessories (continued)		
Note: 16" & 18" spoke wheels do not include handrims.		
18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	NXT-31537	n/c
Note: Rim, hub & spokes are anodized silver.		
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Silver Aluminum Handrim	NXT-31034	\$ 77
18" Black Plastic Coated Handrim	NXT-31036	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
18" Black Pneumatic Tire & Black Spoke Aluminum Wheel	NXT-31538	n/c
Note: Rim, hub & spokes are anodized black.		
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Black Aluminum Handrim	NXT-31035	\$ 77
18" Black Plastic Coated Handrim	NXT-31036	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
20" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31540	\$ 204
20" Airless Inserts (E2213)	NXT-31022	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
22" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31541	\$ 204
22" Airless Inserts (E2213)	NXT-31023	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
24" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31542	\$ 204
24" Airless Inserts (E2213)	NXT-31024	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
Wheel Locks			Accessories		
Push to Lock - Standard	NXT-31501	n/c	Medical Necessities Pouch	NXT-23207	\$ 68
Pull to Lock	NXT-31506	n/c	Note: Black nylon.		
Extension Handles (E0961)	NXT-31801	\$ 61	Footrest Hanger Pads - Pair	SFDI-9HP	\$ 92
Foot Activated Wheel Lock (E2228)	NXT-31508	n/c	Note: Wrap-around pads made from black neoprene.		
Note: Includes manual wheel locks to be used with public transportation.			Frame Pads - Pair	10911-2	\$ 82
Note: Not available with 12" mag wheels with low pro poly tires.			Note: Wrap-around pads made from black neoprene.		
Note: Not available with ALL vent and battery trays.			Note: Typically used to pad exposed back post tubes.		
Hand Activated Brake and Wheel Lock System (E2228)	NXT-30275	n/c	Flip-Down Rigidizing Bar	NXT-23210	\$ 214
Note: Includes manual wheel locks to be used with public transportation.			Note: Mounts to back post tubes to add rigidity to frame.		
Note: Not available with 12" mag wheels with low pro poly tires.			Note: Does NOT adjust in width to grow with frame.		
Note: Used to control direction, speed, brake and lock wheels from strolling position.			Freedom Backpack	NXT-10900	\$ 56
Combination Dual Hand Activated Tilt, Brake and Wheel Lock System (E2228)	NXT-30277	n/c	Note: Black Nylon only with Freedom Designs logo.		
Note: Includes manual wheel locks to be used with public transportation.			Note: Monogramming is not available on Freedom Backpack.		
Note: Not available with 12" mag wheels with low pro poly tires.			Freedom Canopy - Flip-Back	NXT-10905	\$ 128
Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.			Note: Canopy assembly includes removable mounting hardware.		
			Note: Canopy is available in black only with Freedom Designs logo.		
			Note: Canopy includes rear flap with zippered storage pocket.		
			Note: Canopy includes clear-view window.		
			Anterior Support Surfaces – (Trays)		
			3/8" Lexan Tray - Channel Bracket Hardware	(E0950) 08-2000	\$ 271
			3/8" Lexan Tray - Adjustable Knob Hardware & Ribs	(E0950) 08-6000	\$ 219
			1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs	(E0950) 08-13000	\$ 183
			1/4" Lexan Tray - Hook & Loop Straps	(E0950) 08-16000	\$ 106
Rear Anti-Tippers					
Curb Assist Anti-Tip Assembly Flip-Up, Adjustable	(E0971) NXT-31600	\$ 97			
Note: Pair includes one anti-tip tube with curb assist feature.					
Standard Anti-Tip Assembly Flip-Up, Adjustable	(E0971) NXT-31603	\$ 97			
Seat Pan Options					
Fixed Solid Seat Pan - Bolteron	(E2231) NXT-39010	\$ 150			
Note: Not available with folding frame.					
Fixed Solid Seat Pan - Aluminum	(E2231) NXT-39030	\$ 170			
Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".					
Note: Not available with folding frame.					
Removable Solid Seat Pan - Bolteron	(E2231) NXT-39012	\$ 170			
Removable Solid Seat Pan - Aluminum	(E2231) NXT-39032	\$ 255			
Note: Removable Seat Pan raises the effective seat to floor height by 5/8".					
Note: Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.					
Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.					
Note: All seat pans include 2" x 2" notches for growth adjustment.					
Note: All seat pans made to fit frame dimensions unless specified below.					
Seat Width _____" X	Seat Depth _____"				

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

NXT Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Choose NXT Seat to Floor Height _____"

Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height.

Note: Caster sizes vary slightly in height by type.

Note: Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".

NXT MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.
Note: Without extensions bottom hanger hole not accessible for hardware.
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

NXT Wheel / Caster Seat to Floor Heights

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5-13"	20"	3"	Tall	N/A
12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"
12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"
12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"
12"	5"	Short	15.5"	20"	5"	Tall	16-17"
12"	5"	Tall	17"	20"	6"	Short	15.5-16"
12"	6"	Short	16"	20"	6"	Tall	17-17.5"
12"	6"	Tall	17.5"	20"	7"	Short	16.5"
16"	3"	Short	12.5-13"	20"	7"	Tall	18"
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A
16"	5"	Short	14.5"	22"	4"	Tall	16.5"
16"	5"	Tall	16"	22"	5"	Short	N/A
16"	6"	Short	16"	22"	5"	Tall	16.5-17"
16"	6"	Tall	17.5"	22"	6"	Short	N/A
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"
16"	7"	Tall	18"	22"	7"	Short	N/A
18"	3"	Short	12.5-13"	22"	7"	Tall	18"
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A
18"	5"	Tall	16-17"	24"	5"	Short	N/A
18"	6"	Short	15.5"	24"	5"	Tall	N/A
18"	6"	Tall	17"	24"	6"	Short	N/A
18"	7"	Short	N/A	24"	6"	Tall	17.5"
18"	7"	Tall	N/A	24"	7"	Short	N/A
20"	3"	Short	N/A	24"	7"	Tall	18"