



Our kids are our passion.<sup>™</sup>  
2241 N. Madera Rd

**NXT**  
Generation NEXT

Simi Valley, CA 93065

# Adult Manual Wheelchair Transport Package Included

Phone: (800) 331-8551 Fax: (888) 582-1509

www.freedomdesigns.com



Date: \_\_\_\_\_ Buyer: \_\_\_\_\_  
Company: \_\_\_\_\_  
Marked For: \_\_\_\_\_  
Email: \_\_\_\_\_

Quote  Order PO# \_\_\_\_\_  
Location: \_\_\_\_\_  
Phone: \_\_\_\_\_  
Fax: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail			
<b>NXT Adult Frame Options</b>			<b>NXT Frame Colors</b>					
<input type="radio"/> <b>NXT Adult Folding Transport Included</b>	(E1161) NXT-30303	<b>\$3175</b>	<input type="radio"/> Sugar Cookie	n/c	<input type="radio"/> Pacific Blue	n/c		
<i>Note: Includes Freedom Future Fit Program.</i>			<input type="radio"/> Neon Yellow	n/c	<input type="radio"/> Blue Blast	n/c		
<input type="radio"/> <b>NXT Adult Rigid Transport Included</b>	(E1161) NXT-30313	<b>\$3075</b>	<input type="radio"/> Yellow	n/c	<input checked="" type="radio"/> Electric Blue	n/c		
<i>Note: Includes Freedom Future Fit Program.</i>			<input type="radio"/> Bright Pink	n/c	<input type="radio"/> Galaxy Blue	n/c		
<b>Transport Option</b>			<input type="radio"/> Pink Parfait	n/c	<input type="radio"/> Passion Purple	n/c		
<input checked="" type="radio"/> <b>Transport Package Included</b>	NXT-99999-1	n/c	<input type="radio"/> Sugar Plum	n/c	<input type="radio"/> Octane Orange	n/c		
<i>Note: Must be installed by Freedom Designs manufacturing.</i>			<input type="radio"/> Vigorous Violet	n/c	<input type="radio"/> Red Sunset	n/c		
<i>Note: Not as parts.</i>			<input type="radio"/> Candy Violet	n/c	<input type="radio"/> Really Red	n/c		
<i>Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.</i>			<input type="radio"/> Candy Raspberry	n/c	<input checked="" type="radio"/> Fiery Red	n/c		
<i>Note: Most Freedom seating, hardware, and accessories has been WC20 tested and approved for transit.</i>			<input type="radio"/> Extreme Green	n/c	<input type="radio"/> Twilight	n/c		
<input type="radio"/> <b>NXT LEAF Seat Frame - Pair</b>	(K0108) NXT-30324	<b>\$700</b>	<input type="radio"/> Sparkle Green	n/c	<input checked="" type="radio"/> Stormy Gray	n/c		
<input type="radio"/> <b>NXT LEAF Seat Frame - Left</b>	(K0108) NXT-30324-L	<b>\$350</b>	<input type="radio"/> Mystic Teal	n/c	<input type="radio"/> Black	n/c		
<input type="radio"/> <b>NXT LEAF Seat Frame - Right</b>	(K0108) NXT-30324-R	<b>\$350</b>	<input type="radio"/> Blizzard Blue	n/c				
<i>Note: LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.</i>			<b>Frame Dimensions</b>					
<i>Note: Freedom customized seats can be ordered to fit preconfigured LEAF frame models.</i>			<b>Specify Frame Width Setting</b>					
<i>Note: No seating system is provided with NXT frames. A seating system must be purchased separately.</i>			<input type="radio"/> 14"	<input type="radio"/> 15"	<input type="radio"/> 16"	Adjusts 14" to 16"	NXT-30240	n/c
<input type="radio"/> <b>NXT SRO - Seat Rail Only - Pair</b>	(K0108) NXT-30325	<b>\$300</b>	<input type="radio"/> 15"	<input type="radio"/> 16"	<input type="radio"/> 17"	Adjusts 15" to 17"	NXT-30290	n/c
<input type="radio"/> <b>NXT SRO - Seat Rail Only - Left</b>	(K0108) NXT-30325-L	<b>\$150</b>	<input type="radio"/> 16"	<input type="radio"/> 17"	<input type="radio"/> 18"	Adjusts 16" to 18"	NXT-30260	n/c
<input type="radio"/> <b>NXT SRO - Seat Rail Only - Right</b>	(K0108) NXT-30325-R	<b>\$150</b>	<input type="radio"/> 17"	<input type="radio"/> 18"		Adjusts 17" to 18"	NXT-30261	n/c
<i>Note: SRO is manufactured without standard footrest hanger attachments.</i>			<input type="radio"/> 18"			Not adjustable	NXT-30262	n/c
<i>Note: SRO used when customer attaches leg/footrest(s) directly to underside of seat.</i>			<input type="radio"/> 19"			Not adjustable	NXT-30263	n/c
<b>Tilt Activations</b>			<input type="radio"/> 20"			Not adjustable (E2201)	NXT-30264	<b>\$300</b>
<input type="radio"/> <b>Dual Lever Hand Activated Tilt</b>	NXT-30272	n/c	<i>Note: Frame widths measured across and from outside of seat frame tubes.</i>					
<input type="radio"/> <b>Foot Activated Tilt</b>	NXT-30273	n/c	<b>Back Post Locations</b>					
<input checked="" type="radio"/> <b>Combination Hand Activated Tilt, Activated Tilt and Wheel Lock System</b>	NXT-30277	n/c	<b>Specify Frame Type and Frame Length</b>					
<i>Note: Includes manual wheel locks to be used with public transportation.</i>			<i>Note: Extensions may limit tilt range with 22" and 24" wheels and may affect stability when in tilt positions. Adjusting seat rail forward compensates for tilt range limitations and instability.</i>					
<b>Vent &amp; Battery Trays</b>			<input type="radio"/> <b>Medium Frame &amp; Short Tilt Bracket</b>	NXT-30512	n/c			
<i>Note: Vent and battery trays are not available with transport option.</i>			<input type="radio"/> 14"	NXT-30506	n/c			
			<input type="radio"/> 15"	NXT-30521	n/c			
			<input type="radio"/> 16"	NXT-30522	n/c			
			<input type="radio"/> 17"	NXT-30523	n/c			
			<input type="radio"/> 18"	NXT-30524	n/c			
			* <input type="radio"/> 19"	Must select 2" extension		NXT-30527	n/c	
			* <input type="radio"/> 20"	Must select 2" extension	(E2203)	NXT-30528	<b>\$300</b>	
			* <input type="radio"/> 21"	Must select 3" extension	(E2203)	NXT-30529	<b>\$300</b>	
			<input type="radio"/> 2" extension - Pair	(K0108)	NXT-30525	<b>\$ 25</b>		
			<input type="radio"/> 3" extension - Pair	(K0108)	NXT-30530	<b>\$ 25</b>		
			<input type="radio"/> <b>Specify location</b> (if other than specified above) _____"					
			<i>Note: Floor to seat frame height range is 12.5" to 16.5".</i>					
			<i>* Note: Must choose extension for 19", 20" and 21" frame lengths.</i>					
			<i>Note: 22" &amp; 24" wheels not with medium frame and short tilt bracket.</i>					

(Back Post Locations continued on page 2)

Date: \_\_\_\_\_ P.O.# \_\_\_\_\_

Marked for: \_\_\_\_\_

Description	Part No.	Retail
<b>Back Post Locations</b> Specify Frame Type and Frame Length		
<input type="radio"/> <b>Medium Frame &amp; Tall Tilt Bracket</b>	NXT-30540	n/c
<input type="radio"/> 14"	NXT-30506TP	n/c
<input type="radio"/> 15"	NXT-30521TP	n/c
<input type="radio"/> 16"	NXT-30522TP	n/c
<input type="radio"/> 17"	NXT-30523TP	n/c
<input type="radio"/> 18"	NXT-30524TP	n/c
* <input type="radio"/> 19" Must select 2" extension	NXT-30527TP	n/c
* <input type="radio"/> 20" Must select 2" extension (E2203)	NXT-30528TP	\$300
* <input type="radio"/> 21" Must select 3" extension (E2203)	NXT-30529TP	\$300
* <input type="radio"/> 22" Must select 4" extension (E2204)	NXT-30535TP	\$500
<input type="radio"/> 2" extension - Pair (K0108)	NXT-30525	\$ 25
<input type="radio"/> 3" extension - Pair (K0108)	NXT-30530	\$ 25
<input type="radio"/> 4" extension - Pair (K0108)	NXT-30531	\$ 25
<input type="radio"/> <b>Specify location</b> (if other than specified above) _____"		
<b>Note:</b> Tall tilt bracket increases floor to seat frame height by 1.5".		
<b>Note:</b> Floor to seat frame height range is 14" to 18".		
<b>*Note:</b> Must choose extension for 19", 20", 21" and 22" frame lengths.		
<input type="radio"/> <b>Long Frame &amp; Tall Tilt Bracket</b>	NXT-30520	n/c
<input type="radio"/> 16"	NXT-30522TP	n/c
<input type="radio"/> 17"	NXT-30523TP	n/c
<input type="radio"/> 18"	NXT-30524TP	n/c
<input type="radio"/> 19"	NXT-30527TP	n/c
<input type="radio"/> 20" (E2203)	NXT-30528TP	\$300
* <input type="radio"/> 21" Must select 2" extension (E2203)	NXT-30529TP	\$300
* <input type="radio"/> 22" Must select 2" extension (E2204)	NXT-30535TP	\$500
* <input type="radio"/> 23" Must select 3" extension (E2204)	NXT-30537TP	\$500
* <input type="radio"/> 24" Must select 4" extension (E2204)	NXT-30538TP	\$500
<input type="radio"/> 2" extension - Pair (K0108)	NXT-30525	\$ 25
<input type="radio"/> 3" extension - Pair (K0108)	NXT-30530	\$ 25
<input type="radio"/> 4" extension - Pair (K0108)	NXT-30531	\$ 25
<input type="radio"/> <b>Specify location</b> (if other than specified above) _____"		
<b>Note:</b> Tall tilt bracket increases floor to seat frame height by 1.5".		
<b>Note:</b> Floor to seat frame height range is 14" to 18".		
<b>*Note:</b> Must choose extension for 21", 22", 23" and 24" frame lengths.		
<input type="radio"/> <b>5° Anterior Tilt Seat Frame</b> (K0108)	NXT-30276	n/c
<b>Note:</b> Limits posterior tilt to 40°, includes split collars for adjustment.		
<b>Note:</b> Tall tilt bracket only.		
NXT has 4" of depth adjustment for growth. Depth measured from front of back post to end of seat frame. For back post location add seating / hardware thickness (typically 2") to clients seat depth. Example: Seat depth of 12" + seat / hardware thickness of 2" = back post tube location set at 14".		
<b>Back Post Tubes</b> Specify Push Handle, Back Post Height, and Extensions		
<b>Note:</b> All stroller handle ext. add approximately 6.5" to push handle height.		
<b>Note:</b> Back posts measured from top of seat frame to top of handle bend.		
<input type="radio"/> <b>Standard 90° Bend Push Handles and Back Post Tubes</b>	NXT-30518	n/c
<input type="radio"/> <b>Ergo Bend Push Handles and Back Post Tubes</b>	NXT-30519	n/c
<b>Note:</b> All push / stroller handle ext. not recommended with Ergo Bend.		
<input type="radio"/> <b>Seating Dynamics Rocker Back</b> (K0108)	NXT-30532	\$1700
<b>Note:</b> Angle adjustable from 72° to 108° in 6° increments.		
<b>Note:</b> Supplied with a complete set of 4 different durometers of elastomers to allow spring adjustment in the field.		
<b>Note:</b> Locks into a non dynamic position when necessary.		
<b>Note:</b> Does not allow the back post tubes to fold forward.		
<input type="radio"/> <b>Flip Down Rigidizing Bar</b> (K0108)	NXT-23210	\$ 210
<b>Note:</b> Not width adjustable.		
<b>Back Post Height</b>		
<input type="radio"/> <b>Short Push Handle - 19" High</b>	NXT-30515	n/c
<input type="radio"/> <b>Medium Push Handle - 23" High</b>	NXT-30514	n/c
<input type="radio"/> <b>Tall Push Handle - 26" High</b>	NXT-30516	n/c

Description	Part No.	Retail
<b>Stroller Handle Extensions</b>		
<input type="radio"/> <b>Push Handle Extension Width Adjustable Removable</b>	NXT-30513	n/c
<input type="radio"/> <b>Stroller Handle Extension Width Adjustable Removable</b>	NXT-30526	\$125
<input type="radio"/> <b>Push Button Angle Adjustable Stroller Handle Extensions (PB)</b>	NXT-30536	\$175
<b>Note:</b> PB stroller handles not with frames less than 14" wide.		
<b>Note:</b> PB stroller handle not removable and not with folding frame.		
<input type="radio"/> <b>Individual Stroller Handle Extensions Removable - Pair</b>	NXT-30517	\$100
<b>Footrest Hangers and Extensions</b>		
<input type="radio"/> <b>Omit Footrest Hangers &amp; Footplates</b>	NXT-33000	(\$30)
<input type="radio"/> <b>75° Pop-Up Swing Away Hangers</b>	NXT-30280	n/c
<b>Note:</b> Not with high mount bracket.		
<input type="radio"/> <b>90° Pop-Up Swing Away Hangers</b>	NXT-30281	n/c
<b>Note:</b> Not with 8" adjustable extensions.		
<input type="radio"/> <b>75° Latch Style Swing Away Hangers</b>	NXT-30282	n/c
<b>Note:</b> Not with high mount bracket.		
<input type="radio"/> <b>90° Latch Style Swing Away Hangers</b>	NXT-30283	n/c
<b>Note:</b> Not with 8" adjustable extensions.		
<input type="radio"/> <b>4" Short Length Adjustable Extensions</b>	NXT-30285	n/c
<input type="radio"/> <b>6" Medium Length Adjustable Extensions</b>	NXT-30286	n/c
<input type="radio"/> <b>8" Long Length Adjustable Extensions</b>	NXT-30287	n/c
<b>Note:</b> Not with all 90° hangers.		
<input type="radio"/> <b>4.5" Heavy Duty Extensions</b>	NXT-36402	\$ 55
<b>Note:</b> Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.		
<input type="radio"/> <b>6.5" Heavy Duty Extensions</b>	NXT-36403	\$ 65
<b>Note:</b> Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.		
<input type="radio"/> <b>High Mount Brackets</b> (K0037)	NXT-36407	\$100
<b>Note:</b> Includes flip up mounting hardware for growth.		
<b>Note:</b> Not with 75° hangers.		
<input type="radio"/> <b>Pop-Up Swing Away Elevating Legrests (ELR)</b> (E0990)	NXT-30279	\$235
<b>Note:</b> Not with all length adjustable extensions.		
<b>Note:</b> Not with 6.5" heavy duty extensions.		
<input type="radio"/> <b>4.5" Heavy Duty Extensions (ELR)</b>	NXT-36401	\$ 55
<input type="radio"/> <b>Specify Knee to Heel Measurement</b> _____" (see chart below)		
(See Matrix on Page 5)		



NXT Rider Weight Capacity is 225 lbs.

NXT Adult Transport Included



Date: \_\_\_\_\_ P.O.# \_\_\_\_\_

Marked for: \_\_\_\_\_

Description	Part No.	Retail
<b>Foot Plates</b>		
<input type="radio"/> Omit Footrest Plates	NXT-36404	(\$15)
<input type="radio"/> Composite Flip Up	NXT-17070	n/c
<input type="radio"/> Kids Aluminum Angle Adj. Flip Up (K0040)	NXT-36020	\$150
Note: For frames 14" wide and under.		
<input type="radio"/> Heel Loops (E0951)	NXT-36022	\$ 50
<input type="radio"/> Adult Aluminum Angle Adj. Flip Up Footrest Plates w/Heel Loops (K0040 / E0951)	NXT-36021	\$150
Note: For frames 15" wide and greater.		
<input type="radio"/> Omit Heel Loops	NXT-36023	n/c
<input type="radio"/> Aluminum Adjustable Angle Footboard Removable	NXT-36029	\$150
<input type="radio"/> MAPS (Multi Axis Position System) - Pair	NXT-36040	\$300
<input type="radio"/> MAPS (Multi Axis Position System) - Left	NXT-36040-L	\$150
<input type="radio"/> MAPS (Multi Axis Position System) - Right	NXT-36040-R	\$150
Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.		
Note: Must choose aluminum footplates.		
Note: MAPS not with composite footrest plates or footboard.		
Note: MAPS not available on 10" & 11" wide frames.		
<b>Footrest Hangers and Extensions Specialized Footrest Systems</b>		
 <b>Super Contracture Footrest Hanger (K0108)</b>		
<input type="radio"/> Short Left - 11" to 15"	NXT-30291-L	\$1750
<input type="radio"/> Short Right - 11" to 15"	NXT-30291-R	\$1750
<input type="radio"/> Short Pair - 11" to 15"	NXT-30291-2	\$3500
<input type="radio"/> Long Left - 13" to 18"	NXT-30292-L	\$1750
<input type="radio"/> Long Right - 13" to 18"	NXT-30292-R	\$1750
<input type="radio"/> Long Pair - 13" to 18"	NXT-30292-2	\$3500
Note: SC for NXT frames 15" wide or wider.		
Note: SC is both swing away and removable.		
Note: SC Has adjustable knee joints that have 48 positions that adjust in 7.5° increments.		
Note: SC has short or long length adjustable extension tubes that have adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar/dorsiflexion positioning.		
Note: SC has 5"x 7" flip up aluminum foot plates that have multiple holes for location adjustment and secondary positioning accessory attachment.		
Note: NXT frame and/or front caster choices may limit SC knee joint adjustment range and/or cause foot plate/front caster interference.		
Note: SC position (in contracture position) may limit ankle joint/foot plate adjustment range.		
 <b>2D Dynamic Footrest System (K0108)</b>		
<input type="radio"/> Short Left - 14" to 16"	NXT-30293-L	\$2000
<input type="radio"/> Short Right - 14" to 16"	NXT-30293-R	\$2000
<input type="radio"/> Short Pair - 14" to 16"	NXT-30293-2	\$4000
<input type="radio"/> Long Left - 16" to 18"	NXT-30294-L	\$2000
<input type="radio"/> Long Right - 16" to 18"	NXT-30294-R	\$2000
<input type="radio"/> Long Pair - 16" to 18"	NXT-30294-2	\$4000
Note: 2D for NXT frames 15" wide or wider.		
Note: 2D System is both swing away and removable.		
Note: 2D has adjustable knee joints that have 13 positions that adjust in 7.5° increments ranging from 0° to 100°.		
Note: 2D has short or long length adjustable extension tubes with 1.5" of dynamic telescoping travel.		
Note: 2D has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.		
Note: 2D has 5"x 7" angle adjustable aluminum foot plates that have multiple holes for location adjustment and secondary positioning accessory attachment.		
Note: NXT frame and/or front caster choices may limit 2D knee adjustment range and/or cause foot plate/front caster interference.		

Description	Part No.	Retail
 <b>3D Dynamic Footrest System (K0108)</b>		
<input type="radio"/> Short Left - 13" to 15"	NXT-30295-L	\$2250
<input type="radio"/> Short Right 13" to 15"	NXT-30295-R	\$2250
<input type="radio"/> Short Pair - 13" to 15"	NXT-30295-2	\$4500
<input type="radio"/> Long Left - 15" to 17"	NXT-30296-L	\$2250
<input type="radio"/> Long Right - 15" to 17"	NXT-30296-R	\$2250
<input type="radio"/> Long Pair - 15" to 17"	NXT-30296-2	\$4500
Note: 3D for NXT frames 15" wide or wider.		
Note: 3D is both swing away and removable.		
Note: 3D has adjustable dynamic elevating knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.		
Note: 3D has short or long length adjustable extension tubes with 1.5" of dynamic telescoping travel.		
Note: 3D has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.		
Note: 3D has 5"x 7" angle adjustable aluminum foot plates that have multiple holes for location adjustment and secondary positioning accessory attachment.		
Note: NXT frame and/or front caster choices may limit 3D knee joint adjustment range and/or cause foot plate/caster interference.		
<b>Armrests and Pads</b>		
<input type="radio"/> Omit Flip Back Armrests	NXT-38209	(\$20)
<input type="radio"/> 10.5" Long Flip Back	NXT-38204	n/c
<input type="radio"/> 13.5" Long Flip Back	NXT-38202	n/c
<input type="radio"/> 10.5" to 13.5" Long Flip Back Telescoping (E0973)	NXT-38208	\$100
<input type="radio"/> 10.5" Long Flip Back Locking (E0973)	NXT-38201	\$235
<input type="radio"/> 13.5" Long Flip Back Locking (E0973)	NXT-38200	\$235
<input type="radio"/> 10.5" to 13.5" Long Flip Back Telescoping Locking (E0973)	NXT-38210	\$260
<input type="radio"/> 17" to 20" Long Flip Back Telescoping Locking (E0973)	NXT-38215	\$260
Note: Full length pad on 13.5" or longer armrests only.		
Note: Flip back armrest length measured from front of backpost to end of armrest with foam wrap installed.		
<input type="radio"/> Foam Wraps - Standard	NXT-38216	n/c
<input type="radio"/> Desk Length Pads	NXT-38213	n/c
<input type="radio"/> Full Length Pads	NXT-38214	n/c
Note: Full length on 13.5" or longer armrests only.		
<input type="radio"/> Single Post Height Adj. Removable Short (E0973)	NXT-38217	\$235
<input type="radio"/> Single Post Height Adj. Removable Tall (E0973)	NXT-38203	\$235
<input type="radio"/> Dual Post Fixed Height Removable	NXT-38205	n/c
<input type="radio"/> Dual Post Flip Back Height Adj. Removable (E0973)	NXT-38206	\$245
Note: All dual post armrests not with all 22" and 24" wheels.		
<input type="radio"/> Desk Length Pads	NXT-38213	n/c
<input type="radio"/> Full Length Pads	NXT-38214	n/c
Note: Full length on 13.5" or longer armrests only.		
<input type="radio"/> Specify Seat Rail to Top of Armrest Measurement _____"		
(See Matrix on Page 5)		

Date: \_\_\_\_\_ P.O.# \_\_\_\_\_ Marked for: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
<b>Front Tires - Caster Wheels</b>			<b>Rear Tires - Wheels - Accessories</b>		
<input type="radio"/> Omit Caster Wheels	NXT-34333	n/c	<input type="radio"/> 20" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31527	n/c
<input type="radio"/> 3" Poly	NXT-34888	\$ 55	<input type="radio"/> 20" Airless Insert (E2213)	NXT-31022	\$ 75
<input type="radio"/> 3" Light Up Poly	NXT-34111	\$ 55	<input type="radio"/> 20" Plastic Coated Handrim	NXT-31016	\$117
<small>Note: 3" casters recommended for low seat to floor heights.</small>			<input type="radio"/> 22" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31528	n/c
<input type="radio"/> 4" Poly	NXT-40000	\$ 65	<input type="radio"/> 22" Airless Insert (E2213)	NXT-31023	\$ 75
<input type="radio"/> 4" Light Up Poly	NXT-41000	\$ 99	<input type="radio"/> 22" Plastic Coated Handrim	NXT-31017	\$117
<input type="radio"/> 4" x 1.4" Poly Aluminum Mag	NXT-43000	\$ 90	<input type="radio"/> 24" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31529	n/c
<input type="radio"/> 5" x 1" Poly Plastic Mag	NXT-34555	n/c	<input type="radio"/> 24" Airless Inserts (E2213)	NXT-31024	\$ 75
<input type="radio"/> 5" Light Up Poly	NXT-35050	\$ 99	<input type="radio"/> 24" Plastic Coated Handrim	NXT-31018	\$117
<input type="radio"/> 5" x 1.5" Poly Plastic Mag	NXT-34222	\$ 85	<small>Note: 22" &amp; 24" wheels not with short tilt bracket. Note: 22" and 24" wheels not with all dual post armrests.</small>		
<input type="radio"/> 5" x 1.4" Poly Aluminum Mag	NXT-42000	\$ 90	<input type="radio"/> Omit Handrims from 20", 22" or 24"	NXT-31006	(\$25)
<input type="radio"/> 6" x 1" Poly Plastic Mag	NXT-34666	n/c	<input type="radio"/> 20" Pneumatic Spoke with Aluminum Handrim	NXT-31540	\$200
<input type="radio"/> 6" x 1.5" Poly Plastic Mag	NXT-34999	\$ 85	<input type="radio"/> 20" Airless Insert (E2213)	NXT-31022	\$ 75
<input type="radio"/> 6" x 1.4" Poly Aluminum Mag	NXT-44000	\$ 90	<input type="radio"/> 22" Pneumatic Spoke with Aluminum Handrim	NXT-31541	\$200
<input type="radio"/> 6" x 2" Poly with Tread	NXT-34444	\$ 90	<input type="radio"/> 22" Airless Insert (E2213)	NXT-31023	\$ 75
<input type="radio"/> 6" Pneumatic Plastic Mag (E2214 / E2215)	NXT-34667	\$ 75	<input type="radio"/> 24" Pneumatic Spoke with Aluminum Handrim	NXT-31542	\$200
<input type="radio"/> 7" Poly Plastic Mag	NXT-34777	\$ 75	<input type="radio"/> 24" Airless Insert (E2213)	NXT-31024	\$ 75
<input type="radio"/> 7" x 1.5" Poly Plastic Mag	NXT-34000	\$116	<small>Note: 22" and 24" wheels not with short tilt bracket. Note: 22" and 24" wheels not with all dual post armrests. Note: 20", 22" and 24" spoke wheels not with hand/foot brake. Note: Plastic coated handrims not with all spoke wheels.</small>		
<input type="radio"/> 7" x 2" Poly	NXT-34200	\$116	<input type="radio"/> Omit Handrims from 20", 22" or 24"	NXT-31006	(\$25)
<small>Note: All 7" casters not with all 12" rear wheels.</small>			<b>Wheel Locks</b>		
<b>Rear Tires - Wheels - Accessories</b>			<input type="radio"/> Push to Lock - Standard	NXT-31501	n/c
<small>Note: All Wheel Selections Include Patented Double Locking Safety Axles.</small>			<input type="radio"/> Pull to Lock	NXT-31506	n/c
<input type="radio"/> Omit Wheels and Axles	NXT-31524	(\$55)	<input type="radio"/> Extension Handles (E0961)	NXT-31801	\$ 60
<input type="radio"/> 12" Low Profile Poly Mag (E2218)	NXT-31526	\$ 35	<input type="radio"/> Foot Activated Wheel Lock (E2228)	NXT-31508	n/c
<small>Note: Not with any foot/hand activated brake/wheel lock. Note: Not with airless inserts.</small>			<small>Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" poly plastic mag.</small>		
<input type="radio"/> 12" Semi Pneumatic Plastic Mag	NXT-31523	\$165	<input type="radio"/> Hand Activated Brake and Wheel Lock System (E2228)	NXT-30275	n/c
<input type="radio"/> 12" Pneumatic Plastic Mag	NXT-31525	n/c	<small>Note: Used to control direction, speed, brake and lock wheels from strolling position. Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" poly plastic mag.</small>		
<small>Note: All 12" wheels not with all handrims. Note: All 12" wheels not with all 7" casters.</small>			<input type="radio"/> Combination Dual Lever Hand Activated Tilt and Wheel Lock System (E2228)	NXT-30277	n/c
<input type="radio"/> 12" Airless Inserts (E2213)	NXT-31019	\$ 75	<small>Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position. Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" poly plastic mag.</small>		
<input type="radio"/> 16" Pneumatic Plastic Mag	NXT-31530	\$ 40	<b>Rear Anti-Tippers</b>		
<small>Note: All 16" wheels not with plastic coated handrims.</small>			<input type="radio"/> Curb Assist Anti-Tip Assembly (E0971)	NXT-31600	\$ 95
<input type="radio"/> 16" Airless Inserts (E2213)	NXT-31020	\$ 75	<small>Note: Pair includes one anti-tip tube with curb assist feature.</small>		
<input type="radio"/> 16" Aluminum Handrims	NXT-31014	\$166	<input type="radio"/> Standard Anti-Tip Assembly (E0971)	NXT-31603	\$ 95
<input type="radio"/> 18" Pneumatic Plastic Mag	NXT-31533	\$ 40	<small>Note: Flip Up, Adjustable</small>		
<small>Note: All 18" wheels not with plastic coated handrims.</small>			<b>Solid Seat Options</b>		
<input type="radio"/> 18" Airless Inserts (E2213)	NXT-31021	\$ 75	<input type="radio"/> Fixed Solid Seat Pan	NXT-39010	\$100
<input type="radio"/> 18" Aluminum Handrims	NXT-31015	\$166	<small>Note: Not with folding frame.</small>		
			<input type="radio"/> Removable Solid Seat Pan	NXT-39012	\$150
			<small>Note: All seat pans with 2" x 2" notches for growth. Note: All seat pans fit to frame unless otherwise specified below.</small>		
			<input type="radio"/> Seat Width _____" x Seat Depth _____"		

NXT Rider Weight Capacity is 225 lbs.



Date: \_\_\_\_\_ P.O.# \_\_\_\_\_

Marked for: \_\_\_\_\_

Description Part No. Retail

Description Part No. Retail

**Accessories**

- **Medical Necessities Pouch** NXT-23207 \$ 67  
Note: Black nylon.
- **Footrest Hanger Pads - Pair** SFDI-9HP \$ 90  
Note: Wrap-around pads made from black neoprene.
- **Frame Pads - Pair** 10911-2 \$ 80  
Note: Wrap-around pads made from black neoprene.  
Note: Typically used to pad exposed back post tubes.
- **Freedom Backpack**  NXT-10900 \$ 55  
Note: Black Nylon only with Freedom Designs logo.  
Note: Monogramming is not available on Freedom Backpack.
- **Freedom Canopy - Flip-Back**  NXT-10905 \$ 125  
Note: Canopy assembly includes removable mounting hardware.  
Note: Canopy is available in black only with Freedom Designs logo.  
Note: Canopy includes rear flap with zippered storage pocket.  
Note: Canopy includes clear-view window.  
Note: See offer for **FREE** Backpack and Canopy on last page of order form.

**HCPCS CODES** are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

- **Choose NXT Seat to Floor Height** \_\_\_\_\_"  
Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels and front casters size and tilt pivot height.  
Note: Caster size and types vary slightly in height.  
Note: Seat to floor heights are +/- 3/8".

**NXT Footrest / Extension Lengths**

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

**NXT MAPS Footrest / Extension Lengths**

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.  
Note: Without extensions bottom hanger hole not accessible for hardware.  
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

**Freedom Designs Armrest Heights**

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

**NXT Wheel / Caster Seat to Floor Heights**

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5"	18"	5"	Tall	16.5"
12"	3"	Tall	14"	18"	5"	Short	15.5"
12"	3"	Short	13"	18"	5"	Tall	17"
12"	3"	Tall	14.5"	18"	6"	Short	15.5"
12"	4"	Short	13.5"	18"	6"	Tall	17"
12"	4"	Tall	15"	20"	4"	Short	14"
12"	5"	Short	15.5"	20"	4"	Tall	15.5"
12"	5"	Tall	17"	20"	4"	Short	14.5"
12"	6"	Short	16"	20"	4"	Tall	16"
12"	6"	Tall	17.5"	20"	4"	Short	15"
16"	3"	Short	12.5"	20"	4"	Tall	16.5"
16"	3"	Tall	14"	20"	5"	Short	14.5"
16"	3"	Short	13"	20"	5"	Tall	16"
16"	3"	Tall	14.5"	20"	5"	Short	15"
16"	4"	Short	13.5"	20"	5"	Tall	16.5"
16"	4"	Tall	15"	20"	5"	Short	15.5"
16"	4"	Short	14"	20"	5"	Tall	17"
16"	4"	Tall	15.5"	20"	6"	Short	15.5"
16"	4"	Short	14.5"	20"	6"	Tall	17"
16"	4"	Tall	16"	20"	6"	Short	16"
16"	5"	Short	14.5"	20"	6"	Tall	17.5"
16"	5"	Tall	16"	20"	7"	Short	16.5"
16"	7"	Short	16.5"	20"	7"	Tall	18"
16"	7"	Tall	18"				
18"	3"	Short	12.5"	22"	4"	Tall	16.5"
18"	3"	Tall	14"				
18"	3"	Short	13"	22"	5"	Tall	16.5"
18"	3"	Tall	14.5"				
18"	4"	Short	13.5"	22"	5"	Tall	17"
18"	4"	Tall	15"				
18"	4"	Short	14"	22"	6"	Tall	17"
18"	4"	Tall	15.5"				
18"	4"	Short	14.5"	22"	6"	Tall	17.5"
18"	4"	Tall	16"				
18"	4"	Short	15"	22"	7"	Tall	18"
18"	4"	Tall	16.5"				
18"	5"	Short	14.5"	24"	6"	Tall	17.5"
18"	5"	Tall	16"				
18"	5"	Short	15"	24"	7"	Tall	18"

Note: All 22" and 24" wheels not with short tilt pivot.

**NXT Rider Weight Capacity is 225 lbs.**

**NXT Adult Transport Included**



Date: \_\_\_\_\_ P.O.# \_\_\_\_\_

Marked for: \_\_\_\_\_

**\* FREE \* FREE \* FREE \* FREE \* FREE \* FREE \***

# Freedom Designs & Pindot

Are Pleased To Announce  
**“A Sizzling Summer Sensation”**  
 And For a Limited Time It's...

**FREE \* FREE \* FREE**

When ordering a Freedom  
 ...Wheelchair and a  
 complete Freedom or Pindot  
 seating system including seat,  
 back, and attachment hardware,  
 simply provide the seating Quote  
 or seating PO number  
 below and receive a  
 Freedom Canopy and  
 Freedom Backpack

**\* FREE \***

- Freedom Backpack with Complete Freedom Seating System**  NXT-SYS-10900 n/c

**Note:** Black Nylon only with Freedom Designs logo.

**Note:** Monogramming is not available on Freedom backpack.

**Note:** Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

- Provide Freedom/PinDot Seating Quote/PO # \_\_\_\_\_

- Freedom Canopy - Flip-Back with Complete Freedom Seating System**  NXT-SYS-10905 n/c

**Note:** Canopy assembly includes removable mounting hardware.

**Note:** Canopy is available in black only with Freedom Designs logo.

**Note:** Canopy includes rear flap with zippered storage pocket.

**Note:** Canopy includes clear-view window.

**Note:** Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

- Provide Freedom/PinDot Seating Quote/PO # \_\_\_\_\_

**\* FREE \* FREE \* FREE \* FREE \* FREE \* FREE \***