



Our kids are our passion.™

NXT

Generation NEXT

Pediatric Manual Wheelchair



Yes, you can.™

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com

Date: _____ Buyer: _____
Company: _____
Marked For: _____
Email: _____

Quote Order PO# _____
Location: _____
Phone: _____
Fax: _____

Description	Part No.	Retail
Transport Optional Base		
NXT Pediatric Folding Transport Optional	(E1234) NXT-30300	\$2800
Note: Transportation Option selected separately below. Note: Includes Freedom Future Fit Program.		
NXT Pediatric Rigid Transport Optional	(E1233) NXT-30310	\$2700
Note: Transportation Option selected separately below. Note: Includes Freedom Future Fit Program.		
Transportation Option		
Note: Used when Transport Optional frame is selected above.		
Transportation Package	NXT-99998	\$ 250
Note: Must be installed by Freedom Designs manufacturing. Note: Not available with all vent and battery trays. Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved. Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit		

Transport Included Base		
NXT Pediatric Folding Transport Included	(E1234) NXT-30301	\$2950
Note: Includes Freedom Future Fit Program. Note: Not available with all vent and battery trays.		
NXT Pediatric Rigid Transport Included	(E1233) NXT-30311	\$2850
Note: Includes Freedom Future Fit Program. Note: Not available with all vent and battery trays.		

Seat Rail Options		
Note: Not available with ALL Contracture, 2D & 3D Dynamic footrest systems.		
NXT LEAF Seat Frame- Pair	(K0108) NXT-30324	\$ 714
NXT LEAF Seat Frame- Left	(K0108) NXT-30324-L	\$ 357
NXT LEAF Seat Frame- Right	(K0108) NXT-30324-R	\$ 357
Note: LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.		
Note: Freedom customized seats can be ordered to fit pre-configured LEAF frame models		
Note: No seating system is provided with NXT frames. A seating system must be purchased separately.		
NXT SRO Seat Frame- Pair	(K0108) NXT-30325	\$ 306
NXT SRO Seat Frame- Left	(K0108) NXT-30325-L	\$ 153
NXT SRO Seat Frame- Right	(K0108) NXT-30325-R	\$ 153
Note: SRO is manufactured without standard footrest hanger attachments.		
Note: SRO is used when customer attaches leg/footrest(s) directly to the underside of the solid seat.		

Tilt Activations		
Dual Lever Hand Activated Tilt	NXT-30272	n/c
Note: Must be selected for ALL vent and battery trays.		
Foot Activated Tilt	NXT-30273	n/c
Note: Foot pedal installed on right rear side frame. Note: Not with ALL vent and battery trays.		
Combination Dual Hand Activated Tilt, Brake and Wheel Lock System	NXT-30277	n/c
Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" Low Profile Poly Mag Wheels		

Description	Part No.	Retail
Vent & Battery Trays		
Note: All vent and battery trays not with transport.		
Note: All vent and battery trays not with short frames.		
Note: All vent trays on NXT only with Dual Lever Hand Activated Tilt.		
Note: All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.		
Vent Tray - Fixed	(E1029) 11-7000	\$ 439
Battery Tray - Fixed	(K0108) 11-8000	\$ 429
Quick Release Tray Hardware	(E1028) 563	\$ 272
Note: For LP10 style vent and battery tray on frames 16" wide and greater.		
Note: Vent tray is 4" narrower than frame width and 16" deep.		
Note: Battery tray is 4" narrower than frame width and 10" deep.		

Frame Colors			
Sugar Cookie	n/c	Pacific Blue	n/c
Neon Yellow	n/c	Blue Blast	n/c
Yellow	n/c	Electric Blue	n/c
Bright Pink	n/c	Galaxy Blue	n/c
Pink Parfait	n/c	Passion Purple	n/c
Sugar Plum	n/c	Octane Orange	n/c
Candy Violet	n/c	Red Sunset	n/c
Candy Raspberry	n/c	Really Red	n/c
Extreme Green	n/c	Fiery Red	n/c
Sparkle Green	n/c	Twilight	n/c
Mystic Teal	n/c	Stormy Gray	n/c
		Black	n/c

Frame Dimensions				
Specify Frame Width Setting				
10"	11"	12"	Adjusts 10" to 12"	NXT-30250 n/c
11"	12"	13"	Adjusts 11" to 13"	NXT-30210 n/c
12"	13"	14"	Adjusts 12" to 14"	NXT-30220 n/c
13"	14"	15"	Adjusts 13" to 15"	NXT-30230 n/c
14"	15"	16"	Adjusts 14" to 16"	NXT-30240 n/c
15"	16"	17"	Adjusts 15" to 17"	NXT-30290 n/c
16"	17"	18"	Adjusts 16" to 18"	NXT-30260 n/c
17"	18"		Adjusts 17" to 18"	NXT-30261 n/c
18"			Not width adjustable	NXT-30262 n/c
19"			Not adjustable, not w/ short frame	NXT-30263 n/c
20"			Not adjustable, not w/ short frame (E2201)	NXT-30264 \$ 400
Note: Frame widths measured across and from outside of seat frame tubes.				
Note: Seat frame widths of 19" and 20" not available with short frame.				
Note: Pediatric Long frame only with 15" wide and greater frames				
Note: Vent and battery trays only with 14" and wider frames.				

NXT Rider Weight Capacity is 225 lbs.

Page 1 of 7

NXT Pediatric

CE

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
<p align="center">Back Post Locations Specify Frame Type and Frame Length</p> <p>Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.</p> <p>Short Frame & Short Tilt Bracket NXT-30511 n/c</p> <p>12" NXT-30504 n/c</p> <p>13" NXT-30505 n/c</p> <p>14" NXT-30506 n/c</p> <p>15" NXT-30521 n/c</p> <p>16" NXT-30522 n/c</p> <p>* 17" Must select 2" extension NXT-30523 n/c</p> <p>* 18" Must select 2" extension NXT-30524 n/c</p> <p>2" Extension - Pair (K1008) NXT-30525 \$ 26</p> <p>Specify Back Post Location (if other than specified above) _____ "</p> <p>Note: Floor to seat frame height range is 12.5 to 16.5".</p> <p>Note: Seat frame widths of 19" and 20" not available with short frame.</p> <p>* Note: Must choose an extension for 17" and 18" frame lengths.</p> <p>Note: 22" & 24" wheels not available with short frame and short tilt bracket.</p> <p>Note: All vent and battery trays not available with short frame.</p> <p>Medium Frame & Short Tilt Bracket NXT-30512 n/c</p> <p>14" NXT-30506 n/c</p> <p>15" NXT-30521 n/c</p> <p>16" NXT-30522 n/c</p> <p>17" NXT-30523 n/c</p> <p>18" NXT-30524 n/c</p> <p>* 19" Must select 2" extension NXT-30527 n/c</p> <p>* 20" Must select 2" extension (E2203) NXT-30528 \$ 500</p> <p>* 21" Must select 3" extension (E2203) NXT-30529 \$ 500</p> <p>2" Extension - Pair (K1008) NXT-30525 \$ 26</p> <p>3" Extension - Pair (K1008) NXT-30530 \$ 26</p> <p>Specify Back Post Location (if other than specified above) _____ "</p> <p>Note: Floor to seat frame height range is 12.5" to 16.5".</p> <p>* Note: Must choose an extension for 19", 20", and 21" frame lengths.</p> <p>Note: 22" & 24" wheels not available with medium frame and short tilt bracket.</p> <p>Medium Frame & Tall Tilt Bracket NXT-30540 n/c</p> <p>14" NXT-30506TP n/c</p> <p>15" NXT-30521TP n/c</p> <p>16" NXT-30522TP n/c</p> <p>17" NXT-30523TP n/c</p> <p>18" NXT-30524TP n/c</p> <p>* 19" Must select 2" extension NXT-30527TP n/c</p> <p>* 20" Must select 2" extension (E2203) NXT-30528TP \$ 500</p> <p>* 21" Must select 3" extension (E2203) NXT-30529TP \$ 500</p> <p>* 22" Must select 4" extension (E2204) NXT-30535TP \$ 900</p> <p>2" Extension - Pair (K1008) NXT-30525 \$ 26</p> <p>3" Extension - Pair (K1008) NXT-30530 \$ 26</p> <p>4" Extension - Pair (K1008) NXT-30531 \$ 26</p> <p>Specify Back Post Location (if other than specified above) _____ "</p> <p>Note: Tall tilt bracket increases floor to seat rail height 1.5"</p> <p>Note: Floor to seat frame height range is 14" to 18".</p> <p>* Note: Must choose an extension for 19", 20", 21", and 22" frame lengths.</p> <p>5° Anterior Tilt Seat Frame NXT-30276 n/c</p> <p>Note: Limits posterior tilt to 40°, includes split collars for adjustment.</p> <p>Note: Medium with Tall Tilt and Long frames only.</p> <p>Note: NXT has 4" of depth adjustment for growth. Depth is measured from front of back post to end of the seat frame. For back post location add seating / hardware thickness (typically 2") to client's actual seat depth. Example: seat depth of 12" + back / hardware thickness of 2" = back post location set at 14"</p>			<p>Notes:</p>		
			<p align="center">Back Post Tubes</p> <p align="center">Specify Push Handle, Back Post and Extensions</p> <p>Standard 90° Bend Push Handles and Back Post Tubes NXT-30518 n/c</p> <p>Ergo Bend Push Handles and Back Post Tubes NXT-30519 n/c</p> <p>Note: All push/stroller handle extensions not recommended with Ergo Bend.</p> <p>Seating Dynamics Rocker Back (K0108) NXT-30532 \$1734</p> <p>Note: Angle adjustable from 72° to 108° in 6° increments.</p> <p>Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.</p> <p>Note: Locks into a non-dynamic position when necessary.</p> <p>Note: Does not allow the back post tubes to fold forward.</p> <p align="center">Back Post Height</p> <p>Note: Back Posts are measured from top of seat rail to top of handle bend.</p> <p>Short Push Handle - 20" high NXT-30515 n/c</p> <p>Medium Push Handle - 23" high NXT-30514 n/c</p> <p>Tall Push Handle - 26" high NXT-30516 n/c</p>		

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
Back Post Tubes (continued)			Legrests, Footrests and Extensions (continued)		
Stroller Handle Extensions			Footplates & Footboards		
<p>Note: All stroller handle ext. add approximately 6.5" to push handle height. Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.</p>			<p>Omit Footrest Plates NXT-36404</p> <p>Composite Flip-Up NXT-17070 n/c</p> <p>Kids Aluminum Angle Adjustable, Flip-up (K0040) NXT-36020 \$ 160 Note: For frames 14" wide and narrower.</p> <p>Heel Loops (E0951) NXT-36022 \$ 41</p> <p>Calf Strap - Webbing (K0038) NXT-36024 \$ 41</p> <p>Adult Aluminum Angle Adj. Flip-up w/ Heel Loops (K0040/E0951) NXT-36021 \$ 180 Note: For frames 15" wide and greater.</p> <p>Omit Heel Loops NXT-36023 n/c</p> <p>Aluminum Angle Adjustable Footboard Removable (K0040) NXT-36029 \$ 160</p>		
			<p style="text-align: center;">MAPS Footplate Mounting Hardware</p> <p>MAPS (Multi Axis Positioning System) - Pair NXT-36040 \$ 306</p> <p>MAPS (Multi Axis Positioning System) - Left NXT-36040-L \$ 153</p> <p>MAPS (Multi Axis Positioning System) - Right NXT-36040-R \$ 153 Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning. Note: Must choose aluminum footplates. Note: Not available on 10" & 11" wide frames Note: MAPS not available with composite footplates or aluminum footboard.</p>		
<p>Push Handle Extension - Width Adjustable NXT-30560 n/c Note: Not removable, but allows frame to fold. Note: Mounts direct to back post push handle</p> <p>Push Button, Angle-Adjustable Folding Stroller Handle Extension NXT-30570 \$ 112 Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below. Note: Not removable, but allows frame to fold. Note: Mounts direct to back post push handle</p> <p>Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty NXT-30575 \$ 179 Note: Available 10W-22w wheelchair frames. Note: Not removable, but allows frame to fold. Note: Same as above but with metal hinge mechanism. Note: Mounts direct to back post push handle.</p> <p>Push Button, Angle-Adjustable, Folding Push Handle Extension NXT-30580 \$ 112 Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below. Note: Not removable, but allows frame to fold. Note: Mounts direct to straight back post tubes.</p> <p>Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty NXT-30585 \$ 179 Note: Available 10W-22w wheelchair frames. Note: Not removable, but allows frame to fold. Note: Same as above but with metal hinge mechanism. Note: Mounts direct to straight back post tubes.</p> <p>Individual Stroller Handle Extensions Removable - Pair NXT-30517 \$ 102</p>			<p style="text-align: center;">Legrest Extensions</p> <p>4" Short Length Adjustable Extensions NXT-30285 n/c</p> <p>6" Medium Length Adjustable Extensions NXT-30286 n/c</p> <p>8" Long Length Adjustable Extensions NXT-30287 n/c Note: Not available with all 90° hangers.</p> <p>4.5" Heavy Duty ELR Extensions NXT-36401 \$ 56</p> <p>4.5" Heavy Duty Extensions NXT-36402 \$ 56 Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p> <p>6.5" Heavy Duty Extensions NXT-36403 \$ 66 Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</p> <p>High Mount Brackets NXT-36407 \$ 102 Note: Includes flip-up mounting hardware for future growth adjustment. Note: Not available with all 75° hangers.</p> <p>Specify Knee to Heel Measurement _____" (see matrix page)</p>		
Legrests, Footrests and Extensions			MAPS Footplate Mounting Hardware		
Legrests			MAPS Footplate Mounting Hardware		
<p>Omit Legrest Hangers & Footrests NXT-33000</p> <p>75° Pop-Up Swing Away Hangers NXT-30280 n/c Note: Not available with high mount bracket.</p> <p>90° Pop-Up Swing Away Hangers NXT-30281 n/c Note: Not available with 8" adjustable extensions.</p> <p>90° Pop-Up Swing Away Elevating Legrests (ELRs) (E0990) NXT-30279 \$ 255 Note: Not available with all length adjustable extensions. Note: Not available with 6.5" heavy duty extensions.</p> <p>75° Latch Style Swing Away Hangers NXT-30282 n/c Note: Not available with high mount bracket.</p> <p>90° Latch Style Swing Away Hangers NXT-30283 n/c Note: Not available with 8" adjustable extensions.</p>			<p style="text-align: center;">MAPS Footplate Mounting Hardware</p> <p>MAPS (Multi Axis Positioning System) - Pair NXT-36040 \$ 306</p> <p>MAPS (Multi Axis Positioning System) - Left NXT-36040-L \$ 153</p> <p>MAPS (Multi Axis Positioning System) - Right NXT-36040-R \$ 153 Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning. Note: Must choose aluminum footplates. Note: Not available on 10" & 11" wide frames Note: MAPS not available with composite footplates or aluminum footboard.</p>		
NXT Rider Weight Capacity is 225 lbs.			NXT Pediatric		

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Legrests, Footrests and Extensions Specialized Footrest Systems

Note: ALL legrest measurements are from the seat rail to footplate.
Note: Not available with LEAF or SRO seat rail options.

Super Contracture Footrest Hanger (K0108)

Short Pair - 11" to 15"	NXT-30291-2	\$3570
Long Pair - 13" to 18"	NXT-30292-2	\$3570

Note: Super Contracture for frames 15" wide and greater.
Note: Super Contracture is both swing-away and removable.
Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.
Note: Super Contracture has short or long length-adjustable extension tubes with adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar / dorsiflexion positioning.
Note: Super Contracture includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.
Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.
Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

2D Dynamic Footrest Hanger (K0108)

Short Pair - 14" to 16"	NXT-30293-2	\$4080
Long Pair - 16" to 18"	NXT-30294-2	\$4080

Note: 2D Dynamic for frames 15" wide and greater.
Note: 2D Dynamic is both swing-away and removable.
Note: 2D Dynamic has adjustable knee joints that have 48 positions that adjust in 7.5° increments.
Note: 2D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.
Note: 2D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.
Note: 2D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.
Note: Frame and/or front caster choice may limit 2D Dynamic adjustment range and/or cause footplate/front caster interference.

3D Dynamic Footrest Hanger (K0108)

Short Pair - 13" to 15"	NXT-30295-2	\$4590
Long Pair - 15" to 17"	NXT-30296-2	\$4590

Note: 3D Dynamic for frames 15" wide and greater.
Note: 3D Dynamic is both swing-away and removable.
Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.
Note: 3D Dynamic has short or long length-adjustable extension tubes with 1.5" of dynamic telescoping travel.
Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.
Note: 3D Dynamic includes a 5"x 7" aluminum footplate with multiple holes for location adjustment as well as secondary positioning accessory attachment.
Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

Armrests and Armrest Pads

Omit Armrests NXT-38209

Back Post Mounted Armrests

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

10.5" Long Flip-Back	NXT-38204	n/c
Note: Full length pad on 13.5" or longer armrests only.		
13.5" Long Flip-Back	NXT-38202	n/c
Note: Full length pad on 13.5" or longer armrests only.		
10.5" Long Flip-Back Locking (E0973)	NXT-38201	\$ 245
Note: Full length pad on 13.5" or longer armrests only.		
13.5" Long Flip-Back Locking (E0973)	NXT-38200	\$ 245
Note: Full length pad on 13.5" or longer armrests only.		

Armrests and Armrest Pads (continued)

10.5" to 13.5" Long Flip-Back, Telescoping (E0973)	NXT-38208	\$ 240
Note: Full length pad on 13.5" or longer armrests only.		
10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973)	NXT-38210	\$ 265
Note: Full length pad on 13.5" or longer armrests only.		
17" to 20" Long Flip-Back Locking, Telescoping (E0973)	NXT-38215	\$ 265

Foam Wraps - Standard	NXT-38216	n/c
Desk Length Pads	NXT-38213	n/c
Full Length Pads	NXT-38214	n/c
Note: 8.5" armrests not available with desk and full length armrest pads.		
Note: Full length pad on 13.5" or longer armrests only.		

Frame Mounted Armrests

Single Post Fixed Height, Removable - Short (K0015)	NXT-38222	\$ 225
Single Post Fixed Height, Removable - Tall (K0015)	NXT-38224	\$ 225
Note: Single Post Fixed Height arms include a removable side guard.		
Note: Adjustable in height with tools only. See matrix for height ranges.		

Single Post Height Adjustable, Removable - Short (E0973)	NXT-38217	\$ 265
Single Post Height Adjustable, Removable - Tall (E0973)	NXT-38203	\$ 265

Dual Post Fixed Height Removable	NXT-38205	n/c
Note: All dual post armrests not available with short frame.		
Note: All dual post armrests not available with 22" or 24" wheels.		
Dual Post Flip-Back Height Adjustable, Removable (E0973)	NXT-38206	\$ 275
Note: All dual post armrests not available with short frame.		
Note: All dual post armrests not available with 22" or 24" wheels.		

Desk Length Pads	NXT-38213	n/c
Full Length Pads	NXT-38214	n/c

Specify Seat Rail to Top of Armrest Measurement _____ "

Front Wheels - Casters

Omit Caster Wheels	NXT-34333	
3" Polyurethane	NXT-34888	\$ 56
3" Polyurethane Light-Up	NXT-34111	\$ 56
Note: 3" casters recommended for low seat to floor heights.		
4" Polyurethane	NXT-40000	\$ 66
4" Polyurethane Light-Up	NXT-41000	\$ 101
4" x 1.4" Poly, Aluminum Mag	NXT-43000	\$ 92
5" x 1" Poly, Plastic Mag	NXT-34555	n/c
5" Polyurethane Light-Up	NXT-35050	\$ 101
5" x 1.5" Poly, Plastic Mag	NXT-34222	\$ 87
5" x 1.4" Poly, Aluminum Mag	NXT-42000	\$ 92
6" x 1" Poly, Plastic Mag	NXT-34666	n/c
6" x 1.5" Poly, Plastic Mag	NXT-34999	\$ 87
6" x 1.4" Poly, Aluminum Mag	NXT-44000	\$ 92
6" x 2" Poly with Tread	NXT-34444	\$ 92
6" x 1" Pneumatic Plastic Mag (E2214/E2215)	NXT-34667	\$ 77
7" x 1" Poly, Plastic Mag	NXT-34777	\$ 77
7" x 2" Polyurethane	NXT-34200	\$ 118
Note: All 7" casters not available with all 12" rear wheels.		

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
Rear Wheels - Tires - Accessories		
Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.		
Omit Rear Wheels and Axles	NXT-31524	
Mag Wheels		
12" Low Profile Poly Mag (E2218)	NXT-31526	\$ 36
Note: Not available with any foot/hand activated brake/wheel lock system.		
Note: Not available with airless inserts.		
Note: All 12" wheels not available with all 7" casters.		
12" Semi Pneumatic Plastic Mag (K0069)	NXT-31523	\$ 168
Note: Not available with airless inserts.		
Note: All 12" wheels not available with all 7" casters.		
12" Pneumatic Plastic Mag	NXT-31525	n/c
12" Airless Inserts (E2213)	NXT-31019	\$ 77
Note: All 12" wheels not available with all handrims.		
Note: All 12" wheels not available with all 7" casters.		
16" Pneumatic Plastic Mag	NXT-31530	\$ 41
16" Airless Inserts (E2213)	NXT-31020	\$ 77
16" Aluminum Handrims	NXT-31014	\$ 169
Note: 16" mag wheels not available with plastic coated handrims.		
18" Pneumatic Plastic Mag	NXT-31533	\$ 41
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Aluminum Handrims	NXT-31015	\$ 169
Note: 18" mag wheels not available with plastic coated handrims.		
20" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31527	n/c
20" Airless Inserts (E2213)	NXT-31022	\$ 77
20" Plastic Coated Handrims	NXT-31016	\$ 119
22" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31528	n/c
22" Airless Inserts (E2213)	NXT-31023	\$ 77
22" Plastic Coated Handrims	NXT-31017	\$ 119
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
24" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31529	n/c
24" Airless Inserts (E2213)	NXT-31024	\$ 77
24" Plastic Coated Handrims	NXT-31018	\$ 119
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	
Spoke Wheels		
16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	NXT-31534	n/c
Note: Rim, hub & spokes are anodized silver.		
16" Airless Inserts (E2213)	NXT-31020	\$ 77
16" Silver Aluminum Handrim	NXT-31030	\$ 77
16" Black Plastic Coated Handrim	NXT-31032	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
16" Black Pneumatic Tire & Black Spoke Aluminum Wheel	NXT-31535	n/c
Note: Rim, hub & spokes are anodized black.		
16" Airless Inserts (E2213)	NXT-31020	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
16" Black Aluminum Handrim	NXT-31031	\$ 77
16" Black Plastic Coated Handrim	NXT-31032	\$ 77

Description	Part No.	Retail
Rear Wheels - Tires - Accessories (continued)		
Note: 16" & 18" spoke wheels do not include handrims.		
18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel	NXT-31537	n/c
Note: Rim, hub & spokes are anodized silver.		
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Silver Aluminum Handrim	NXT-31034	\$ 77
18" Black Plastic Coated Handrim	NXT-31036	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
18" Black Pneumatic Tire & Black Spoke Aluminum Wheel	NXT-31538	n/c
Note: Rim, hub & spokes are anodized black.		
18" Airless Inserts (E2213)	NXT-31021	\$ 77
18" Black Aluminum Handrim	NXT-31035	\$ 77
18" Black Plastic Coated Handrim	NXT-31036	\$ 77
Note: 16" & 18" spoke wheels do not include handrims.		
20" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31540	\$ 204
20" Airless Inserts (E2213)	NXT-31022	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
22" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31541	\$ 204
22" Airless Inserts (E2213)	NXT-31023	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
24" Pneumatic Spoke Wheel with Aluminum Handrim	NXT-31542	\$ 204
24" Airless Inserts (E2213)	NXT-31024	\$ 77
Note: Not available with foot/hand activated wheel locks.		
Note: All 22" & 24" wheels not available with short tilt bracket.		
Note: All 22" & 24" wheels not available with dual post armrests.		
Note: 20", 22" & 24" spoke wheels not available plastic coated handrims.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	

Continue To Next Page

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Wheel Locks

Push to Lock - Standard	NXT-31501	n/c
Pull to Lock	NXT-31506	n/c
Extension Handles (E0961)	NXT-31801	\$ 61

Foot Activated Wheel Lock (E2228)	NXT-31508	n/c
Note: Includes manual wheel locks to be used with public transportation.		
Note: Not available with 12" poly plastic mag wheels.		
Note: Not available with ALL vent and battery trays.		

Hand Activated Brake and Wheel Lock System (E2228)	NXT-30275	n/c
Note: Includes manual wheel locks to be used with public transportation.		
Note: Not available with 12" mag wheels with low pro poly tires.		
Note: Used to control direction, speed, brake and lock wheels from strolling position.		
Note: Not available with all vent and battery trays.		

Combination Dual Hand Activated Tilt, Brake and Wheel Lock System (E2228)	NXT-30277	n/c
Note: Includes manual wheel locks to be used with public transportation.		
Note: Not available with 12" poly plastic mag wheels.		
Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position.		
Note: Not available with all vent and battery trays.		

Rear Anti-Tippers

Curb Assist Anti-Tip Assembly Flip-Up, Adjustable (E0971)	NXT-31600	\$ 97
Note: Pair includes one anti-tip tube with curb assist feature.		

Standard Anti-Tip Assembly Flip-Up, Adjustable (E0971)	NXT-31603	\$ 97
---	-----------	-------

Seat Pan Options

Fixed Solid Seat Pan - Bolteron (E2231)	NXT-39010	\$ 150
Note: Not available with folding frame.		

Fixed Solid Seat Pan - Aluminum (E2231)	NXT-39030	\$ 170
Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".		
Note: Not available with folding frame.		

Removable Solid Seat Pan - Bolteron (E2231)	NXT-39012	\$ 170
--	-----------	--------

Removable Solid Seat Pan - Aluminum (E2231)	NXT-39032	\$ 255
Note: Removable Seat Pan raises the effective seat to floor height by 5/8".		
Note: Aluminum Removable Solid Seat Pan adjusts 2" in width by mounting hardware to grow with NXT frame.		
Note: Removable solid seat pan includes 4-points of EZ-Lite mounting hardware.		
Note: All seat pans include 2" x 2" notches for growth adjustment.		
Note: All seat pans made to fit frame dimensions unless specified below.		

Seat Width _____" X Seat Depth _____"

Accessories

Medical Necessities Pouch	NXT-23207	\$ 68
Note: Black nylon.		

Footrest Hanger Pads - Pair	SFDI-9HP	\$ 92
Note: Wrap-around pads made from black neoprene.		

Frame Pads - Pair	10911-2	\$ 82
Note: Wrap-around pads made from black neoprene.		
Note: Typically used to pad exposed back post tubes.		

Flip-Down Rigidizing Bar	NXT-23210	\$ 214
Note: Mounts to back post tubes to add rigidity to frame.		
Note: Does NOT adjust in width to grow with frame.		

Freedom Backpack	NXT-10900	\$ 56
Note: Black Nylon only with Freedom Designs logo.		
Note: Monogramming is not available on Freedom Backpack.		

Freedom Canopy - Flip-Back	NXT-10905	\$ 128
Note: Canopy assembly includes removable mounting hardware.		
Note: Canopy is available in black only with Freedom Designs logo.		
Note: Canopy includes rear flap with zippered storage pocket.		
Note: Canopy includes clear-view window.		

Anterior Support Surfaces – (Trays)

3/8" Lexan Tray - Channel Bracket Hardware (E0950)	08-2000	\$ 271
---	---------	--------

3/8" Lexan Tray - Adjustable Knob Hardware & Ribs (E0950)	08-6000	\$ 219
--	---------	--------

1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs (E0950)	08-13000	\$ 183
--	----------	--------

1/4" Lexan Tray - Hook & Loop Straps (E0950)	08-16000	\$ 106
---	----------	--------

Continue To Next Page

Date: _____ PO# _____

Marked For: _____

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

NXT Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Choose NXT Seat to Floor Height _____"

Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height.
Note: Caster sizes vary slightly in height by type.
Note: Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".
Note: Seat pan choice increases seat rail to floor height 3/8" to 3/4".

NXT MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.
Note: Without extensions bottom hanger hole not accessible for hardware.
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

NXT Wheel / Caster Seat to Floor Heights

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5-13"	20"	3"	Tall	N/A
12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"
12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"
12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"
12"	5"	Short	15.5"	20"	5"	Tall	16-17"
12"	5"	Tall	17"	20"	6"	Short	15.5-16"
12"	6"	Short	16"	20"	6"	Tall	17-17.5"
12"	6"	Tall	17.5"	20"	7"	Short	16.5"
16"	3"	Short	12.5-13"	20"	7"	Tall	18"
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A
16"	5"	Short	14.5"	22"	4"	Tall	16.5"
16"	5"	Tall	16"	22"	5"	Short	N/A
16"	6"	Short	16"	22"	5"	Tall	16.5-17"
16"	6"	Tall	17.5"	22"	6"	Short	N/A
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"
16"	7"	Tall	18"	22"	7"	Short	N/A
18"	3"	Short	12.5-13"	22"	7"	Tall	18"
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A
18"	5"	Tall	16-17"	24"	5"	Short	N/A
18"	6"	Short	15.5"	24"	5"	Tall	N/A
18"	6"	Tall	17"	24"	6"	Short	N/A
18"	7"	Short	N/A	24"	6"	Tall	17.5"
18"	7"	Tall	N/A	24"	7"	Short	N/A
20"	3"	Short	N/A	24"	7"	Tall	18"