



Our kids are our passion.
2241 N. Madera Rd

NXT
Generation **NXT**

Simi Valley, CA 93065

Pediatric Manual Wheelchair Transport Package Optional

Phone: (800) 331-8551 Fax: (888) 582-1509

www.freedomdesigns.com



Date: _____ Buyer: _____
Company: _____
Marked For: _____
Email: _____

Quote Order PO# _____
Location: _____
Phone: _____
Fax: _____

Description	Part No.	Retail
NXT Pediatric Frame Options		
<input type="radio"/> NXT Pediatric Folding Transport Optional Note: Includes Freedom Future Fit Program.	(E1234) NXT-30300	\$2800
<input type="radio"/> NXT Pediatric Rigid Transport Optional Note: Includes Freedom Future Fit Program.	(E1233) NXT-30310	\$2700
Transport Option		
<input type="radio"/> Transport Package Note: Must be installed by Freedom Designs manufacturing. Note: Not as parts. Note: Not with all vent and battery trays. Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved.	NXT-99998	\$150
<input type="radio"/> NXT LEAF Seat Frame - Pair	(K0108) NXT-30324	\$700
<input type="radio"/> NXT LEAF Seat Frame - Left	(K0108) NXT-30324-L	\$350
<input type="radio"/> NXT LEAF Seat Frame - Right Note: LEAF (Lower Extremity Accommodation Frame) front frame angle adjusts 22° from midline for abduction, adduction and windswept lower extremity positioning. LEAF seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy. Note: Freedom customized seats can be ordered to fit preconfigured LEAF frame models. Note: No seating system is provided with NXT frames. A seating system must be purchased separately.	(K0108) NXT-30324-R	\$350
<input type="radio"/> NXT SRO - Seat Rail Only - Pair	(K0108) NXT-30325	\$300
<input type="radio"/> NXT SRO - Seat Rail Only - Left	(K0108) NXT-30325-L	\$150
<input type="radio"/> NXT SRO - Seat Rail Only - Right Note: SRO is manufactured without standard footrest hanger attachments. Note: SRO used when customer attaches leg/footrest(s) directly to underside of seat.	(K0108) NXT-30325-R	\$150
Tilt Activations		
<input type="radio"/> Dual Lever Hand Activated Tilt Note: Must be selected for all vent and battery trays.	NXT-30272	n/c
<input type="radio"/> Foot Activated Tilt Note: Foot pedal installed on right rear side frame. Note: Not with all vent and battery trays.	NXT-30273	n/c
<input checked="" type="radio"/> Combination Dual Lever Hand Activated Tilt and Wheel Lock System Note: Includes manual wheel locks to be used with public transportation. Note: Not with all vent and battery trays.	(E2228) NXT-30277	n/c
Vent & Battery Trays		
<input type="radio"/> Vent Tray - Fixed	(E1029) 11-7000	\$430
<input type="radio"/> Battery Tray - Fixed	(K0108) 11-8000	\$421
<input type="radio"/> Quick Release Tray Hardware Note: For LP10 vent and battery tray with frames 16" wide and greater. Note: Vent tray fixed is 14" wide x 16" deep. Note: Battery tray fixed is 12" wide x 10" deep.	(E1028) 563	\$267
<input type="radio"/> Custom Vent Tray - Fixed	(E1029) 11-1000	\$430
<input type="radio"/> Custom Battery Tray - Fixed	(K0108) 11-5000	\$421
<input type="radio"/> Quick Release Tray Hardware Note: Custom vent / battery trays for 14" and 15" frame sizes. Note: All vent and battery trays not with transport. Note: All vent and battery trays not with short frames. Note: All vent and battery trays only with dual lever hand activated tilt. Note: All vent and battery trays only with standard manual wheel locks. Note: All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.	(E1028) 563	\$267

Description	Part No.	Retail
NXT Frame Colors		
<input type="radio"/> Sugar Cookie	n/c	<input type="radio"/> Pacific Blue n/c
<input type="radio"/> Neon Yellow	n/c	<input type="radio"/> Blue Blast n/c
<input type="radio"/> Yellow	n/c	<input checked="" type="radio"/> Electric Blue n/c
<input type="radio"/> Bright Pink	n/c	<input type="radio"/> Galaxy Blue n/c
<input type="radio"/> Pink Parfait	n/c	<input type="radio"/> Passion Purple n/c
<input type="radio"/> Sugar Plum	n/c	<input type="radio"/> Octane Orange n/c
<input type="radio"/> Vigorous Violet	n/c	<input type="radio"/> Red Sunset n/c
<input type="radio"/> Candy Violet	n/c	<input type="radio"/> Really Red n/c
<input type="radio"/> Candy Raspberry	n/c	<input checked="" type="radio"/> Fiery Red n/c
<input type="radio"/> Extreme Green	n/c	<input type="radio"/> Twilight n/c
<input type="radio"/> Sparkle Green	n/c	<input checked="" type="radio"/> Stormy Gray n/c
<input type="radio"/> Mystic Teal	n/c	<input type="radio"/> Black n/c
<input type="radio"/> Blizzard Blue	n/c	
Frame Dimensions		
Specify Frame Width Setting		
<input type="radio"/> 10"	<input type="radio"/> 11"	<input type="radio"/> 12" Adjusts 10" to 12" NXT-30250 n/c
<input type="radio"/> 11"	<input type="radio"/> 12"	<input type="radio"/> 13" Adjusts 11" to 13" NXT-30210 n/c
<input type="radio"/> 12"	<input type="radio"/> 13"	<input type="radio"/> 14" Adjusts 12" to 14" NXT-30220 n/c
<input type="radio"/> 13"	<input type="radio"/> 14"	<input type="radio"/> 15" Adjusts 13" to 15" NXT-30230 n/c
<input type="radio"/> 14"	<input type="radio"/> 15"	<input type="radio"/> 16" Adjusts 14" to 16" NXT-30240 n/c
<input type="radio"/> 15"	<input type="radio"/> 16"	<input type="radio"/> 17" Adjusts 15" to 17" NXT-30290 n/c
<input type="radio"/> 16"	<input type="radio"/> 17"	<input type="radio"/> 18" Adjusts 16" to 18" NXT-30260 n/c
<input type="radio"/> 17"	<input type="radio"/> 18"	<input type="radio"/> Adjusts 17" to 18" NXT-30261 n/c
<input type="radio"/> 18"		Not adjustable NXT-30262 n/c
<input type="radio"/> 19"		Not adjustable, not w/short frame NXT-30263 n/c
<input type="radio"/> 20"		Not adjustable, not w/short frame (E2201) NXT-30264 \$300
Note: Frame widths measured across and from outside of seat frame tubes. Note: Seat frame widths of 19" and 20" not with short frame. Note: Vent and battery trays only with 14" and wider frames.		
Back Post Locations		
Specify Frame Type and Frame Length		
Note: Extensions may limit tilt range with 22" and 24" wheels and may affect stability when in tilt positions. Adjusting seat rail forward compensates for tilt range limitations and instability.		
<input type="radio"/> Short Frame & Short Tilt Bracket	NXT-30511	n/c
<input type="radio"/> 12"	NXT-30504	n/c
<input type="radio"/> 13"	NXT-30505	n/c
<input type="radio"/> 14"	NXT-30506	n/c
<input type="radio"/> 15"	NXT-30521	n/c
<input type="radio"/> 16"	NXT-30522	n/c
* <input type="radio"/> 17" Must select 2" extension	NXT-30523	n/c
* <input type="radio"/> 18" Must select 2" extension	NXT-30524	n/c
<input type="radio"/> 2" extension - Pair (K0108)	NXT-30525	\$ 25
<input type="radio"/> Specify location (if other than specified above) _____"		
Note: Floor to seat frame height range is 12.5" to 16.5". Note: Seat frame widths of 19" and 20" not with short frame. * Note: Must choose 2" extension for 17" and 18" frame lengths. Note: 22" & 24" wheels not with short frame and short tilt bracket. Note: All vent and battery trays not with short frame.		

(Back Post Locations continued on page 2)

NXT Rider Weight Capacity is 225 lbs.



Date: _____ P.O.# _____ Marked for: _____

Description	Part No.	Retail	Description	Part No.	Retail
Back Post Locations Specify Frame Type and Frame Length			Stroller Handle Extensions		
<input type="radio"/> Medium Frame & Short Tilt Bracket	NXT-30512	n/c	<input type="radio"/> Push Handle Extension Width Adjustable Removable	NXT-30513	n/c
<input type="radio"/> 14"	NXT-30506	n/c	<input type="radio"/> Stroller Handle Extension Width Adjustable Removable	NXT-30526	\$125
<input type="radio"/> 15"	NXT-30521	n/c	<input type="radio"/> Push Button Angle Adjustable Stroller Handle Extensions (PB)	NXT-30536	\$175
<input type="radio"/> 16"	NXT-30522	n/c	Note: PB stroller handles not with frames less than 14" wide. Note: PB stroller handle not removable and not with folding frame.		
<input type="radio"/> 17"	NXT-30523	n/c	<input type="radio"/> Individual Stroller Handle Extensions Removable - Pair	NXT-30517	\$100
<input type="radio"/> 18"	NXT-30524	n/c	Footrest Hangers and Extensions		
<input type="radio"/> 19" Must select 2" extension	NXT-30527	n/c	<input type="radio"/> Omit Footrest Hangers & Footplates	NXT-33000	(\$30)
<input type="radio"/> 20" Must select 2" extension (E2203)	NXT-30528	\$300	<input type="radio"/> 75° Pop-Up Swing Away Hangers	NXT-30280	n/c
<input type="radio"/> 21" Must select 3" extension (E2203)	NXT-30529	\$300	Note: Not with high mount bracket.		
<input type="radio"/> 2" extension - Pair (K0108)	NXT-30525	\$ 25	<input type="radio"/> 90° Pop-Up Swing Away Hangers	NXT-30281	n/c
<input type="radio"/> 3" extension - Pair (K0108)	NXT-30530	\$ 25	Note: Not with 8" adjustable extensions.		
<input type="radio"/> Specify location (if other than specified above) _____			<input type="radio"/> 75° Latch Style Swing Away Hangers	NXT-30282	n/c
Note: Floor to seat frame height range is 12.5" to 16.5".			Note: Not with high mount bracket.		
* Note: Must choose extension for 19", 20" and 21" frame lengths.			<input type="radio"/> 90° Latch Style Swing Away Hangers	NXT-30283	n/c
Note: 22" & 24" wheels not with medium frame and short tilt bracket.			Note: Not with 8" adjustable extensions.		
<input type="radio"/> Medium Frame & Tall Tilt Bracket	NXT-30540	n/c	<input type="radio"/> 4" Short Length Adjustable Extensions	NXT-30285	n/c
<input type="radio"/> 14"	NXT-30506TP	n/c	<input type="radio"/> 6" Medium Length Adjustable Extensions	NXT-30286	n/c
<input type="radio"/> 15"	NXT-30521TP	n/c	<input type="radio"/> 8" Long Length Adjustable Extensions	NXT-30287	n/c
<input type="radio"/> 16"	NXT-30522TP	n/c	Note: Not with all 90° hangers.		
<input type="radio"/> 17"	NXT-30523TP	n/c	<input type="radio"/> 4.5" Heavy Duty Extensions	NXT-36402	\$ 55
<input type="radio"/> 18"	NXT-30524TP	n/c	Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.		
<input type="radio"/> 19" Must select 2" extension	NXT-30527TP	n/c	<input type="radio"/> 6.5" Heavy Duty Extensions	NXT-36403	\$ 65
<input type="radio"/> 20" Must select 2" extension (E2203)	NXT-30528TP	\$300	Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.		
<input type="radio"/> 21" Must select 3" extension (E2203)	NXT-30529TP	\$300	<input type="radio"/> High Mount Brackets (K0037)	NXT-36407	\$100
<input type="radio"/> 22" Must select 4" extension (E2204)	NXT-30535TP	\$500	Note: Includes flip up mounting hardware for growth. Note: Not with 75° hangers.		
<input type="radio"/> 2" extension - Pair (K0108)	NXT-30525	\$ 25	<input type="radio"/> Pop-Up Swing Away Elevating (E0990)	NXT-30279	\$235
<input type="radio"/> 3" extension - Pair (K0108)	NXT-30530	\$ 25	Legrests (ELR)		
<input type="radio"/> 4" extension - Pair (K0108)	NXT-30531	\$ 25	Note: Not with all length adjustable extensions. Note: Not with 6.5" heavy duty extensions.		
<input type="radio"/> Specify location (if other than specified above) _____			<input type="radio"/> 4.5" Heavy Duty Extensions (ELR)	NXT-36401	\$ 55
Note: Tall tilt bracket increases floor to seat frame height by 1.5".			<input type="radio"/> Specify Knee to Heel Measurement _____ (see chart below)		
Note: Floor to seat frame height range is 14" to 18".			(See Matrix on Page 5)		
* Note: Must choose extension for 19", 20", 21" and 22" frame lengths.					
<input type="radio"/> 5° Anterior Tilt Seat Frame	NXT-30276	n/c			
Note: Limits posterior tilt to 40°, includes split collars for adjustment.					
Note: Tall tilt bracket only.					
NXT has 4" of depth adjustment for growth. Depth measured from front of back post to end of seat frame. For back post location add seating / hardware thickness (typically 2") to clients seat depth. Example: Seat depth of 12" + seat / hardware thickness of 2" = back post tube location set at 14".					
Back Post Tubes					
Specify Push Handle, Back Post Height, and Extensions					
Note: All stroller handle ext. add approximately 6.5" to push handle height.					
Note: Back posts measured from top of seat frame to top of handle bend.					
<input type="radio"/> Standard 90° Bend Push Handles and Back Post Tubes	NXT-30518	n/c			
<input type="radio"/> Ergo Bend Push Handles and Back Post Tubes	NXT-30519	n/c			
Note: All push / stroller handle ext. not recommended with Ergo Bend.					
<input type="radio"/> Seating Dynamics Rocker Back (K0108)	NXT-30532	\$1700			
Note: Angle adjustable from 72° to 108° in 6° increments.					
Note: Supplied with a complete set of 4 different durometers of elastomers to allow spring adjustment in the field.					
Note: Locks into a non dynamic position when necessary.					
Note: Does not allow the back post tubes to fold forward.					
<input type="radio"/> Flip Down Rigidizing Bar	(K0108) NXT-23210	\$ 210			
Note: Not width adjustable.					
Back Post Height					
<input type="radio"/> Short Push Handle - 19" High	NXT-30515	n/c			
<input type="radio"/> Medium Push Handle - 23" High	NXT-30514	n/c			
<input type="radio"/> Tall Push Handle - 26" High	NXT-30516	n/c			

NXT Rider Weight Capacity is 225 lbs.

Date: _____ P.O.# _____

Marked for: _____

Description	Part No.	Retail
-------------	----------	--------

Description	Part No.	Retail
-------------	----------	--------

Foot Plates

- Omit Footrest Plates** NXT-36404 **(\$15)**
- Composite Flip Up** NXT-17070 **n/c**
- Kids Aluminum Angle Adj. Flip Up** (K0040) NXT-36020 **\$150**
 Note: For frames 14" wide and under.
- Heel Loops** (E0951) NXT-36022 **\$ 50**
- Adult Aluminum Angle Adj. Flip Up Footrest Plates w/Heel Loops** (K0040 / E0951) NXT-36021 **\$150**
 Note: For frames 15" wide and greater.
- Omit Heel Loops** NXT-36023 **n/c**
- Aluminum Adjustable Angle Footboard Removable** NXT-36029 **\$150**
- MAPS (Multi Axis Position System) - Pair** NXT-36040 **\$300**
- MAPS (Multi Axis Position System) - Left** NXT-36040-L **\$150**
- MAPS (Multi Axis Position System) - Right** NXT-36040-R **\$150**
 Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.
 Note: Must choose aluminum footplates.
 Note: MAPS not with composite footrest plates or footboard.
 Note: MAPS not available on 10" & 11" wide frames.

**Footrest Hangers and Extensions
Specialized Footrest Systems**

 **Super Contracture Footrest Hanger (K0108)**

- Short Left - 11" to 15"** NXT-30291-L **\$1750**
- Short Right - 11" to 15"** NXT-30291-R **\$1750**
- Short Pair - 11" to 15"** NXT-30291-2 **\$3500**
- Long Left - 13" to 18"** NXT-30292-L **\$1750**
- Long Right - 13" to 18"** NXT-30292-R **\$1750**
- Long Pair - 13" to 18"** NXT-30292-2 **\$3500**
 Note: SC for NXT frames 15" wide or wider.
 Note: SC is both swing away and removable.
 Note: SC Has adjustable knee joints that have 48 positions that adjust in 7.5° increments.
 Note: SC has short or long length adjustable extension tubes that have adjustable ankle joints that have 48 positions that adjust in 7.5° increments for plantar/dorsiflexion positioning.
 Note: SC has 5"x 7" flip up aluminum foot plates that have multiple holes for location adjustment and secondary positioning accessory attachment.
 Note: NXT frame and/or front caster choices may limit SC knee joint adjustment range and/or cause foot plate/front caster interference.
 Note: SC position (in contracture position) may limit ankle joint/foot plate adjustment range.

 **2D Dynamic Footrest System (K0108)**

- Short Left - 14" to 16"** NXT-30293-L **\$2000**
- Short Right - 14" to 16"** NXT-30293-R **\$2000**
- Short Pair - 14" to 16"** NXT-30293-2 **\$4000**
- Long Left - 16" to 18"** NXT-30294-L **\$2000**
- Long Right - 16" to 18"** NXT-30294-R **\$2000**
- Long Pair - 16" to 18"** NXT-30294-2 **\$4000**
 Note: 2D for NXT frames 15" wide or wider.
 Note: 2D System is both swing away and removable.
 Note: 2D has adjustable knee joints that have 13 positions that adjust in 7.5° increments ranging from 0° to 100°.
 Note: 2D has short or long length adjustable extension tubes with 1.5" of dynamic telescoping travel.
 Note: 2D has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.
 Note: 2D has 5"x 7" angle adjustable aluminum foot plates that have multiple holes for location adjustment and secondary positioning accessory attachment.
 Note: NXT frame and/or front caster choices may limit 2D knee adjustment range and/or cause foot plate/front caster interference.

 **3D Dynamic Footrest System (K0108)**

- Short Left - 13" to 15"** NXT-30295-L **\$2250**
- Short Right 13" to 15"** NXT-30295-R **\$2250**
- Short Pair - 13" to 15"** NXT-30295-2 **\$4500**
- Long Left - 15" to 17"** NXT-30296-L **\$2250**
- Long Right - 15" to 17"** NXT-30296-R **\$2250**
- Long Pair - 15" to 17"** NXT-30296-2 **\$4500**
 Note: 3D for NXT frames 15" wide or wider.
 Note: 3D is both swing away and removable.
 Note: 3D has adjustable dynamic elevating knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.
 Note: 3D has short or long length adjustable extension tubes with 1.5" of dynamic telescoping travel.
 Note: 3D has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.
 Note: 3D has 5"x 7" angle adjustable aluminum foot plates that have multiple holes for location adjustment and secondary positioning accessory attachment.
 Note: NXT frame and/or front caster choices may limit 3D knee joint adjustment range and/or cause foot plate/caster interference.

Armrests and Pads

- Omit Flip Back Armrests** NXT-38209 **(\$20)**
- 10.5" Long Flip Back** NXT-38204 **n/c**
- 13.5" Long Flip Back** NXT-38202 **n/c**
- 10.5" to 13.5" Long Flip Back Telescoping** (E0973) NXT-38208 **\$100**
- 10.5" Long Flip Back Locking** (E0973) NXT-38201 **\$235**
- 13.5" Long Flip Back Locking** (E0973) NXT-38200 **\$235**
- 10.5" to 13.5" Long Flip Back Telescoping Locking** (E0973) NXT-38210 **\$260**
- 17" to 20" Long Flip Back Telescoping Locking** (E0973) NXT-38215 **\$260**
 Note: Full length pad on 13.5" or longer armrests only.
 Note: Flip back armrest length measured from front of backpost to end of armrest with foam wrap installed.
- Foam Wraps - Standard** NXT-38216 **n/c**
- Desk Length Pads** NXT-38213 **n/c**
- Full Length Pads** NXT-38214 **n/c**
 Note: Full length on 13.5" or longer armrests only.
- Single Post Height Adj. Removable Short** (E0973) NXT-38217 **\$235**
- Single Post Height Adj. Removable Tall** (E0973) NXT-38203 **\$235**
- Dual Post Fixed Height Removable** NXT-38205 **n/c**
- Dual Post Flip Back Height Adj. Removable** (E0973) NXT-38206 **\$245**
 Note: All dual post armrests not with short frame.
 Note: All dual post armrests not with all 22" and 24" wheels.
- Desk Length Pads** NXT-38213 **n/c**
- Full Length Pads** NXT-38214 **n/c**
 Note: Full length on 13.5" or longer armrests only.
- Specify Seat Rail to Top of Armrest Measurement _____"**

(See Matrix on Page 5)

Date: _____ P.O.# _____

Marked for: _____

Description	Part No.	Retail
Front Tires - Caster Wheels		
<input type="radio"/> Omit Caster Wheels	NXT-34333	n/c
<input type="radio"/> 3" Poly	NXT-34888	\$ 55
<input type="radio"/> 3" Light Up Poly	NXT-34111	\$ 55
<small>Note: 3" casters recommended for low seat to floor heights.</small>		
<input type="radio"/> 4" Poly	NXT-40000	\$ 65
<input type="radio"/> 4" Light Up Poly	NXT-41000	\$ 99
<input type="radio"/> 4" x 1.4" Poly Aluminum Mag	NXT-43000	\$ 90
<input type="radio"/> 5" x 1" Poly Plastic Mag	NXT-34555	n/c
<input type="radio"/> 5" Light Up Poly	NXT-35050	\$ 99
<input type="radio"/> 5" x 1.5" Poly Plastic Mag	NXT-34222	\$ 85
<input type="radio"/> 5" x 1.4" Poly Aluminum Mag	NXT-42000	\$ 90
<input type="radio"/> 6" x 1" Poly Plastic Mag	NXT-34666	n/c
<input type="radio"/> 6" x 1.5" Poly Plastic Mag	NXT-34999	\$ 85
<input type="radio"/> 6" x 1.4" Poly Aluminum Mag	NXT-44000	\$ 90
<input type="radio"/> 6" x 2" Poly with Tread	NXT-34444	\$ 90
<input type="radio"/> 6" Pneumatic Plastic Mag (E2214 / E2215)	NXT-34667	\$ 75
<input type="radio"/> 7" Poly Plastic Mag	NXT-34777	\$ 75
<input type="radio"/> 7" x 1.5" Poly Plastic Mag	NXT-34000	\$116
<input type="radio"/> 7" x 2" Poly	NXT-34200	\$116
<small>Note: All 7" casters not with all 12" rear wheels.</small>		
Rear Tires - Wheels - Accessories		
<small>Note: All Wheel Selections Include Patented Double Locking Safety Axles.</small>		
<input type="radio"/> Omit Wheels and Axles	NXT-31524	(\$55)
<input type="radio"/> 12" Low Profile Poly Mag (E2218)	NXT-31526	\$ 35
<small>Note: Not with any foot/hand activated brake/wheel lock. Note: Not with airless inserts.</small>		
<input type="radio"/> 12" Semi Pneumatic Plastic Mag	NXT-31523	\$165
<input type="radio"/> 12" Pneumatic Plastic Mag	NXT-31525	n/c
<small>Note: All 12" wheels not with all handrims. Note: All 12" wheels not with all 7" casters.</small>		
<input type="radio"/> 12" Airless Inserts (E2213)	NXT-31019	\$ 75
<input type="radio"/> 16" Pneumatic Plastic Mag	NXT-31530	\$ 40
<small>Note: All 16" wheels not with plastic coated handrims.</small>		
<input type="radio"/> 16" Airless Inserts (E2213)	NXT-31020	\$ 75
<input type="radio"/> 16" Aluminum Handrims	NXT-31014	\$166
<input type="radio"/> 18" Pneumatic Plastic Mag	NXT-31533	\$ 40
<small>Note: All 18" wheels not with plastic coated handrims.</small>		
<input type="radio"/> 18" Airless Inserts (E2213)	NXT-31021	\$ 75
<input type="radio"/> 18" Aluminum Handrims	NXT-31015	\$166

Description	Part No.	Retail
Rear Tires - Wheels - Accessories		
<input type="radio"/> 20" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31527	n/c
<input type="radio"/> 20" Airless Insert (E2213)	NXT-31022	\$ 75
<input type="radio"/> 20" Plastic Coated Handrim	NXT-31016	\$117
<input type="radio"/> 22" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31528	n/c
<input type="radio"/> 22" Airless Insert (E2213)	NXT-31023	\$ 75
<input type="radio"/> 22" Plastic Coated Handrim	NXT-31017	\$117
<input type="radio"/> 24" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31529	n/c
<input type="radio"/> 24" Airless Inserts (E2213)	NXT-31024	\$ 75
<input type="radio"/> 24" Plastic Coated Handrim	NXT-31018	\$117
<small>Note: 22" & 24" wheels not with short tilt bracket. Note: 22" and 24" wheels not with all dual post armrests.</small>		
<input type="radio"/> Omit Handrims from 20", 22" or 24"	NXT-31006	(\$25)
<input type="radio"/> 20" Pneumatic Spoke with Aluminum Handrim	NXT-31540	\$200
<input type="radio"/> 20" Airless Insert (E2213)	NXT-31022	\$ 75
<input type="radio"/> 22" Pneumatic Spoke with Aluminum Handrim	NXT-31541	\$200
<input type="radio"/> 22" Airless Insert (E2213)	NXT-31023	\$ 75
<input type="radio"/> 24" Pneumatic Spoke with Aluminum Handrim	NXT-31542	\$200
<input type="radio"/> 24" Airless Insert (E2213)	NXT-31024	\$ 75
<small>Note: 22" and 24" wheels not with short tilt bracket. Note: 22" and 24" wheels not with all dual post armrests. Note: 20", 22" and 24" spoke wheels not with hand/foot brake. Note: Plastic coated handrims not with all spoke wheels.</small>		
<input type="radio"/> Omit Handrims from 20", 22" or 24"	NXT-31006	(\$25)
Wheel Locks		
<input type="radio"/> Push to Lock - Standard	NXT-31501	n/c
<input type="radio"/> Pull to Lock	NXT-31506	n/c
<input type="radio"/> Extension Handles (E0961)	NXT-31801	\$ 60
<input type="radio"/> Foot Activated Wheel Lock (E2228)	NXT-31508	n/c
<small>Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" poly plastic mag.</small>		
<input type="radio"/> Hand Activated Brake and Wheel Lock System (E2228)	NXT-30275	n/c
<small>Note: Used to control direction, speed, brake and lock wheels from strolling position. Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" poly plastic mag.</small>		
<input type="radio"/> Combination Dual Lever Hand Activated Tilt and Wheel Lock System (E2228)	NXT-30277	n/c
<small>Note: Used to control direction, speed, brake and lock wheels and tilt frame from strolling position. Note: Includes manual wheel locks to be used with public transportation. Note: Not with 12" poly plastic mag. Note: Not with all vent and battery trays.</small>		
Rear Anti-Tippers		
<input type="radio"/> Curb Assist Anti-Tip Assembly (E0971)	NXT-31600	\$ 95
<small>Note: Pair includes one anti-tip tube with curb assist feature.</small>		
<input type="radio"/> Standard Anti-Tip Assembly (E0971)	NXT-31603	\$ 95
Solid Seat Options		
<input type="radio"/> Fixed Solid Seat Pan	NXT-39010	\$100
<small>Note: Not with folding frame.</small>		
<input type="radio"/> Removable Solid Seat Pan	NXT-39012	\$150
<small>Note: All seat pans with 2" x 2" notches for growth. Note: All seat pans fit to frame unless otherwise specified below.</small>		
<input type="radio"/> Seat Width _____" x Seat Depth _____"		

NXT Rider Weight Capacity is 225 lbs.

Date: _____ P.O.# _____

Marked for: _____

Description	Part No.	Retail
Accessories		
<input type="radio"/> Medical Necessities Pouch Note: Black nylon.	NXT-23207	\$ 67
<input type="radio"/> Footrest Hanger Pads - Pair Note: Wrap-around pads made from black neoprene.	SFDI-9HP	\$ 90
<input type="radio"/> Frame Pads - Pair Note: Wrap-around pads made from black neoprene. Note: Typically used to pad exposed back post tubes.	10911-2	\$ 80
<input type="radio"/> Freedom Backpack  Note: Black Nylon only with Freedom Designs logo. Note: Monogramming is not available on Freedom Backpack.	NXT-10900	\$ 55
<input type="radio"/> Freedom Canopy - Flip-Back  Note: Canopy assembly includes removable mounting hardware. Note: Canopy is available in black only with Freedom Designs logo. Note: Canopy includes rear flap with zippered storage pocket. Note: Canopy includes clear-view window. Note: See offer for FREE Backpack and Canopy on last page of order form.	NXT-10905	\$ 125

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Choose NXT Seat to Floor Height _____"
Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels and front casters size and tilt pivot height.
Note: Caster size and types vary slightly in height.
Note: Seat to floor heights are +/- 3/8".

NXT Footrest / Extension Lengths					
Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A

NXT MAPS Footrest / Extension Lengths					
Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.
Note: Without extensions bottom hanger hole not accessible for hardware.
Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Freedom Designs Armrest Heights			
Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short	N/A	8" - 14"	8" - 14"
Single Post Height Adj. Tall	N/A	10" - 16"	10" - 16"
Dual Post Fixed Height	N/A	9.5"	9.5"
Dual Post Flip Back Height Adj.	N/A	9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

NXT Wheel / Caster Seat to Floor Heights							
Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5"	18"	5"	Tall	16.5"
12"	3"	Tall	14"	18"	5"	Short	15.5"
12"	3"	Short	13"	18"	5"	Tall	17"
12"	3"	Tall	14.5"	18"	6"	Short	15.5"
12"	4"	Short	13.5"	18"	6"	Tall	17"
12"	4"	Tall	15"	20"	4"	Short	14"
12"	5"	Short	15.5"	20"	4"	Tall	15.5"
12"	5"	Tall	17"	20"	4"	Short	14.5"
12"	6"	Short	16"	20"	4"	Tall	16"
12"	6"	Tall	17.5"	20"	4"	Short	15"
16"	3"	Short	12.5"	20"	4"	Tall	16.5"
16"	3"	Tall	14"	20"	5"	Short	14.5"
16"	3"	Short	13"	20"	5"	Tall	16"
16"	3"	Tall	14.5"	20"	5"	Short	15"
16"	4"	Short	13.5"	20"	5"	Tall	16.5"
16"	4"	Tall	15"	20"	5"	Short	15.5"
16"	4"	Short	14"	20"	5"	Tall	17"
16"	4"	Tall	15.5"	20"	6"	Short	15.5"
16"	4"	Short	14.5"	20"	6"	Tall	17"
16"	4"	Tall	16"	20"	6"	Short	16"
16"	5"	Short	14.5"	20"	6"	Tall	17.5"
16"	5"	Tall	16"	20"	7"	Short	16.5"
16"	7"	Short	16.5"	20"	7"	Tall	18"
16"	7"	Tall	18"				
18"	3"	Short	12.5"	22"	4"	Tall	16.5"
18"	3"	Tall	14"				
18"	3"	Short	13"	22"	5"	Tall	16.5"
18"	3"	Tall	14.5"				
18"	4"	Short	13.5"	22"	5"	Tall	17"
18"	4"	Tall	15"				
18"	4"	Short	14"	22"	6"	Tall	17"
18"	4"	Tall	15.5"				
18"	4"	Short	14.5"	22"	6"	Tall	17.5"
18"	4"	Tall	16"				
18"	4"	Short	15"	22"	7"	Tall	18"
18"	4"	Tall	16.5"				
18"	5"	Short	14.5"	24"	6"	Tall	17.5"
18"	5"	Tall	16"				
18"	5"	Short	15"	24"	7"	Tall	18"

Note: All 22" and 24" wheels not with short tilt pivot.

Date: _____ P.O.# _____

Marked for: _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***

Freedom Designs & Pindot

Are Pleased To Announce
"A Sizzling Summer Sensation"
And For a Limited Time It's...

FREE * FREE * FREE

When ordering a Freedom
Wheelchair and a
complete Freedom or Pindot
seating system including seat,
back, and attachment hardware,
simply provide the seating Quote
or seating PO number
below and receive a
Freedom Canopy and
Freedom Backpack

*** FREE ***

- Freedom Backpack with Complete Freedom Seating System**  NXT-SYS-10900 n/c

Note: Black Nylon only with Freedom Designs logo.

Note: Monogramming is not available on Freedom backpack.

Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

- Provide Freedom/PinDot Seating Quote/PO # _____

- Freedom Canopy - Flip-Back with Complete Freedom Seating System**  NXT-SYS-10905 n/c

Note: Canopy assembly includes removable mounting hardware.

Note: Canopy is available in black only with Freedom Designs logo.

Note: Canopy includes rear flap with zippered storage pocket.

Note: Canopy includes clear-view window.

Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

- Provide Freedom/PinDot Seating Quote/PO # _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***