



Freedom Designs, Inc.

TriPod LT

Adult Manual Wheelchair Transport Package Optional

Our kids are our passion.™ 2241 N. Madera Road Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

Bill to: Date: PO#: Company: Buyer: Marked For:

Ship to: Company: Address: City: State: Zip: Phone:

Table with 6 columns: Description, Part No., Retail, Description, Part No., Retail

Frame Type

- TriPod LT (Rigid) L ATP-30270 \$2315 (K0005)

Front Frame Configuration

- 90 Degree Straight L ATP-30330 n/c



Transport Option

- Transport Package L ATP-99999 \$150 Available with FDI WC19 approved seating or on frame only. No Transport Package L ATP-99999-1 n/c

Frame & Axle Colors

- Sugar Cookie Neon Yellow Yellow Bright Pink Pink Parfait Sugar Plum Vigorous Violet Candy Violet Candy Raspberry Extreme Green Sparkle Green Mystic Teal Blizzard Blue Pacific Blue Blue Blast Electric Blue Galaxy Blue Passion Purple Octane Orange Red Sunset Really Red Fiery Red Twilight Stormy Gray Black

Axle Colors

- Black Yellow Same as Frame

Frame Width

- Specify Frame Width: 12" Wide L ATP-30220 n/c 14" Wide L ATP-30240 n/c 16" Wide L ATP-30260 n/c 18" Wide L ATP-30280 n/c

Frame Length

Specify Frame Length:

- Small (11" - 13") L ATP-30507 n/c Set-up Depth: 11" Deep L ATP-11 n/c 12" Deep L ATP-12 n/c 13" Deep L ATP-13 n/c Medium (12" - 15") L ATP-30508 n/c Set-up Depth: 12" Deep L ATP-12 n/c 13" Deep L ATP-13 n/c 14" Deep L ATP-14 n/c 15" Deep L ATP-15 n/c Long (15" - 18") L ATP-30510 n/c Set-up Depth: 15" Deep L ATP-15 n/c 16" Deep L ATP-16 n/c 17" Deep L ATP-17 n/c 18" Deep L ATP-18 n/c

Shorter seat depth may be achieved through seating system. Seat depth indicates back cane placement.

Seating Options

- Omit Seat Pan or Sling L ATP-39011 n/c Seat Pan - Fixed L ATP-39010 \$100 Seat Pan w/EZ Lite L ATP-39013 \$150 Seat Sling L ATP-39012 n/c

Specify Sling Length: Front of seat rail to back upholstery.

- Seat Sling Length: 11" Deep L ATP-30501 n/c 12" Deep L ATP-30502 n/c 13" Deep L ATP-30503 n/c 14" Deep L ATP-30504 n/c 15" Deep L ATP-30505 n/c 16" Deep L ATP-30506 n/c 17" Deep L ATP-30537 n/c 18" Deep L ATP-30538 n/c

Cushions

- 1" Cushion L ATP-39003 \$85 2" Cushion L ATP-39007 \$95 Includes 4# foam & black knit cover.

Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Posts / Push Handles
Hand Grip Bend Options**

- Straight** - No push handles
With fold over sling back
LATP-30600 n/c
- Specify height from rear seat rail to top of straight cane:
- 10" High** LATP-30610 n/c
 - 11" High** LATP-30611 n/c
 - 12" High** LATP-30612 n/c
 - 13" High** LATP-30613 n/c
 - 14" High** LATP-30614 n/c
 - 15" High** LATP-30615 n/c
- 3/4" Telescoping Push Handles** LATP-30616 **\$175**
(10" - 16") These are push handles only.
Back upholstery does not attach to telescoping handles.
For use with straight back posts only.

- Traditional Bend Push Handles** LATP-30518 n/c
- Ergo Bend Push Handles** LATP-30519 n/c



Traditional Push Handle



Ergo Push Handle

- Specify height from rear seat rail to top of bend on cane:
- Short (19" High)** LATP-30515 n/c
 - Medium (23" High)** LATP-30514 n/c
 - Tall (26" High)** LATP-30516 n/c
 - Push Handle Extension Quick Release** LATP-30517 **\$100**
 - 1 Piece Push Handle Ext. Quick Release** LATP-30526 **\$125**
 - One Piece Push Button Angle Adjustable Stroller Handle** LATP-30527 **\$175**
Available on chairs 14"W and up.

Back canes can achieve an additional height of 6.5" when either quick release push handle extensions are ordered.

Sling Backs

- Omit Sling Back** LATP-39005 n/c
- Sling Back** LATP-39015 n/c
- Tension Adjustable Back** LATP-30008 **\$231**
Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery.

Specify Sling Back Height:
Top of seat to top of back upholstery.

- 10" High** LATP-30910 n/c
- 11" High** LATP-30911 n/c
- 12" High** LATP-30912 n/c
- 13" High** LATP-30913 n/c
- 14" High** LATP-30914 n/c
- 15" High** LATP-30915 n/c
- 16" High** LATP-30916 n/c
- 17" High** LATP-30917 n/c
- 18" High** LATP-30918 n/c


Footrest Options

- Angle Adjustable, Flip-Up Footboard** LATP-36028 **\$125**
(K0040)
Adjustable leg length from 3 1/2" - 12".
Flip up footboard available with straight (non taper) frames only.
- Individual Flip Up/Angle Adjustable Footplates** (K0040) LATP-36030 **\$105**
- Composite Individual Flip Up Footplates** LATP-36031 n/c
Adjustable leg length from 4.5" - 6.5".
All measurements are without cushions.



Angle Adjustable Flip Up Footplate

Flip Up Footboard



- Knee to Heel Measurement** _____ Inches

Casters



- Omit Casters Wheels** LATP-34333 n/c
- 3" Poly Casters** LATP-34888 **\$ 55**
- 3" Light Up Casters** LATP-34111 **\$ 55**
- 4" Poly Casters** LATP-40000 **\$ 65**
- 4" Light Up Casters** LATP-41000 **\$ 99**
- 5" Poly Casters** LATP-34555 n/c
- 5" x 1.5" Poly Casters** LATP-34222 **\$ 65**
- 6" Poly Casters** LATP-34666 n/c
- 6" x 1.5" Poly Casters** LATP-34999 **\$ 85**
- 6" Pneuomatic Casters** LATP-34667 **\$ 75**

Center Wheel Camber Angle

- 0 Degrees** LATP-31010 n/c
- 3 Degrees** LATP-31011 n/c
- 6 Degrees** LATP-31012 n/c



Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Rear Wheels & Wheel Accessories

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** L ATP-31524 **(\$55)**



Mag Wheels

- 20" Pneumatic Mag** L ATP-31527 **n/c**
Only available with 3" casters.
- 22" Pneumatic Mag** L ATP-31528 **n/c**
- 24" Pneumatic Mag** L ATP-31529 **n/c**

Aluminum handrims are standard on 20", 22" & 24" wheels.

To omit handrims, please check below.

24" wheels are not available with small frame small lengths.



Spoke Wheels

- 20" Pneumatic Spoke** L ATP-31540 **\$200**
Only available with 3" casters.
- 22" Pneumatic Spoke** L ATP-31541 **\$200**
- 24" Pneumatic Spoke** L ATP-31542 **\$200**

Aluminum handrims are included on 20", 22" & 24" wheels.

To omit handrims, please check below.

Plastic coated handrims are not available with spoke wheels.

24" wheels are not available with small frame lengths.

Wheel Accessories

- Omit Handrims** L ATP-31006 **(\$15)**
- Plastic Coated Handrims (pr.)** L ATP-31005 **\$117**
Plastic coated handrims are not available with spoke wheels
- Natural Fit Handrims**
Please contact your customer service representative for details.
- Airless Inserts (pr.)** (E2213) L ATP-31000 **\$ 75**
For pneumatic tires only.

Wheel Locks



Specify lever type:

- Push to Lock** L ATP-31501 **n/c**
- Pull to Lock** L ATP-31506 **n/c**
- Wheel Lock Handle Ext.** (E0961) L ATP-31801 **\$ 60**

Armrests—Back Cane Mount

- Omit Flip Back Armrest** L ATP-38209 **(\$20)**

Back Cane Mount

- Flip Back Armrest—13" Long** L ATP-38204 **n/c**
- Flip Back Armrest—16" Long** L ATP-38202 **n/c**
- Locking Flip Back Armrest—13" Long** (E0973) L ATP-38201 **\$235**
- Locking Flip Back Armrest—16" Long** (E0973) L ATP-38200 **\$235**

- Flip Back Armrests—Telescoping** L ATP-38208 **\$100**
- Locking Flip Back Armrests—Telescoping** L ATP-38207 **\$260**

Frame Mount

- Height Adjustable Padded Armrests** (E0973) L ATP-38203 **\$235**
Not available with sling or seat pan.



Flip Back Armrest



Height Adjustable Padded Armrest



Locking Flip Back Armrest

- Desk Length Flat Armrest Pad** L ATP-38213 **n/c**
- Full Length Flat Armrest Pad** L ATP-38214 **n/c**
Only available on 16" armrests.
- Seat Rail to Top of Armrest** _____ **Inches**

Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Accessories

- Clothing Guards** (Fabric) L ATP-CGF \$ 70
- Clothing Guards** (Quick Release - Rigid) L ATP-CGR \$120
- Medical Necessities Pouch** L ATP-23207 \$ 67
Note: Black nylon.
- Footrest Hanger Pads** - Pair SFDI-9HP \$ 90
Note: Wrap around pads made from black neoprene.
- Frame Pads** - Pair 10911-2 \$ 80
Note: Wrap around pads made from black neoprene.
Note: Typically used to pad exposed back post tubes.
- Freedom Backpack** L ATP-10900 \$ 55
Note: Back Nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom Backpack.
- Freedom Canopy - Flip Back** L ATP-10905 \$125
Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.
Note: See offer for ~~FREE~~ Backpack and Canopy on last page of order form.

Wheel & Caster Matrix

Wheel Size	Caster Size	Floor to Front of Seat Rail	Arc Over Rail
20"	4"	14.75"	7"
22"	4"	15.75"	7"
24"	4"	16"	8"
20"	5"	16.25"	7"
22"	5"	17.25"	7"
24"	5"	17.5"	8"
20"	6"	16.5"	7"
22"	6"	17.75"	7"
24"	6"	18"	8"

Maximum Rider Weight—200 lbs.
Frame weight approximately 16.5 lbs.

Notes:

Notes:

Date: _____ PO#: _____

Marked For: _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***

Freedom Designs & Pindot

Are Pleased To Announce
"A Sizzling Summer Sensation"
And For a Limited Time It's...

FREE * FREE * FREE

When ordering a Freedom
Wheelchair and a
complete Freedom or Pindot
seating system including seat,
back, and attachment hardware,
simply provide the seating Quote
or seating PO number
below and receive a
Freedom Canopy and
Freedom Backpack

*** FREE ***



Freedom Backpack with Complete Freedom Seating System LATP-SYS-10900 n/c

Note: Black Nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom backpack.
Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

Provide Freedom/PinDot Seating Quote/PO # _____



Freedom Canopy - Flip-Back with Complete Freedom Seating System LATP-SYS-10905 n/c

Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.
Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

Provide Freedom/PinDot Seating Quote/PO # _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***