



# Pediatric Manual Wheelchair Transport Package Optional

Our kids are our passion.™  
2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

**Bill to:**  
 Date: \_\_\_\_\_ PO#: \_\_\_\_\_  
 Company: \_\_\_\_\_  
 Buyer: \_\_\_\_\_  
 Marked For: \_\_\_\_\_

**Ship to:**  
 Company: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_  
 Phone: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type		
<input type="radio"/> TriPod LT Pediatric (Rigid) (E1237)	LPTP-30270	\$2315
<b>Front Frame Configuration</b>		
<input type="radio"/> 90 Degree Straight	LPTP-30330	n/c



Transport Option		
<input type="radio"/> Transport Package	LPTP-99999	\$150
<small>Available with FDI WC19 approved seating or on frame only.</small>		
<input type="radio"/> No Transport Package	LPTP-99999-1	n/c

Frame & Axle Colors		
<input type="radio"/> Sugar Cookie	<input type="radio"/> Pacific Blue	
<input type="radio"/> Neon Yellow	<input type="radio"/> Blue Blast	
<input type="radio"/> Yellow	<input type="radio"/> Electric Blue	
<input type="radio"/> Bright Pink	<input type="radio"/> Galaxy Blue	
<input type="radio"/> Pink Parfait	<input type="radio"/> Passion Purple	
<input type="radio"/> Sugar Plum	<input type="radio"/> Octane Orange	
<input type="radio"/> Vigorous Violet	<input type="radio"/> Red Sunset	
<input type="radio"/> Candy Violet	<input type="radio"/> Really Red	
<input type="radio"/> Candy Raspberry	<input type="radio"/> Fiery Red	
<input type="radio"/> Extreme Green	<input type="radio"/> Twilight	
<input type="radio"/> Sparkle Green	<input type="radio"/> Stormy Gray	
<input type="radio"/> Mystic Teal	<input type="radio"/> Black	
<input type="radio"/> Blizzard Blue		

Axle Colors		
<input type="radio"/> Black	<input type="radio"/> Yellow	<input type="radio"/> Same as Frame

Frame Widths		
<u>Specify Frame Width:</u>		
<input type="radio"/> 10" Wide	LPTP-30200	n/c
<input type="radio"/> 12" Wide	LPTP-30220	n/c
<input type="radio"/> 14" Wide	LPTP-30240	n/c
<input type="radio"/> 16" Wide	LPTP-30260	n/c
<input type="radio"/> 18" Wide	LPTP-30280	n/c

Frame Length			
<u>Specify Frame Length:</u>			
<input type="radio"/> Small (11" - 13")	LPTP-30507	n/c	
<input type="radio"/> 11" Deep	LPTP-11	n/c	
<input type="radio"/> 12" Deep	LPTP-12	n/c	
<input type="radio"/> 13" Deep	LPTP-13	n/c	

Shorter seat depth may be achieved through seating system. Seat depth indicates back cane placement.

Seating Options			
<input type="radio"/> Omit Seat Pan or Sling	LPTP-39011	n/c	
<input type="radio"/> Seat Pan - Fixed	LPTP-39010	\$100	
<input type="radio"/> Seat Pan w/EZ Lite	LPTP-39013	\$150	
<input type="radio"/> Seat Sling	LPTP-39012	n/c	

<u>Specify Sling Length:</u> <small>Front of seat rail to back upholstery.</small>			
<b>Seat Sling Length:</b>			
<input type="radio"/> 11" Deep	LPTP-30501	n/c	
<input type="radio"/> 12" Deep	LPTP-30502	n/c	
<input type="radio"/> 13" Deep	LPTP-30503	n/c	
<input type="radio"/> 14" Deep	LPTP-30504	n/c	

Cushions			
<input type="radio"/> 1" Cushion	LPTP-39003	\$ 85	
<input type="radio"/> 2" Cushion	LPTP-39007	\$ 95	
<small>Includes 4# foam &amp; black knit cover.</small>			

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles  
Hand Grip Bend Options**

- Straight** - No push handles LPTP-30600 n/c  
With fold over sling back
- Specify height from rear seat rail to top of straight cane:
- 10" High** LPTP-30610 n/c
- 11" High** LPTP-30611 n/c
- 12" High** LPTP-30612 n/c
- 13" High** LPTP-30613 n/c
- 14" High** LPTP-30614 n/c
- 15" High** LPTP-30615 n/c
- 3/4" Telescoping Push Handles** LPTP-30616 **\$175**  
(10" - 15") These are push handles only.  
Back upholstery does not attach to telescoping handles.  
For use with straight back posts only.
- Traditional Bend Push Handles** LPTP-30518 n/c
- Ergo Bend Push Handles\*** LPTP-30519 n/c



- Specify height from rear seat rail to top of bend on cane:
  - Short (19" High)** LPTP-30515 n/c
  - Medium (23" High)** LPTP-30514 n/c
  - Tall (26" High)** LPTP-30516 n/c
  - Push Handle Extension-Quick Release** LPTP-30517 **\$100**
  - One Piece Push Handle Extension - Quick Release** LPTP-30526 **\$125**
  - One Piece Push Button Angle Adjustable Stroller Handle** LPTP-30527 **\$175**  
Available on chairs 14"W and up.
- Back canes can achieve an additional height of 6.5" when either quick release push handle extensions are ordered.

**Back Upholstery Height**

- Omit Sling Back** LPTP-39005 n/c
  - Sling Back\*** LPTP-39015 n/c
  - Tension Adjustable Back\*** LPTP-30008 **\$231**
- Specify Sling Back Height:  
(Measured from top of seat to top of back upholstery)
- 10" High** LPTP-30910 n/c
  - 11" High** LPTP-30911 n/c
  - 12" High** LPTP-30912 n/c
  - 13" High** LPTP-30913 n/c
  - 14" High** LPTP-30914 n/c
  - 15" High** LPTP-30915 n/c
  - 16" High** LPTP-30916 n/c
  - 17" High** LPTP-30917 n/c
  - 18" High** LPTP-30918 n/c

\*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery.

**Footrest Options**

- Flip-Up Footboard** (K0040) LPTP-36028 **\$125**  
Adjustable leg length from 3 1/2" - 12".  
Flip up footboard available with straight (non taper) frames only.
- Individual Flip Up/Angle Adjustable Footplates** (K0040) LPTP-36030 **\$105**
- Composite Individual Flip Up Footplates** LPTP-36031 n/c  
Adjustable leg length from 4 .5" - 6 .5".  
All measurements are without a cushion.



**Knee to Heel Measurement** \_\_\_\_\_ Inches

**Casters**



- Omit Casters Wheels** LPTP-34333 n/c
- 3" Poly Casters** LPTP-34888 **\$ 55**
- 3" Light Up Poly Casters** LPTP-34111 **\$ 55**
- 4" Poly Casters** LPTP-40000 n/c
- 4" Light Up Poly Casters** LPTP-41000 **\$ 99**

**Center Wheel Camber Angle**

- Specify:
- 0 Degrees** LPTP-31010 n/c
  - 3 Degrees** LPTP-31011 n/c
  - 6 Degrees** LPTP-31012 n/c

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Wheels & Wheel Accessories**

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** LPTP-31524 **(\$55)**

**Mag Wheels**



- 16" Pneumatic Mag** LPTP-31530 n/c
- 18" Pneumatic Mag** LPTP-31533 n/c

**Wheel Accessories**

- Aluminum Handrims 16" & 18" (pr.)** LPTP-31007 n/c
- Natural Fit Handrims**  
Please contact your customer service representative for details.
- Airless Inserts (pr.) (E2213)** LPTP-31000 **\$ 75**

**Wheel Locks**



Specify lever type:

- Push to Lock** LPTP-31501 n/c
- Pull to Lock** LPTP-31506 n/c
- Wheel Lock Handle Ext. (E0961)** LPTP-31801 **\$ 60**

**Armrests**

- Omit Flip Back Armrests** LPTP-38209 **(\$20)**
- Back Cane Mount**
- Flip Back Armrests—13" Long** LPTP-38204 n/c
  - Flip Back Armrests—16" Long** LPTP-38202 n/c
  - Locking Flip Back Armrests—13" Long (E0973)** LPTP-38201 **\$235**
  - Locking Flip Back Armrests—16" Long (E0973)** LPTP-38200 **\$235**
  - Flip Back Armrests—Telescoping** LPTP-38208 **\$100**
  - Locking Flip Back Armrests—Telescoping** LPTP-38207 **\$260**

**Frame Mount**

- Height Adjustable Padded Armrests (E0973)** LPTP-38203 **\$235**  
Not available with sling or seat pan.



Flip Back Armrest



Locking Flip Back Armrest



Height Adjustable Padded Armrest

- Desk Length Flat Armrest Pad** LPTP-38213 n/c
- Full Length Flat Armrest Pad** LPTP38214 n/c  
Only available on 16" armrests.
- Seat Rail to Top of Armrest** \_\_\_\_\_ Inches

**Accessories**

- Clothing Guards (Fabric)** LPTP-CGF **\$ 70**
- Clothing Guards (Quick Release - Rigid)** LPTP-CGR **\$120**
- Medical Necessities Pouch** LPTP-23207 **\$ 67**  
Note: Black nylon.
- Footrest Hanger Pads - Pair** SFDI-9HP **\$ 90**  
Note: Wrap around pads made from black neoprene.
- Frame Pads - Pair** 10911-2 **\$ 80**  
Note: Wrap around pads made from black neoprene.  
Note: Typically used to pad exposed back post tubes.
- Freedom Backpack** LPTP-10900 **\$ 55**  
Note: Back nylon only with Freedom Designs logo.  
Note: Monogramming is not available on Freedom Backpack.
- Freedom Canopy - Flip Back** LPTP-10905 **\$125**  
Note: Canopy assembly includes removable mounting hardware.  
Note: Canopy is available in black only with Freedom Designs logo.  
Note: Canopy includes rear flap with zippered storage pocket.  
Note: Canopy includes clear-view window.  
Note: See offer for **FREE** Backpack and Canopy on last page of order form.

Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Wheel & Caster Matrix**

Wheel Size	Caster Size	Floor to Front of Seat Rail	Arc Over Rail
16"	3"	11"	4.5"
18"	3"	11.5"	5.5"
16"	4"	11.5"	4.5"
18"	4"	12"	5.5"

**Maximum Rider Weight—200 lbs.**

Frame weight approximately 16.5 lbs.

**Notes:**

**Notes:**



Date: \_\_\_\_\_ PO#: \_\_\_\_\_

Marked For: \_\_\_\_\_

**\* FREE \* FREE \* FREE \* FREE \* FREE \* FREE \***

# Freedom Designs & Pindot

Are Pleased To Announce  
*"A Sizzling Summer Sensation"*  
And For a Limited Time It's...

**FREE \* FREE \* FREE**

When ordering a Freedom  
Wheelchair and a  
complete Freedom or Pindot  
seating system including seat,  
back, and attachment hardware,  
simply provide the seating Quote  
or seating PO number  
below and receive a  
Freedom Canopy and  
Freedom Backpack

**\* FREE \***



**Freedom Backpack with Complete LPTP-SYS-10900 n/c  
Freedom Seating System**

**Note:** Black Nylon only with Freedom Designs logo.

**Note:** Monogramming is not available on Freedom backpack.

**Note:** Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

Provide Freedom/PinDot Seating Quote/PO # \_\_\_\_\_



**Freedom Canopy - Flip-Back with LPTP-SYS-10905 n/c  
Complete Freedom Seating System**

**Note:** Canopy assembly includes removable mounting hardware.

**Note:** Canopy is available in black only with Freedom Designs logo.

**Note:** Canopy includes rear flap with zippered storage pocket.

**Note:** Canopy includes clear-view window.

**Note:** Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

Provide Freedom/PinDot Seating Quote/PO # \_\_\_\_\_

**\* FREE \* FREE \* FREE \* FREE \* FREE \* FREE \***