



TriPod™ LT

Pediatric Manual Wheelchair Transport Package Optional

Our kids are our passion.™
2241 N. Madera Rd. Simi Valley, CA 93065 Phone: (800) 331-8551 Fax (888) 582-1509 www.freedomdesigns.com

Bill to:
Date: _____ PO#: _____
Company: _____
Buyer: _____
Marked For: _____

Ship to:
Company: _____
Address: _____
City: _____ State: _____ Zip: _____
Phone: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type		
<input type="radio"/> TriPod LT Pediatric (Rigid) (E1237)	LPTP-30270	\$2315
Front Frame Configuration		
<input type="radio"/> 90 Degree Straight	LPTP-30330	n/c



Transport Option		
<input type="radio"/> Transport Package	LPTP-99999	\$150
<small>Available with FDI WC19 approved seating or on frame only.</small>		
<input type="radio"/> No Transport Package	LPTP-99999-1	n/c

Frame & Axle Colors		
<input type="radio"/> Sugar Cookie	<input type="radio"/> Pacific Blue	
<input type="radio"/> Neon Yellow	<input type="radio"/> Blue Blast	
<input type="radio"/> Yellow	<input type="radio"/> Electric Blue	
<input type="radio"/> Bright Pink	<input type="radio"/> Galaxy Blue	
<input type="radio"/> Pink Parfait	<input type="radio"/> Passion Purple	
<input type="radio"/> Sugar Plum	<input type="radio"/> Octane Orange	
<input type="radio"/> Vigorous Violet	<input type="radio"/> Red Sunset	
<input type="radio"/> Candy Violet	<input type="radio"/> Really Red	
<input type="radio"/> Candy Raspberry	<input type="radio"/> Fiery Red	
<input type="radio"/> Extreme Green	<input type="radio"/> Twilight	
<input type="radio"/> Sparkle Green	<input type="radio"/> Stormy Gray	
<input type="radio"/> Mystic Teal	<input type="radio"/> Black	
<input type="radio"/> Blizzard Blue		

Axle Colors		
<input type="radio"/> Black	<input type="radio"/> Yellow	<input type="radio"/> Same as Frame

Frame Widths		
<u>Specify Frame Width:</u>		
<input type="radio"/> 10" Wide	LPTP-30200	n/c
<input type="radio"/> 12" Wide	LPTP-30220	n/c
<input type="radio"/> 14" Wide	LPTP-30240	n/c
<input type="radio"/> 16" Wide	LPTP-30260	n/c
<input type="radio"/> 18" Wide	LPTP-30280	n/c

Frame Length			
<u>Specify Frame Length:</u>			
<input type="radio"/> Small (11" - 13")	LPTP-30507	n/c	
<input type="radio"/> 11" Deep	LPTP-11	n/c	
<input type="radio"/> 12" Deep	LPTP-12	n/c	
<input type="radio"/> 13" Deep	LPTP-13	n/c	

Shorter seat depth may be achieved through seating system. Seat depth indicates back cane placement.

Seating Options			
<input type="radio"/> Omit Seat Pan or Sling	LPTP-39011	n/c	
<input type="radio"/> Seat Pan - Fixed	LPTP-39010	\$100	
<input type="radio"/> Seat Pan w/EZ Lite	LPTP-39013	\$150	
<input type="radio"/> Seat Sling	LPTP-39012	n/c	

<u>Specify Sling Length:</u> <small>Front of seat rail to back upholstery.</small>			
<input type="radio"/> Seat Sling Length:			
<input type="radio"/> 11" Deep	LPTP-30501	n/c	
<input type="radio"/> 12" Deep	LPTP-30502	n/c	
<input type="radio"/> 13" Deep	LPTP-30503	n/c	
<input type="radio"/> 14" Deep	LPTP-30504	n/c	

Cushions			
<input type="radio"/> 1" Cushion	LPTP-39003	\$ 85	
<input type="radio"/> 2" Cushion	LPTP-39007	\$ 95	
<small>Includes 4# foam & black knit cover.</small>			

Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Back Canes / Push Handles
Hand Grip Bend Options**

- Straight** - No push handles LPTP-30600 n/c
With fold over sling back
- Specify height from rear seat rail to top of straight cane:
 - 10" High** LPTP-30610 n/c
 - 11" High** LPTP-30611 n/c
 - 12" High** LPTP-30612 n/c
 - 13" High** LPTP-30613 n/c
 - 14" High** LPTP-30614 n/c
 - 15" High** LPTP-30615 n/c
- 3/4" Telescoping Push Handles** LPTP-30616 **\$175**
(10" - 15") These are push handles only.
Back upholstery does not attach to telescoping handles.
For use with straight back posts only.
- Traditional Bend Push Handles** LPTP-30518 n/c
- Ergo Bend Push Handles*** LPTP-30519 n/c



Traditional Push Handle



Ergo Push Handle

- Specify height from rear seat rail to top of bend on cane:
 - Short (19" High)** LPTP-30515 n/c
 - Medium (23" High)** LPTP-30514 n/c
 - Tall (26" High)** LPTP-30516 n/c
 - Push Handle Extension-Quick Release** LPTP-30517 **\$100**
 - One Piece Push Handle Extension - Quick Release** LPTP-30526 **\$125**
 - One Piece Push Button Angle Adjustable Stroller Handle** LPTP-30527 **\$175**
Available on chairs 14"W and up.
- Back canes can achieve an additional height of 6.5" when either quick release push handle extensions are ordered.

Back Upholstery Height

- Omit Sling Back** LPTP-39005 n/c
 - Sling Back*** LPTP-39015 n/c
 - Tension Adjustable Back*** LPTP-30008 **\$231**
- Specify Sling Back Height:
(Measured from top of seat to top of back upholstery)
- 10" High** LPTP-30910 n/c
 - 11" High** LPTP-30911 n/c
 - 12" High** LPTP-30912 n/c
 - 13" High** LPTP-30913 n/c
 - 14" High** LPTP-30914 n/c
 - 15" High** LPTP-30915 n/c
 - 16" High** LPTP-30916 n/c
 - 17" High** LPTP-30917 n/c
 - 18" High** LPTP-30918 n/c

*Back upholstery height must be chosen when ordering sling or tension adjustable back upholstery.


Footrest Options

- Flip-Up Footboard** (K0040) LPTP-36028 **\$125**
Adjustable leg length from 3 1/2" - 12".
Flip up footboard available with straight (non taper) frames only.
- Individual Flip Up/Angle Adjustable Footplates** (K0040) LPTP-36030 **\$105**
- Composite Individual Flip Up Footplates** LPTP-36031 n/c
Adjustable leg length from 4.5" - 6.5".
All measurements are without a cushion.



Flip Up Footplate

Flip Up Footboard



- Knee to Heel Measurement** _____ Inches

Casters



- Omit Casters Wheels** LPTP-34333 n/c
- 3" Poly Casters** LPTP-34888 **\$ 55**
- 3" Light Up Poly Casters** LPTP-34111 **\$ 55**
- 4" Poly Casters** LPTP-40000 n/c
- 4" Light Up Poly Casters** LPTP-41000 **\$ 99**

Center Wheel Camber Angle

- Specify:
- 0 Degrees** LPTP-31010 n/c
 - 3 Degrees** LPTP-31011 n/c
 - 6 Degrees** LPTP-31012 n/c



Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Wheels & Wheel Accessories

Patented double lock axle receiver included at no charge.

- Omit Wheels & Axles** LPTP-31524 (\$55)

Mag Wheels



- 16" Pneumatic Mag** LPTP-31530 n/c
- 18" Pneumatic Mag** LPTP-31533 n/c

Wheel Accessories

- Aluminum Handrims 16" & 18" (pr.)** LPTP-31007 n/c
- Natural Fit Handrims**
Please contact your customer service representative for details.
- Airless Inserts (pr.) (E2213)** LPTP-31000 \$ 75

Wheel Locks



Specify lever type:

- Push to Lock** LPTP-31501 n/c
- Pull to Lock** LPTP-31506 n/c
- Wheel Lock Handle Ext. (E0961)** LPTP-31801 \$ 60

Armrests

- Omit Flip Back Armrests** LPTP-38209 (\$20)
- Back Cane Mount**
- Flip Back Armrests—13" Long** LPTP-38204 n/c
 - Flip Back Armrests—16" Long** LPTP-38202 n/c
 - Locking Flip Back Armrests—13" Long (E0973)** LPTP-38201 \$235
 - Locking Flip Back Armrests—16" Long (E0973)** LPTP-38200 \$235
 - Flip Back Armrests—Telescoping** LPTP-38208 \$100
 - Locking Flip Back Armrests—Telescoping** LPTP-38207 \$260

Frame Mount

- Height Adjustable Padded Armrests (E0973)** LPTP-38203 \$235
Not available with sling or seat pan.



Flip Back Armrest



Locking Flip Back Armrest



Height Adjustable Padded Armrest

- Desk Length Flat Armrest Pad** LPTP-38213 n/c
- Full Length Flat Armrest Pad** LPTP38214 n/c
Only available on 16" armrests.
- Seat Rail to Top of Armrest** _____ Inches

Accessories

- Clothing Guards (Fabric)** LPTP-CGF \$ 70
- Clothing Guards (Quick Release - Rigid)** LPTP-CGR \$120
- Medical Necessities Pouch** LPTP-23207 \$ 67
Note: Black nylon.
- Footrest Hanger Pads - Pair** SFDI-9HP \$ 90
Note: Wrap around pads made from black neoprene.
- Frame Pads - Pair** 10911-2 \$ 80
Note: Wrap around pads made from black neoprene.
Note: Typically used to pad exposed back post tubes.
- Freedom Backpack** LPTP-10900 \$ 55
Note: Back nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom Backpack.
- Freedom Canopy - Flip Back** LPTP-10905 \$125
Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.
Note: See offer for **FREE** Backpack and Canopy on last page of order form.



Date: _____ PO#: _____

Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Wheel & Caster Matrix

Wheel Size	Caster Size	Floor to Front of Seat Rail	Arc Over Rail
16"	3"	11"	4.5"
18"	3"	11.5"	5.5"
16"	4"	11.5"	4.5"
18"	4"	12"	5.5"

Maximum Rider Weight—200 lbs.

Frame weight approximately 16.5 lbs.

Notes:

Notes:



Date: _____ PO#: _____ Marked For: _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***

Freedom Designs & Pindot

Are Pleased To Announce
"A Sizzling Summer Sensation"
And For a Limited Time It's...


FREE * FREE * FREE

When ordering a Freedom
Wheelchair and a
complete Freedom or Pindot
seating system including seat,
back, and attachment hardware,
simply provide the seating Quote
or seating PO number
below and receive a
Freedom Canopy and
Freedom Backpack

*** FREE ***

-  **Freedom Backpack with Complete LPTP-SYS-10900 n/c
Freedom Seating System**
Note: Black Nylon only with Freedom Designs logo.
Note: Monogramming is not available on Freedom backpack.
Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for backpack offer.

Provide Freedom/PinDot Seating Quote/PO # _____

-  **Freedom Canopy - Flip-Back with LPTP-SYS-10905 n/c
Complete Freedom Seating System**
Note: Canopy assembly includes removable mounting hardware.
Note: Canopy is available in black only with Freedom Designs logo.
Note: Canopy includes rear flap with zippered storage pocket.
Note: Canopy includes clear-view window.
Note: Complete Freedom Designs seating system including solid seat, solid back, and mounting hardware must be ordered on Freedom frame to be eligible for canopy offer.

Provide Freedom/PinDot Seating Quote/PO # _____

*** FREE * FREE * FREE * FREE * FREE * FREE ***