




NXT Mini Pediatric Manual Wheelchair



Thank you for using Freedom Designs / Invacare eForms

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.freedomdesigns.com



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-331-8551



SEAT RAIL EXTENSIONS

NXT-30525 2" Extension - Pair\$45 (K0108)

Specify Back Post Location (if other than specified above) _____ "

Vent Options

11-1100 Back Post Mounted Vent Bracket\$633 (E1029)

For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

Auxillary Device Mounts

11-3100 O2 Holder - Aluminum w/ QR Hdwe\$631 (K0108)

11-4000 O2 Holder -
IV Pole w/ Mounting Hardware \$311 (K0108)

Tilt Activations

NXT-30272 Dual Lever Hand Activated Tilt n/c

NXT-30273 Foot Activated Tilt n/c

Frame Colors

Sugar Cookien/c	Candy Violet.....n/c	Electric Blue.....n/c	Fiery Redn/c
Neon Yellown/c	Candy Raspberry.....n/c	Galaxy Blue.....n/c	Twilightn/c
Yellown/c	Extreme Greenn/c	Passion Purplen/c	Stormy Gray.....n/c
Bright Pink.....n/c	Sparkle Green.....n/c	Octane Orange.....n/c	Blackn/c
Pink Parfait.....n/c	Mystic Teal.....n/c	Red Sunset.....n/c	Color to be
Sugar Plum.....n/c	Pacific Blue.....n/c	Really Red.....n/c	Determined.....n/c

Back Posts

BACK POST OPTIONS

NXT-30545 Fold-Down Back Posts.....n/c

PUSH HANDLE OPTIONS

NXT-30518 90° Bend Back Post Tubes.....n/c

NXT-30544 Straight Back Post Tubes.....n/c

Requires selection of Push Button Angle Adjustable Stroller Handle (NXT-30570 or NXT-30575).

PUSH HANDLE HEIGHT

NXT-30515 Short - 20" Highn/c

NXT-30514 Medium - 23" Highn/c

NXT-30516 Tall - 26" Highn/c

Push Handle Extension Options



All push handle extensions except NXT-30560 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

NXT-30560 Push Handle Extension -
Width Adjustablen/c
Not available with Straight Back Posts.

NXT-30570 Push Button,
Angle-Adjustable Stroller Handle Extension\$121
Available on wheelchair frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. NXT-30570 or NXT-30575 required with Straight Back Post Tubes (NXT-30544)
Not removable, but does allow frame to fold.

NXT-30517 Individual Stroller
Handle Extensions Removable - Pair\$119
Not available with Straight Back Posts.

NXT-30590 Stroller Handle Extension -
Width Adjustable, Removable.....\$103
Not available with Straight Back Posts.

NXT-30517 Individual Stroller Handle
Extensions Removable - Pair\$193
Not available with Straight Back Posts.

Footrests

HANGERS

NXT-30370 75° Fixed Legrest.....n/c

NXT-30372 75° Fixed Legrest - Extended..... \$49
4" longer than standard. Seat rail must be moved 1.5" forward for base frame clearance

NXT-30327 15° Abduction Seat Rail Insert with 75°
Fixed Legrest \$172
Not available with extended legrest.

Mid-Line Mounted Footplates

Abd Mounted Footplates

NXT-30371 85° Fixed Legrest..... n/c

NXT-30373 85° Fixed Legrest - Extended..... \$49
4" longer than standard. Seat rail must be moved 1.5" forward for base frame clearance

Specify Seat Rail to Foot Plate Setting

3"	3.5"	4"	4.5"	5"	5.5"	6"
6.5"	7"	7.5"	Requires Extended Fixed Legrest			
8"	8.5"	9"	Requires Extended Fixed Legrest			
9.5"	10"		Requires Extended Fixed Legrest			

FOOTPLATES & FOOTBOARDS

NXT-36033 Aluminum Angle
Adj. Flip-up Footplates \$205 (K0040)
Does not include Heel Loops, select below.)

NXT-17070 Composite Flip-Upn/c

NXT-36030 Bolteron Angle
Adjustable Flip-Up Footboard n/c (K0040)

NXT-36404 Omit Footrest Plates

FOOTPLATE OPTIONS

NXT-36022 Heel Loops\$49 (E0951)

NXT-36024 Calf Strap - Webbing..... \$49 (K0038)

Armrests

BACK POST MOUNTED ARMRESTS

Omit Armrests

Foam wrap pad STD on all back post mounted armrests.

NXT-38218 8.5" Flip-Back.....n/c

NXT-38204 10.5" Flip-Back.....n/c

NXT-38210 10.5" to 13.5" Flip-Back

Locking, Telescoping\$307 (E0973)

Full length pad on 13.5" or longer armrests only.

NXT-38201 10.5" Flip-Back Locking\$285 (E0973)

FRAME MOUNTED ARMRESTS

NXT-38217 Single Post Height Adjustable,

Removable - Low\$307 (E0973)

NXT-38203 Single Post Height Adjustable,

Removable - Standard\$307 (E0973)

ARM PAD OPTIONS

NXT-38213 Desk Length Padsn/c

NXT-38214 Full Length Padsn/c

13.5" or longer armrests only.

Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"



Casters

NXT-34888 3" x 1" Polyurethanen/c

NXT-34111 3" Polyurethane Light-Up..... \$94

NXT-40000 4" x 1" Polyurethane.....n/c

NXT-43000 4" x 1.4" Poly, Aluminum Mag \$180

NXT-34555 5" x 1" Polyurethane.....n/c

NXT-34222 5" x 1.5" Polyurethane..... \$94

NXT-35050 5" Polyurethane Light-Up..... \$110

NXT-42000 5" x 1.4" Poly, Aluminum Mag \$180

NXT-34666 6" x 1" Polyurethanen/c

NXT-34999 6" x 1.5" Polyurethane..... \$94

NXT-44000 6" x 1.4" Poly, Aluminum Mag \$180

NXT-34444 6" x 2" Polyurethane, Treaded\$105

NXT-34777 7" x 1" Polyurethane.....n/c

NXT-34200 7" x 2" Foam Filled\$215

7" casters not available with 12" mag wheels.



Rear Wheels

MAGS

NXT-31525 12" Mag Wheel.....n/c
12" mag wheels not available with 7" casters.

NXT-31530 16" Mag Wheel.....\$45

NXT-31533 18" Mag Wheel\$45

SPOKE

NXT-31535 16" Spoke Wheel - Black w/ Black Tire..... n/c

NXT-31538 18" Spoke Wheel - Black w/ Black Tire.....n/c

TIRE OPTIONS

NXT-31040 Pneumatic Tire \$85

NXT-31000 Pneumatic Tire
with Airless Insert \$85 (E2213)

NXT-31045 Low Profile Polyurethane Tiren/c (E2218)
12" Mag wheel only, not available with drum brakes

Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Caster Wheel				
	3"	4"	5"	6"	7"
12"	15"-15.5"	16"	18"	18.5"	
16"	15"-15.5"	16.5"-17"	17"	18"	19"
18"	15"-15.5"	17.5"-19"	18"-19"	18"	
Clear Selection					

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".

Wheel Locks

STANDARD

NXT-31501 Push to Lock.....n/c

NXT-31506 Pull to Lockn/c

WHEEL LOCK OPTIONS

NXT-31801 Extension Handles \$81 (E0961)

DRUM LOCK

Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires.
Not Available with 20",22" & 24" Spoke Wheels.

NXT-31508 Foot Activated Wheel Lock.....n/c
Not available with ALL vent and battery trays.

NXT-30275 Hand Activated Brake
and Wheel Lock System\$1000 (E2228)

Anti-Tippers

NXT-31605 Swing-Away Anti-Tip
Assembly Foot Activated \$194 (E0971)
Not available with 16" or 18" wheels and 7" x 2" casters

NXT-31606 Fixed Anti-Tip Assembly \$135 (E0971)



Seat Pans

Only choose one Seat Pan option. ALL seat pans include 2" x 2" notches for growth adjustment.

NON-ADJUSTABLE

NXT-39010 Fixed Solid
Seat Pan - Composite.....\$183 (E2231)

ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width _____" X Seat Depth _____"

NXT-39030 Fixed Solid Seat Pan - Aluminum\$215 (E2231)
Fixed Seat Pan raises the effective seat to floor height by 1/2".

Available in the following sizes ONLY:

- 14" wide x 16" deep
- 18" wide x 18" deep
- 16" wide x 18" deep

ADJUSTABLE

NXT-39040 Fixed, Width Adjustable
Solid Seat Pan - Aluminum\$312 (E2231)

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with NXT frame.

Available in the following sizes ONLY:

- 10"-12" wide x 12" deep
- 12"-14" wide x 14" deep
- 14"-16" wide x 16" deep

MOUNTING HARDWARE OPTIONS

Use to convert fixed seat pan to removable.

NXT-EZLITE Removable Pan Hardware Kit.....\$223 (K0108)

Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

Accessories

NXT-23210 Flip-Down Rigidizing Bar.....\$189 (K0108)

Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.

NXT-10900 Freedom Backpack.....\$65 (K0108)

Black Nylon only with Freedom Designs logo. Monogramming is not available on Freedom Backpack.

Specifications

Frame Style	Tilt-In-Space
Base Frame Widths	14" - 20"
Overall Width	Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel)
Seat Widths	14" - 20"
Seat Depths	14" - 26"
Back Heights	20", 23", 26"
Adjustable Angle Back	90°, 100°, 110°, 120°
Tilt Range	0° to 55°
Recline Back	N/A
Weight Limit	225 LBS.
Transport Weight	Approximately 22 lbs
Warranty	5 year warranty on side frames and quick release axles

Special Configurations Through Customer Service

Asymmetrical Seat Depth	Seat rails may be configured up to 4" offset
-------------------------	--

HCPCS Codes

are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's

Special Instructions:



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-331-8551