



NXT Mini Pediatric Manual Wheelchair



Thank you for using Freedom Designs / Invacare eForms

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.freedomdesigns.com



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-331-8551



Account Information

Request Type: Quote Order

Date: _____

Purchase Order #: _____

Account #: _____

CONTACT

Name: _____

Company: _____

Back Up Contact: _____

SHIP TO

Name: _____

Phone: _____

Address: _____

Email: _____

City: _____ State: _____ Zip: _____

Comments: _____

Country: _____

MF: Please do not provide with a customer name

Base Options

A seating system must be purchased separately. No seating system is provided with the NXT frame. No seat or back upholstery is available as an option.

TRANSPORT OPTIONAL

NXT-30330 NXT-Mini Pediatric Folding.....\$2850 (E1234)

NXT-30340 NXT-Mini Pediatric Rigid..... \$2725 (E1233)

Transportation Option selected separately below.
Includes Freedom Future Fit Program

NXT-99998 Transport Brackets.....\$300

Must be installed by Freedom Designs manufacturing.
Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.
Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

TRANSPORT INCLUDED

NXT-30331 NXT-Mini Pediatric Folding.....\$2995 (E1234)

NXT-30341 NXT-Mini Pediatric Rigid \$2875 (E1233)

Includes Freedom Future Fit Program

Frame Dimensions

FRAME WIDTH SETTING & CROSS BRACE

NXT-30250n/c

NXT-30220n/c

Adjusts 10" to 12"

Adjusts 12" to 14"

10" 11" 12"

12" 13" 14"

FRAME DEPTH & BACK POST SETTING

Mini Frame, X-Tall Tiltn/c

Requires 2" Extension

Seat rail to floor height range is 15" to 19". All vent and battery trays not available with Mini frame. All dual post armrests not available with short frame.

8" 11"

14"

9" 12"

Requires 2" Extension

10" 13"



SEAT RAIL EXTENSIONS

NXT-30525 2" Extension - Pair.....\$41 (K0108)

Specify Back Post Location (if other than specified above) _____”

Vent Options

11-1100 Back Post Mounted Vent Bracket..... \$550 (E1029)

For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

Auxillary Device Mounts

11-3100 O2 Holder - Aluminum w/ QR Hdwe....\$548 (K0108)

11-4000 O2 Holder - IV Pole w/ Mounting Hardware..... \$270 (K0108)

Tilt Activations

NXT-30272 Dual Lever Hand Activated Tilt.....n/c

NXT-30273 Foot Activated Tilt.....n/c

Frame Colors

- | | | | |
|-----------------------|------------------------|-------------------------|---------------------------------|
| Sugar Cookie..... n/c | Candy Violet..... n/c | Electric Blue..... n/c | Fiery Red..... n/c |
| Neon Yellow..... n/c | Candy Raspberry n/c | Galaxy Blue..... n/c | Twilight..... n/c |
| Yellow..... n/c | Extreme Green..... n/c | Passion Purple..... n/c | Stormy Gray..... n/c |
| Bright Pink..... n/c | Sparkle Green..... n/c | Octane Orange..... n/c | Black..... n/c |
| Pink Parfait..... n/c | Pacific Blue..... n/c | Red Sunset..... n/c | Color to be Determined..... n/c |
| Sugar Plum..... n/c | Blue Blast..... n/c | Really Red..... n/c | |

Back Posts

BACK POST OPTIONS

NXT-30545 Fold-Down Back Posts.....n/c

PUSH HANDLE OPTIONS

NXT-30518 90° Bend Back Post Tubes.....n/c

NXT-30544 Straight Back Post Tubes.....n/c

Requires selection of Push Button Angle Adjustable Stroller Handle (NXT-30570 or NXT-30575).

PUSH HANDLE HEIGHT

NXT-30515 Short - 20" High.....n/c

NXT-30514 Medium - 23" High.....n/c

NXT-30516 Tall - 26" High.....n/c

Push Handle Extension Options



All push handle extensions except NXT-30560 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

NXT-30560 Push Handle Extension - Width Adjustablen/c
Not available with Straight Back Posts.

NXT-30570 Push Button, Angle-Adjustable Stroller Handle Extension \$112
Available on wheelchair frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. NXT-30570 or NXT-30575 required with Straight Back Post Tubes (NXT-30544)
Not removable, but does allow frame to fold.

NXT-30517 Individual Stroller Handle Extensions Removable - Pair \$110
Not available with Straight Back Posts.

NXT-30590 Stroller Handle Extension - Width Adjustable, Removable \$95
Not available with Straight Back Posts.

NXT-30517 Individual Stroller Handle Extensions Removable - Pair \$110
Not available with Straight Back Posts.

Footrests

HANGERS

NXT-30370 75° Fixed Legrest.....n/c

NXT-30372 75° Fixed Legrest - Extended..... \$45
4" longer than standard. Seat rail must be moved 1.5" forward for base frame clearance

NXT-30327 15° Abduction Seat Rail Insert with 75° Fixed Legrest.....\$160
Not available with extended legrest.

Mid-Line Mounted Footplates

Abd Mounted Footplates

NXT-30371 85° Fixed Legrest.....n/c

NXT-30373 85° Fixed Legrest - Extended.....\$45
4" longer than standard. Seat rail must be moved 1.5" forward for base frame clearance

Specify Seat Rail to Foot Plate Setting

3"	3.5"	4"	4.5"	5"	5.5"	6"
6.5"	7"	7.5"	Requires Extended Fixed Legrest			
8"	8.5"	9"	Requires Extended Fixed Legrest			
9.5"	10"		Requires Extended Fixed Legrest			

FOOTPLATES & FOOTBOARDS

NXT-36033 Aluminum Angle Adj. Flip-up Footplates..... \$190 (K0040)
Does not include Heel Loops, select below.)

NXT-17070 Composite Flip-Up.....n/c

NXT-36030 Bolteron Angle Adjustable Flip-Up Footboard n/c (K0040)

NXT-36404 Omit Footrest Plates

FOOTPLATE OPTIONS

NXT-36022 Heel Loops.....\$45 (E0951)

NXT-36024 Calf Strap - Webbing.....\$45 (K0038)

Armrests

BACK POST MOUNTED ARMRESTS

Omit Armrests

Foam wrap pad STD on all back post mounted armrests.

NXT-38218 8.5" Flip-Back.....n/c

NXT-38204 10.5" Flip-Back.....n/c

NXT-38210 10.5" to 13.5" Flip-Back

Locking, Telescoping\$285 (E0973)

Full length pad on 13.5" or longer armrests only.

NXT-38201 10.5" Flip-Back Locking\$265 (E0973)

FRAME MOUNTED ARMRESTS

NXT-38217 Single Post Height Adjustable,

Removable - Low\$285 (E0973)

NXT-38203 Single Post Height Adjustable,

Removable - Standard\$285 (E0973)

ARM PAD OPTIONS


NXT-38213 Desk Length Padsn/c

NXT-38214 Full Length Pads.....n/c

13.5" or longer armrests only.

Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"



Casters

NXT-34888 3" x 1" Polyurethane.....n/c

NXT-43111 3" Polyurethane Light-Up..... \$87

NXT-40000 4" x 1" Polyurethane.....n/c

NXT-43000 4" x 1.4" Poly, Aluminum Mag\$125

NXT-34555 5" x 1" Polyurethane.....n/c

NXT-34222 5" x 1.5" Polyurethane..... \$87

NXT-35050 5" Polyurethane Light-Up..... \$101

NXT-42000 5" x 1.4" Poly, Aluminum Mag\$125

NXT-34666 6" x 1" Polyurethane.....n/c

NXT-34999 6" x 1.5" Polyurethane..... \$87

NXT-44000 6" x 1.4" Poly, Aluminum Mag\$125

NXT-34444 6" x 2" Polyurethane, Treaded\$97

NXT-34777 7" x 1" Polyurethane.....n/c

NXT-34200 7" x 2" Foam Filled.....\$200

7" casters not available with 12" mag wheels.



NXT Mini Pediatric Manual Wheelchair

Rear Wheels

MAGS

NXT-31525 12" Mag Wheel.....n/c
12" mag wheels not available with 7" casters.

NXT-31533 18" Mag Wheel.....\$41

NXT-31530 16" Mag Wheel..... \$41

SPOKE

NXT-31535 16" Spoke Wheel - Black w/ Black Tiren/c

NXT-31538 18" Spoke Wheel - Black w/ Black Tiren/c

TIRE OPTIONS

NXT-31040 Pneumatic Tire..... \$85

NXT-31000 Pneumatic Tire with Airless Insert\$85 (E2213)

NXT-31045 Low Profile Polyurethane Tiren/c (E2218)
12" Mag wheel only, not available with drum brakes

Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Caster Wheel				
	3"	4"	5"	6"	7"
12"	15"-15.5"	16"	18"	18.5"	
16"	15"-15.5"	16.5"-17"	17"	18"	19"
18"	15"-15.5"	17.5"-19"	18"-19"	18"	
Clear Selection					

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".

Wheel Locks

STANDARD

NXT-31501 Push to Lockn/c

NXT-31506 Pull to Lock.....n/c

WHEEL LOCK OPTIONS

NXT-31801 Extension Handles \$75 (E0961)

DRUM LOCK

Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires. Not Available with 20",22" & 24" Spoke Wheels.

NXT-31508 Foot Activated Wheel Lock.....\$200
Not available with ALL vent and battery trays.

NXT-30275 Hand Activated Brake and Wheel Lock Systemn/c (E2228)

Anti-Tippers

NXT-31605 Swing-Away Anti-Tip Assembly Foot Activated\$180 (E0971)
Not available with 16" or 18" wheels and 7" x 2" casters

NXT-31606 Fixed Anti-Tip Assembly \$125 (E0971)



Yes, you can.



NXT Mini Pediatric Manual Wheelchair

Seat Pans

ALL seat pans include 2" x 2" notches for growth adjustment.

NON-ADJUSTABLE

NXT-39010 Fixed Solid
Seat Pan - Composite.....\$170 (E2231)

ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width _____" X Seat Depth _____"

NXT-39030 Fixed Solid Seat Pan - Aluminum..... \$200 (E2231)
Fixed Seat Pan raises the effective seat to floor height by 1/2".

Available in the following sizes ONLY:

- 14" wide x 16" deep
- 18" wide x 18" deep
- 16" wide x 18" deep

ADJUSTABLE

NXT-39040 Fixed, Width Adjustable
Solid Seat Pan - Aluminum..... \$290 (E2231)

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with NXT frame.

Available in the following sizes ONLY:

- 10"-12" wide x 12" deep
- 12"-14" wide x 14" deep
- 14"-16" wide x 16" deep

MOUNTING HARDWARE OPTIONS

Use to convert fixed seat pan to removable.

NXT-EZLITE Removable Pan Hardware Kit.....\$207 (K0108)
Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".
Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

Accessories

NXT-23210 Flip-Down Rigidizing Bar \$175 (K0108)
Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.

NXT-10900 Freedom Backpack \$60 (K0108)
Black Nylon only with Freedom Designs logo. Monogramming is not available on Freedom Backpack.

NXT-10905 Freedom Canopy - Flip-Back.....n/c (K0108)
Canopy assembly includes removable mounting hardware. Available in black only with Freedom Designs logo. Includes rear flap with zippered storage pocket and clear-view window.

Specifications

Frame Style	Tilt-In-Space
Base Frame Widths	14" - 20"
Overall Width	Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel)
Seat Widths	14" - 20"
Seat Depths	14" - 26"
Back Heights	20", 23", 26"
Adjustable Angle Back	90°, 100°, 110°, 120°
Tilt Range	0° to 55°
Recline Back	N/A
Weight Limit	225 LBS.
Transport Weight	
Warranty	Lifetime on frame

Special Configurations Through Customer Service

Asymmetrical Seat Depth	Seat rails may be configured up to 4" offset
-------------------------	--

HCPCS Codes

are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-331-8551