



# NXT<sup>TM</sup> Mini Pediatric Manual Wheelchair

Generation NEXT

Clear Form

Email Form to Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: \_\_\_\_\_ Buyer: \_\_\_\_\_ Quote Order Install Seating? PinDot Matrix  
 Company: \_\_\_\_\_ PO#: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Location: \_\_\_\_\_  
 Email: \_\_\_\_\_ Phone: \_\_\_\_\_

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
-------------	-------	----------	--------	-------------	-------	----------	--------

## Base Options

A seating system must be purchased separately. No seating system is provided with the NXT frame. No seat or back upholstery is available as an option.

### Transport Optional

NXT-Mini Pediatric Folding	(E1234)	NXT-30330	\$2850
NXT-Mini Pediatric Rigid	(E1233)	NXT-30340	\$2725

Transportation Option selected separately below.  
 Includes Freedom Future Fit Program

Transport Brackets		NXT-99998	\$ 300
--------------------	--	-----------	--------

Must be installed by Freedom Designs manufacturing.

Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.

Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

### Transport Included

NXT-Mini Pediatric Folding	(E1234)	NXT-30331	\$2995
NXT-Mini Pediatric Rigid	(E1233)	NXT-30341	\$2875

Includes Freedom Future Fit Program

## Frame Dimensions

### Frame Width Setting & Cross Brace

10"	11"	12"	Adjusts 10" to 12"	NXT-30250	n/c
12"	13"	14"	Adjusts 12" to 14"	NXT-30220	n/c

### Frame Depth & Back Post Setting

Mini Frame, X-Tall Tilt Seat rail to floor height range is 15" to 19". All vent and battery trays not available with Mini frame. All dual post armrests not available with short frame.

8" n/c	9" n/c	10" n/c	11" n/c	12" n/c	13" n/c	14" n/c
					Requires 2" Extension	Requires 2" Extension

## Frame Dimensions - Continued

### Seat Rail Extensions

2" Extension - Pair	(K0108)	NXT-30525	\$ 41	Specify Back Post Location (if other than specified above) _____"
---------------------	---------	-----------	-------	---

## Vent Options

Back Post Mounted Vent Bracket	(E1029)	11-1100	\$ 550
--------------------------------	---------	---------	--------

For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

## Auxillary Device Mounts

O2 Holder - Aluminum w/ QR Hdwe	(K0108)	11-3100	\$ 548	IV Pole w/ Mounting Hardware	(K0108)	11-4000	\$ 270
---------------------------------	---------	---------	--------	------------------------------	---------	---------	--------

## Tilt Activations

Dual Lever Hand Activated Tilt	NXT-30272	n/c	Foot Activated Tilt	NXT-30273	n/c
--------------------------------	-----------	-----	---------------------	-----------	-----

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
<b>Frame Colors</b>							
Sugar Cookie			n/c	Blue Blast			n/c
Neon Yellow			n/c	Electric Blue			n/c
Yellow			n/c	Galaxy Blue			n/c
Bright Pink			n/c	Passion Purple			n/c
Pink Parfait			n/c	Octane Orange			n/c
Sugar Plum			n/c	Red Sunset			n/c
Candy Violet			n/c	Really Red			n/c
Candy Raspberry			n/c	Fiery Red			n/c
Extreme Green			n/c	Twilight			n/c
Sparkle Green			n/c	Stormy Gray			n/c
Mystic Teal			n/c	Black			n/c
Pacific Blue			n/c	Color to be Determined			

**Back Posts**

<b>Back Post Options</b>							
Fold-Down Back Posts		NXT-30545	n/c				

<b>Push Handle Options</b>							
90° Bend Back Post Tubes		NXT-30518	n/c	Straight Back Post Tubes		NXT-30544	n/c
				Requires selection of Push Button Angle Adjustable Stroller Handle (NXT-30570 or NXT-30575).			

<b>Push Handle Height</b>							
Short - 20" High		NXT-30515	n/c	Tall - 26" High		NXT-30516	n/c
Medium - 23" High		NXT-30514	n/c				

**Push Handle Extension Options**



All push handle extensions except NXT-30560 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

Push Handle Extension - Width Adjustable		NXT-30560	n/c	Stroller Handle Extension - Width Adjustable, Removable		NXT-30590	\$ 95
Not available with Straight Back Posts.				Not available with Straight Back Posts.			
Push Button, Angle-Adjustable Stroller Handle Extension		NXT-30570	\$ 112	Individual Stroller Handle Extensions Removable - Pair		NXT-30517	\$ 110
Available on wheelchair frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. NXT-30570 or NXT-30575 required with Straight Back Post Tubes (NXT-30544). Not removable, but does allow frame to fold.				Not available with Straight Back Posts.			

**Footrests**

<b>Hangers</b>							
75° Fixed Legrest		NXT-30370	n/c	85° Fixed Legrest		NXT-30371	n/c
75° Fixed Legrest - Extended		NXT-30372	\$ 45	85° Fixed Legrest - Extended		NXT-30373	\$ 45

4" longer than standard. Seat rail must be moved 1.5" forward for base frame clearance

4" longer than standard. Seat rail must be moved 1.5" forward for base frame clearance

15° Abduction Seat Rail Insert with 75° Fixed Legrest		NXT-30327	\$ 160
---	--	-----------	--------

Not available with extended legrest.

<b>Specify Seat Rail to Foot Plate Setting</b>									
	3"	3.5"	4"	4.5"	5"	5.5"	6"		
	6.5"	7"	7.5"	Requires Extended Fixed Legrest					
	8"	8.5"	9"	Requires Extended Fixed Legrest					
	9.5"	10"		Requires Extended Fixed Legrest					

Mid-Line Mounted Footplates Abd Mounted Footplates

**Footplates & Footboards**

Aluminum Angle Adj. Flip-up Footplates	(K0040)	NXT-36033	\$ 190	Bolteron Angle Adjustable Flip-Up Footboard	(K0040)	NXT-36030	\$ 190
Does not include Heel Loops, select below.							
Composite Flip-Up		NXT-17070	n/c	Omit Footrest Plates		NXT-36404	

<b>Footplate Options</b>							
Heel Loops	(E0951)	NXT-36022	\$ 45	Calf Strap - Webbing	(K0038)	NXT-36024	\$ 45

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
-------------	-------	----------	--------	-------------	-------	----------	--------

### Armrests

#### Back Post Mounted Armrests

Omit Armrests

Foam wrap pad STD on all back post mounted armrests..

8.5" Flip-Back		NXT-38218	n/c				
10.5" Flip-Back		NXT-38204	n/c	10.5" Flip-Back Locking	(E0973)	NXT-38201	\$ 265
10.5" to 13.5" Flip-Back Locking, Telescoping	(E0973)	NXT-38210	\$ 285				

Full length pad on 13.5" or longer armrests only.

#### Frame Mounted Armrests

Single Post Height Adjustable, Removable - Low	(E0973)	NXT-38217	\$ 285	Single Post Height Adjustable, Removable - Standard	(E0973)	NXT-38203	\$ 285
--	---------	-----------	--------	---	---------	-----------	--------

#### Arm Pad Options


Desk Length Pads		NXT-38213	n/c	Full Length Pads		NXT-38214	n/c
------------------	--	-----------	-----	------------------	--	-----------	-----

13.5" or longer armrests only.

#### Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"

All heights from top of seat frame to top of armrest pad.



#### Casters

3" x 1" Polyurethane		NXT-34888	n/c	3" Polyurethane Light-Up		NXT-34111	\$ 87
4" x 1" Polyurethane		NXT-40000	n/c	4" x 1.4" Poly, Aluminum Mag		NXT-43000	\$ 125
5" x 1" Polyurethane		NXT-34555	n/c	5" Polyurethane Light-Up		NXT-35050	\$ 101
5" x 1.5" Polyurethane		NXT-34222	\$ 87	5" x 1.4" Poly, Aluminum Mag		NXT-42000	\$ 125
6" x 1" Polyurethane		NXT-34666	n/c	6" x 1.4" Poly, Aluminum Mag		NXT-44000	\$ 125
6" x 1.5" Polyurethane		NXT-34999	\$ 87	6" x 2" Polyurethane, Treaded		NXT-34444	\$ 97
7" x 1" Polyurethane		NXT-34777	n/c	7" x 2" Foam Filled		NXT-34200	\$ 200

7" casters not available with 12" mag wheels.

#### Rear Wheels

##### Mags

12" Mag Wheel		NXT-31525	n/c
16" Mag Wheel		NXT-31530	\$ 41
18" Mag Wheel		NXT-31533	\$ 41

12" mag wheels not available with 7" casters.

##### Spoke

16" Spoke Wheel - Black w/ Black Tire	NXT-31535	n/c
18" Spoke Wheel - Black w/ Black Tire	NXT-31538	n/c

#### Tire Options

Pneumatic Tire		NXT-31040	\$ 85
Pneumatic Tire with Airless Insert	(E2213)	NXT-31000	\$ 85
Low Profile Polyurethane Tire	(E2218)	NXT-31045	n/c

12" Mag wheel only, not available with drum brakes

#### Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Caster Wheel				
	3"	4"	5"	6"	7"
12"	15"-15.5"	16"	18"	18.5"	
16"	15.5"-16"	16.5"-17"	17"	18"	19"
18"	15"-15.5"	17.5"-19"	18"-19"	18"	

Clear selection

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".

#### Wheel Locks

##### Standard

Push to Lock		NXT-31501	n/c
Pull to Lock		NXT-31506	n/c

##### Wheel Lock Options

Extension Handles	(E0961)	NXT-31801	\$ 75
-------------------	---------	-----------	-------

##### Drum Lock

Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires.			
Foot Activated Wheel Lock		NXT-31508	\$ 200
Hand Activated Brake and Wheel Lock System	(E2228)	NXT-30275	\$ 1000

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
<b>Anti-Tippers</b>							
Swing-Away Anti-Tip Assembly Foot Activated	(E0971)	NXT-31605	\$ 180	Fixed Anti-Tip Assembly	(E0971)	NXT-31606	\$ 125
Not available with 16" or 18" wheels and 7" x 2" casters							

**Seat Pans**

ALL seat pans include 2" x 2" notches for growth adjustment.

**Non-Adjustable**

Fixed Solid Seat Pan - Composite (E2231) NXT-39010 \$ 170

ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width \_\_\_\_\_" X Seat Depth \_\_\_\_\_"

Fixed Solid Seat Pan - Aluminum (E2231) NXT-39030 \$ 200

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Available in the following sizes ONLY:

10" wide x 12" deep                      14" wide x 16" deep

12" wide x 14" deep

**Adjustable**

Fixed, Width Adjustable Solid Seat Pan - Aluminum (E2231) NXT-39040 \$ 290

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with NXT frame.

Available in the following sizes ONLY:

10"-12" wide x 12" deep                      14"-16" wide x 16" deep

12"-14" wide x 14" deep

**Mounting Hardware Options**

Use to convert fixed seat pan to removable.

Removable Pan Hardware Kit (K0108) NXT-EZLITE \$ 207

Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

**Accessories**

Flip-Down Rigidizing Bar (K0108) NXT-23210 \$ 175

Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.

Freedom Backpack (K0108) NXT-10900 \$ 60

Black Nylon only with Freedom Designs logo. Monogramming is not available on Freedom Backpack.

Freedom Canopy - Flip-Back (K0108) NXT-10905 n/c

Canopy assembly includes removable mounting hardware. Available in black only with Freedom Designs logo. Includes rear flap with zippered storage pocket and clear-view window.

**Notes:**

**Specifications:**

Frame Style	Tilt-In-Space
Base Frame Widths	14' - 20"
Overall Width	Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel)
Seat Widths	14" - 20"
Seat Depths	14" - 26"
Back Heights	20", 23", 26"
Adjustable Angle Back	90°, 100°, 110°, 120°
Tilt Range	0° to 55°
Recline Back	N/A
Weight Limit	225 LBS.
Transport Weight	
Warranty	Lifetime on frame

**Special Configurations Through Customer Service:**

Asymmetrical Seat Depth      Seat rails may be configured up to 4" offset

**HCPCS Codes:**

are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.