



Yes, you can.®

# NXT-SE Pediatric Manual Wheelchair



## Thank you for using Freedom Designs / Invacare eForms

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email\* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit [www.freedomdesigns.com](http://www.freedomdesigns.com)



### Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**800-331-8551**



# Account Information

Request Type:      Quote                      Order

Date: \_\_\_\_\_

Account #: \_\_\_\_\_

Company: \_\_\_\_\_

SHIP TO  
Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_ State: \_\_\_\_\_ Zip: \_\_\_\_\_

Country: \_\_\_\_\_

Purchase Order #: \_\_\_\_\_

CONTACT  
Name: \_\_\_\_\_

Back Up Contact: \_\_\_\_\_

Phone: \_\_\_\_\_

Email: \_\_\_\_\_

Comments: \_\_\_\_\_

MF: Please do not provide with a customer name

# Base Options

A seating system must be purchased separately. No seating system is provided with the NXT frame. No seat or back upholstery is available as an option.

## TRANSPORT OPTIONAL

- NXT-30304 NXT-SE Pediatric Folding ..... \$2700 (E1234)
  - NXT-30314 NXT-SE Pediatric Rigid ..... \$2550 (E1233)
- Transportation Option selected separately below. SE Frames Do Not include Freedom Future Fit Program

- NXT-99997 Transport Brackets - SE ..... \$300

Must be installed by Freedom Designs manufacturing.  
Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.  
Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

# Frame Dimensions

## FRAME WIDTH SETTING & CROSS BRACE

NXT-30250 ..... n/c	Adjusts 10" to 12"		
	10"	11"	12"
NXT-30210 ..... n/c	Adjusts 11" to 13"		
	11"	12"	13"
NXT-30220 ..... n/c	Adjusts 12" to 14"		
	12"	13"	14"
NXT-30230 ..... n/c	Adjusts 13" to 15"		
	13"	14"	15"
NXT-30240 ..... n/c	Adjusts 14" to 16"		
	14"	15"	16"
NXT-30290 ..... n/c	Adjusts 15" to 17"		
	15"	16"	17"

NXT-30260 ..... n/c	Adjusts 16" to 18"		
	16"	17"	18"
NXT-30261 ..... n/c	Adjusts 17" to 18"		
	17"	18"	
NXT-30262 ..... n/c	Not width adjustable		
	18"		
NXT-30263 ..... n/c	Not with short frame		
	19"		
NXT-30264 ..... \$500 (E2201)	Not with short frame		
	20"		

Frame widths measures across and from outside of seat frame tubes. Vent and battery trays only with 14" wide and greater frame.



FRAME DEPTH & BACK POST SETTING

Short Frame, Low Tilt

Seat rail to floor height range is 12.5" to 16.5". Seat frame widths of 19" and 20" not available with short frame. 22" & 24" wheels not available with low tilt bracket. All vent and battery trays not available with short frame. All dual post armrests not available with short frame.

- 12" n/c
13" n/c
14" n/c
15" n/c
16" n/c
17" n/c Requires 2" Extension
18" n/c Requires 2" Extension

Medium Frame, Low Tilt

Seat rail to floor height range is 12.5" to 16.5". 22" & 24" wheels not available with low tilt

- 14" n/c
15" n/c
16" n/c
17" n/c
18" n/c
19" n/c Requires 2" Extension
20" \$550 (E2203) Includes 2" Extension
21" \$550 (E2203) Includes 4" Extension
22" \$990 (E2204) Includes 4" Extension

Medium Frame, Tall Tilt

Seat rail to floor height range is 14" to 18".

- 14" n/c
15" n/c
16" n/c
17" n/c
18" n/c
19" n/c Requires 2" Extension
20" \$550 (E2203) Includes 2" Extension
21" \$550 (E2203) Includes 4" Extension
22" \$990 (E2204) Includes 4" Extension

Pediatric Long Frame, Tall Tilt \$153 (E1298)

Seat rail to floor height range is 14" to 18".

- 16" n/c
17" n/c
18" n/c
19" n/c
20" \$550 (E2203)
21" \$550 (E2203)
22" \$990 (E2204)
23" \$990 (E2204) Includes 2" Extension
24" \$990 (E2204) Includes 2" Extension
25" \$990 (E2204) Includes 4" Extension
26" \$990 (E2204) Includes 4" Extension

SEAT RAIL EXTENSIONS

NXT-30525 2" Extension - Pair \$41 (K0108)

NXT-30531 4" Extension - Pair \$41 (K0108)

Specify Back Post Location (if other than specified above) "



## Seat Rail Options

### LEAF - LOWER EXTREMITY ACCOMMODATION FRAME

Not available with ALL Contracture & 3D Dynamic footrest systems.

NXT-30324 LEAF Seat Frame- Pair .....\$775 (K0108)

Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

NXT-30276 5° Anterior Tilt Seat Frame.....n/c

Limits posterior tilt to 40°. Includes split collars for adjustment. Medium with Tall Tilt and Long frames only.

### SRO - SEAT RAIL ONLY

NXT-30325 SRO Seat Frame- Pair .....\$325 (K0108)

SRO is manufactured without standard footrest hanger attachments.

## Vent Options

### LEAF - LOWER EXTREMITY ACCOMMODATION FRAME

Not available with ALL Contracture & 3D Dynamic footrest systems.

11-7000 Vent Tray - Fixed .....\$461 (E1029)

Not with Foot activated tilt. For LP10 style vent, base must be 16" or greater. Tray measures 4" narrower than the width of the base frame and 16" deep. All vent and battery trays not available with short frame.

11-1100 Back Post Mounted Vent Bracket .....\$550 (E1029)

For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

11-8000 Battery Tray - Fixed .....\$450 (K0108)

Tray measures 4" narrower than the width of the base frame and 10" deep. All vent and battery trays not available with short frame.

563 Quick Release Tray Hardware .....\$286 (K0108)

## Auxillary Device Mounts

11-3100 O2 Holder - Aluminum w/ QR Hdwe.....\$548 (E0108)

11-4000 IV Pole w/ Mounting Hardware .....\$270 (K0108)

## Tilt Activations

NXT-30272 Dual Lever Hand Activated Tilt.....n/c

NXT-30273 Foot Activated Tilt .....n/c

Not available with vent tray.

## Frame Colors

Sugar Cookie .....n/c

Candy Violet.....n/c

Electric Blue.....n/c

Fiery Red .....n/c

Neon Yellow .....n/c

Candy Raspberry.....n/c

Galaxy Blue.....n/c

Twilight .....n/c

Yellow .....n/c

Extreme Green .....n/c

Passion Purple .....n/c

Stormy Gray.....n/c

Bright Pink.....n/c

Sparkle Green.....n/c

Octane Orange.....n/c

Black .....n/c

Pink Parfait.....n/c

Mystic Teal.....n/c

Red Sunset.....n/c

Color to be

Sugar Plum.....n/c

Pacific Blue.....n/c

Really Red.....n/c

Determined.....n/c

## Back Posts

### BACK POST OPTIONS

NXT-30545 Fold-Down Back Posts.....n/c

NXT-30532 Dynamic Rocker Back - Seating Dynamics..... \$1985 (K0108)  
Angle adjustable from 72° to 108° in 6° increments. Can be locked-out when needed. Includes a set of 4 elastomers in different durometers to allow spring adjustment in the field.

### PUSH HANDLE OPTIONS

NXT-30518 90° Bend Back Post Tubes.....n/c

NXT-30544 Straight Back Post Tubes.....n/c  
Requires selection of Push Button Angle Adjustable Stroller Handle (NXT-30570 or NXT-30575).

### PUSH HANDLE HEIGHT

NXT-30515 Short - 20" High .....n/c

NXT-30516 Tall - 26" High .....n/c

NXT-30514 Medium - 23" High.....n/c

## Push Handle Extension Options



All push handle extensions except NXT-30560 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

NXT-30560 Push Handle Extension - Width Adjustable.....n/c  
Not available with Straight Back Posts.

NXT-30590 Stroller Handle Extension - Width Adjustable, Removable.....\$95  
Not available with Straight Back Posts.

NXT-30570 Push Button, Angle-Adjustable Stroller Handle Extension.....\$112  
Available on wheelchair frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. NXT-30570 or NXT-30575 required with Straight Back Post Tubes (NXT-30544)  
Not removable, but does allow frame to fold.

NXT-30575 Heavy Duty Push Button, Angle-Adjustable Stroller Handle Extension..... \$179  
Available on any size frame. Metal hinge up to 180° range of adjustment. NXT-30570 or NXT-30575 required with Straight Back Post Tubes (NXT-30544)  
Not removable, but does allow frame to fold.

NXT-30517 Individual Stroller Handle Extensions Removable - Pair.....\$110  
Not available with Straight Back Posts.

## Footrests

### HANGERS

NXT-30280 75° Pop-Up Swing Away Hangers.....n/c

NXT-30283 90° Latch Style Swing Away Hangers .....n/c  
Not available with 8" or 10" Telescoping Extensions.

NXT-30281 90° Pop-Up Swing Away Hangers.....n/c  
Not available with 8" or 10" Telescoping Extensions.

NXT-33000 Omit Legrest Hangers & Footrests.....n/c

NXT-30279 90° Pop-Up Swing Away Elevating Legrests (ELRs).....\$295 (E099)  
Not available with all Telescoping extensions.

NXT-30282 75° Latch Style Swing Away Hangers ..... n/c

## FOOTPLATES & FOOTBOARDS

NXT-36033 Aluminum Angle  
Adj. Flip-up Footplates ..... \$190 (K0040)  
Does not include Heel Loops, select below.

NXT-17070 Composite Flip-Up .....n/c

NXT-36029 Aluminum Angle  
Adjustable Footboard Removable..... \$190 (K0040)

NXT-36404 Omit Footrest Plates.....n/c

### Footrest / Extension Lengths

Extension/ Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-12"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	13"-17"
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	13"-15"
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	13"-18"
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	14"-20"
10" Adjustable	16"-19"	14.5"-18.5"	N/A	N/A	16"-22"



### MAPS Footrest / Extension Lengths

Extension/ Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-14.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	12.5"-16.5"
4" Adjust- able	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	12.5"-14.5"
6" Adjust- able	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	12.5"-17.5"
8" Adjust- able	12.5"-15.5"	12"-16"	N/A	N/A	13.5"-19.5"
10" Adjust- able	14.5"-17.5"	14"-18"	N/A	N/A	15.5"-21.5"

Length is measured from top of footrest plate to top of seat rail. Without extensions bottom hanger hole not accessible for hardware. Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.



## FOOTPLATE OPTIONS

NXT-36040 MAPS (Multi Axis Positioning System) ..... \$325  
Requires selection of Aluminum Angle Adjustable Footplates. Not compatible with Composite footplates or Aluminum Footboard. Not compatible with Seating Dynamics hangers.

NXT-36022 Heel Loops .....\$45 (E0951)

NXT-36024 Calf Strap - Webbing.....\$45 (K0038)

## HANGER EXTENSION OPTIONS

NXT-30285 4" Telescoping Extensions.....n/c

NXT-36402 4.5" Heavy Duty Extensions .....\$58

NXT-30286 6" Telescoping Extensions.....n/c

NXT-36403 6.5" Heavy Duty Extensions .....\$68

NXT-30287 8" Telescoping Extensions.....n/c  
Not available with all 90° hangers.

Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

NXT-30288 10" Telescoping Extensions.....n/c  
Not available with all 90° hangers.

NXT-36410 High Mount Clamp & Tubes..... \$102

Clear Selection

## Specialized Footrests

ALL legrest measurements are from the seat rail to footplate. Not available with LEAF or SRO seat rail options. Frame and/or front caster choice may limit adjustment range and/or cause interference.

### CONTRACTURE FOOTREST HANGERS

Super Contracture for frames 14" wide and greater only.

NXT-30291-2 Contracture Short - 11" to 15" ..... \$3570  
Includes short hanger extension.

NXT-30292-2 Contracture Long - 15" to 19" ..... \$3570  
Includes long hanger extension.

### 3D DYNAMIC FOOTREST HANGERS

NXT-30297-2 3D Dynamic Footrest - Hi Pivot .....\$6215

NXT-30298-2 3D Dynamic Footrest - STD ..... \$5775

### DYNAMIC FOOTREST EXTENSION OPTIONS

NXT-6405-01 Dynamic Adj Plantar/Dorsi Ext, 0-2" .....n/c

NXT-6409-01 Dynamic Adj Plantar/Dorsi Ext, 1-3" .....n/c

NXT-6406-01 Dynamic Adj Plantar/Dorsi Ext, 2-4" .....n/c

NXT-6410-01 Dynamic Adj Plantar/Dorsi Ext, 3-5" .....n/c

NXT-6407-01 Dynamic Adj Plantar/Dorsi Ext, 4-6" .....n/c

NXT-6411-01 Dynamic Adj Plantar/Dorsi Ext, 5-7" .....n/c

NXT-6408-01 Dynamic Adj Plantar/Dorsi Ext, 6-8" .....n/c

### SPECIALIZED FOOTPLATES

Must pick a footplate and extension

NXT-6503 4" x 6" HD Aluminum Footplates .....n/c

NXT-6504 5" x 8" HD Aluminum Footplates .....n/c  
For frames 16" wide and greater only.

NXT-6505 6" x 10" HD Aluminum Footplates.....n/c  
For frames 18" wide and greater only. Not available with Contracture Hangers

## Specialized Footrest Seat Rail to Footplate Measurements

Part Number	Extension	Super Contracture	3D Dynamic Hi Pivot	3D Dynamic STD
6402-01	Short	10.5" - 13.5"		
6402-02	Long	13.5" - 18.5"		
6405-01	0"-2" Extension		6.25" - 10.75"	12.75-15.25"
6409-01	1"-3" Extension		7.25-11.75"	13.75-16.25"
6406-01	2"-4" Extension		8.25-12.75"	14.75-17.25"
6410-01	3"-5" Extension		9.25-13.75"	15.75-18.25"
6407-01	4"-6" Extension		10.25-14.75"	16.75-19.25"
6411-01	5"-7" Extension		11.25-15.75"	17.75-20.25"
6508-01	6"-8" Extension		12.25-16.75"	18.75-21.25"



Length measured from top of footrest plate to top of seat rail. 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.

## Armrests

### BACK POST MOUNTED ARMRESTS

Omit Armrests

Foam wrap pad STD on all back post mounted armrests.

NXT-38204 10.5" Flip-Back.....n/c

NXT-38202 13.5" Flip-Back.....n/c

NXT-38210 10.5" to 13.5" Flip-Back  
Locking, Telescoping.....\$285 (E0973)

Full length pad on 13.5" or longer armrests only.

NXT-38201 10.5" Flip-Back Locking.....\$265 (E0973)

NXT-38200 13.5" Flip-Back Locking.....\$265 (E0973)

NXT-38215 17" to 20" Flip-Back  
Locking, Telescoping.....\$265 (E0973)

### FRAME MOUNTED ARMRESTS

NXT-38217 Single Post Height Adjustable,  
Removable - Low .....\$285 (E0973)

NXT-38203 Single Post Height Adjustable,  
Removable - Standard .....\$285 (E0973)

NXT-38206 Dual Post Flip-Back Height  
Adjustable, Removable .....\$295 (E0973)

NXT-38205 Dual Post Fixed Height Removable.....n/c

Not available with short frame. 22" & 24" Wheels not available with dual post armrests

### ARM PAD OPTIONS

NXT-38213 Desk Length Pads.....n/c

NXT-38214 Full Length Pads.....n/c  
13.5" or longer armrests only.

## Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

All heights from top of seat frame to top of armrest pad.







## Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Tilt Pivot	Caster Wheel				
		3"	4"	5"	6"	7"
12"	Low	12.5"-13"	13.5	15.5"	16"	
	Tall	14"-14.5"	15"	17"	17.5"	
16"	Low	12.5"-13"	13.5"-14"	14.5"	16"	16.5"
	Tall	14"-14.5"	15"-15.5"	16"	17.5"	18"
18"	Low	12.5"-13"	13.5"-14.5"	14.5"-15.5"	15.5"	
	Tall	14"-14.5"	15"-16.5"	16"-17"	17"	
20"	Low		14"-15"	14.5"-15.5"	15.5"-16"	16.5"
	Tall		15.5"-16.5"	16"-17"	17"-17.5"	18"
22"	Low					
	Tall		16.5"	16.5-17"	17"-17.5"	18"
24"	Low					
	Tall				17.5"	18"
Clear Selection						

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".

## Casters

NXT-34888 3" x 1" Polyurethane .....n/c	NXT-42000 5" x 1.4" Poly, Aluminum Mag ..... \$125
NXT-43111 3" Polyurethane Light-Up..... \$87	NXT-34666 6" x 1" Polyurethane .....n/c
NXT-40000 4" x 1" Polyurethane.....n/c	NXT-34999 6" x 1.5" Polyurethane..... \$87
NXT-43000 4" x 1.4" Poly, Aluminum Mag ..... \$125	NXT-44000 6" x 1.4" Poly, Aluminum Mag..... \$125
NXT-34555 5" x 1" Polyurethane.....n/c	NXT-34444 6" x 2" Polyurethane, Treaded .....\$97
NXT-34222 5" x 1.5" Polyurethane..... \$87	NXT-34777 7" x 1" Polyurethane.....n/c
NXT-35050 5" Polyurethane Light-Up..... \$101	NXT-34200 7" x 2" Foam Filled .....\$200

7" casters not available with 12" mag wheels.



## Rear Wheels

### MAGS

- NXT-31525 12" Mag Wheel.....n/c  
12" mag wheels not available with 7" casters.
- NXT-31530 16" Mag Wheel.....\$41
- NXT-31533 18" Mag Wheel .....\$41

- NXT-31527 20" Mag Wheel - Aluminum Handrim..... n/c
- NXT-31528 22" Mag Wheel - Aluminum Handrim..... n/c
- NXT-31529 24" Mag Wheel - Aluminum Handrim .....n/c  
22" & 24" Wheels not available with dual post armrests.

### SPOKE

- NXT-31535 16" Spoke Wheel - Black w/ Black Tire..... n/c
- NXT-31538 18" Spoke Wheel - Black w/ Black Tire.....n/c
- NXT-31540 20" Spoke Wheel - Aluminum Handrim..... \$204

- NXT-31541 22" Spoke Wheel - Aluminum Handrim..... \$204
- NXT-31542 24" Spoke Wheel - Aluminum Handrim..... \$204  
20", 22" & 24" Spoke Wheels not available with drum brakes. 22" & 24" Wheels not available with dual post armrests. Spoke wheels not available with plastic coated handrims.

### TIRE OPTIONS

- NXT-31040 Pneumatic Tire ..... \$85
- NXT-31000 Pneumatic Tire  
with Airless Insert ..... \$85 (E2213)

- NXT-31045 Low Profile Polyurethane Tire .....n/c (E2218)  
12" Mag wheel only, not available with drum brakes

### HANDRIM OPTIONS

- NXT-31005 Plastic Coated Handrims .....\$119  
20", 22" & 24" Mag wheels only, not available with drum brakes

- NXT-31006 Omit Handrims

## Wheel Locks

### STANDARD

- NXT-31501 Push to Lock.....n/c

- NXT-31506 Pull to Lock .....n/c

### WHEEL LOCK OPTIONS

- NXT-31801 Extension Handles ..... \$75 (E0961)

### DRUM LOCK

Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires. Not Available with 20", 22" & 24" Spoke Wheels.

- NXT-31508 Foot Activated Wheel Lock.....n/c  
Not available with ALL vent and battery trays.

- NXT-30275 Hand Activated Brake  
and Wheel Lock System .....\$1,000 (E2228)

## Anti-Tippers

- NXT-31600 Curb Assist Anti-Tip  
Assembly Flip-Up, Adjustable ..... \$125 (E0971)  
Pair includes one anti-tip tube with curb assist feature.

- NXT-31603 Standard Anti-Tip  
Assembly Flip-Up, Adjustable ..... \$125 (E0971)



## Seat Pans

Only choose one Seat Pan option. ALL seat pans include 2" x 2" notches for growth adjustment.

### NON-ADJUSTABLE

NXT-39010 Fixed Solid Seat Pan - Composite .....\$170 (E2231)

ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width \_\_\_\_\_" X Seat Depth \_\_\_\_\_"

NXT-39030 Fixed Solid Seat Pan - Aluminum .....\$200 (E2231)

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Available in the following sizes ONLY:

14" wide x 16" deep	16" wide x 18" deep	18" wide x 18" deep
------------------------	------------------------	------------------------

### ADJUSTABLE

NXT-39040 Fixed, Width Adjustable  
Solid Seat Pan - Aluminum .....\$290 (E2231)

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with NXT frame.

Available in the following sizes ONLY:

14"-16" wide x 16" deep

16"-18" wide x 18" deep

## MOUNTING HARDWARE OPTIONS

Use to convert fixed seat pan to removable.

NXT-EZLITE Removable Pan Hardware Kit.....\$207 (K0108)

Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

## Accessories

SFDI-9HP Footrest Hanger Pads - Pair ..... \$97 (K0108)

Wrap-around pads made from black neoprene.

NXT-23210 Flip-Down Rigidizing Bar..... \$175 (K0108)

Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.

NXT-10900 Freedom Backpack..... \$60 (K0108)

Black Nylon only with Freedom Designs logo. Monogramming is not available on Freedom Backpack.

NXT-10905 Freedom Canopy - Flip-Back..... n/c (K0108)

Canopy assembly includes removable mounting hardware. Available in black only with Freedom Designs logo. Includes rear flap with zippered storage pocket and clear-view window.

## Specifications

Frame Style	Tilt-In-Space
Base Frame Widths	14' - 20"
Overall Width	Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel)
Seat Widths	14" - 20"
Seat Depths	14" - 26"
Back Heights	20", 23", 26"
Adjustable Angle Back	90°, 100°, 110°, 120°
Tilt Range	0° to 55°
Recline Back	N/A
Weight Limit	225 LBS.
Transport Weight	Approximately 22 lbs
Warranty	5 year warranty on side frames and quick release axles

## Special Configurations Through Customer Service

Asymmetrical Seat Depth	Seat rails may be configured up to 4" offset
-------------------------	--

## HCPCS Codes

are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Special Instructions:



## Adobe Acrobat Reader DC

\* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

---

### Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form.

Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

### Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

### Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

**800-331-8551**