



# Pediatric Manual Wheelchair



Yes, you can!

Clear Form

Email Form to Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: \_\_\_\_\_ Buyer: \_\_\_\_\_ Quote Order PO# \_\_\_\_\_  
 Company: \_\_\_\_\_ Location: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Email: \_\_\_\_\_ Fax: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Transport Optional Base

<b>NXT Pediatric Folding</b>	(E1234) NXT-30300	<b>\$2800</b>
<i>Transport Optional</i>		
<b>NXT Pediatric Rigid</b>	(E1233) NXT-30310	<b>\$2700</b>
<i>Transport Optional</i>		
<b>Note:</b> Transportation Option selected separately below.		
<b>Note:</b> Includes Freedom Future Fit Program		

<b>Transport Package</b>	NXT-99998	<b>\$ 250</b>
<b>Note:</b> Must be installed by Freedom Designs manufacturing.		
<b>Note:</b> Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.		
<b>Note:</b> Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.		

### Vent Tray Options

**Note:** All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.

<b>Vent Tray - Fixed</b>	(E1029) 11-7000	<b>\$ 439</b>
<b>Note:</b> Not with Foot Activated Tilt.		
<b>Battery Tray - Fixed</b>	(K0108) 11-8000	<b>\$ 429</b>
<b>Quick Release Tray Hardware</b>	(E1028) 563	<b>\$ 272</b>
<b>Note:</b> For LP10 style vent and battery tray on frames 16" wide and wider.		
<b>Note:</b> Vent tray is 4" narrower than frame width and 16" deep.		
<b>Note:</b> Battery tray is 4" narrower than frame width and 10" deep.		
<b>Back Post Mounted Vent Bracket</b>	(E1029) 11-1100	<b>\$ 550</b>
<b>Note:</b> For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.		

### Transport Included Base

<b>NXT Pediatric Folding</b>	(E1234) NXT-30301	<b>\$2950</b>
<i>Transport Included</i>		
<b>NXT Pediatric Rigid</b>	(E1233) NXT-30311	<b>\$2850</b>
<i>Transport Included</i>		
<b>Note:</b> Includes Freedom Future Fit Program		
<b>Note:</b> Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.		

### Frame Color Options

Sugar Cookie	n/c	Blue Blast	n/c
Neon Yellow	n/c	Electric Blue	n/c
Yellow	n/c	Galaxy Blue	n/c
Bright Pink	n/c	Passion Purple	n/c
Pink Parfait	n/c	Octane Orange	n/c
Sugar Plum	n/c	Red Sunset	n/c
Candy Violet	n/c	Really Red	n/c
Candy Raspberry	n/c	Fiery Red	n/c
Extreme Green	n/c	Twilight	n/c
Sparkle Green	n/c	Stormy Gray	n/c
Mystic Teal	n/c	Black	n/c
Pacific Blue	n/c		

### Seat Rail Options

**Note:** Not available with ALL Contracture & 3D Dynamic footrest systems.

**LEAF - Lower Extremity Accommodation Frame**

<b>LEAF Seat Frame- Pair</b>	(K0108) NXT-30324	<b>\$ 714</b>
<b>LEAF Seat Frame- Left</b>	(K0108) NXT-30324-L	<b>\$ 357</b>
<b>LEAF Seat Frame- Right</b>	(K0108) NXT-30324-R	<b>\$ 357</b>

**Note:** Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

### SRO - Seat Rail Only

<b>SRO Seat Frame- Pair</b>	(K0108) NXT-30325	<b>\$ 306</b>
<b>SRO Seat Frame- Left</b>	(K0108) NXT-30325-L	<b>\$ 153</b>
<b>SRO Seat Frame- Right</b>	(K0108) NXT-30325-R	<b>\$ 153</b>

**Note:** SRO is manufactured without standard footrest hanger attachments.

### Tilt Activation Options

<b>Dual Lever Hand Activated Tilt</b>	NXT-30272	n/c
<b>Foot Activated Tilt</b>	NXT-30273	n/c
<b>Note:</b> Not with vent & battery trays.		
<b>Combination Dual Hand Activated Tilt, Brake and Wheel Lock System</b>	NXT-30277	n/c

**Note:** Includes manual wheel locks to be used with public transportation.  
**Note:** Not available with 12" mag wheels with low pro poly tires.

### Frame Dimensions

#### Specify Frame Width Setting

10"	11"	12"	Adjusts 10" to 12"	NXT-30250	n/c
11"	12"	13"	Adjusts 11" to 13"	NXT-30210	n/c
12"	13"	14"	Adjusts 12" to 14"	NXT-30220	n/c
13"	14"	15"	Adjusts 13" to 15"	NXT-30230	n/c
14"	15"	16"	Adjusts 14" to 16"	NXT-30240	n/c
15"	16"	17"	Adjusts 15" to 17"	NXT-30290	n/c
16"	17"	18"	Adjusts 16" to 18"	NXT-30260	n/c
17"	18"		Adjusts 17" to 18"	NXT-30261	n/c
18"			Not width adjustable	NXT-30262	n/c
19"			Not width adjustable, not with short frame	NXT-30263	n/c
20"			Not width adjustable, not with short frame	(E2201) NXT-30264	<b>\$ 400</b>

**Note:** Frame widths measures across and from outside of seat frame tubes.

**Note:** Seat frame widths of 19" and 20" not available with short frame.

**Note:** Vent and battery trays only with 14" wide and greater frame.

**NXT Rider Weight Capacity is 225 lbs.**  
All prices subject to change without notice.

**NXT Pediatric**  
Note: All Specifications and dimensions are approximate.

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Frame Type Options

#### Specify Frame Type and Frame Length

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

**Short Frame and Short Tilt Bracket** NXT-30511 n/c

#### Specify Back Post Depth Setting

12"	NXT-30504	n/c	
13"	NXT-30505	n/c	
14"	NXT-30506	n/c	
15"	NXT-30521	n/c	
16"	NXT-30522	n/c	
* 17" Must select 2" extension	NXT-30523	n/c	
* 18" Must select 2" extension	NXT-30524	n/c	
2" Extension - Pair	(K0108) NXT-30525	\$ 26	

#### Specify Back Post Location (if other than specified above) \_\_\_\_\_"

- Note: Seat rail to floor height range is 12.5" to 16.5".  
 Note: Seat frame widths of 19" and 20" not available with short frame.  
 \* Note: Must choose extension for 17" and 18" frame lengths.  
 Note: 22" & 24" wheels not available with short tilt bracket.  
 Note: All vent and battery trays not available with short frame.

**Medium Frame and Short Tilt Bracket** NXT-30512 n/c

#### Specify Back Post Depth Setting

14"	NXT-30506	n/c	
15"	NXT-30521	n/c	
16"	NXT-30522	n/c	
17"	NXT-30523	n/c	
18"	NXT-30524	n/c	
* 19" Must select 2" extension	NXT-30527	n/c	
* 20" Must select 2" extension	(E2203) NXT-30528	\$ 500	
* 21" Must select 3" extension	(E2203) NXT-30529	\$ 500	
2" Extension - Pair	(K0108) NXT-30525	\$ 26	
3" Extension - Pair	(K0108) NXT-30530	\$ 26	

#### Specify Back Post Location (if other than specified above) \_\_\_\_\_"

- Note: Seat rail to floor height range is 12.5" to 16.5".  
 \* Note: Must choose extension for 19" to 21" frame lengths.  
 Note: 22" & 24" wheels not available with short tilt bracket.

**Medium Frame and Tall Tilt Bracket** NXT-30540 n/c

#### Specify Back Post Depth Setting

14"	NXT-30506TP	n/c	
15"	NXT-30521TP	n/c	
16"	NXT-30522TP	n/c	
17"	NXT-30523TP	n/c	
18"	NXT-30524TP	n/c	
* 19" Must select 2" extension	NXT-30527TP	n/c	
* 20" Must select 2" extension	(E2203) NXT-30528TP	\$ 500	
* 21" Must select 3" extension	(E2203) NXT-30529TP	\$ 500	
* 22" Must select 4" extension	(E2204) NXT-30535TP	\$ 900	
2" Extension - Pair	(K0108) NXT-30525	\$ 26	
3" Extension - Pair	(K0108) NXT-30530	\$ 26	
4" Extension - Pair	(K0108) NXT-30531	\$ 26	

#### Specify Back Post Location (if other than specified above) \_\_\_\_\_"

- Note: Tall tilt bracket increases seat rail to floor 1.5".  
 Note: Seat rail to floor height range is 14" to 18".  
 \* Note: Must choose extension for 19" to 22" frame lengths.

### Frame Type Options - Continued

**Pediatric Long Frame & Tall Tilt Bracket** (E1298) NXT-30520P \$ 153

#### Specify Back Post Depth Setting

16"	NXT-30522TP	n/c	
17"	NXT-30523TP	n/c	
18"	NXT-30524TP	n/c	
19"	NXT-30527TP	n/c	
20"	(E2203) NXT-30528TP	\$ 500	
21"	(E2203) NXT-30529TP	\$ 500	
22"	(E2204) NXT-30535TP	\$ 900	
* 23" Must select 2" extension	(E2204) NXT-30537TP	\$ 900	
* 24" Must select 2" extension	(E2204) NXT-30538TP	\$ 900	
* 25" Must select 3" extension	(E2204) NXT-30542TP	\$ 900	
* 26" Must select 4" extension	(E2204) NXT-30543TP	\$ 900	
2" Extension - Pair	(K0108) NXT-30525	\$ 26	
3" Extension - Pair	(K0108) NXT-30530	\$ 26	
4" Extension - Pair	(K0108) NXT-30531	\$ 26	

#### Specify Back Post Location (if other than specified above) \_\_\_\_\_"

- Note: Tall tilt bracket increases seat rail to floor 1.5".  
 Note: Seat rail to floor height range is 14" to 18".  
 \* Note: Must choose extension for 23" to 26" frame lengths.

**5° Anterior Tilt Seat Frame** NXT-30276 n/c

- Note: Limits posterior tilt to 40°. Includes split collars for adjustment.  
 Note: Medium with Tall Tilt and Long frames only.

Note: NXT has 4" of depth adjustment for growth. Depth is measured from front of back-post to the front end of the seat frame. For back post location add seating & hardware thickness (typically 2") to client's actual seat depth. Example: 12" client seat depth + 2" back & hardware thickness = 14" back post location.

### Back Post Options

#### Specify Push Handle, Back Post and Extensions

**Standard 90° Bend Push Handle Back Post Tubes** NXT-30518 n/c

**Ergo Bend Push Handle Back Post Tubes** NXT-30519 n/c

Note: All push/stroller handle extensions not recommended with Ergo Bend.

**Seating Dynamics Rocker Back** (K0108) NXT-30532 \$1985

- Note: Angle adjustable from 72° to 108° in 6° increments.  
 Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.  
 Note: Locks into a non-dynamic position when necessary.  
 Note: Does not allow the back post tubes to fold forward.

#### Back Post Height Setting

Note: Back Posts are measured from top of seat rail to top of handle bend.

Short Push Handle - 20" high	NXT-30515	n/c
Medium Push Handle - 23" high	NXT-30514	n/c
Tall Push Handle - 26" high	NXT-30516	n/c

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Stroller Handle Extension Options

Note: All stroller handle ext. add approximately 6.5" to push handle height.  
 Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.



**Push Handle Extension - Width Adjustable - Mounts direct to back post push handle** NXT-30560 n/c

**Stroller Handle Extension - Width Adjustable, Removable - Mounts direct to back post push handle** **NEW** NXT-30590 \$ 95

**Push Button, Angle-Adjustable Stroller Handle Extension - Mounts direct to back post push handle** NXT-30570 \$ 112

**Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty - Mounts direct to back post tubes** NXT-30575 \$ 179

Note: Available 10W-20w wheelchair frames.  
 Note: Not removable, but allows frame to fold.

**Push Button, Angle-Adjustable, Folding Push Handle Extension - Mounts direct to straight back post tubes** NXT-30580 \$ 112

**Push Button, Angle-Adjustable, Folding Push Handle Extension - Heavy Duty - Mounts direct to straight back post tubes** NXT-30585 \$ 179

**Individual Stroller Handle Extensions Removable - Pair** NXT-30517 \$ 102

### Legrests, Footrests and Extension Options

#### Legrests

**Omit Legrest Hangers & Footrests** NXT-33000

**75° Pop-Up Swing Away Hangers** NXT-30280 n/c

**90° Pop-Up Swing Away Hangers** NXT-30281 n/c

Note: Not available with 8" adjustable extensions.  
**90° Pop-Up Swing Away Elevating Legrests (ELRs)** (E0990) NXT-30279 \$ 255

Note: Not available with all length adjustable extensions.  
 Note: Not available with 6.5" heavy duty extensions.  
**75° Latch Style Swing Away Hangers** NXT-30282 n/c

**90° Latch Style Swing Away Hangers** NXT-30283 n/c

Note: Not available with 8" adjustable extensions.

#### Footplates & Footboards

**Omit Footrest Plates** NXT-36404

**Composite Flip-Up** NXT-17070 n/c

**Kids Aluminum Angle Adjustable, Flip-Up** (K0040) NXT-36020 \$ 160

Note: For frames 14" wide and narrower.  
**Heel Loops** (E0951) NXT-36022 \$ 41

**Calf Strap - Webbing** (K0038) NXT-36024 \$ 41

**Adult Aluminum Angle Adj. Flip-up w/ Heel Loops** (K0040/E0951) NXT-36021 \$ 180

Note: For frames 15" wide and greater.  
**Omit Heel Loops** NXT-36023 n/c

**Aluminum Angle Adjustable Footboard Removable** (K0040) NXT-36029 \$ 160

#### MAPS Footplate Mounting Hardware

**MAPS (Multi Axis Positioning System) - Pair** NXT-36040 \$ 306

**MAPS (Multi Axis Positioning System) - Left** NXT-36040-L \$ 153

**MAPS (Multi Axis Positioning System) - Right** NXT-36040-R \$ 153

Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.

Note: Must choose aluminum footplates.  
 Note: MAPS not available with composite footplates or aluminum footboard.

### Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	1.5"-3"	1.5"-3"	1.5"-3"	1.5"-3"	
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	
8" Adjustable	14"-17"	12.5"-16.5"			

### MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount					
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	
8" Adjustable	12.5"-15.5"	12"-16"			

Note: Length from top of footrest plate to top of seat rail.  
 Note: Without extensions bottom hanger hole not accessible for hardware.  
 Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Legrests, Footrests and Extension Options

#### Legrest Extensions

<b>4" Length Adjustable Extensions</b>	NXT-30285	n/c
<b>6" Length Adjustable Extensions</b>	NXT-30286	n/c
<b>8" Length Adjustable Extensions</b>	NXT-30287	n/c

Note: 8" & 10" Not available with all 90° hangers.

<b>4.5" Heavy Duty ELR Extensions</b>	NXT-36401	\$ 56
<b>4.5" Heavy Duty Extensions</b>	NXT-36402	\$ 56
<b>6.5" Heavy Duty Extensions</b>	NXT-36403	\$ 66

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

<b>High Mount Brackets</b>	NXT-36407	\$ 102
----------------------------	-----------	--------

Note: Includes flip-up mounting hardware for future growth adjustment.

Note: 75° hangers use a tube & clamp style hardware.

Specify Knee to Heel Measurement \_\_\_\_\_" (see matrix)

**Notes:**

### Specialized Footrest Systems

Note: ALL legrest measurements are from the seat rail to footplate.  
 Note: Not available with LEAF or SRO seat rail options.

#### Super Contracture Footrest Hanger (K0108)

Note: Super Contracture for frames 14" wide and greater.

<b>Contracture Short - 11" to 15"</b>	NXT-30291-2	\$3570
---------------------------------------	-------------	--------

Note: See Matrix for seat to rail footplate measurements.

#### Includes:

<b>Adj Contracture Upper</b>	NXT-6302-03	Incl.
<b>Footrest Adapter HD</b>	NXT-6220-01	Incl.
<b>Adj Plantar/Dorsi Ext. Short</b>	NXT-6401-02	Incl.

#### Footplate Options

<b>4" x 6" HD Alu Footplates</b>	NXT-6503	n/c
<b>5" x 7" HD Alu Footplates</b>	NXT-6504	n/c

Note: for frames 16" wide and greater only.

<b>Contracture Long - 15" to 19"</b>	NXT-30292-2	\$3570
--------------------------------------	-------------	--------

Note: See Matrix for seat to rail footplate measurements.

#### Includes:

<b>Adj Contracture Upper</b>	NXT-6302-03	Incl.
<b>Footrest Adapter HD</b>	NXT-6220-01	Incl.
<b>Adj Plantar/Dorsi Ext. Long</b>	NXT-6402-02	Incl.

#### Footplate Options

<b>4" x 6" HD Alu Footplates</b>	NXT-6503	n/c
<b>5" x 7" HD Alu Footplates</b>	NXT-6504	n/c

Note: Super Contracture is both swing-away and removable.

Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.

Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.

Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

#### 3D Dynamic Footrest Hanger (K0108)

Note: 3D Dynamic for frames 14" wide and greater.

<b>3D Dynamic Footrest - Hi Pivot</b>	NXT-30297-2	\$6215
---------------------------------------	-------------	--------

#### Includes:

<b>HiPivot Dyn Adj Knee &amp; Dyn Tel Up</b>	NXT-6304-02	Incl.
<b>Footrest Adapter HD</b>	NXT-6220-01	Incl.

#### Extension Options

Note: See Matrix for seat to rail footplate measurements.

<b>Dynamic Adj Plantar/Dorsi Ext, 0-2"</b>	NXT-6405-01	n/c
<b>Dynamic Adj Plantar/Dorsi Ext, 1-3"</b>	NXT-6409-01	n/c
<b>Dynamic Adj Plantar/Dorsi Ext, 2-4"</b>	NXT-6406-01	n/c
<b>Dynamic Adj Plantar/Dorsi Ext, 3-5"</b>	NXT-6410-01	n/c
<b>Dynamic Adj Plantar/Dorsi Ext, 4-6"</b>	NXT-6407-01	n/c
<b>Dynamic Adj Plantar/Dorsi Ext, 5-7"</b>	NXT-6411-01	n/c
<b>Dynamic Adj Plantar/Dorsi Ext, 6-8"</b>	NXT-6408-01	n/c

#### Footplate Options

<b>4" x 6" HD Alu Footplates</b>	NXT-6503	n/c
<b>5" x 7" HD Alu Footplates</b>	NXT-6504	n/c

Note: for frames 16" wide and greater only.

<b>6" x 10" HD Alu Footplates</b>	NXT-6505	n/c
-----------------------------------	----------	-----

Note: for frames 18" wide and greater only.

### Specialized Footrest Seat Rail to Footplate Measurements

Part Number	Extension	Super Contracture	3D Dynamic Hi Pivot	3D Dynamic STD
6402-01	Short	10.5" - 13.5"		
6402-02	Long	13.5" - 18.5"		
6405-01	0"-2" Extension		6.25" - 10.75"	12.75-15.25"
6409-01	1"-3" Extension		7.25-11.75"	13.75-16.25"
6406-01	2"-4" Extension		8.25-12.75"	14.75-17.25"
6410-01	3"-5" Extension		9.25-13.75"	15.75-18.25"
6407-01	4"-6" Extension		10.25-14.75"	16.75-19.25"
6411-01	5"-7" Extension		11.25-15.75"	17.75-20.25"
6508-01	6"-8" Extension		12.25-16.75"	18.75-21.25"

Note: Length from top of footrest plate to top of seat rail.

Note: 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Specialized Footrest Systems - Continued

**3D Dynamic Footrest - STD** NXT-30298-2 \$5775

**Includes:**

**Dynamic Adj Knee & Dyn Tel Up** NXT-6304-01 Incl.

**Footrest Adapter HD** NXT-6220-01 Incl.

**Extension Options**

Note: See Matrix for seat to rail footplate measurements.

**Dynamic Adj Plantar/Dorsi Ext, 0-2"** NXT-6405-01 n/c

**Dynamic Adj Plantar/Dorsi Ext, 1-3"** NXT-6409-01 n/c

**Dynamic Adj Plantar/Dorsi Ext, 2-4"** NXT-6406-01 n/c

**Dynamic Adj Plantar/Dorsi Ext, 3-5"** NXT-6410-01 n/c

**Dynamic Adj Plantar/Dorsi Ext, 4-6"** NXT-6407-01 n/c

**Dynamic Adj Plantar/Dorsi Ext, 5-7"** NXT-6411-01 n/c

**Dynamic Adj Plantar/Dorsi Ext, 6-8"** NXT-6408-01 n/c

**Footplate Options**

**4" x 6" HD Alu Footplates** NXT-6503 n/c

**5" x 7" HD Alu Footplates** NXT-6504 n/c

Note: for frames 16" wide and greater only.

**6" x 10" HD Alu Footplates** NXT-6505 n/c

Note: for frames 18" wide and greater only.

Note: 3D Dynamic is both swing-away and removable.

Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.

Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.

Note: Angle adjustable footplate included.

Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

### Armrests and Armrest Pad Options

**Omit Armrests** NXT-38209

**Back Post Mounted Armrests**

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

**10.5" Long Flip-Back** NXT-38204 n/c

Note: Not available with full length arm pads

**13.5" Long Flip-Back** NXT-38202 n/c

**10.5" Long Flip-Back Locking** (E0973) NXT-38201 \$ 245

Note: Not available with full length arm pads

**13.5" Long Flip-Back Locking** (E0973) NXT-38200 \$ 245

**10.5" to 13.5" Long Flip-Back, Telescoping** (E0973) NXT-38208 \$ 240

Note: Full length pad on 13.5" or longer armrests only.

**10.5" to 13.5" Long Flip-Back Locking, Telescoping** (E0973) NXT-38210 \$ 265

Note: Full length pad on 13.5" or longer armrests only.

**17" to 20" Long Flip-Back Locking, Telescoping** (E0973) NXT-38215 \$ 265

**Foam Wraps - Standard** NXT-30216 n/c

**Desk Length Pads** NXT-38213 n/c

**Full Length Pads** NXT-38214 n/c

Note: Full length pad on 13.5" or longer armrests only.

**Frame Mounted Armrests**

**Single Post Fixed Height, Removable - Short** (K0015) NXT-38222 \$ 225

**Single Post Fixed Height, Removable - Tall** (K0015) NXT-38224 \$ 225

Note: Single Post Fixed Height arms include a removable side guard.

Note: Adjustable in height with tools only. See matrix for height ranges.

**Single Post Height Adjustable, Removable - Short** (E0973) NXT-38217 \$ 265

**Single Post Height Adjustable, Removable - Tall** (E0973) NXT-38203 \$ 265

**Dual Post Fixed Height Removable** NXT-38205 n/c

**Dual Post Flip-Back Height Adjustable, Removable** (E0973) NXT-38206 \$ 275

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

**Desk Length Pads** NXT-38213 n/c

**Full Length Pads** NXT-38214 n/c

Specify Seat Rail to Top of Armrest Measurement \_\_\_\_\_"

### Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Single Post Fxd Hgt Short		8" - 10.5"	8" - 10.5"
Single Post Fxd Hgt Short		9" - 11.5"	9" - 11.5"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

Continue To Next Page

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Choose NXT Seat Rail to Floor Height Below:

**Note:** Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height.  
**Note:** Caster sizes vary slightly in height by type.  
**Note:** Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".  
**Note:** Seat pan choice increases seat rail to floor height 3/8" to 3/4".

### NXT Wheel / Caster Seat to Floor Heights

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5-13"	20"	3"	Tall	N/A
12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"
12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"
12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"
12"	5"	Short	15.5"	20"	5"	Tall	16-17"
12"	5"	Tall	17"	20"	6"	Short	15.5-16"
12"	6"	Short	16"	20"	6"	Tall	17-17.5"
12"	6"	Tall	17.5"	20"	7"	Short	16.5"
16"	3"	Short	12.5-13"	20"	7"	Tall	18"
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A
16"	5"	Short	14.5"	22"	4"	Tall	16.5"
16"	5"	Tall	16"	22"	5"	Short	N/A
16"	6"	Short	16	22"	5"	Tall	16.5-17"
16"	6"	Tall	17.5	22"	6"	Short	N/A
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"
16"	7"	Tall	18"	22"	7"	Short	N/A
18"	3"	Short	12.5-13"	22"	7"	Tall	18"
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A
18"	5"	Tall	16-17"	24"	5"	Short	N/A
18"	6"	Short	15.5"	24"	5"	Tall	N/A
18"	6"	Tall	17"	24"	6"	Short	N/A
18"	7"	Short	N/A	24"	6"	Tall	17.5"
18"	7"	Tall	N/A	24"	7"	Short	N/A
20"	3"	Short	N/A	24"	7"	Tall	18"

### Front Wheels - Caster Options

<b>Omit Caster Wheels</b>	NXT-34333	
<b>3" Polyurethane</b>	NXT-34888	n/c
<b>3" Polyurethane Light-Up</b>	NXT-34111	\$ 101
<b>4" Polyurethane</b>	NXT-40000	n/c
<b>4" Polyurethane Light-Up</b>	NXT-41000	\$ 101
<b>4" x 1.4" Poly, Aluminum Mag</b>	NXT-43000	\$ 175
<b>5" x 1" Poly, Plastic Mag</b>	NXT-34555	n/c
<b>5" Polyurethane Light-Up</b>	NXT-35050	n/c
<b>5" x 1.5" Poly, Plastic Mag</b>	NXT-34222	\$ 87
<b>5" x 1.4" Poly, Aluminum Mag</b>	NXT-42000	\$ 175

### Front Wheels - Caster Options

<b>6" x 1" Poly, Plastic Mag</b>	NXT-34666	n/c
<b>6" x 1.5" Poly, Plastic Mag</b>	NXT-34999	\$ 87
<b>6" x 1.4" Poly, Aluminum Mag</b>	NXT-44000	\$ 175
<b>6" x 2" Poly with Tread</b>	NXT-34444	\$ 92
<b>6" x 1" Pneumatic Plastic Mag</b>	(E2214/ E2215) NXT-34667	\$ 77
<b>7" x 1" Poly, Plastic Mag</b>	NXT-34777	n/c
<b>7" x 2" Polyurethane</b>	NXT-34200	\$ 118

**Note:** All 7" casters not available with all 12" rear wheels.

### Rear Wheels - Tires - Accessories

**Note:** All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.

<b>Omit Rear Wheels and Axles</b>	NXT-31524	
<b>Mag Wheels</b>		
<b>12" Low Profile Poly Mag</b>	(E2218) NXT-31526	\$ 36
<b>Note:</b> Not available with any foot/hand activated brake/wheel lock system.		
<b>12" Semi Pneumatic Plastic Mag</b>	NXT-31523	\$ 168
<b>Note:</b> Not available with airless inserts.		
<b>12" Pneumatic Plastic Mag</b>	NXT-31525	n/c
<b>12" Airless Inserts</b>	(E2213) NXT-31019	\$ 77
<b>Note:</b> All 12" wheels not available with all 7" casters.		
<b>16" Pneumatic Plastic Mag</b>	NXT-31530	\$ 41
<b>16" Airless Inserts</b>	(E2213) NXT-31020	\$ 77
<b>16" Aluminum Handrims</b>	NXT-31014	\$ 169
<b>18" Pneumatic Plastic Mag</b>	NXT-31533	\$ 41
<b>18" Airless Inserts</b>	(E2213) NXT-31021	\$ 77
<b>18" Aluminum Handrims</b>	NXT-31015	\$ 169
<b>20" Pneumatic Plastic Mag with Aluminum Handrim</b>	NXT-31527	n/c
<b>20" Airless Inserts</b>	(E2213) NXT-31022	\$ 77
<b>20" Plastic Coated Handrims</b>	NXT-31016	\$ 119
<b>Note:</b> Plastic Coated Handrims NOT available with Drum Braking Systems.		
<b>22" Pneumatic Plastic Mag with Aluminum Handrim</b>	NXT-31528	n/c
<b>22" Airless Inserts</b>	(E2213) NXT-31023	\$ 77
<b>22" Plastic Coated Handrims</b>	NXT-31017	\$ 119
<b>Note:</b> Not available with Short Tilt		
<b>Note:</b> Not available with dual post armrests.		
<b>Note:</b> Plastic Coated Handrims NOT available with Drum Braking Systems.		
<b>24" Pneumatic Plastic Mag with Aluminum Handrim</b>	NXT-31529	n/c
<b>24" Airless Inserts</b>	(E2213) NXT-31024	\$ 77
<b>24" Plastic Coated Handrims</b>	NXT-31018	\$ 119
<b>Note:</b> Not available with Short Tilt		
<b>Note:</b> Not available with dual post armrests.		
<b>Note:</b> Plastic Coated Handrims NOT available with Drum Braking Systems.		
<b>Omit Handrims from 20", 22" or 24" Wheels</b>	NXT-31006	

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

**Rear Wheels - Tires - Accessories - Cont.**

**Spoke Wheels**

Note: All spoke wheels not available plastic coated handrims.

**16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel** NXT-31534 n/c

Note: Rim, hub & spokes are anodized silver.

**16" Airless Inserts** (E2213) NXT-31020 \$ 77

**16" Silver Aluminum Handrim** NXT-31030 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**16" Black Pneumatic Tire & Black Spoke Aluminum Wheel** NXT-31535 n/c

Note: Rim, hub & spokes are anodized black.

**16" Airless Inserts** (E2213) NXT-31020 \$ 77

**16" Black Aluminum Handrim** NXT-31031 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel** NXT-31537 n/c

Note: Rim, hub & spokes are anodized silver.

**18" Airless Inserts** (E2213) NXT-31021 \$ 77

**18" Silver Aluminum Handrim** NXT-31034 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**18" Black Pneumatic Tire & Black Spoke Aluminum Wheel** NXT-31538 n/c

Note: Rim, hub & spokes are anodized black.

**18" Airless Inserts** (E2213) NXT-31021 \$ 77

**18" Black Aluminum Handrim** NXT-31035 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**20" Pneumatic Spoke Wheel with Aluminum Handrim** NXT-31540 \$ 204

**20" Airless Inserts** (E2213) NXT-31022 \$ 77

**22" Pneumatic Spoke Wheel with Aluminum Handrim** NXT-31541 \$ 204

**22" Airless Inserts** (E2213) NXT-31023 \$ 77

Note: Not available with Short Tilt

Note: Not available with dual post armrests.

**24" Pneumatic Spoke Wheel with Aluminum Handrim** NXT-31542 \$ 204

**24" Airless Inserts** (E2213) NXT-31024 \$ 77

Note: Not available with Short Tilt

Note: Not available with dual post armrests.

**Omit Handrims from 20", 22" or 24" Wheels** NXT-31006

**Wheel Lock Options**

**Push to Lock - Standard** NXT-31501 n/c

**Pull to Lock** NXT-31506 n/c

**Extension Handles** (E0961) NXT-31801 \$ 61

**Foot Activated Wheel Lock** (E2228) NXT-31508 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Not available with ALL vent and battery trays.

**Hand Activated Brake and Wheel Lock System** (E2228) NXT-30275 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Used to control direction, speed, brake and lock wheels from strolling position.

**Combination Dual Hand Activated Tilt, Brake and Wheel Lock System** (E2228) NXT-30277 n/c

Note: Includes manual wheel locks to be used with public transportation.

Note: Not available with 12" mag wheels with low pro poly tires.

**Rear Anti-Tipper Options**

**Curb Assist Anti-Tip Assembly Flip-Up, Adjustable** (E0971) NXT-31600 \$ 97

Note: Pair includes one anti-tip tube with curb assist feature.

**Standard Anti-Tip Assembly Flip-Up, Adjustable** (E0971) NXT-31603 \$ 97

**Notes:**

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Seat Pan Options

Note: ALL seat pans include 2" x 2" notches for growth adjustment.

**Fixed Solid Seat Pan - Bolteron** (E2231) NXT-39010 \$ 170

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Note: ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width \_\_\_\_\_" X Seat Depth \_\_\_\_\_"

**Fixed Solid Seat Pan - Alumi-num** (E2231) NXT-39030 \$ 200

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

10" wide x 12" deep 16" wide x 18" deep

12" wide x 14" deep 18" wide x 18" deep

14" wide x 16" deep

**Fixed, Width Adjustable Solid Seat Pan - Aluminum** (E2231) NXT-39040 \$ 290

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with NXT frame.

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

10"-12" wide x 12" deep 14"-16" wide x 16" deep

12"-14" wide x 14" deep 16"-18" wide x 18" deep

### Mounting Hardware Options

Note: Use to convert fixed seat pan to removable.

**Removable Pan Hardware Kit** (K0108) NXT-EZLITE \$ 198

Note: Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Note: Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

### Accessories

**Medical Necessities Pouch** (K0108) NXT-23207 \$ 68

Note: Black nylon.

**Footrest Hanger Pads - Pair** (K0108) SFDI-9HP \$ 92

Note: Wrap-around pads made from black neoprene.

**Frame Pads - Pair** (K0108) 10911-2 \$ 82

Note: Wrap-around pads made from black neoprene.

Note: Typically used to pad exposed back post tubes.

**Flip-Down Rigidizing Bar** (K0108) NXT-23210 \$ 194

Note: Mounts to back post tubes to add rigidity to frame.

Note: Does NOT adjust in width to grow with frame.

**Freedom Backpack** (K0108) NXT-10900 \$ 56

Note: Black Nylon only with Freedom Designs logo.

Note: Monogramming is not available on Freedom Backpack.

**Freedom Canopy - Flip-Back** (K0108) NXT-10905 \$ 128

Note: Canopy assembly includes removable mounting hardware.

Note: Canopy is available in black only with Freedom Designs logo.

Note: Canopy includes rear flap with zippered storage pocket.

Note: Canopy includes clear-view window.

### Anterior Support Surfaces – (Trays)

**3/8" Lexan Tray - Channel Bracket Hardware** (E0950) 08-2000 \$ 271

**3/8" Lexan Tray - Adjustable Knob Hardware & Ribs** (E0950) 08-6000 \$ 219

**1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs** (E0950) 08-13000 \$ 183

**1/4" Lexan Tray - Hook & Loop Straps** (E0950) 08-16000 \$ 106

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

### Notes: