



NXT
Generation NEXT

Pediatric Manual Wheelchair



Yes, you can!

Clear Form

Email Form to
Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: _____	Buyer: _____	Quote _____	Order _____	PO# _____
Company: _____		Location: _____		
Marked For: _____		Phone: _____		
Email: _____		Fax: _____		

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Transport Optional Base

NXT Pediatric Folding	(E1234) NXT-30300	\$2800
<i>Transport Optional</i>		
NXT Pediatric Rigid	(E1233) NXT-30310	\$2700
<i>Transport Optional</i>		
Note: Transportation Option selected separately below.		
Note: Includes Freedom Future Fit Program		

Transport Package	NXT-99998	\$ 250
Note: Must be installed by Freedom Designs manufacturing.		
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.		
Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.		

Vent Tray Options

Note: All vent and battery trays best with standard 90° bend push handles and push button angle adjustable stroller handle extensions.

Vent Tray - Fixed	(E1029) 11-7000	\$ 439
Note: Not with Foot Activated Tilt.		
Battery Tray - Fixed	(K0108) 11-8000	\$ 429
Quick Release Tray Hardware	(E1028) 563	\$ 272
Note: For LP10 style vent and battery tray on frames 16" wide and wider.		
Note: Vent tray is 4" narrower than frame width and 16" deep.		
Note: Battery tray is 4" narrower than frame width and 10" deep.		
Back Post Mounted Vent Bracket	(E1029) 11-1100	\$ 550
Note: For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.		

Transport Included Base

NXT Pediatric Folding	(E1234) NXT-30301	\$2950
<i>Transport Included</i>		
NXT Pediatric Rigid	(E1233) NXT-30311	\$2850
<i>Transport Included</i>		
Note: Includes Freedom Future Fit Program		
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.		

Frame Color Options

Sugar Cookie	n/c	Blue Blast	n/c
Neon Yellow	n/c	Electric Blue	n/c
Yellow	n/c	Galaxy Blue	n/c
Bright Pink	n/c	Passion Purple	n/c
Pink Parfait	n/c	Octane Orange	n/c
Sugar Plum	n/c	Red Sunset	n/c
Candy Violet	n/c	Really Red	n/c
Candy Raspberry	n/c	Fiery Red	n/c
Extreme Green	n/c	Twilight	n/c
Sparkle Green	n/c	Stormy Gray	n/c
Mystic Teal	n/c	Black	n/c
Pacific Blue	n/c		

Seat Rail Options

Note: Not available with ALL Contracture & 3D Dynamic footrest systems.

LEAF - Lower Extremity Accommodation Frame

LEAF Seat Frame- Pair	(K0108) NXT-30324	\$ 714
LEAF Seat Frame- Left	(K0108) NXT-30324-L	\$ 357
LEAF Seat Frame- Right	(K0108) NXT-30324-R	\$ 357

Note: Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

SRO - Seat Rail Only

SRO Seat Frame- Pair	(K0108) NXT-30325	\$ 306
SRO Seat Frame- Left	(K0108) NXT-30325-L	\$ 153
SRO Seat Frame- Right	(K0108) NXT-30325-R	\$ 153

Note: SRO is manufactured without standard footrest hanger attachments.

Tilt Activation Options

Dual Lever Hand Activated Tilt	NXT-30272	n/c
Foot Activated Tilt	NXT-30273	n/c
Note: Not with vent & battery trays.		
Combination Dual Hand Activated Tilt, Brake and Wheel Lock System	NXT-30277	n/c

Note: Includes manual wheel locks to be used with public transportation.
Note: Not available with 12" mag wheels with low pro poly tires.

Frame Dimensions

Specify Frame Width Setting

10"	11"	12"	Adjusts 10" to 12"	NXT-30250	n/c
11"	12"	13"	Adjusts 11" to 13"	NXT-30210	n/c
12"	13"	14"	Adjusts 12" to 14"	NXT-30220	n/c
13"	14"	15"	Adjusts 13" to 15"	NXT-30230	n/c
14"	15"	16"	Adjusts 14" to 16"	NXT-30240	n/c
15"	16"	17"	Adjusts 15" to 17"	NXT-30290	n/c
16"	17"	18"	Adjusts 16" to 18"	NXT-30260	n/c
17"	18"		Adjusts 17" to 18"	NXT-30261	n/c
18"	Not width adjustable			NXT-30262	n/c
19"	Not width adjustable, not with short frame			NXT-30263	n/c
20"	Not width adjustable, not with short frame		(E2201)	NXT-30264	\$ 400

Note: Frame widths measures across and from outside of seat frame tubes.

Note: Seat frame widths of 19" and 20" not available with short frame.

Note: Vent and battery trays only with 14" wide and greater frame.

NXT Rider Weight Capacity is 225 lbs.
All prices subject to change without notice.

NXT Pediatric
Note: All Specifications and dimensions are approximate.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Type Options

Specify Frame Type and Frame Length

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

Short Frame and Short Tilt Bracket	NXT-30511	n/c
---	-----------	-----

Specify Back Post Depth Setting

12"	NXT-30504	n/c
13"	NXT-30505	n/c
14"	NXT-30506	n/c
15"	NXT-30521	n/c
16"	NXT-30522	n/c
* 17" Must select 2" extension	NXT-30523	n/c
* 18" Must select 2" extension	NXT-30524	n/c
2" Extension - Pair (K0108)	NXT-30525	\$ 26

Specify Back Post Location (if other than specified above) _____"

Note: Seat rail to floor height range is 12.5" to 16.5".

Note: Seat frame widths of 19" and 20" not available with short frame.

* Note: Must choose extension for 17" and 18" frame lengths.

Note: 22" & 24" wheels not available with short tilt bracket.

Note: All vent and battery trays not available with short frame.

Medium Frame and Short Tilt Bracket	NXT-30512	n/c
--	-----------	-----

Specify Back Post Depth Setting

14"	NXT-30506	n/c
15"	NXT-30521	n/c
16"	NXT-30522	n/c
17"	NXT-30523	n/c
18"	NXT-30524	n/c
* 19" Must select 2" extension	NXT-30527	n/c
* 20" Must select 2" extension (E2203)	NXT-30528	\$ 500
* 21" Must select 3" extension (E2203)	NXT-30529	\$ 500
2" Extension - Pair (K0108)	NXT-30525	\$ 26
3" Extension - Pair (K0108)	NXT-30530	\$ 26

Specify Back Post Location (if other than specified above) _____"

Note: Seat rail to floor height range is 12.5" to 16.5".

* Note: Must choose extension for 19" to 21" frame lengths.

Note: 22" & 24" wheels not available with short tilt bracket.

Medium Frame and Tall Tilt Bracket	NXT-30540	n/c
---	-----------	-----

Specify Back Post Depth Setting

14"	NXT-30506TP	n/c
15"	NXT-30521TP	n/c
16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
* 19" Must select 2" extension	NXT-30527TP	n/c
* 20" Must select 2" extension (E2203)	NXT-30528TP	\$ 500
* 21" Must select 3" extension (E2203)	NXT-30529TP	\$ 500
* 22" Must select 4" extension (E2204)	NXT-30535TP	\$ 900
2" Extension - Pair (K0108)	NXT-30525	\$ 26
3" Extension - Pair (K0108)	NXT-30530	\$ 26
4" Extension - Pair (K0108)	NXT-30531	\$ 26

Specify Back Post Location (if other than specified above) _____"

Note: Tall tilt bracket increases seat rail to floor 1.5".

Note: Seat rail to floor height range is 14" to 18".

* Note: Must choose extension for 19" to 22" frame lengths.

Frame Type Options - Continued

Pediatric Long Frame & Tall Tilt Bracket	(E1298) NXT-30520P	\$ 153
---	--------------------	--------

Specify Back Post Depth Setting

16"	NXT-30522TP	n/c
17"	NXT-30523TP	n/c
18"	NXT-30524TP	n/c
19"	NXT-30527TP	n/c
20"	(E2203) NXT-30528TP	\$ 500
21"	(E2203) NXT-30529TP	\$ 500
22"	(E2204) NXT-30535TP	\$ 900
* 23" Must select 2" extension (E2204)	NXT-30537TP	\$ 900
* 24" Must select 2" extension (E2204)	NXT-30538TP	\$ 900
* 25" Must select 3" extension (E2204)	NXT-30542TP	\$ 900
* 26" Must select 4" extension (E2204)	NXT-30543TP	\$ 900
2" Extension - Pair (K0108)	NXT-30525	\$ 26
3" Extension - Pair (K0108)	NXT-30530	\$ 26
4" Extension - Pair (K0108)	NXT-30531	\$ 26

Specify Back Post Location (if other than specified above) _____"

Note: Tall tilt bracket increases seat rail to floor 1.5".

Note: Seat rail to floor height range is 14" to 18".

* Note: Must choose extension for 23" to 26" frame lengths.

5° Anterior Tilt Seat Frame	NXT-30276	n/c
------------------------------------	-----------	-----

Note: Limits posterior tilt to 40°. Includes split collars for adjustment.

Note: Medium with Tall Tilt and Long frames only.

Note: NXT has 4" of depth adjustment for growth. Depth is measured from front of back-post to the front end of the seat frame. For back post location add seating & hardware thickness (typically 2") to client's actual seat depth. Example: 12" client seat depth + 2" back & hardware thickness = 14" back post location.

Back Post Options

Specify Push Handle, Back Post and Extensions

Standard 90° Bend Push Handle Back Post Tubes	NXT-30518	n/c
Ergo Bend Push Handle Back Post Tubes	NXT-30519	n/c

Note: All push/stroller handle extensions not recommended with Ergo Bend.

Seating Dynamics Rocker Back	(K0108) NXT-30532	\$1985
-------------------------------------	-------------------	--------

Note: Angle adjustable from 72° to 108° in 6° increments.

Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.

Note: Locks into a non-dynamic position when necessary.

Note: Does not allow the back post tubes to fold forward.

Back Post Height Setting

Note: Back Posts are measured from top of seat rail to top of handle bend.

Short Push Handle - 20" high	NXT-30515	n/c
Medium Push Handle - 23" high	NXT-30514	n/c
Tall Push Handle - 26" high	NXT-30516	n/c

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Stroller Handle Extension Options

Note: All stroller handle ext. add approximately 6.5" to push handle height.
 Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.



Push Handle Extension - Width Adjustable - Mounts direct to back post push handle NXT-30560 n/c

Stroller Handle Extension - Width Adjustable, Removable - Mounts direct to back post push handle **NEW** NXT-30590 \$ 95
 Note: Not available with Super Reclining Back.

Push Button, Angle-Adjustable Stroller Handle Extension - Mounts direct to back post push handle NXT-30570 \$ 112
 Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.
 Note: Not removable, but allows frame to fold.

Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty - Mounts direct to back post tubes NXT-30575 \$ 179
 Note: Available 10W-20w wheelchair frames.
 Note: Not removable, but allows frame to fold.

Push Button, Angle-Adjustable, Folding Push Handle Extension - Mounts direct to straight back post tubes NXT-30580 \$ 112
 Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.
 Note: Not removable, but allows frame to fold.

Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty - Mounts direct to straight back post tubes NXT-30585 \$ 179
 Note: Available 10W-20w wheelchair frames.
 Note: Not removable, but allows frame to fold.

Individual Stroller Handle Extensions Removable - Pair NXT-30517 \$ 102

Legrests, Footrests and Extension Options

Legrests

Omit Legrest Hangers & Footrests	NXT-33000	
75° Pop-Up Swing Away Hangers	NXT-30280	n/c
90° Pop-Up Swing Away Hangers	NXT-30281	n/c
90° Pop-Up Swing Away Elevating Legrests (ELRs)	(E0990) NXT-30279	\$ 255
75° Latch Style Swing Away Hangers	NXT-30282	n/c
90° Latch Style Swing Away Hangers	NXT-30283	n/c

Note: Not available with 8" adjustable extensions.
 Note: Not available with all length adjustable extensions.
 Note: Not available with 6.5" heavy duty extensions.
 Note: Not available with 8" adjustable extensions.

Footplates & Footboards

Omit Footrest Plates	NXT-36404	
Composite Flip-Up	NXT-17070	n/c
Kids Aluminum Angle Adjustable, Flip-Up	(K0040) NXT-36020	\$ 160
Heel Loops	(E0951) NXT-36022	\$ 41
Calf Strap - Webbing	(K0038) NXT-36024	\$ 41
Adult Aluminum Angle Adj. Flip-up w/ Heel Loops	(K0040/ E0951) NXT-36021	\$ 180
Omit Heel Loops	NXT-36023	n/c
Aluminum Angle Adjustable Footboard Removable	(K0040) NXT-36029	\$ 160

Note: For frames 14" wide and narrower.
 Note: For frames 15" wide and greater.

MAPS Footplate Mounting Hardware

MAPS (Multi Axis Positioning System) - Pair	NXT-36040	\$ 306
MAPS (Multi Axis Positioning System) - Left	NXT-36040-L	\$ 153
MAPS (Multi Axis Positioning System) - Right	NXT-36040-R	\$ 153

Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.

Note: Must choose aluminum footplates.

Note: MAPS not available with composite footplates or aluminum footboard.

Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	1.5"-3"	1.5"-3"	1.5"-3"	1.5"-3"	
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	
8" Adjustable	14"-17"	12.5"-16.5"			

MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount					
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	
8" Adjustable	12.5"-15.5"	12"-16"			

Note: Length from top of footrest plate to top of seat rail.

Note: Without extensions bottom hanger hole not accessible for hardware.

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Continue To Next Page

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Legrests, Footrests and Extension Options

Legrest Extensions

4" Length Adjustable Extensions	NXT-30285	n/c
6" Length Adjustable Extensions	NXT-30286	n/c
8" Length Adjustable Extensions	NXT-30287	n/c

Note: 8" & 10" Not available with all 90° hangers.

4.5" Heavy Duty ELR Extensions	NXT-36401	\$ 56
4.5" Heavy Duty Extensions	NXT-36402	\$ 56
6.5" Heavy Duty Extensions	NXT-36403	\$ 66

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

High Mount Brackets	NXT-36407	\$ 102
----------------------------	-----------	--------

Note: Includes flip-up mounting hardware for future growth adjustment.

Note: 75° hangers use a tube & clamp style hardware.

Specify Knee to Heel Measurement _____" (see matrix)

Notes:

Specialized Footrest Systems

Note: ALL legrest measurements are from the seat rail to footplate.
Note: Not available with LEAF or SRO seat rail options.

Super Contracture Footrest Hanger (K0108)

Note: Super Contracture for frames 14" wide and greater.

Contracture Short - 11" to 15"	NXT-30291-2	\$3570
---------------------------------------	-------------	--------

Note: See Matrix for seat to rail footplate measurements.

Includes:

Adj Contracture Upper	NXT-6302-03	Incl.
Footrest Adapter HD	NXT-6220-01	Incl.
Adj Plantar/Dorsi Ext. Short	NXT-6401-02	Incl.

Footplate Options

4" x 6" HD Alu Footplates	NXT-6503	n/c
5" x 7" HD Alu Footplates	NXT-6504	n/c

Note: for frames 16" wide and greater only.

Contracture Long - 15" to 19"	NXT-30292-2	\$3570
--------------------------------------	-------------	--------

Note: See Matrix for seat to rail footplate measurements.

Includes:

Adj Contracture Upper	NXT-6302-03	Incl.
Footrest Adapter HD	NXT-6220-01	Incl.
Adj Plantar/Dorsi Ext. Long	NXT-6402-02	Incl.

Footplate Options

4" x 6" HD Alu Footplates	NXT-6503	n/c
5" x 7" HD Alu Footplates	NXT-6504	n/c

Note: Super Contracture is both swing-away and removable.

Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.

Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.

Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

3D Dynamic Footrest Hanger (K0108)

Note: 3D Dynamic for frames 14" wide and greater.

3D Dynamic Footrest - Hi Pivot	NXT-30297-2	\$6215
---------------------------------------	-------------	--------

Includes:

HiPivot Dyn Adj Knee & Dyn Tel Up	NXT-6304-02	Incl.
Footrest Adapter HD	NXT-6220-01	Incl.

Extension Options

Note: See Matrix for seat to rail footplate measurements.

Dynamic Adj Plantar/Dorsi Ext, 0-2"	NXT-6405-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 1-3"	NXT-6409-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 2-4"	NXT-6406-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 3-5"	NXT-6410-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 4-6"	NXT-6407-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 5-7"	NXT-6411-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 6-8"	NXT-6408-01	n/c

Footplate Options

4" x 6" HD Alu Footplates	NXT-6503	n/c
5" x 7" HD Alu Footplates	NXT-6504	n/c

Note: for frames 16" wide and greater only.

6" x 10" HD Alu Footplates	NXT-6505	n/c
-----------------------------------	----------	-----

Note: for frames 18" wide and greater only.

Specialized Footrest Seat Rail to Footplate Measurements

Part Number	Extension	Super Contracture	3D Dynamic Hi Pivot	3D Dynamic STD
6402-01	Short	10.5" - 13.5"		
6402-02	Long	13.5" - 18.5"		
6405-01	0"-2" Extension		6.25" - 10.75"	12.75-15.25"
6409-01	1"-3" Extension		7.25-11.75"	13.75-16.25"
6406-01	2"-4" Extension		8.25-12.75"	14.75-17.25"
6410-01	3"-5" Extension		9.25-13.75"	15.75-18.25"
6407-01	4"-6" Extension		10.25-14.75"	16.75-19.25"
6411-01	5"-7" Extension		11.25-15.75"	17.75-20.25"
6508-01	6"-8" Extension		12.25-16.75"	18.75-21.25"

Note: Length from top of footrest plate to top of seat rail.

Note: 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Specialized Footrest Systems - Continued

3D Dynamic Footrest - STD NXT-30298-2 \$5775

Includes:

Dynamic Adj Knee & Dyn Tel Up NXT-6304-01 Incl.

Footrest Adapter HD NXT-6220-01 Incl.

Extension Options

Note: See Matrix for seat to rail footplate measurements.

Dynamic Adj Plantar/Dorsi Ext, 0-2" NXT-6405-01 n/c

Dynamic Adj Plantar/Dorsi Ext, 1-3" NXT-6409-01 n/c

Dynamic Adj Plantar/Dorsi Ext, 2-4" NXT-6406-01 n/c

Dynamic Adj Plantar/Dorsi Ext, 3-5" NXT-6410-01 n/c

Dynamic Adj Plantar/Dorsi Ext, 4-6" NXT-6407-01 n/c

Dynamic Adj Plantar/Dorsi Ext, 5-7" NXT-6411-01 n/c

Dynamic Adj Plantar/Dorsi Ext, 6-8" NXT-6408-01 n/c

Footplate Options

4" x 6" HD Alu Footplates NXT-6503 n/c

5" x 7" HD Alu Footplates NXT-6504 n/c

Note: for frames 16" wide and greater only.

6" x 10" HD Alu Footplates NXT-6505 n/c

Note: for frames 18" wide and greater only.

Note: 3D Dynamic is both swing-away and removable.

Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.

Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.

Note: Angle adjustable footplate included.

Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

Armrests and Armrest Pad Options

Omit Armrests NXT-38209

Back Post Mounted Armrests

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

10.5" Long Flip-Back NXT-38204 n/c

Note: Not available with full length arm pads

13.5" Long Flip-Back NXT-38202 n/c

10.5" Long Flip-Back Locking (E0973) NXT-38201 \$ 245

Note: Not available with full length arm pads

13.5" Long Flip-Back Locking (E0973) NXT-38200 \$ 245

10.5" to 13.5" Long Flip-Back, Telescoping (E0973) NXT-38208 \$ 240

Note: Full length pad on 13.5" or longer armrests only.

10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973) NXT-38210 \$ 265

Note: Full length pad on 13.5" or longer armrests only.

17" to 20" Long Flip-Back Locking, Telescoping (E0973) NXT-38215 \$ 265

Foam Wraps - Standard NXT-30216 n/c

Desk Length Pads NXT-38213 n/c

Full Length Pads NXT-38214 n/c

Note: Full length pad on 13.5" or longer armrests only.

Frame Mounted Armrests

Single Post Fixed Height, Removable - Short (K0015) NXT-38222 \$ 225

Single Post Fixed Height, Removable - Tall (K0015) NXT-38224 \$ 225

Note: Single Post Fixed Height arms include a removable side guard.

Note: Adjustable in height with tools only. See matrix for height ranges.

Single Post Height Adjustable, Removable - Short (E0973) NXT-38217 \$ 265

Single Post Height Adjustable, Removable - Tall (E0973) NXT-38203 \$ 265

Dual Post Fixed Height Removable NXT-38205 n/c

Dual Post Flip-Back Height Adjustable, Removable (E0973) NXT-38206 \$ 275

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

Desk Length Pads NXT-38213 n/c

Full Length Pads NXT-38214 n/c

Specify Seat Rail to Top of Armrest Measurement _____"

Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Single Post Fxd Hgt Short		8" - 10.5"	8" - 10.5"
Single Post Fxd Hgt Short		9" - 11.5"	9" - 11.5"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

Continue To Next Page

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Choose NXT Seat Rail to Floor Height Below:

Note: Use chart to determine seat to floor height. Seat to floor height determined by frame type, rear wheels, front caster size and tilt pivot height.
Note: Caster sizes vary slightly in height by type.
Note: Seat to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 3/8".
Note: Seat pan choice increases seat rail to floor height 3/8" to 3/4".

NXT Wheel / Caster Seat to Floor Heights

Wheel	Caster	Pivot	Height	Wheel	Caster	Pivot	Height
12"	3"	Short	12.5-13"	20"	3"	Tall	N/A
12"	3"	Tall	14-14.5"	20"	4"	Short	14-15"
12"	4"	Short	13.5"	20"	4"	Tall	15.5-16.5"
12"	4"	Tall	15"	20"	5"	Short	14.5-15.5"
12"	5"	Short	15.5"	20"	5"	Tall	16-17"
12"	5"	Tall	17"	20"	6"	Short	15.5-16"
12"	6"	Short	16"	20"	6"	Tall	17-17.5"
12"	6"	Tall	17.5"	20"	7"	Short	16.5"
16"	3"	Short	12.5-13"	20"	7"	Tall	18"
16"	3"	Tall	14-14.5"	22"	3"	Short	N/A
16"	4"	Short	13.5-14"	22"	3"	Tall	N/A
16"	4"	Tall	15-15.5"	22"	4"	Short	N/A
16"	5"	Short	14.5"	22"	4"	Tall	16.5"
16"	5"	Tall	16"	22"	5"	Short	N/A
16"	6"	Short	16	22"	5"	Tall	16.5-17"
16"	6"	Tall	17.5	22"	6"	Short	N/A
16"	7"	Short	16.5"	22"	6"	Tall	17-17.5"
16"	7"	Tall	18"	22"	7"	Short	N/A
18"	3"	Short	12.5-13"	22"	7"	Tall	18"
18"	3"	Tall	14-14.5"	24"	3"	Short	N/A
18"	4"	Short	13.5-14.5"	24"	3"	Tall	N/A
18"	4"	Tall	15-16.5"	24"	4"	Short	N/A
18"	5"	Short	14.5-15.5"	24"	4"	Tall	N/A
18"	5"	Tall	16-17"	24"	5"	Short	N/A
18"	6"	Short	15.5"	24"	5"	Tall	N/A
18"	6"	Tall	17"	24"	6"	Short	N/A
18"	7"	Short	N/A	24"	6"	Tall	17.5"
18"	7"	Tall	N/A	24"	7"	Short	N/A
20"	3"	Short	N/A	24"	7"	Tall	18"

Front Wheels - Caster Options

Omit Caster Wheels	NXT-34333	
3" Polyurethane	NXT-34888	n/c
3" Polyurethane Light-Up	NXT-34111	\$ 101
4" Polyurethane	NXT-40000	n/c
4" Polyurethane Light-Up	NXT-41000	\$ 101
4" x 1.4" Poly, Aluminum Mag	NXT-43000	\$ 175
5" x 1" Poly, Plastic Mag	NXT-34555	n/c
5" Polyurethane Light-Up	NXT-35050	\$ 101
5" x 1.5" Poly, Plastic Mag	NXT-34222	\$ 87
5" x 1.4" Poly, Aluminum Mag	NXT-42000	\$ 175

Front Wheels - Caster Options

6" x 1" Poly, Plastic Mag	NXT-34666	n/c
6" x 1.5" Poly, Plastic Mag	NXT-34999	\$ 87
6" x 1.4" Poly, Aluminum Mag	NXT-44000	\$ 175
6" x 2" Poly with Tread	NXT-34444	\$ 92
6" x 1" Pneumatic Plastic Mag	(E2214/ E2215) NXT-34667	\$ 77
7" x 1" Poly, Plastic Mag	NXT-34777	n/c
7" x 2" Polyurethane	NXT-34200	\$ 118

Note: All 7" casters not available with all 12" rear wheels.

Rear Wheels - Tires - Accessories

Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.

Omit Rear Wheels and Axles	NXT-31524	
Mag Wheels		
12" Low Profile Poly Mag	(E2218) NXT-31526	\$ 36
Note: Not available with any foot/hand activated brake/wheel lock system.		
12" Semi Pneumatic Plastic Mag	NXT-31523	\$ 168
Note: Not available with airless inserts.		
12" Pneumatic Plastic Mag	NXT-31525	n/c
12" Airless Inserts	(E2213) NXT-31019	\$ 77
Note: All 12" wheels not available with all 7" casters.		
16" Pneumatic Plastic Mag	NXT-31530	\$ 41
16" Airless Inserts	(E2213) NXT-31020	\$ 77
16" Aluminum Handrims	NXT-31014	\$ 169
18" Pneumatic Plastic Mag	NXT-31533	\$ 41
18" Airless Inserts	(E2213) NXT-31021	\$ 77
18" Aluminum Handrims	NXT-31015	\$ 169
20" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31527	n/c
20" Airless Inserts	(E2213) NXT-31022	\$ 77
20" Plastic Coated Handrims	NXT-31016	\$ 119
Note: Plastic Coated Handrims NOT available with Drum Braking Systems.		
22" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31528	n/c
22" Airless Inserts	(E2213) NXT-31023	\$ 77
22" Plastic Coated Handrims	NXT-31017	\$ 119
Note: Not available with Short Tilt		
Note: Not available with dual post armrests.		
Note: Plastic Coated Handrims NOT available with Drum Braking Systems.		
24" Pneumatic Plastic Mag with Aluminum Handrim	NXT-31529	n/c
24" Airless Inserts	(E2213) NXT-31024	\$ 77
24" Plastic Coated Handrims	NXT-31018	\$ 119
Note: Not available with Short Tilt		
Note: Not available with dual post armrests.		
Note: Plastic Coated Handrims NOT available with Drum Braking Systems.		
Omit Handrims from 20", 22" or 24" Wheels	NXT-31006	

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Rear Wheels - Tires - Accessories - Cont.

Spoke Wheels

Note: All spoke wheels not available plastic coated handrims.

16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel NXT-31534 n/c

Note: Rim, hub & spokes are anodized silver.

16" Airless Inserts (E2213) NXT-31020 \$ 77

16" Silver Aluminum Handrim NXT-31030 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

16" Black Pneumatic Tire & Black Spoke Aluminum Wheel NXT-31535 n/c

Note: Rim, hub & spokes are anodized black.

16" Airless Inserts (E2213) NXT-31020 \$ 77

16" Black Aluminum Handrim NXT-31031 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel NXT-31537 n/c

Note: Rim, hub & spokes are anodized silver.

18" Airless Inserts (E2213) NXT-31021 \$ 77

18" Silver Aluminum Handrim NXT-31034 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

18" Black Pneumatic Tire & Black Spoke Aluminum Wheel NXT-31538 n/c

Note: Rim, hub & spokes are anodized black.

18" Airless Inserts (E2213) NXT-31021 \$ 77

18" Black Aluminum Handrim NXT-31035 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

20" Pneumatic Spoke Wheel with Aluminum Handrim NXT-31540 \$ 204

20" Airless Inserts (E2213) NXT-31022 \$ 77

22" Pneumatic Spoke Wheel with Aluminum Handrim NXT-31541 \$ 204

22" Airless Inserts (E2213) NXT-31023 \$ 77

Note: Not available with Short Tilt

Note: Not available with dual post armrests.

24" Pneumatic Spoke Wheel with Aluminum Handrim NXT-31542 \$ 204

24" Airless Inserts (E2213) NXT-31024 \$ 77

Note: Not available with Short Tilt

Note: Not available with dual post armrests.

Omit Handrims from 20", 22" or 24" Wheels NXT-31006

Wheel Lock Options

Push to Lock - Standard NXT-31501 n/c

Pull to Lock NXT-31506 n/c

Extension Handles (E0961) NXT-31801 \$ 61

Foot Activated Wheel Lock (E2228) NXT-31508 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Not available with ALL vent and battery trays.

Hand Activated Brake and Wheel Lock System (E2228) NXT-30275 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Used to control direction, speed, brake and lock wheels from strolling position.

Combination Dual Hand Activated Tilt, Brake and Wheel Lock System (E2228) NXT-30277 n/c

Note: Includes manual wheel locks to be used with public transportation.

Note: Not available with 12" mag wheels with low pro poly tires.

Rear Anti-Tipper Options

Curb Assist Anti-Tip Assembly Flip-Up, Adjustable (E0971) NXT-31600 \$ 97

Note: Pair includes one anti-tip tube with curb assist feature.

Standard Anti-Tip Assembly Flip-Up, Adjustable (E0971) NXT-31603 \$ 97

Notes:

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Seat Pan Options

Note: ALL seat pans include 2" x 2" notches for growth adjustment.

Fixed Solid Seat Pan - Bolteron (E2231) NXT-39010 \$ 170

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Note: ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width _____" X Seat Depth _____"

Fixed Solid Seat Pan - Alumi-num (E2231) NXT-39030 \$ 200

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

10" wide x 12" deep 16" wide x 18" deep

12" wide x 14" deep 18" wide x 18" deep

14" wide x 16" deep

Fixed, Width Adjustable Solid Seat Pan - Aluminum (E2231) NXT-39040 \$ 290

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with NXT frame.

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

10"-12" wide x 12" deep 14"-16" wide x 16" deep

12"-14" wide x 14" deep 16"-18" wide x 18" deep

Mounting Hardware Options

Note: Use to convert fixed seat pan to removable.

Removable Pan Hardware Kit (K0108) NXT-EZLITE \$ 198

Note: Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Note: Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

Accessories

Medical Necessities Pouch (K0108) NXT-23207 \$ 68

Note: Black nylon.

Footrest Hanger Pads - Pair (K0108) SFDI-9HP \$ 92

Note: Wrap-around pads made from black neoprene.

Frame Pads - Pair (K0108) 10911-2 \$ 82

Note: Wrap-around pads made from black neoprene.

Note: Typically used to pad exposed back post tubes.

Flip-Down Rigidizing Bar (K0108) NXT-23210 \$ 194

Note: Mounts to back post tubes to add rigidity to frame.

Note: Does NOT adjust in width to grow with frame.

Freedom Backpack (K0108) NXT-10900 \$ 56

Note: Black Nylon only with Freedom Designs logo.

Note: Monogramming is not available on Freedom Backpack.

Freedom Canopy - Flip-Back (K0108) NXT-10905 \$ 128

Note: Canopy assembly includes removable mounting hardware.

Note: Canopy is available in black only with Freedom Designs logo.

Note: Canopy includes rear flap with zippered storage pocket.

Note: Canopy includes clear-view window.

Anterior Support Surfaces – (Trays)

3/8" Lexan Tray - Channel Bracket Hardware (E0950) 08-2000 \$ 271

3/8" Lexan Tray - Adjustable Knob Hardware & Ribs (E0950) 08-6000 \$ 219

1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs (E0950) 08-13000 \$ 183

1/4" Lexan Tray - Hook & Loop Straps (E0950) 08-16000 \$ 106

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Notes: