



Yes, you can.®



PRO-CG™-SE Adult Manual Wheelchair

Clear Form

Email Form to
Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: _____ Buyer: _____ Quote Order Install Seating? PinDot Matrix
 Company: _____ PO#: _____
 Marked For: _____ Location: _____
 Email: _____ Phone: _____

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
-------------	-------	----------	--------	-------------	-------	----------	--------

Base Options

A seating system must be purchased separately. No seating system is provided with the PRO-CG frame. No seat or back upholstery is available as an option.

Transport Optional

PRO-CG Adult SE (E1161) PRO-30312 \$2775

Transportation Option selected separately below.

SE frame does not include Freedom Future Fit Program

Transport Brackets – SE PRO-99997 \$ 300

Must be installed by Freedom Designs manufacturing.

Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.

Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

Frame Dimensions

16" Base Width PRO-30240

Specify Seat Width Setting

14"	PRO-14W	n/c
15"	PRO-15W	n/c
16"	PRO-16W	n/c

Seat frame adjusts 14" to 16" wide in 1" increments.

Specify Seat Depth Setting

14"	PRO-30506	n/c
15"	PRO-30521	n/c
16"	PRO-30522	n/c
17"	PRO-30523	n/c
18"	PRO-30524	n/c
19" - 2" extension required	PRO-30527	n/c
20" - 2" extension included	(E2203) PRO-30528	\$ 550
21" - 4" extension included	(E2203) PRO-30529	\$ 550
22" - 4" extension included	(E2204) PRO-30535	\$ 990

18" Base Width PRO-30260

Specify Seat Width Setting

16"	PRO-16W	n/c
17"	PRO-17W	n/c
18"	PRO-18W	n/c

Seat frame adjusts 16" to 18" wide in 1" increments.

Specify Seat Depth Setting

16"	PRO-30522	n/c
17"	PRO-30523	n/c
18"	PRO-30524	n/c
19"	PRO-30527	n/c
20"	(E2203) PRO-30528	\$ 550
21" - 2" extension included	(E2203) PRO-30529	\$ 550
22" - 2" extension included	(E2204) PRO-30535	\$ 990
23" - 4" extension included	(E2204) PRO-30536	\$ 990
24" - 4" extension included	(E2204) PRO-30537	\$ 990

Notes:

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
Frame Dimensions (continued)							
20" Base Width				PRO-30261			
Specify Seat Width Setting				Specify Seat Depth Setting			
18"		PRO-18W	n/c	18"		PRO-30524	n/c
19"		PRO-19W	n/c	19"		PRO-30527	n/c
20"	(E2201)	PRO-20W	\$ 500	20"	(E2203)	PRO-30528	\$ 550
Seat frame adjusts 18" to 20" wide in 1" increments.				21"	(E2203)	PRO-30529	\$ 550
				22"	(E2204)	PRO-30535	\$ 990
				23" - 2" extension included	(E2204)	PRO-30536	\$ 990
				24" - 2" extension included	(E2204)	PRO-30537	\$ 990
				25" - 4" extension included	(E2204)	PRO-30538	\$ 990
				26" - 4" extension included	(E2204)	PRO-30539	\$ 990

Seat Rail Extensions			
2" Extension - Pair	(K0108)	PRO-30525	\$ 41
4" Extension - Pair	(K0108)	PRO-30531	\$ 41
Specify Back Post Location (if other than specified above) _____"			

Seat Rail Options

Not available with ALL Contracture & 3D Dynamic footrest systems.

LEAF - Lower Extremity Accommodation Frame				SRO - Seat Rail Only			
LEAF Seat Frame- Pair	(K0108)	PRO-30324	\$ 775	SRO Seat Frame- Pair	(K0108)	PRO-30325	\$ 325
Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.				SRO is manufactured without standard footrest hanger attachments.			

Vent Options

Vent Tray Base W/ Tray	(E1029)	PRO-30335	\$ 589	Back Post Mounted Vent Bracket	(E1029)	11-1100	\$ 550
Extended wheelchair base including vent tray. For LP10 style vent, base must be 16" or greater. Tray measures full width of the base frame and 12" deep				For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.			

Auxillary Device Mounts

O2 Holder - Aluminum w/ QR Hdwe	(K0108)	11-3100	\$ 548	IV Pole w/ Mounting Hardware	(K0108)	11-4000	\$ 270
---------------------------------	---------	---------	--------	------------------------------	---------	---------	--------

Tilt Activations

Dual Lever Hand Activated Tilt		PRO-30272	n/c	Foot Activated Tilt		PRO-30273	n/c
				Not available with vent tray.			

Frame Colors

Sugar Cookie		n/c	Blue Blast		n/c
Neon Yellow		n/c	Electric Blue		n/c
Yellow		n/c	Galaxy Blue		n/c
Bright Pink		n/c	Passion Purple		n/c
Pink Parfait		n/c	Octane Orange		n/c
Sugar Plum		n/c	Red Sunset		n/c
Candy Violet		n/c	Really Red		n/c
Candy Raspberry		n/c	Fiery Red		n/c
Extreme Green		n/c	Twilight		n/c
Sparkle Green		n/c	Stormy Gray		n/c
Mystic Teal		n/c	Black		n/c
Pacific Blue		n/c	Color to be Determined		

Notes:

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
-------------	-------	----------	--------	-------------	-------	----------	--------

Back Posts

Back Post Options

Fold-Down Back Posts		PRO-30545	n/c	Dynamic Rocker Back - Seating Dynamics	(E2398)	PRO-30542	\$1985
				Angle adjustable from 72° to 108° in 6° increments. Can be locked-out when needed.			
				Includes a set of 4 elastomers in different durometers to allow spring adjustment in the field.			
Reduced Shear Recline Back	(E1226)	PRO-30562	\$1138	Reduced Shear Recline Back with Angle-Adjustable Stroller Handle Extension		PRO-30564	\$1335

Recline range is 65° from 97° to 162°. Pivot point is 2.5" above the seat rail. Solid Back Insert installation requires 4-points of attachment hardware. Seat-to-Back hinges are not compatible with any Super Reclining Reduced Shear Back.

Push Handle Options

90° Bend Back Post Tubes		PRO-30518	n/c	Straight Back Post Tubes		PRO-30544	n/c
				Requires selection of Push Button Angle Adjustable Stroller Handle (PRO-30571 or PRO-30576). Not available with Reduced Shear Recline Back.			

Push Handle Height

Short - 20" High		PRO-30515	n/c	Tall - 26" High		PRO-30516	n/c
Medium - 23" High		PRO-30514	n/c				

Push Handle Extension Options



All push handle extensions except PRO-30561 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

Push Handle Extension - Width Adjustable		PRO-30561	n/c	Stroller Handle Extension - Width Adjustable, Removable		PRO-30590	\$ 95
				Not available with Straight Back Posts. Not available with Recline.			
Push Button, Angle-Adjustable Stroller Handle Extension		PRO-30571	\$ 112	Heavy Duty Push Button, Angle-Adjustable Stroller Handle Extension		PRO-30576	\$ 179
				Available on wheelchairs frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. PRO-30571 or PRO-30576 required with Straight Back Post Tubes (PRO-30544)			
Individual Stroller Handle Extensions Removable - Pair		PRO-30517	\$ 110				
				Available on any size frame. Metal hinge up to 180° range of adjustment. PRO-30571 or PRO-30576 required with Straight Back Post Tubes (PRO-30544)			
				Not available with Straight Back Posts. Not available with Recline.			

Footrests


Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-12"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	13"-17"
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	13"-15"
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	13"-18"
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	14"-20"
10" Adjustable	16"-19"	14.5"-18.5"	N/A	N/A	16"-22"

MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90 °Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-14.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	12.5"-16.5"
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	12.5"-14.5"
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	12.5"-17.5"
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	13.5"-19.5"
10" Adjustable	14.5"-17.5"	14"-18"	N/A	N/A	15.5"-21.5"

Length is measured from top of footrest plate to top of seat rail. Without extensions bottom hanger hole not accessible for hardware. Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.



Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
Footrests - continued							
Hangers							
75° Pop-Up Swing Away Hangers		PRO-30280	n/c	75° Latch Style Swing Away Hangers		PRO-30282	n/c
90° Pop-Up Swing Away Hangers		PRO-30281	n/c	90° Latch Style Swing Away Hangers		PRO-30283	n/c
<i>Not available with 8" or 10" Telescoping Extensions.</i>				<i>Not available with 8" or 10" Telescoping Extensions.</i>			
90° Pop-Up Swing Away Elevating Legrests (ELRs)	(E0990)	PRO-30279	\$ 295	Omit Legrest Hangers & Footrests		PRO-33000	
<i>Not available with all Telescoping extensions.</i>							

Footplates & Footboards							
Aluminum Angle Adj. Flip-up Footplates	(K0040)	PRO-36033	\$ 190	Aluminum Angle Adjustable Foot-board Removable	(K0040)	PRO-36029	\$ 190
<i>Does not include Heel Loops, select below.</i>							
Composite Flip-Up		PRO-17070	n/c	Omit Footrest Plates		PRO-36404	

Footplate Options							
MAPS (Multi Axis Positioning System)		PRO-36040	\$ 325	Heel Loops	(E0951)	PRO-36022	\$ 45
<i>Requires selection of Aluminum Angle Adjustable Footplates. Not compatible with Composite footplates or Aluminum Footboard. Not compatible with Seating Dynamics hangers.</i>				Calf Strap - Webbing	(K0038)	PRO-36024	\$ 45


Hanger Extension Options							
4" Telescoping Extensions		PRO-30285	n/c	6" Telescoping Extensions		PRO-30286	n/c
8" Telescoping Extensions		PRO-30287	n/c	10" Telescoping Extensions		PRO-30288	n/c
<i>Not available with all 90° hangers.</i>				<i>Not available with all 90° hangers.</i>			
4.5" Heavy Duty Extensions		PRO-36402	\$ 58	6.5" Heavy Duty Extensions		PRO-36403	\$ 68
<i>Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</i>							
High Mount Clamp & Tubes		PRO-36410	\$ 102	Clear Selection			

Specialized Footrests
 ALL legrest measurements are from the seat rail to footplate. Not available with LEAF or SRO seat rail options. Frame and/or front caster choice may limit adjustment range and/or cause interference.

Specialized Footrest Seat Rail to Footplate Measurements

Part Number	Extension	Super Contracture	3D Dynamic Hi Pivot	3D Dynamic STD
6402-01	Short	10.5" - 13.5"		
6402-02	Long	13.5" - 18.5"		
6405-01	0"-2" Extension		6.25" - 10.75"	12.75-15.25"
6409-01	1"-3" Extension		7.25-11.75"	13.75-16.25"
6406-01	2"-4" Extension		8.25-12.75"	14.75-17.25"
6410-01	3"-5" Extension		9.25-13.75"	15.75-18.25"
6407-01	4"-6" Extension		10.25-14.75"	16.75-19.25"
6411-01	5"-7" Extension		11.25-15.75"	17.75-20.25"
6508-01	6"-8" Extension		12.25-16.75"	18.75-21.25"

Length measured from top of footrest plate to top of seat rail. 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.



Contracture Footrest Hangers							
<i>Super Contracture for frames 14" wide and greater only.</i>							
Contracture Short - 11" to 15"		PRO-30291-2	\$3570	Contracture Long - 15" to 19"		PRO-30292-2	\$3570
<i>Includes short hanger extension.</i>				<i>Includes long hanger extension.</i>			

3D Dynamic Footrest Hangers							
3D Dynamic Footrest - Hi Pivot		PRO-30297-2	\$6215	3D Dynamic Footrest - STD		PRO-30298-2	\$5775

Dynamic Footrest Extension Options							
Dynamic Adj Plantar/Dorsi Ext, 0-2"		PRO-6405-01	n/c	Dynamic Adj Plantar/Dorsi Ext, 1-3"		PRO-6409-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 2-4"		PRO-6406-01	n/c	Dynamic Adj Plantar/Dorsi Ext, 3-5"		PRO-6410-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 4-6"		PRO-6407-01	n/c	Dynamic Adj Plantar/Dorsi Ext, 5-7"		PRO-6411-01	n/c
Dynamic Adj Plantar/Dorsi Ext, 6-8"		PRO-6408-01	n/c				

Specialized Footplates							
4" x 6" HD Aluminum Footplates		PRO-6503	n/c	5" x 8" HD Aluminum Footplates		PRO-6504	n/c
6" x 10" HD Aluminum Footplates		PRO-6505	n/c	<i>For frames 16" wide and greater only.</i>			
<i>For frames 18" wide and greater only. Not available with Contracture Hangers</i>							

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
-------------	-------	----------	--------	-------------	-------	----------	--------

Armrests

Back Post Mounted Armrests

Omit Armrests

Foam wrap pad STD on all back post mounted armrests..

10.5" Flip-Back		PRO-38204	n/c	10.5" Flip-Back Locking	(E0973)	PRO-38201	\$ 265
13.5" Flip-Back		PRO-38202	n/c	13.5" Flip-Back Locking	(E0973)	PRO-38200	\$ 265
10.5" to 13.5" Flip-Back Locking, Telescoping	(E0973)	PRO-38210	\$ 285	17" to 20" Flip-Back Locking, Telescoping	(E0973)	PRO-38215	\$ 285

Full length pad on 13.5" or longer armrests only.

Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

All heights from top of seat frame to top of armrest pad.



Frame Mounted Armrests

Single Post Height Adjustable, Removable - Low	(E0973)	PRO-38217	\$ 285	Single Post Height Adjustable, Removable - Standard	(E0973)	PRO-38203	\$ 285
Dual Post Fixed Height Removable		PRO-38205	n/c	Dual Post Flip-Back Height Adjustable, Removable	(E0973)	PRO-38206	\$ 295

Arm Pad Options

Desk Length Pads		PRO-38213	n/c	Full Length Pads 13.5" or longer armrests only.		PRO-38214	n/c
------------------	--	-----------	-----	--	--	-----------	-----

Casters

4" x 1" Polyurethane		PRO-40000	n/c	4" x 1.4" Poly, Aluminum Mag		PRO-43000	\$ 125
5" x 1" Polyurethane		PRO-34555	n/c	5" Polyurethane Light-Up		PRO-35050	\$ 101
5" x 1.5" Polyurethane		PRO-34222	\$ 87	5" x 1.4" Poly, Aluminum Mag		PRO-42000	\$ 125
6" x 1" Polyurethane		PRO-34666	n/c	6" x 1.4" Poly, Aluminum Mag		PRO-44000	\$ 125
6" x 1.5" Polyurethane		PRO-34999	\$ 87	6" x 2" Polyurethane, Treaded		PRO-34444	\$ 97
7" x 1" Polyurethane		PRO-34777	n/c	7" x 2" Foam Filled		PRO-34200	\$ 200

Rear Wheels

Mags

12" Mag Wheel		PRO-31525	n/c
16" Mag Wheel		PRO-31530	\$ 41
18" Mag Wheel		PRO-31533	\$ 41
20" Mag Wheel - Aluminum Handrim		PRO-31527	n/c
22" Mag Wheel - Aluminum Handrim		PRO-31528	n/c
24" Mag Wheel - Aluminum Handrim		PRO-31529	n/c

Spoke

16" Spoke Wheel - Black w/ Black Tire	PRO-31535	n/c
18" Spoke Wheel - Black w/ Black Tire	PRO-31538	n/c
20" Spoke Wheel - Aluminum Handrim	PRO-31540	\$ 204
22" Spoke Wheel - Aluminum Handrim	PRO-31541	\$ 204
24" Spoke Wheel - Aluminum Handrim	PRO-31542	\$ 204

20", 22" & 24" Spoke Wheels not available with drum brakes

Tire Options

Pneumatic Tire		PRO-31040	\$ 85
Pneumatic Tire with Airless Insert	(E2213)	PRO-31000	\$ 85
Low Profile Polyurethane Tire	(E2218)	PRO-31045	n/c

12" Mag wheel only, not available with drum brakes

Handrim Options

Plastic Coated Handrims	PRO-31005	\$ 119
Omit Handrims	PRO-31006	

20", 22" & 24" Mag wheels only, not available with drum brakes

Wheel/Caster Seat Rail to Floor Heights							Wheel/Caster Seat Rail to Floor Heights						
Rear Wheel	Caster Wheel						Rear Wheel	Caster Wheel					
	4"	5"	6"	6 x 2"	7"	7 x 2"		4"	5"	6"	6 x 2"	7"	7 x 2"
12"	15"	15"	18.5"	18.5"	18.5"	18.5"	20"	17"	17"	17"	17"	17"	17"
	15.5"	15.5"	19"	19"	19"	19"		17.5"	17.5"	17.5"	17.5"	17.5"	17.5"
	16"	16"						18"	18"	18"	18"	18"	18"
	18.5"	18.5"						18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"						19"	19"	19"	19"	19"	19"
16"	15"	15"	16"	16"	17"	17"	22"	18"	18"	18"	18"	18"	18"
	15.5"	15.5"	16.5"	16.5"	17.5"	17.5"		18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	16"	16"	17"	17"	20"	20"		19"	19"	19"	19"	19"	19"
	16.5"	16.5"	17.5"	17.5"	20.5"	20.5"		19.5"	19.5"	19.5"	19.5"	20"	20"
	17"	17"	20"	20"				20"	20"	20"	20"		
18"	15.5"	15.5"	16"	16"	17"	17"	24"	18.5"	19"	19"	19"	19"	19"
	16"	16"	16.5"	16.5"	17.5"	17.5"		19"	19.5"	19.5"	19.5"	20"	20"
	16.5"	16.5"	17"	17"	18"	18"							
	17"	17"	17.5"	17.5"	18.5"	18.5"							
	17.5"	17.5"	18"	18"	18.5"	18.5"							

Clear selection

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".

Wheel Locks

Standard			
Push to Lock		PRO-31501	n/c
Pull to Lock		PRO-31506	n/c
Wheel Lock Options			
Extension Handles	(E0961)	PRO-31801	\$ 75

Drum Lock			
Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires. Not Available with 20", 22" & 24" Spoke Wheels.			
Foot Activated Wheel Lock		PRO-31508	\$ 0
Not available with ALL vent and battery trays.			
Hand Activated Brake and Wheel Lock System	(E2228)	PRO-30275	\$ 1000

Anti-Tippers

Curb Assist Anti-Tip Assembly Flip-Up, Adjustable	(E0971)	PRO-31600	\$ 125	Spring Loaded Swing-Away Anti-Tip Assy - Foot Activated	(E0971)	PRO-31605	\$ 180
---	---------	-----------	--------	---	---------	-----------	--------

Pair includes one anti-tip tube with curb assist feature.

Seat Pans

ALL seat pans include 2" x 2" notches for growth adjustment.

Non-Adjustable			
Fixed Solid Seat Pan - Composite	(E2231)	PRO-39010	\$ 170
ALL Bolteron seat pans are made to fit frame dimensions unless specified below.			
Seat Width _____" X	Seat Depth _____"		
Fixed Solid Seat Pan - Aluminum	(E2231)	PRO-39030	\$ 200
Fixed Seat Pan raises the effective seat to floor height by 1/2".			
Available in the following sizes ONLY:			
14" wide x 16" deep	18" wide x 18" deep		
16" wide x 18" deep			

Adjustable			
Fixed, Width Adjustable Solid Seat Pan - Aluminum	(E2231)	PRO-39040	\$ 290
Fixed Seat Pan raises the effective seat to floor height by 1/2".			
Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with PRO frame.			
Available in the following sizes ONLY:			
14"-16" wide x 16" deep	16"-18" wide x 18" deep		

Mounting Hardware Options

Use to convert fixed seat pan to removable.

Removable Pan Hardware Kit	(K0108)	PRO-EZLITE	\$ 207	1" Drop Pan Hardware Kit	(K0108)	PRO-DROP1	\$ 207			
Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".				2" Drop Pan Hardware Kit				(K0108)	PRO-DROP2	\$ 207
Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.				Includes 4 each 2" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps.				Not available with ANY Aluminum Solid Seat Pan.		

Description	hccps	Part No.	Retail	Description	hccps	Part No.	Retail
Accessories							
Tilt Range Limit Kit <small>Pins install into outside of tilt channel to limit tilt travel.</small>	(K0108)	PRO-10920	\$ 35	Freedom Backpack <small>Black Nylon only with Freedom Designs logo. Monogramming is not available on Freedom Backpack.</small>	(K0108)	PRO-10900	\$ 60
Footrest Hanger Pads - Pair <small>Wrap-around pads made from black neoprene.</small>	(K0108)	SFDI-9HP	\$ 97	Freedom Canopy - Flip-Back <small>Canopy assembly includes removable mounting hardware. Available in black only with Freedom Designs logo. Includes rear flap with zippered storage pocket and clear-view window.</small>	(K0108)	PRO-10905	n/c
Flip-Down Rigidizing Bar <small>Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.</small>	(K0108)	PRO-23210	\$ 175				

Notes:

Specifications:	
Frame Style	Tilt-In-Space
Base Frame Widths	16", 18", 20"
Overall Width	Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel)
Seat Widths	14" - 20"
Seat Depths	14" - 26"
Back Heights	20", 23", 26"
Adjustable Angle Back	90°, 100°, 110°, 120°
Tilt Range	0° to 55°
Recline Back	162°
Weight Limit	250 LBS.
Transport Weight	
Warranty	Lifetime on frame

Special Configurations Through Customer Service:	
Anterior Tilt	As much as -5° depending on configuration. Reduces posterior tilt angle by same amount.
Asymmetrical Seat Depth	Seat rails may be configured up to 4" offset

HCPCS Codes:
are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.