



PRO-CG™ Adult Manual Wheelchair



Yes, you can!™

Clear Form

Email Form to
Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: _____ Buyer: _____ Quote Order PO# _____
 Company: _____ Location: _____
 Marked For: _____ Phone: _____
 Email: _____ Fax: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Notice: A seating system must be purchased separately. No seating system is provided with the PRO-CG frame. No seat or back upholstery is available as an option.

Transport Optional Base

PRO-CG Adult Rigid (E1161) PRO-30310 \$2930

Transport Optional

Note: Transportation Option selected separately below.
Note: Includes Freedom Future Fit Program

Transport Package PRO-99998 \$ 250

Note: Must be installed by Freedom Designs manufacturing.
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.
Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

Transport Included Base

PRO-CG Adult Rigid (E1161) PRO-30311 \$3180

Transport Included

Note: Includes Freedom Future Fit Program
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.

Seat Rail Options

Note: Not available with ALL Contracture & 3D Dynamic footrest systems.

LEAF - Lower Extremity Accommodation Frame

LEAF Seat Frame- Pair (K0108) PRO-30324 \$ 714

LEAF Seat Frame- Left (K0108) PRO-30324-L \$ 357

LEAF Seat Frame- Right (K0108) PRO-30324-R \$ 357

Note: Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

SRO - Seat Rail Only

SRO Seat Frame- Pair (K0108) PRO-30325 \$ 306

SRO Seat Frame- Left (K0108) PRO-30325-L \$ 153

SRO Seat Frame- Right (K0108) PRO-30325-R \$ 153

Note: SRO is manufactured without standard footrest hanger attachments.

Tilt Activation Options

Dual Lever Hand Activated Tilt PRO-30272 n/c

Foot Activated Tilt PRO-30273 n/c

Note: Not with vent & battery trays.

Vent Tray Options

Note: Vent Ready Base required for Vent Tray.

Note: Vent Tray not with Foot Activated Tilt.

Vent Ready Base (K0108) PRO-30331 \$ 459

Note: Includes Push Button Adjustable Stroller Handle. For LP10 style vent and battery tray on frames 16" wide and greater.

Vent Tray - Fixed (E1029) 11-7000 \$ 439

Note: For LP10 style vent tray on 16" base width and greater.

Note: Vent Tray is full width of the base frame and 11" deep.

Quick Release Vent Tray Hardware (E1028) 563 \$ 272

Back Post Mounted Vent (E1029) 11-1100 \$ 550

Bracket

Note: For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

Frame Color Options

Sugar Cookie n/c Blue Blast n/c

Neon Yellow n/c Electric Blue n/c

Yellow n/c Galaxy Blue n/c

Bright Pink n/c Passion Purple n/c

Pink Parfait n/c Octane Orange n/c

Sugar Plum n/c Red Sunset n/c

Candy Violet n/c Really Red n/c

Candy Raspberry n/c Fiery Red n/c

Extreme Green n/c Twilight n/c

Sparkle Green n/c Stormy Gray n/c

Mystic Teal n/c Black n/c

Pacific Blue n/c

PRO-CG Rider Weight Capacity is 250 lbs.
All Prices subject to change without notice.

PRO-CG Adult
Note: All Specifications and dimensions are approximate

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Frame Dimensions


Specify Frame Width & Frame Depth Settings

Note: *PRO-CG* has 4" of depth adjustment for growth. Depth is measured from front of back post to front edge of the seat frame. For back post location add seating/hardware thickness (typically 2") to client's actual seat depth.

Example: seat depth of 12" + back/hardware thickness of 2" = back post location set at 14"

Note: The *PRO-CG* has 6" of center of gravity (CG) adjustment on all adult frames. The default CG setting is 1 position forward of either the client's actual seat depth or the manufacturer's seat depth setting.

Specify Client's Actual Seat Depth: _____"



<i>PRO-CG</i> Base Frame Width	Manufacturing CG location Setting														
	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	25"	26"
16" Base Frame Width			-3	-2	-1	0	1	2	3	3	3				
18" Base Frame Width				-3	-2	-1	0	1	2	3	3	3			
20" Base Frame Width						-3	-2	-1	0	1	2	3	3	3	

Note: Seat rail extensions required for manufacturing seat depth settings with CG setting in gray area.

16" Base Frame Width – Seat PRO-30240
frame adjusts 14" to 16" wide in 1" increments.

Specify Seat Width Setting

14"	PRO-14W	n/c
15"	PRO-15W	n/c
16"	PRO-16W	n/c

14" To 22" Seat Frame Length PRO-30512

Specify Seat Depth Setting

14"	PRO-30506	n/c
15"	PRO-30521	n/c
16"	PRO-30522	n/c
17"	PRO-30523	n/c
18"	PRO-30524	n/c
* 19" - Must select 2" extension	PRO-30527	n/c
* 20" - Must select 2" extension	(E2203) PRO-30528	\$ 500
** 21" - Must select 4" extension	(E2203) PRO-30529	\$ 500
** 22" - Must select 4" extension	(E2204) PRO-30535	\$ 900
2" Extension - Pair	(K0108) PRO-30525	\$ 26
4" Extension - Pair	(K0108) PRO-30531	\$ 26

* **Note:** Must choose 2" extension for 19" and 20" frame lengths.

** **Note:** Must choose 4" extension for 21" and 22" frame lengths.

Specify Back Post Location (if other than specified above) _____"

Frame Dimensions (continued)

Specify Frame Width & Frame Depth Settings

18" Base Frame Width – Seat frame PRO-30260
adjusts 16" to 18" wide in 1" increments.

Specify Seat Width Setting

16"	PRO-16W	n/c
17"	PRO-17W	n/c
18"	PRO-18W	n/c

16" To 24" Seat Frame Length PRO-30520

Specify Seat Depth Setting

16"	PRO-30522	n/c
17"	PRO-30523	n/c
18"	PRO-30524	n/c
19"	PRO-30527	n/c
20"	(E2203) PRO-30528	\$ 500
* 21" - Must select 2" extension	(E2203) PRO-30529	\$ 500
* 22" - Must select 2" extension	(E2204) PRO-30535	\$ 900
** 23" - Must select 4" extension	(E2204) PRO-30536	\$ 900
** 24" - Must select 4" extension	(E2204) PRO-30537	\$ 900
2" Extension - Pair	(K0108) PRO-30525	\$ 26
4" Extension - Pair	(K0108) PRO-30531	\$ 26

* **Note:** Must choose 2" extension for 21" and 22" frame lengths.

** **Note:** Must choose 4" extension for 23" and 24" frame lengths.

Specify Back Post Location (if other than specified above) _____"

20" Base Frame Width – Seat frame PRO-30261
adjusts 18" to 20" wide in 1" increments.

Specify Seat Width Setting

18"	PRO-18W	n/c
19"	PRO-19W	n/c
20"	(E2201) PRO-20W	\$ 400

18" To 26" Seat Frame Length PRO-30540

Specify Seat Depth Setting

18"	PRO-30524	n/c
19"	PRO-30527	n/c
20"	(E2203) PRO-30528	\$ 500
21"	(E2203) PRO-30529	\$ 500
22"	(E2204) PRO-30535	\$ 900
* 23" - Must select 2" extension	(E2204) PRO-30536	\$ 900
* 24" - Must select 2" extension	(E2204) PRO-30537	\$ 900
** 25" - Must select 4" extension	(E2204) PRO-30538	\$ 900
** 26" - Must select 4" extension	(E2204) PRO-30539	\$ 900
2" Extension - Pair	(K0108) PRO-30525	\$ 26
4" Extension - Pair	(K0108) PRO-30531	\$ 26

* **Note:** Must choose 2" extension for 23" and 24" frame lengths.

** **Note:** Must choose 4" extension for 25" and 26" frame lengths.

Specify Back Post Location (if other than specified above) _____"

Note: Extensions may limit tilt range and may affect stability when in tilted positions. Adjusting seat rail location will compensate for tilt range limitations and possible instability.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Back Post Options

Specify Push Handle, Back Post and Extensions

Standard 90° Bend Push Handle Back Post Tubes PRO-30518 n/c

Ergo Bend Push Handle Back Post Tubes PRO-30519 n/c

Note: All push/stroller handle extensions not recommended with Ergo Bend.

Seating Dynamics Rocker Back (K0108) PRO-30542 **\$1985**

Note: Angle adjustable from 72° to 108° in 6° increments.

Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.

Note: Locks into a non-dynamic position when necessary.

Note: Does not allow the back post tubes to fold forward.

Super Reclining Reduced Shear Back (K0108) PRO-30562 **\$1138**

Super Reclining Reduced Shear Back with Push Button, Angle-Adjustable Stroller Handle Extension (K0108) PRO-30564 **\$1335**

Note: Pivot point is 2.5" above the seat rail.

Note: Recline range is 65° from 97° to 162°.

Note: Includes height adjustable back post tubes for 21" - 27" heights.

Note: Available with standard 90° push handles only.

Note: Only compatible with Push Button, Angle-Adjustable Folding Stroller Handle Extension.

Note: Does not allow the back post tubes to fold forward.

Note: Solid Back Insert installation requires 4-points of attachment hardware. Seat-to-Back hinges are not compatible with any Super Reclining Reduced Shear Back.

Back Post Height Setting

Note: Back Posts are measured from top of seat rail to top of handle bend.

Short Push Handle - 20" high PRO-30515 n/c

Medium Push Handle - 23" high PRO-30514 n/c

Tall Push Handle - 26" high PRO-30516 n/c

Stroller Handle Extension Options

Note: All stroller handle ext. add approximately 6.5" to push handle height.

Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.



-30560 Push Handle Extension - Width Adjustable - Mounts direct to back post push handle PRO-30561 n/c

-30571/6 Stroller Handle Extension - Width Adjustable, Removable - Mounts direct to back post push handle PRO-30590 **\$ 95**

-30581/6 Push Button, Angle-Adjustable Stroller Handle Extension - Mounts direct to back post push handle PRO-30571 **\$ 112**

Note: Not available with Super Reclining Back.

-30517 Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty - Mounts direct to back post tubes PRO-30576 **\$ 179**

Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.

Note: Not removable.

Note: Available 10W-20w wheelchair frames.

Note: Not removable.

Stroller Handle Extension Options - Continued

Note: All stroller handle ext. add approximately 6.5" to push handle height.

Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.

Push Button, Angle-Adjustable, Folding Push Handle Extension - Mounts direct to straight back post tubes PRO-30581 **\$ 112**

Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.

Note: Not removable.

Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty - Mounts direct to straight back post tubes PRO-30586 **\$ 179**

Note: Available 10W-20w wheelchair frames.

Note: Not removable.

Individual Stroller Handle Extensions Removable - Pair PRO-30517 **\$ 102**

PRO Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A
10" Adjustable	16"-19"	14.5"-18.5"	N/A	N/A	N/A

PRO MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A
10" Adjustable	14.5"-17.5"	14"-18"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.

Note: Without extensions bottom hanger hole not accessible for hardware.

Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

Legrests, Footrests and Extension Options

Legrests

Omit Legrest Hangers & Footrests PRO-33000

75° Pop-Up Swing Away Hangers PRO-30280 n/c

90° Pop-Up Swing Away Hangers PRO-30281 n/c

Note: Not available with 8" adjustable extensions.

90° Pop-Up Swing Away Elevating Legrests (ELRs) (E0990) PRO-30279 **\$ 255**

Note: Not available with all length adjustable extensions.

Note: Not available with 6.5" heavy duty extensions.

75° Latch Style Swing Away Hangers PRO-30282 n/c

90° Latch Style Swing Away Hangers PRO-30283 n/c

Note: Not available with 8" adjustable extensions.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Legrests, Footrests and Extension Options

Footplates & Footboards

Omit Footrest Plates	PRO-36404	
Composite Flip-Up	PRO-17070	n/c
Kids Aluminum Angle Adjustable, Flip-Up	(K0040) PRO-36020	\$ 160
<i>Note: For frames 14" wide and narrower.</i>		
Heel Loops	(E0951) PRO-36022	\$ 41
Calf Strap - Webbing	(K0038) PRO-36024	\$ 41
Adult Aluminum Angle Adj. Flip-up w/ Heel Loops	(K0040/E0951) PRO-36021	\$ 180
<i>Note: For frames 15" wide and greater.</i>		
Omit Heel Loops	PRO-36023	n/c
Aluminum Angle Adjustable Footboard Removable	(K0040) PRO-36029	\$ 160

MAPS Footplate Mounting Hardware

MAPS (Multi Axis Positioning System) - Pair	PRO-36040	\$ 306
MAPS (Multi Axis Positioning System) - Left	PRO-36040-L	\$ 153
MAPS (Multi Axis Positioning System) - Right	PRO-36040-R	\$ 153
<i>Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.</i>		
<i>Note: Must choose aluminum footplates.</i>		
<i>Note: MAPS not available with composite footplates or aluminum footboard.</i>		
<i>Note: MAPS not available on 10" & 11" wide frames.</i>		

Legrest Extensions

4" Length Adjustable Extensions	PRO-30285	n/c
6" Length Adjustable Extensions	PRO-30286	n/c
8" Length Adjustable Extensions	PRO-30287	n/c
10" Length Adjustable Extensions	PRO-30288	n/c
<i>Note: 8" & 10" Not available with all 90° hangers.</i>		
4.5" Heavy Duty ELR Extensions	PRO-36401	\$ 56
6.5" Heavy Duty ELR Extensions	PRO-36405	\$ 56
4.5" Heavy Duty Extensions	PRO-36402	\$ 56
6.5" Heavy Duty Extensions	PRO-36403	\$ 66
<i>Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</i>		
High Mount Brackets	PRO-36407	\$ 102
<i>Note: Includes flip-up mounting hardware for future growth adjustment.</i>		
<i>Note: 75° hangers use a tube & clamp style hardware.</i>		

Specify Knee to Heel Measurement _____" (see matrix)

Notes:

Specialized Footrest Seat Rail to Footplate Measurements

Part Number	Extension	Super Contracture	3D Dynamic Hi Pivot	3D Dynamic STD
6402-01	Short	10.5" - 13.5"		
6402-02	Long	13.5" - 18.5"		
6405-01	0"-2" Extension		6.25" - 10.75"	12.75-15.25"
6409-01	1"-3" Extension		7.25-11.75"	13.75-16.25"
6406-01	2"-4" Extension		8.25-12.75"	14.75-17.25"
6410-01	3"-5" Extension		9.25-13.75"	15.75-18.25"
6407-01	4"-6" Extension		10.25-14.75"	16.75-19.25"
6411-01	5"-7" Extension		11.25-15.75"	17.75-20.25"
6508-01	6"-8" Extension		12.25-16.75"	18.75-21.25"

Note: Length from top of footrest plate to top of seat rail.
Note: 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.

Specialized Footrest Systems

Note: ALL legrest measurements are from the seat rail to footplate.
Note: Not available with LEAF or SRO seat rail options.

Super Contracture Footrest Hanger (K0108)

Note: Super Contracture for frames 14" wide and greater.

Contracture Short - 11" to 15" PRO-30291-2 \$3570

Note: See Matrix for seat to rail footplate measurements.

Includes:

Adj Contracture Upper	PRO-6302-03	Incl.
Footrest Adapter HD	PRO-6220-01	Incl.
Adj Plantar/Dorsi Ext. Short	PRO-6401-02	Incl.

Footplate Options

4" x 6" HD Alu Footplates	PRO-6503	n/c
5" x 7" HD Alu Footplates	PRO-6504	n/c

Note: for frames 16" wide and greater only.

Contracture Long - 15" to 19" PRO-30292-2 \$3570

Note: See Matrix for seat to rail footplate measurements.

Includes:

Adj Contracture Upper	PRO-6302-03	Incl.
Footrest Adapter HD	PRO-6220-01	Incl.
Adj Plantar/Dorsi Ext. Long	PRO-6402-02	Incl.

Footplate Options

4" x 6" HD Alu Footplates	PRO-6503	n/c
5" x 7" HD Alu Footplates	PRO-6504	n/c

Note: Super Contracture is both swing-away and removable.

Note: Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.

Note: Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.

Note: Super Contracture position may limit ankle joint/ footplate adjustment range.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Specialized Footrest Systems - Continued

3D Dynamic Footrest Hanger (K0108)

Note: 3D Dynamic for frames 14" wide and greater.

3D Dynamic Footrest - Hi Pivot PRO-30297-2 **\$6215**

Includes:

HiPivot Dyn Adj Knee & Dyn Tel Up PRO-6304-02 **Incl.**

Footrest Adapter HD PRO-6220-01 **Incl.**

Extension Options

Note: See Matrix for seat to rail footplate measurements.

Dynamic Adj Plantar/Dorsi Ext, 0-2" PRO-6405-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 1-3" PRO-6409-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 2-4" PRO-6406-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 3-5" PRO-6410-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 4-6" PRO-6407-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 5-7" PRO-6411-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 6-8" PRO-6408-01 **n/c**

Footplate Options

4" x 6" HD Alu Footplates PRO-6503 **n/c**

5" x 7" HD Alu Footplates PRO-6504 **n/c**

Note: for frames 16" wide and greater only.

6" x 10" HD Alu Footplates PRO-6505 **n/c**

Note: for frames 18" wide and greater only.

3D Dynamic Footrest - STD PRO-30298-2 **\$5775**

Includes:

Dynamic Adj Knee & Dyn Tel Up PRO-6304-01 **Incl.**

Footrest Adapter HD PRO-6220-01 **Incl.**

Extension Options

Note: See Matrix for seat to rail footplate measurements.

Dynamic Adj Plantar/Dorsi Ext, 0-2" PRO-6405-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 1-3" PRO-6409-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 2-4" PRO-6406-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 3-5" PRO-6410-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 4-6" PRO-6407-01 **n/c**

Dynamic Adj Plantar/Dorsi Ext, 5-7" PRO-6411-01 **n/c**

Footplate Options

4" x 6" HD Alu Footplates PRO-6503 **n/c**

5" x 7" HD Alu Footplates PRO-6504 **n/c**

Note: for frames 16" wide and greater only.

6" x 10" HD Alu Footplates PRO-6505 **n/c**

Note: for frames 18" wide and greater only.

Note: 3D Dynamic is both swing-away and removable.

Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.

Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.

Note: Angle adjustable footplate included.

Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Single Post Fxd Hgt Short		8" - 10.5"	8" - 10.5"
Single Post Fxd Hgt Short		9" - 11.5"	9" - 11.5"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

Armrests and Armrest Pad Options

Omit Armrests PRO-38209

Back Post Mounted Armrests

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

10.5" Long Flip-Back PRO-38204 **n/c**

Note: Not available with full length arm pads

13.5" Long Flip-Back PRO-38202 **n/c**

10.5" Long Flip-Back Locking (E0973) PRO-38201 **\$ 245**

Note: Not available with full length arm pads

13.5" Long Flip-Back Locking (E0973) PRO-38200 **\$ 245**

10.5" to 13.5" Long Flip-Back, Telescoping (E0973) PRO-38208 **\$ 240**

Note: Full length pad on 13.5" or longer armrests only.

10.5" to 13.5" Long Flip-Back Locking, Telescoping (E0973) PRO-38210 **\$ 265**

Note: Full length pad on 13.5" or longer armrests only.

17" to 20" Long Flip-Back Locking, Telescoping (E0973) PRO-38215 **\$ 265**

Foam Wraps - Standard PRO-30216 **n/c**

Desk Length Pads PRO-38213 **n/c**

Full Length Pads PRO-38214 **n/c**

Note: Full length pad on 13.5" or longer armrests only.

Frame Mounted Armrests

Single Post Fixed Height, Removable - Short (K0015) PRO-38222 **\$ 225**

Single Post Fixed Height, Removable - Tall (K0015) PRO-38224 **\$ 225**

Note: Single Post Fixed Height arms include a removable side guard.

Note: Adjustable in height with tools only. See matrix for height ranges.

Single Post Height Adjustable, Removable - Short (E0973) PRO-38217 **\$ 265**

Single Post Height Adjustable, Removable - Tall (E0973) PRO-38203 **\$ 265**

Dual Post Fixed Height Removable PRO-38205 **n/c**

Dual Post Flip-Back Height Adjustable, Removable (E0973) PRO-38206 **\$ 275**

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

Desk Length Pads PRO-38213 **n/c**

Full Length Pads PRO-38214 **n/c**

Specify Seat Rail to Top of Armrest Measurement _____ "

PRO-CG Rider Weight Capacity is 250 lbs.
All Prices subject to change without notice.

PRO-CG Adult
Note: All Specifications and dimensions are approximate

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Caster Wheel				
	4"	5"	6"	7"	
12"	15"	15"	18.5"	18.5"	
	15.5"	15.5"	19"	19"	
	16"	16"			
	18.5"	18.5"			
	19"	19"			
16"	15"	15"	16"	17"	
	15.5"	15.5"	16.5"	17.5"	
	16"	16"	17"	20"	
	16.5"	16.5"	17.5"	20.5"	
	17"	17"	20"		
18"	15.5"	15.5"	16"	17"	
	16"	16"	16.5"	17.5"	
	16.5"	16.5"	17"	18"	
	17"	17"	17.5"	18.5"	
	17.5"	17.5"	18"	18.5"	
20"	17"	17"	17"	17"	
	17.5"	17.5"	17.5"	17.5"	
	18"	18"	18"	18"	
	18.5"	18.5"	18.5"	18.5"	
	19"	19"	19"	19"	
22"	18"	18"	18"	18"	
	18.5"	18.5"	18.5"	18.5"	
	19"	19"	19"	19"	
		19.5"	19.5"	20"	
			20"		
24"	18.5"	19"	19"	19"	
	19"	19.5"	19.5"	20"	
			20"	20.5"	

Note: Use chart to determine seat rail to floor height. Seat rail to floor height determined by frame type, rear wheels, and front caster size.

Note: Caster sizes vary slightly in height by type.

Note: Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4".

Note: Seat pan choice increases seat rail to floor height 3/8" to 3/4".

Front Wheels - Caster Options - Continued

6" x 1" Poly, Plastic Mag	PRO-34666	n/c
6" x 1.5" Poly, Plastic Mag	PRO-34999	\$ 87
6" x 1.4" Poly, Aluminum Mag	PRO-44000	\$ 175
6" x 2" Poly with Tread	PRO-34444	\$ 92
6" x 1" Pneumatic Plastic Mag (E2214/E2215)	PRO-34667	\$ 77
7" x 1" Poly, Plastic Mag	PRO-34777	n/c
7" x 2" Polyurethane	PRO-34200	\$ 118

Rear Wheels - Tires - Accessories

Note: All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.

Omit Rear Wheels and Axles	PRO-31524	
Mag Wheels		
12" Low Profile Poly Mag (E2218)	PRO-31526	\$ 36
Note: Not available with any foot/hand activated brake/wheel lock system.		
12" Semi Pneumatic Plastic Mag	PRO-31523	\$ 168
Note: Not available with airless inserts.		
12" Pneumatic Plastic Mag	PRO-31525	n/c
12" Airless Inserts (E2213)	PRO-31019	\$ 77
16" Pneumatic Plastic Mag	PRO-31530	\$ 41
16" Airless Inserts (E2213)	PRO-31020	\$ 77
16" Aluminum Handrims	PRO-31014	\$ 169
18" Pneumatic Plastic Mag	PRO-31533	\$ 41
18" Airless Inserts (E2213)	PRO-31021	\$ 77
18" Aluminum Handrims	PRO-31015	\$ 169
20" Pneumatic Plastic Mag with Aluminum Handrim	PRO-31527	n/c
20" Airless Inserts (E2213)	PRO-31022	\$ 77
20" Plastic Coated Handrims	PRO-31016	\$ 119
Note: Plastic Coated Handrims NOT available with Drum Braking Systems.		
22" Pneumatic Plastic Mag with Aluminum Handrim	PRO-31528	n/c
22" Airless Inserts (E2213)	PRO-31023	\$ 77
22" Plastic Coated Handrims	PRO-31017	\$ 119
Note: Plastic Coated Handrims NOT available with Drum Braking Systems.		
24" Pneumatic Plastic Mag with Aluminum Handrim	PRO-31529	n/c
24" Airless Inserts (E2213)	PRO-31024	\$ 77
24" Plastic Coated Handrims	PRO-31018	\$ 119
Note: Plastic Coated Handrims NOT available with Drum Braking Systems.		
Omit Handrims from 20", 22" or 24" Wheels	PRO-31006	

Front Wheels - Caster Options

Omit Caster Wheels	PRO-34333	
3" Polyurethane	PRO-34888	n/c
3" Polyurethane Light-Up	PRO-34111	\$ 101
4" Polyurethane	PRO-40000	n/c
4" Polyurethane Light-Up	PRO-41000	\$ 101
4" x 1.4" Poly, Aluminum Mag	PRO-43000	\$ 175
5" x 1" Poly, Plastic Mag	PRO-34555	n/c
5" Polyurethane Light-Up	PRO-35050	n/c
5" x 1.5" Poly, Plastic Mag	PRO-34222	\$ 87
5" x 1.4" Poly, Aluminum Mag	PRO-42000	\$ 175

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Rear Wheels - Tires - Accessories - Cont.

Spoke Wheels

Note: All spoke wheels not available plastic coated handrims.

16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel PRO-31534 n/c

Note: Rim, hub & spokes are anodized silver.

16" Airless Inserts (E2213) PRO-31020 \$ 77

16" Silver Aluminum Handrim PRO-31030 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

16" Black Pneumatic Tire & Black Spoke Aluminum Wheel PRO-31535 n/c

Note: Rim, hub & spokes are anodized black.

16" Airless Inserts (E2213) PRO-31020 \$ 77

16" Black Aluminum Handrim PRO-31031 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel PRO-31537 n/c

Note: Rim, hub & spokes are anodized silver.

18" Airless Inserts (E2213) PRO-31021 \$ 77

18" Silver Aluminum Handrim PRO-31034 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

18" Black Pneumatic Tire & Black Spoke Aluminum Wheel PRO-31538 n/c

Note: Rim, hub & spokes are anodized black.

18" Airless Inserts (E2213) PRO-31021 \$ 77

18" Black Aluminum Handrim PRO-31035 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

20" Pneumatic Spoke Wheel with Aluminum Handrim PRO-31540 \$ 204

20" Airless Inserts (E2213) PRO-31022 \$ 77

22" Pneumatic Spoke Wheel with Aluminum Handrim PRO-31541 \$ 204

22" Airless Inserts (E2213) PRO-31023 \$ 77

24" Pneumatic Spoke Wheel with Aluminum Handrim PRO-31542 \$ 204

24" Airless Inserts (E2213) PRO-31024 \$ 77

Omit Handrims from 20", 22" or 24" Wheels PRO-31006

Wheel Lock Options

Push to Lock - Standard PRO-31501 n/c

Pull to Lock PRO-31506 n/c

Extension Handles (E0961) PRO-31801 \$ 61

Foot Activated Wheel Lock PRO-31508 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Not available with ALL vent and battery trays.

Hand Activated Brake and Wheel Lock System (E2228) PRO-30275 \$ 900

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Used to control direction, speed, brake and lock wheels from strolling position.

Rear Anti-Tipper Options

Curb Assist Anti-Tip Assembly Flip-Up, Adjustable (E0971) PRO-31600 \$ 97

Note: Pair includes one anti-tip tube with curb assist feature.

Spring Loaded Swing-Away Anti-Tip Assy - Foot Activated (E0971) PRO-31605 \$ 180

Seat Pan Options

Note: ALL seat pans include 2" x 2" notches for growth adjustment.

Fixed Solid Seat Pan - Bolteron (E2231) PRO-39010 \$ 170

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Note: ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width _____" X Seat Depth _____"

Fixed Solid Seat Pan - Aluminum (E2231) PRO-39030 \$ 200

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

14" wide x 16" deep **18" wide x 18" deep**

16" wide x 18" deep

Fixed, Width Adjustable Solid Seat Pan - Aluminum (E2231) PRO-39040 \$ 290

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with PRO frame.

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

14"-16" wide x 16" deep

16"-18" wide x 18" deep

Mounting Hardware Options

Note: Use to convert fixed seat pan to removable.

Removable Pan Hardware Kit (K0108) PRO-EZLITE \$ 198

Note: Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Note: Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

1" Drop Pan Hardware Kit (K0108) PRO-DROP1 \$ 198

Note: Includes 4 each 1" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps.

Note: Not available with ANY Aluminum Solid Seat Pan.

2" Drop Pan Hardware Kit (K0108) PRO-DROP2 \$ 198

Note: Includes 4 each 2" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps.

Note: Not available with ANY Aluminum Solid Seat Pan.

Date: _____ PO# _____ Marked For: _____

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

Accessories

Tilt Range Limit Kit	(K0108)	PRO-10920	\$ 35
Note: Pins install into outside of tilt channel to limit tilt travel.			
Medical Necessities Pouch	(K0108)	PRO-23207	\$ 68
Note: Black nylon.			
Footrest Hanger Pads - Pair	(K0108)	SFDI-9HP	\$ 92
Note: Wrap-around pads made from black neoprene.			
Frame Pads - Pair	(K0108)	10911-2	\$ 82
Note: Wrap-around pads made from black neoprene.			
Note: Typically used to pad exposed back post tubes.			
Flip-Down Rigidizing Bar	(K0108)	PRO-23210	\$ 194
Note: Mounts to back post tubes to add rigidity to frame.			
Note: Does NOT adjust in width to grow with frame.			
Freedom Backpack	(K0108)	PRO-10900	\$ 56
Note: Black Nylon only with Freedom Designs logo.			
Note: Monogramming is not available on Freedom Backpack.			
Freedom Canopy - Flip-Back	(K0108)	PRO-10905	\$ 128
Note: Canopy assembly includes removable mounting hardware.			
Note: Canopy is available in black only with Freedom Designs logo.			
Note: Canopy includes rear flap with zippered storage pocket.			
Note: Canopy includes clear-view window.			

HCPCS CODES are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Notes:

Anterior Support Surfaces – (Trays)

3/8" Lexan Tray - Channel Bracket Hardware	(E0950)	08-2000	\$ 271
3/8" Lexan Tray - Adjustable Knob Hardware & Ribs	(E0950)	08-6000	\$ 219
1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs	(E0950)	08-13000	\$ 183
1/4" Lexan Tray - Hook & Loop Straps	(E0950)	08-16000	\$ 106