



Yes, you can.®



PRO-CG™ Pediatric Manual Wheelchair

Clear Form

Email Form to
Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: _____ Buyer: _____ Quote Order Install Seating? PinDot Matrix
 Company: _____ PO#: _____
 Marked For: _____ Location: _____
 Email: _____ Phone: _____

| Description | hccps | Part No. | Retail | Description | hccps | Part No. | Retail |
|-------------|-------|----------|--------|-------------|-------|----------|--------|
|-------------|-------|----------|--------|-------------|-------|----------|--------|

Base Options

A seating system must be purchased separately. No seating system is provided with the PRO-CG frame. No seat or back upholstery is available as an option.

Transport Optional

PRO-CG Pediatric SE (E1233) PRO-30302 \$2625

Transportation Option selected separately below.

SE frame does not include Freedom Future Fit Program

Transport Brackets – SE PRO-99997 \$ 300

Must be installed by Freedom Designs manufacturing.

Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.

Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

Frame Dimensions

14" Base Width PRO-30220

Specify Seat Width Setting

| | | |
|-----|---------|-----|
| 12" | PRO-12W | n/c |
| 13" | PRO-13W | n/c |
| 14" | PRO-14W | n/c |

Seat frame adjusts 12" to 14" wide in 1" increments.

Specify Seat Depth Setting

| | | |
|-----------------------------|-------------------|--------|
| 12" | PRO-30504 | n/c |
| 13" | PRO-30505 | n/c |
| 14" | PRO-30506 | n/c |
| 15" | PRO-30521 | n/c |
| 16" | PRO-30522 | n/c |
| 17" - 2" extension required | PRO-30523 | n/c |
| 18" - 2" extension required | PRO-30524 | n/c |
| 19" - 4" extension required | PRO-30527 | n/c |
| 20" - 4" extension included | (E2203) PRO-30528 | \$ 550 |

16" Base Width PRO-30240

Specify Seat Width Setting

| | | |
|-----|---------|-----|
| 14" | PRO-14W | n/c |
| 15" | PRO-15W | n/c |
| 16" | PRO-16W | n/c |

Seat frame adjusts 14" to 16" wide in 1" increments.

Specify Seat Depth Setting

| | | |
|-----------------------------|-------------------|--------|
| 14" | PRO-30506 | n/c |
| 15" | PRO-30521 | n/c |
| 16" | PRO-30522 | n/c |
| 17" | PRO-30523 | n/c |
| 18" | PRO-30524 | n/c |
| 19" - 2" extension required | PRO-30527 | n/c |
| 20" - 2" extension included | (E2203) PRO-30528 | \$ 550 |
| 21" - 4" extension included | (E2203) PRO-30529 | \$ 550 |
| 22" - 4" extension included | (E2204) PRO-30535 | \$ 990 |

Notes:

| Description | hccps | Part No. | Retail | Description | hccps | Part No. | Retail |
|--|-------|----------------------------|--------|-----------------------------|---------|----------------------------|--------|
| Frame Dimensions (continued) | | | | | | | |
| 18" Base Width | | PRO-30260 | | | | | |
| | | Specify Seat Width Setting | | | | Specify Seat Depth Setting | |
| 16" | | PRO-16W | n/c | 16" | | PRO-30522 | n/c |
| 17" | | PRO-17W | n/c | 17" | | PRO-30523 | n/c |
| 18" | | PRO-18W | n/c | 18" | | PRO-30524 | n/c |
| Seat frame adjusts 16" to 18" wide in 1" increments. | | | | 19" | | PRO-30527 | n/c |
| | | | | 20" | (E2203) | PRO-30528 | \$ 550 |
| | | | | 21" - 2" extension included | (E2203) | PRO-30529 | \$ 550 |
| | | | | 22" - 2" extension included | (E2204) | PRO-30535 | \$ 990 |
| | | | | 23" - 4" extension included | (E2204) | PRO-30536 | \$ 990 |
| | | | | 24" - 4" extension included | (E2204) | PRO-30537 | \$ 990 |

| Frame Dimensions (continued) | | | | | | | |
|-------------------------------------|---------|-----------|-------|--|--|--|--|
| Seat Rail Extensions | | | | | | | |
| 2" Extension - Pair | (K0108) | PRO-30525 | \$ 41 | | | | |
| 4" Extension - Pair | (K0108) | PRO-30531 | \$ 41 | | | | |
| | | | | Specify Back Post Location (if other than specified above)_____" | | | |

| Seat Rail Options | | | | | | | |
|---|---------|-----------|--------|---|---------|-----------|--------|
| Not available with ALL Contracture & 3D Dynamic footrest systems. | | | | | | | |
| LEAF - Lower Extremity Accommodation Frame | | | | SRO - Seat Rail Only | | | |
| LEAF Seat Frame- Pair | (K0108) | PRO-30324 | \$ 775 | SRO Seat Frame- Pair | (K0108) | PRO-30325 | \$ 325 |
| Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy. | | | | SRO is manufactured without standard footrest hanger attachments. | | | |

| Vent Options | | | | | | | |
|---|---------|-----------|--------|---|---------|---------|--------|
| Vent Tray Base W/ Tray | (E1029) | PRO-30335 | \$ 589 | Back Post Mounted Vent Bracket | (E1029) | 11-1100 | \$ 550 |
| Extended wheelchair base including vent tray. For LP10 style vent, base must be 16" or greater. Tray measures full width of the base frame and 12" deep | | | | For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs. | | | |

| Auxillary Device Mounts | | | | | | | |
|---------------------------------|---------|---------|--------|------------------------------|---------|---------|--------|
| O2 Holder - Aluminum w/ QR Hdwe | (K0108) | 11-3100 | \$ 548 | IV Pole w/ Mounting Hardware | (K0108) | 11-4000 | \$ 270 |

| Tilt Activations | | | | | | | |
|--------------------------------|--|-----------|-----|-------------------------------|--|-----------|-----|
| Dual Lever Hand Activated Tilt | | PRO-30272 | n/c | Foot Activated Tilt | | PRO-30273 | n/c |
| | | | | Not available with vent tray. | | | |

| Frame Colors | | | | | | | |
|---------------------|--|--|-----|------------------------|--|--|-----|
| Sugar Cookie | | | n/c | Blue Blast | | | n/c |
| Neon Yellow | | | n/c | Electric Blue | | | n/c |
| Yellow | | | n/c | Galaxy Blue | | | n/c |
| Bright Pink | | | n/c | Passion Purple | | | n/c |
| Pink Parfait | | | n/c | Octane Orange | | | n/c |
| Sugar Plum | | | n/c | Red Sunset | | | n/c |
| Candy Violet | | | n/c | Really Red | | | n/c |
| Candy Raspberry | | | n/c | Fiery Red | | | n/c |
| Extreme Green | | | n/c | Twilight | | | n/c |
| Sparkle Green | | | n/c | Stormy Gray | | | n/c |
| Mystic Teal | | | n/c | Black | | | n/c |
| Pacific Blue | | | n/c | Color to be Determined | | | |

Notes:

| Description | hccps | Part No. | Retail | Description | hccps | Part No. | Retail |
|--|---------|-----------|--------|---|---------|-----------|--------|
| Back Posts | | | | | | | |
| Back Post Options | | | | | | | |
| Fold-Down Back Posts | | PRO-30545 | n/c | Dynamic Rocker Back - Seating Dynamics | (E2398) | PRO-30542 | \$1985 |
| | | | | Angle adjustable from 72° to 108° in 6° increments. Can be locked-out when needed. | | | |
| | | | | Includes a set of 4 elastomers in different durometers to allow spring adjustment in the field. | | | |
| Reduced Shear Recline Back | (E1226) | PRO-30562 | \$1138 | Reduced Shear Recline Back with Angle-Adjustable Stroller Handle Extension | | PRO-30564 | \$1335 |
| Recline range is 65° from 97° to 162°. Pivot point is 2.5" above the seat rail. Solid Back Insert installation requires 4-points of attachment hardware. Seat-to-Back hinges are not compatible with any Super Reclining Reduced Shear Back. | | | | | | | |

| Push Handle Options | | | | | | | |
|---|--|-----------|-----|--------------------------|--|-----------|-----|
| 90° Bend Back Post Tubes | | PRO-30518 | n/c | Straight Back Post Tubes | | PRO-30544 | n/c |
| Requires selection of Push Button Angle Adjustable Stroller Handle (PRO-30571 or PRO-30576). Not available with Reduced Shear Recline Back. | | | | | | | |

| Push Handle Height | | | | | | | |
|---------------------------|--|-----------|-----|-----------------|--|-----------|-----|
| Short - 20" High | | PRO-30515 | n/c | Tall - 26" High | | PRO-30516 | n/c |
| Medium - 23" High | | PRO-30514 | n/c | | | | |

Push Handle Extension Options



All push handle extensions except PRO-30561 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

| | | | | | | | |
|--|--|-----------|--------|--|--|-----------|--------|
| Push Handle Extension - Width Adjustable | | PRO-30561 | n/c | Stroller Handle Extension - Width Adjustable, Removable | | PRO-30590 | \$ 95 |
| Not available with Straight Back Posts. Not available with Recline. | | | | Not available with Straight Back Posts. Not available with Recline. | | | |
| Push Button, Angle-Adjustable Stroller Handle Extension | | PRO-30571 | \$ 112 | Heavy Duty Push Button, Angle-Adjustable Stroller Handle Extension | | PRO-30576 | \$ 179 |
| Available on wheelchair frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. PRO-30571 or PRO-30576 required with Straight Back Post Tubes (PRO-30544) | | | | Available on any size frame. Metal hinge up to 180° range of adjustment. PRO-30571 or PRO-30576 required with Straight Back Post Tubes (PRO-30544) | | | |
| Individual Stroller Handle Extensions Removable - Pair | | PRO-30517 | \$ 110 | | | | |
| Not available with Straight Back Posts. Not available with Recline. | | | | | | | |

Footrests


Footrest / Extension Lengths

| Extension/Hanger | 75° Pop | 75° Lch | 90° Pop | 90 °Lch | ELR (90°) |
|------------------|---------|-------------|---------|---------|-----------|
| Hanger Only | 6"-8" | 5.5"-8.5" | 4"-8" | 4"-8" | 6"-12" |
| High Mount | N/A | N/A | 1.5"-3" | N/A | N/A |
| 4.5" H D Ext. | 9"-12" | 9.5"-12.5" | 9"-12" | 9"-12" | 13"-15" |
| 6.5" H D Ext. | 9"-15" | 9.5"-15.5" | 9"-15" | 9"-15" | 13"-17" |
| 4" Adjustable | 10"-13" | 9.5"-12.5" | 11"-13" | 10"-13" | 13"-15" |
| 6" Adjustable | 12"-15" | 10.5"-14.5" | 13"-15" | 12"-15" | 13"-18" |
| 8" Adjustable | 14"-17" | 12.5"-16.5" | N/A | N/A | 14"-20" |
| 10" Adjustable | 16"-19" | 14.5"-18.5" | N/A | N/A | 16"-22" |

MAPS Footrest / Extension Lengths

| Extension/Hanger | 75° Pop | 75° Lch | 90° Pop | 90 °Lch | ELR (90°) |
|------------------|-------------|---------|-------------|------------|--------------------|
| Hanger Only | 5.5"-7.5" | 5"-8" | 5.5"-7.5" | 4.5"-7.5" | 4.5"- 8.5" & 11.5" |
| High Mount | N/A | N/A | N/A | N/A | N/A |
| 4.5" H D Ext. | 8.5"-11.5" | 9"-12" | 8.5"-11.5" | 8.5"-11.5" | 12.5"-14.5" |
| 6.5" H D Ext. | 8.5"-13.5" | 9"-14" | 8.5"-13.5" | 8.5"-13.5" | 12.5"-16.5" |
| 4" Adjustable | 8.5"-11.5" | 9"-12" | 9.5"-11.5" | 8.5"-11.5" | 12.5"-14.5" |
| 6" Adjustable | 10.5"-13.5" | 10"-14" | 11.5"-13.5" | 8.5"-13.5" | 12.5"-17.5" |
| 8" Adjustable | 12.5"-15.5" | 12"-16" | N/A | N/A | 13.5"-19.5" |
| 10" Adjustable | 14.5"-17.5" | 14"-18" | N/A | N/A | 15.5"-21.5" |

Length is measured from top of footrest plate to top of seat rail. Without extensions bottom hanger hole not accessible for hardware. Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.



| Description | hccps | Part No. | Retail | Description | hccps | Part No. | Retail |
|---|---------|-----------|--------|---|-------|-----------|--------|
| Footrests - continued | | | | | | | |
| Hangers | | | | | | | |
| 75° Pop-Up Swing Away Hangers | | PRO-30280 | n/c | 75° Latch Style Swing Away Hangers | | PRO-30282 | n/c |
| 90° Pop-Up Swing Away Hangers | | PRO-30281 | n/c | 90° Latch Style Swing Away Hangers | | PRO-30283 | n/c |
| <i>Not available with 8" or 10" Telescoping Extensions.</i> | | | | <i>Not available with 8" or 10" Telescoping Extensions.</i> | | | |
| 90° Pop-Up Swing Away Elevating Legrests (ELRs) | (E0990) | PRO-30279 | \$ 295 | Omit Legrest Hangers & Footrests | | PRO-33000 | |
| <i>Not available with all Telescoping extensions.</i> | | | | | | | |

| | | | | | | | |
|---|---------|-----------|--------|--|---------|-----------|--------|
| Footplates & Footboards | | | | | | | |
| Aluminum Angle Adj. Flip-up Footplates | (K0040) | PRO-36033 | \$ 190 | Aluminum Angle Adjustable Foot-board Removable | (K0040) | PRO-36029 | \$ 190 |
| <i>Does not include Heel Loops, select below.</i> | | | | | | | |
| Composite Flip-Up | | PRO-17070 | n/c | Omit Footrest Plates | | PRO-36404 | |

| | | | | | | | |
|--|--|-----------|--------|----------------------|---------|-----------|-------|
| Footplate Options | | | | | | | |
| MAPS (Multi Axis Positioning System) | | PRO-36040 | \$ 325 | Heel Loops | (E0951) | PRO-36022 | \$ 45 |
| <i>Requires selection of Aluminum Angle Adjustable Footplates. Not compatible with Composite footplates or Aluminum Footboard. Not compatible with Seating Dynamics hangers.</i> | | | | Calf Strap - Webbing | (K0038) | PRO-36024 | \$ 45 |


| | | | | | | | |
|--|--|-----------|--------|--|--|-----------|-------|
| Hanger Extension Options | | | | | | | |
| 4" Telescoping Extensions | | PRO-30285 | n/c | 6" Telescoping Extensions | | PRO-30286 | n/c |
| 8" Telescoping Extensions | | PRO-30287 | n/c | 10" Telescoping Extensions | | PRO-30288 | n/c |
| <i>Not available with all 90° hangers.</i> | | | | <i>Not available with all 90° hangers.</i> | | | |
| 4.5" Heavy Duty Extensions | | PRO-36402 | \$ 58 | 6.5" Heavy Duty Extensions | | PRO-36403 | \$ 68 |
| <i>Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</i> | | | | | | | |
| High Mount Clamp & Tubes | | PRO-36410 | \$ 102 | Clear Selection | | | |

Specialized Footrests
 ALL legrest measurements are from the seat rail to footplate. Not available with LEAF or SRO seat rail options. Frame and/or front caster choice may limit adjustment range and/or cause interference.

Specialized Footrest Seat Rail to Footplate Measurements

| Part Number | Extension | Super Contracture | 3D Dynamic Hi Pivot | 3D Dynamic STD |
|-------------|-----------------|-------------------|---------------------|----------------|
| 6402-01 | Short | 10.5" - 13.5" | | |
| 6402-02 | Long | 13.5" - 18.5" | | |
| 6405-01 | 0"-2" Extension | | 6.25" - 10.75" | 12.75-15.25" |
| 6409-01 | 1"-3" Extension | | 7.25-11.75" | 13.75-16.25" |
| 6406-01 | 2"-4" Extension | | 8.25-12.75" | 14.75-17.25" |
| 6410-01 | 3"-5" Extension | | 9.25-13.75" | 15.75-18.25" |
| 6407-01 | 4"-6" Extension | | 10.25-14.75" | 16.75-19.25" |
| 6411-01 | 5"-7" Extension | | 11.25-15.75" | 17.75-20.25" |
| 6508-01 | 6"-8" Extension | | 12.25-16.75" | 18.75-21.25" |

Length measured from top of footrest plate to top of seat rail. 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.



| | | | | | | | |
|--|--|-------------|--------|--|--|-------------|--------|
| Contracture Footrest Hangers | | | | | | | |
| <i>Super Contracture for frames 14" wide and greater only.</i> | | | | | | | |
| Contracture Short - 11" to 15" | | PRO-30291-2 | \$3570 | Contracture Long - 15" to 19" | | PRO-30292-2 | \$3570 |
| <i>Includes short hanger extension.</i> | | | | <i>Includes long hanger extension.</i> | | | |

| | | | | | | | |
|------------------------------------|--|-------------|--------|---------------------------|--|-------------|--------|
| 3D Dynamic Footrest Hangers | | | | | | | |
| 3D Dynamic Footrest - Hi Pivot | | PRO-30297-2 | \$6215 | 3D Dynamic Footrest - STD | | PRO-30298-2 | \$5775 |

| | | | | | | | |
|---|--|-------------|-----|-------------------------------------|--|-------------|-----|
| Dynamic Footrest Extension Options | | | | | | | |
| Dynamic Adj Plantar/Dorsi Ext, 0-2" | | PRO-6405-01 | n/c | Dynamic Adj Plantar/Dorsi Ext, 1-3" | | PRO-6409-01 | n/c |
| Dynamic Adj Plantar/Dorsi Ext, 2-4" | | PRO-6406-01 | n/c | Dynamic Adj Plantar/Dorsi Ext, 3-5" | | PRO-6410-01 | n/c |
| Dynamic Adj Plantar/Dorsi Ext, 4-6" | | PRO-6407-01 | n/c | Dynamic Adj Plantar/Dorsi Ext, 5-7" | | PRO-6411-01 | n/c |
| Dynamic Adj Plantar/Dorsi Ext, 6-8" | | PRO-6408-01 | n/c | | | | |

| | | | | | | | |
|---|--|----------|-----|--|--|----------|-----|
| Specialized Footplates | | | | | | | |
| 4" x 6" HD Aluminum Footplates | | PRO-6503 | n/c | 5" x 8" HD Aluminum Footplates | | PRO-6504 | n/c |
| 6" x 10" HD Aluminum Footplates | | PRO-6505 | n/c | <i>For frames 16" wide and greater only.</i> | | | |
| <i>For frames 18" wide and greater only. Not available with Contracture Hangers</i> | | | | | | | |

| Description | hccps | Part No. | Retail | Description | hccps | Part No. | Retail |
|-------------|-------|----------|--------|-------------|-------|----------|--------|
|-------------|-------|----------|--------|-------------|-------|----------|--------|

Armrests

Back Post Mounted Armrests

Omit Armrests

Foam wrap pad STD on all back post mounted armrests..

| | | | | | | | |
|---|---------|-----------|--------|---|---------|-----------|--------|
| 10.5" Flip-Back | | PRO-38204 | n/c | 10.5" Flip-Back Locking | (E0973) | PRO-38201 | \$ 265 |
| 13.5" Flip-Back | | PRO-38202 | n/c | 13.5" Flip-Back Locking | (E0973) | PRO-38200 | \$ 265 |
| 10.5" to 13.5" Flip-Back Locking, Telescoping | (E0973) | PRO-38210 | \$ 285 | 17" to 20" Flip-Back Locking, Telescoping | (E0973) | PRO-38215 | \$ 285 |

Full length pad on 13.5" or longer armrests only.

Armrest Heights

| Armrest and Pad | Foam | Desk | Full |
|---------------------------------|--------------|--------------|--------------|
| Flip Back | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Flip Back Locking | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Telescoping Flip Back | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Flip Back Locking Telescoping | 6.5" - 11.5" | 7.5" - 12.5" | 7.5" - 12.5" |
| Single Post Height Adj. Short | | 8" - 14" | 8" - 14" |
| Single Post Height Adj. Tall | | 10" - 16" | 10" - 16" |
| Dual Post Fixed Height | | 9.5" | 9.5" |
| Dual Post Flip Back Height Adj. | | 9" - 13" | 9" - 13" |

All heights from top of seat frame to top of armrest pad.



Frame Mounted Armrests

| | | | | | | | |
|--|---------|-----------|--------|---|---------|-----------|--------|
| Single Post Height Adjustable, Removable - Low | (E0973) | PRO-38217 | \$ 285 | Single Post Height Adjustable, Removable - Standard | (E0973) | PRO-38203 | \$ 285 |
| Dual Post Fixed Height Removable | | PRO-38205 | n/c | Dual Post Flip-Back Height Adjustable, Removable | (E0973) | PRO-38206 | \$ 295 |

Arm Pad Options

| | | | | | | | |
|------------------|--|-----------|-----|------------------|--|-----------|-----|
| Desk Length Pads | | PRO-38213 | n/c | Full Length Pads | | PRO-38214 | n/c |
|------------------|--|-----------|-----|------------------|--|-----------|-----|

13.5" or longer armrests only.

Casters

| | | | | | | | |
|------------------------|--|-----------|-------|-------------------------------|--|-----------|--------|
| 4" x 1" Polyurethane | | PRO-40000 | n/c | 4" x 1.4" Poly, Aluminum Mag | | PRO-43000 | \$ 125 |
| 5" x 1" Polyurethane | | PRO-34555 | n/c | 5" Polyurethane Light-Up | | PRO-35050 | \$ 101 |
| 5" x 1.5" Polyurethane | | PRO-34222 | \$ 87 | 5" x 1.4" Poly, Aluminum Mag | | PRO-42000 | \$ 125 |
| 6" x 1" Polyurethane | | PRO-34666 | n/c | 6" x 1.4" Poly, Aluminum Mag | | PRO-44000 | \$ 125 |
| 6" x 1.5" Polyurethane | | PRO-34999 | \$ 87 | 6" x 2" Polyurethane, Treaded | | PRO-34444 | \$ 97 |
| 7" x 1" Polyurethane | | PRO-34777 | n/c | 7" x 2" Foam Filled | | PRO-34200 | \$ 200 |

Rear Wheels

Mags

| | | | |
|----------------------------------|--|-----------|-------|
| 12" Mag Wheel | | PRO-31525 | n/c |
| 16" Mag Wheel | | PRO-31530 | \$ 41 |
| 18" Mag Wheel | | PRO-31533 | \$ 41 |
| 20" Mag Wheel - Aluminum Handrim | | PRO-31527 | n/c |
| 22" Mag Wheel - Aluminum Handrim | | PRO-31528 | n/c |
| 24" Mag Wheel - Aluminum Handrim | | PRO-31529 | n/c |

Spoke

| | | | |
|---------------------------------------|--|-----------|--------|
| 16" Spoke Wheel - Black w/ Black Tire | | PRO-31535 | n/c |
| 18" Spoke Wheel - Black w/ Black Tire | | PRO-31538 | n/c |
| 20" Spoke Wheel - Aluminum Handrim | | PRO-31540 | \$ 204 |
| 22" Spoke Wheel - Aluminum Handrim | | PRO-31541 | \$ 204 |
| 24" Spoke Wheel - Aluminum Handrim | | PRO-31542 | \$ 204 |

20", 22" & 24" Spoke Wheels not available with drum brakes

Tire Options

| | | | |
|------------------------------------|---------|-----------|-------|
| Pneumatic Tire | | PRO-31040 | \$ 85 |
| Pneumatic Tire with Airless Insert | (E2213) | PRO-31000 | \$ 85 |
| Low Profile Polyurethane Tire | (E2218) | PRO-31045 | n/c |

12" Mag wheel only, not available with drum brakes

Handrim Options

| | | | |
|-------------------------|--|-----------|--------|
| Plastic Coated Handrims | | PRO-31005 | \$ 119 |
| Omit Handrims | | PRO-31006 | |

20", 22" & 24" Mag wheels only, not available with drum brakes

| Wheel/Caster Seat Rail to Floor Heights | | | | | | | Wheel/Caster Seat Rail to Floor Heights | | | | | | |
|---|--------------|-------|-------|--------|-------|--------|---|--------------|-------|-------|--------|-------|--------|
| Rear Wheel | Caster Wheel | | | | | | Rear Wheel | Caster Wheel | | | | | |
| | 4" | 5" | 6" | 6 x 2" | 7" | 7 x 2" | | 4" | 5" | 6" | 6 x 2" | 7" | 7 x 2" |
| 12" | 15" | 15" | 18.5" | 18.5" | 18.5" | 18.5" | 20" | 17" | 17" | 17" | 17" | 17" | 17" |
| | 15.5" | 15.5" | 19" | 19" | 19" | 19" | | 17.5" | 17.5" | 17.5" | 17.5" | 17.5" | 17.5" |
| | 16" | 16" | | | | | | 18" | 18" | 18" | 18" | 18" | 18" |
| | 18.5" | 18.5" | | | | | | 18.5" | 18.5" | 18.5" | 18.5" | 18.5" | 18.5" |
| | 19" | 19" | | | | | | 19" | 19" | 19" | 19" | 19" | 19" |
| 16" | 15" | 15" | 16" | 16" | 17" | 17" | 22" | 18" | 18" | 18" | 18" | 18" | 18" |
| | 15.5" | 15.5" | 16.5" | 16.5" | 17.5" | 17.5" | | 18.5" | 18.5" | 18.5" | 18.5" | 18.5" | 18.5" |
| | 16" | 16" | 17" | 17" | 20" | 20" | | 19" | 19" | 19" | 19" | 19" | 19" |
| | 16.5" | 16.5" | 17.5" | 17.5" | 20.5" | 20.5" | | 19.5" | 19.5" | 19.5" | 19.5" | 20" | 20" |
| | 17" | 17" | 20" | 20" | | | | 20" | 20" | 20" | 20" | | |
| 18" | 15.5" | 15.5" | 16" | 16" | 17" | 17" | 24" | 18.5" | 19" | 19" | 19" | 19" | 19" |
| | 16" | 16" | 16.5" | 16.5" | 17.5" | 17.5" | | 19" | 19.5" | 19.5" | 19.5" | 20" | 20" |
| | 16.5" | 16.5" | 17" | 17" | 18" | 18" | | | | | | | |
| | 17" | 17" | 17.5" | 17.5" | 18.5" | 18.5" | | | | | | | |
| | 17.5" | 17.5" | 18" | 18" | 18.5" | 18.5" | | | | | | | |

Clear selection

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".

| Wheel Locks | | | | Drum Lock | | | |
|---------------------------|---------|-----------|-------|---|---------|-----------|---------|
| Standard | | | | Drum Lock | | | |
| Push to Lock | | PRO-31501 | n/c | Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires. Not Available with 20", 22" & 24" Spoke Wheels. | | | |
| Pull to Lock | | PRO-31506 | n/c | Foot Activated Wheel Lock | | PRO-31508 | \$ 0 |
| Wheel Lock Options | | | | Not available with ALL vent and battery trays. | | | |
| Extension Handles | (E0961) | PRO-31801 | \$ 75 | Hand Activated Brake and Wheel Lock System | (E2228) | PRO-30275 | \$ 1000 |

| Anti-Tippers | | | |
|---|---------|-----------|--------|
| Curb Assist Anti-Tip Assembly Flip-Up, Adjustable | (E0971) | PRO-31600 | \$ 125 |
| Pair includes one anti-tip tube with curb assist feature. | | | |
| Spring Loaded Swing-Away Anti-Tip Assy - Foot Activated | (E0971) | PRO-31605 | \$ 180 |

| Seat Pans | | | |
|---|-------------------------|-----------|--------|
| ALL seat pans include 2" x 2" notches for growth adjustment. | | | |
| Non-Adjustable | | | |
| Fixed Solid Seat Pan - Composite | (E2231) | PRO-39010 | \$ 170 |
| ALL Bolteron seat pans are made to fit frame dimensions unless specified below. | | | |
| Seat Width _____" X | Seat Depth _____" | | |
| Fixed Solid Seat Pan - Aluminum | (E2231) | PRO-39030 | \$ 200 |
| Fixed Seat Pan raises the effective seat to floor height by 1/2". | | | |
| Available in the following sizes ONLY: | | | |
| 10" wide x 12" deep | 16" wide x 18" deep | | |
| 12" wide x 14" deep | 18" wide x 18" deep | | |
| 14" wide x 16" deep | | | |
| Adjustable | | | |
| Fixed, Width Adjustable Solid Seat Pan - Aluminum | (E2231) | PRO-39040 | \$ 290 |
| Fixed Seat Pan raises the effective seat to floor height by 1/2". | | | |
| Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with PRO frame. | | | |
| Available in the following sizes ONLY: | | | |
| 10"-12" wide x 12" deep | 14"-16" wide x 16" deep | | |
| 12"-14" wide x 14" deep | 16"-18" wide x 18" deep | | |

| Mounting Hardware Options | | | |
|--|---------|------------|--------|
| Use to convert fixed seat pan to removable. | | | |
| Removable Pan Hardware Kit | (K0108) | PRO-EZLITE | \$ 207 |
| Removable Seat Pan Hardware raises the effective seat to floor height by 5/8". | | | |
| Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps. | | | |
| 1" Drop Pan Hardware Kit | (K0108) | PRO-DROP1 | \$ 207 |
| 2" Drop Pan Hardware Kit | (K0108) | PRO-DROP2 | \$ 207 |
| Includes 4 each 2" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps. | | | |
| Not available with ANY Aluminum Solid Seat Pan. | | | |

| Description | hccps | Part No. | Retail | Description | hccps | Part No. | Retail |
|---|---------|-----------|--------|--|---------|-----------|--------|
| Accessories | | | | | | | |
| Tilt Range Limit Kit <small>Pins install into outside of tilt channel to limit tilt travel.</small> | (K0108) | PRO-10920 | \$ 35 | Freedom Backpack <small>Black Nylon only with Freedom Designs logo. Monogramming is not available on Freedom Backpack.</small> | (K0108) | PRO-10900 | \$ 60 |
| Footrest Hanger Pads - Pair <small>Wrap-around pads made from black neoprene.</small> | (K0108) | SFDI-9HP | \$ 97 | Freedom Canopy - Flip-Back <small>Canopy assembly includes removable mounting hardware. Available in black only with Freedom Designs logo. Includes rear flap with zippered storage pocket and clear-view window.</small> | (K0108) | PRO-10905 | n/c |
| Flip-Down Rigidizing Bar <small>Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.</small> | (K0108) | PRO-23210 | \$ 175 | | | | |

Notes:

| Specifications: | |
|------------------------|---|
| Frame Style | Tilt-In-Space |
| Base Frame Widths | 16", 18", 20" |
| Overall Width | Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel) |
| Seat Widths | 14" - 20" |
| Seat Depths | 14" - 26" |
| Back Heights | 20", 23", 26" |
| Adjustable Angle Back | 90°, 100°, 110°, 120° |
| Tilt Range | 0° to 55° |
| Recline Back | 162° |
| Weight Limit | 250 LBS. |
| Transport Weight | |
| Warranty | Lifetime on frame |

| Special Configurations Through Customer Service: | |
|---|---|
| Anterior Tilt | As much as -5° depending on configuration. Reduces posterior tilt angle by same amount. |
| Asymmetrical Seat Depth | Seat rails may be configured up to 4" offset |

HCPCS Codes:
are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.