



PRO-CG™

Pediatric Manual Wheelchair



Yes, you can!

Clear Form

Email Form to Customer Service

2241 N Madera Rd Simi Valley, CA 93065 Phone: (800) 331-8551 Fax: (888) 582-1509 www.freedomdesigns.com email: customerservice@freedomdesigns.com

Date: \_\_\_\_\_ Buyer: \_\_\_\_\_ Quote Order PO# \_\_\_\_\_
Company: \_\_\_\_\_ Location: \_\_\_\_\_
Marked For: \_\_\_\_\_ Phone: \_\_\_\_\_
Email: \_\_\_\_\_ Fax: \_\_\_\_\_

Table with 6 columns: Description, Part No., Retail, Description, Part No., Retail

Notice: A seating system must be purchased separately. No seating system is provided with the PRO-CG frame. No seat or back upholstery is available as an option.

Transport Optional Base

PRO-CG Pediatric Rigid (E1233) PRO-30300 \$2800

Transport Optional

Note: Transportation Option selected separately below.
Note: Includes Freedom Future Fit Program

Transport Package PRO-99998 \$ 250

Note: Must be installed by Freedom Designs manufacturing.
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.
Note: Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

Transport Included Base

PRO-CG Pediatric Rigid (E1233) PRO-30301 \$3050

Transport Included

Note: Includes Freedom Future Fit Program
Note: Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D.

Seat Rail Options

Note: Not available with ALL Contracture & 3D Dynamic footrest systems.

LEAF - Lower Extremity Accommodation Frame

LEAF Seat Frame- Pair (K0108) PRO-30324 \$ 714

LEAF Seat Frame- Left (K0108) PRO-30324-L \$ 357

LEAF Seat Frame- Right (K0108) PRO-30324-R \$ 357

Note: Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

SRO - Seat Rail Only

SRO Seat Frame- Pair (K0108) PRO-30325 \$ 306

SRO Seat Frame- Left (K0108) PRO-30325-L \$ 153

SRO Seat Frame- Right (K0108) PRO-30325-R \$ 153

Note: SRO is manufactured without standard footrest hanger attachments.

Tilt Activation Options

Dual Lever Hand Activated Tilt PRO-30272 n/c

Foot Activated Tilt PRO-30273 n/c

Note: Not with vent & battery trays.

Vent Tray Options

Note: Vent Ready Base required for Vent Tray.

Note: Vent Tray not with Foot Activated Tilt.

Vent Ready Base (K0108) PRO-30331 \$ 459

Note: Includes Push Button Adjustable Stroller Handle. For LP10 style vent and battery tray on frames 16" wide and greater.

Vent Tray - Fixed (E1029) 11-7000 \$ 439

Note: For LP10 style vent tray on 16" base width and greater.

Note: Vent Tray is full width of the base frame and 11" deep.

Quick Release Vent Tray Hardware (E1028) 563 \$ 272

Back Post Mounted Vent Bracket (E1029) 11-1100 \$ 550

Note: For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

Frame Color Options

Sugar Cookie n/c Blue Blast n/c

Neon Yellow n/c Electric Blue n/c

Yellow n/c Galaxy Blue n/c

Bright Pink n/c Passion Purple n/c

Pink Parfait n/c Octane Orange n/c

Sugar Plum n/c Red Sunset n/c

Candy Violet n/c Really Red n/c

Candy Raspberry n/c Fiery Red n/c

Extreme Green n/c Twilight n/c

Sparkle Green n/c Stormy Gray n/c

Mystic Teal n/c Black n/c

Pacific Blue n/c

PRO-CG Rider Weight Capacity is 250 lbs. All Prices subject to change without notice.

PRO-CG Pediatric Note: All Specifications and dimensions are approximate

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Frame Dimensions


#### Specify Frame Width & Frame Depth Settings

**Note:** *PRO-CG* has 4" of depth adjustment for growth. Depth is measured from front of back post to front edge of the seat frame. For back post location add seating/hardware thickness (typically 2") to client's actual seat depth.

**Example:** seat depth of 12" + back/hardware thickness of 2" = back post location set at 14"

**Note:** The *PRO-CG* has 6" of center of gravity (CG) adjustment on all adult frames. The default CG setting is 1 position forward of either the client's actual seat depth or the manufacturer's seat depth setting.

Specify Client's Actual Seat Depth: \_\_\_\_\_"



<i>PRO-CG</i> Base Frame Width	Actual Client Seat Depth Specification or Specified Manufacturer Seat Depth Setting															
	10"	11"	12"	13"	14"	15"	16"	17"	18"	19"	20"	21"	22"	23"	24"	
14" Base Frame Width			-2*	-2	-1	0	1	2	3	3	3					
16" Base Frame Width					-3	-2	-1	0	1	2	3	3	3			
18" Base Frame Width								-3	-2	-1	0	1	2	3	3	

**Note:** Seat rail extensions required for manufacturing seat depth settings with CG setting in gray area.

\* **Note:** Maximum negative CG setting on 14" Base Frame Width is -2

**14" Base Frame Width – Seat frame** PRO-30220  
adjusts 12" to 14" wide in 1" increments.

#### Specify Seat Width Setting

12"	PRO-12W	n/c
13"	PRO-13W	n/c
14"	PRO-14W	n/c

**12" To 20" Seat Frame Length** PRO-30511

#### Specify Seat Depth Setting

12"	PRO-30504	n/c
13"	PRO-30505	n/c
14"	PRO-30506	n/c
15"	PRO-30521	n/c
16"	PRO-30522	n/c
* 17" - Must select 2" extension	PRO-30523	n/c
* 18" - Must select 2" extension	PRO-30524	n/c
** 19" - Must select 4" extension	PRO-30527	n/c
** 20" - Must select 4" extension	(E2203) PRO-30528	\$ 500
2" Extension - Pair	(K0108) PRO-30525	\$ 26
4" Extension - Pair	(K0108) PRO-30531	\$ 26

\* **Note:** Must choose 2" extension for 17" and 18" frame lengths.

\*\* **Note:** Must choose 4" extension for 19" and 20" frame lengths.

Specify Back Post Location (if other than specified above) \_\_\_\_\_"

### Frame Dimensions (continued)

#### Specify Frame Width & Frame Depth Settings

**16" Base Frame Width – Seat** PRO-30240  
frame adjusts 14" to 16" wide in 1" increments.

#### Specify Seat Width Setting

14"	PRO-14W	n/c
15"	PRO-15W	n/c
16"	PRO-16W	n/c

**14" To 22" Seat Frame Length** PRO-30512

#### Specify Seat Depth Setting

14"	PRO-30506	n/c
15"	PRO-30521	n/c
16"	PRO-30522	n/c
17"	PRO-30523	n/c
18"	PRO-30524	n/c
* 19" - Must select 2" extension	PRO-30527	n/c
* 20" - Must select 2" extension	(E2203) PRO-30528	\$ 500
** 21" - Must select 4" extension	(E2203) PRO-30529	\$ 500
** 22" - Must select 4" extension	(E2204) PRO-30535	\$ 900
2" Extension - Pair	(K0108) PRO-30525	\$ 26
4" Extension - Pair	(K0108) PRO-30531	\$ 26

\* **Note:** Must choose 2" extension for 19" and 20" frame lengths.

\*\* **Note:** Must choose 4" extension for 21" and 22" frame lengths.

Specify Back Post Location (if other than specified above) \_\_\_\_\_"

**18" Base Frame Width – Seat frame** PRO-30260  
adjusts 16" to 18" wide in 1" increments.

#### Specify Seat Width Setting

16"	PRO-16W	n/c
17"	PRO-17W	n/c
18"	PRO-18W	n/c

**16" To 24" Seat Frame Length** PRO-30520

#### Specify Seat Depth Setting

16"	PRO-30522	n/c
17"	PRO-30523	n/c
18"	PRO-30524	n/c
19"	PRO-30527	n/c
20"	(E2203) PRO-30528	\$ 500
* 21" - Must select 2" extension	(E2203) PRO-30529	\$ 500
* 22" - Must select 2" extension	(E2204) PRO-30535	\$ 900
** 23" - Must select 4" extension	(E2204) PRO-30536	\$ 900
** 24" - Must select 4" extension	(E2204) PRO-30537	\$ 900
2" Extension - Pair	(K0108) PRO-30525	\$ 26
4" Extension - Pair	(K0108) PRO-30531	\$ 26

\* **Note:** Must choose 2" extension for 21" and 22" frame lengths.

\*\* **Note:** Must choose 4" extension for 23" and 24" frame lengths.

Specify Back Post Location (if other than specified above) \_\_\_\_\_"

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Back Post Options

#### Specify Push Handle, Back Post and Extensions

**Standard 90° Bend Push Handle Back Post Tubes** PRO-30518 n/c

**Ergo Bend Push Handle Back Post Tubes** PRO-30519 n/c

Note: All push/stroller handle extensions not recommended with Ergo Bend.

**Seating Dynamics Rocker Back** (K0108) PRO-30542 **\$1985**

Note: Angle adjustable from 72° to 108° in 6° increments.  
 Note: Supplied with a complete set of 4 elastomers in different durometers to allow spring adjustment in the field.  
 Note: Locks into a non-dynamic position when necessary.  
 Note: Does not allow the back post tubes to fold forward.

**Super Reclining Reduced Shear Back** (K0108) PRO-30562 **\$1138**

**Super Reclining Reduced Shear Back with Push Button, Angle-Adjustable Stroller Handle Extension** (K0108) PRO-30564 **\$1335**

Note: Pivot point is 2.5" above the seat rail.  
 Note: Recline range is 65° from 97° to 162°.  
 Note: Includes height adjustable back post tubes for 21" - 27" heights.  
 Note: Available with standard 90° push handles only.  
 Note: Only compatible with Push Button, Angle-Adjustable Folding Stroller Handle Extension.  
 Note: Does not allow the back post tubes to fold forward.  
 Note: Solid Back Insert installation requires 4-points of attachment hardware. Seat-to-Back hinges are not compatible with any Super Reclining Reduced Shear Back.

#### Back Post Height Setting

Note: Back Posts are measured from top of seat rail to top of handle bend.

**Short Push Handle - 20" high** PRO-30515 n/c

**Medium Push Handle - 23" high** PRO-30514 n/c

**Tall Push Handle - 26" high** PRO-30516 n/c

### Stroller Handle Extension Options

Note: All stroller handle ext. add approximately 6.5" to push handle height.  
 Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.



**Push Handle Extension - Width Adjustable - Mounts direct to back post push handle** PRO-30561 n/c

**Stroller Handle Extension - Width Adjustable, Removable - Mounts direct to back post push handle** PRO-30590 **\$ 95**

**Push Button, Angle-Adjustable Stroller Handle Extension - Mounts direct to back post push handle** PRO-30571 **\$ 112**

**Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty - Mounts direct to back post tubes** PRO-30576 **\$ 179**

Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.  
 Note: Not removable.

**Push Button, Angle-Adjustable Folding Stroller Handle Extension - Heavy Duty - Mounts direct to back post tubes** PRO-30576 **\$ 179**

Note: Available 10W-20w wheelchair frames.  
 Note: Not removable.

### Stroller Handle Extension Options - Continued

Note: All stroller handle ext. add approximately 6.5" to push handle height.  
 Note: All stroller/push handle extensions except Individual include a black neoprene wrap around pad.

**Push Button, Angle-Adjustable, Folding Push Handle Extension - Mounts direct to straight back post tubes** PRO-30581 **\$ 112**

Note: Available 10W-15W wheelchair frames only. For frames wider than 15W must choose Heavy Duty below.

Note: Not removable.

**Push Button, Angle-Adjustable Folding Push Handle Extension - Heavy Duty - Mounts direct to straight back post tubes** PRO-30586 **\$ 179**

Note: Available 10W-20w wheelchair frames.  
 Note: Not removable.

**Individual Stroller Handle Extensions Removable - Pair** PRO-30517 **\$ 102**

### PRO Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-10"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	N/A
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	N/A
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	N/A
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	N/A
10" Adjustable	16"-19"	14.5"-18.5"	N/A	N/A	N/A

### PRO MAPS Footrest / Extension Lengths

Extension/Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"- 8.5" & 11.5"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-15.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	N/A
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	N/A
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	N/A
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	N/A
10" Adjustable	14.5"-17.5"	14"-18"	N/A	N/A	N/A

Note: Length from top of footrest plate to top of seat rail.  
 Note: Without extensions bottom hanger hole not accessible for hardware.  
 Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.

### Legrests, Footrests and Extension Options

#### Legrests

**Omit Legrest Hangers & Footrests** PRO-33000

**75° Pop-Up Swing Away Hangers** PRO-30280 n/c

**90° Pop-Up Swing Away Hangers** PRO-30281 n/c

Note: Not available with 8" adjustable extensions.  
**90° Pop-Up Swing Away Elevating Legrests (ELRs)** (E0990) PRO-30279 **\$ 255**

Note: Not available with all length adjustable extensions.  
 Note: Not available with 6.5" heavy duty extensions.

**75° Latch Style Swing Away Hangers** PRO-30282 n/c

**90° Latch Style Swing Away Hangers** PRO-30283 n/c

Note: Not available with 8" adjustable extensions.

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Legrests, Footrests and Extension Options

#### Footplates & Footboards

<b>Omit Footrest Plates</b>	PRO-36404	
<b>Composite Flip-Up</b>	PRO-17070	n/c
<b>Kids Aluminum Angle Adjustable, Flip-Up</b>	(K0040) PRO-36020	\$ 160
<i>Note: For frames 14" wide and narrower.</i>		
<b>Heel Loops</b>	(E0951) PRO-36022	\$ 41
<b>Calf Strap - Webbing</b>	(K0038) PRO-36024	\$ 41
<b>Adult Aluminum Angle Adj. Flip-up w/ Heel Loops</b>	(K0040/E0951) PRO-36021	\$ 180
<i>Note: For frames 15" wide and greater.</i>		
<b>Omit Heel Loops</b>	PRO-36023	n/c
<b>Aluminum Angle Adjustable Footboard Removable</b>	(K0040) PRO-36029	\$ 160

#### MAPS Footplate Mounting Hardware

<b>MAPS (Multi Axis Positioning System) - Pair</b>	PRO-36040	\$ 306
<b>MAPS (Multi Axis Positioning System) - Left</b>	PRO-36040-L	\$ 153
<b>MAPS (Multi Axis Positioning System) - Right</b>	PRO-36040-R	\$ 153
<i>Note: MAPS connects footrest hangers to aluminum footrest plates to adjust for fixed foot and ankle positioning. MAPS footrest plates rotate and adjust for plantar and dorsiflexion and inversion and eversion positioning.</i>		
<i>Note: Must choose aluminum footplates.</i>		
<i>Note: MAPS not available with composite footplates or aluminum footboard.</i>		
<i>Note: MAPS not available on 10" &amp; 11" wide frames.</i>		

#### Legrest Extensions

<b>4" Length Adjustable Extensions</b>	PRO-30285	n/c
<b>6" Length Adjustable Extensions</b>	PRO-30286	n/c
<b>8" Length Adjustable Extensions</b>	PRO-30287	n/c
<b>10" Length Adjustable Extensions</b>	PRO-30288	n/c
<i>Note: 8" &amp; 10" Not available with all 90° hangers.</i>		
<b>4.5" Heavy Duty ELR Extensions</b>	PRO-36401	\$ 56
<b>6.5" Heavy Duty ELR Extensions</b>	PRO-36405	\$ 56
<b>4.5" Heavy Duty Extensions</b>	PRO-36402	\$ 56
<b>6.5" Heavy Duty Extensions</b>	PRO-36403	\$ 66
<i>Note: Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.</i>		
<b>High Mount Brackets</b>	PRO-36407	\$ 102
<i>Note: Includes flip-up mounting hardware for future growth adjustment.</i>		
<i>Note: 75° hangers use a tube &amp; clamp style hardware.</i>		

Specify Knee to Heel Measurement \_\_\_\_\_" (see matrix)

#### Notes:

### Specialized Footrest Seat Rail to Footplate Measurements

Part Number	Extension	Super Contracture	3D Dynamic Hi Pivot	3D Dynamic STD
6402-01	Short	10.5" - 13.5"		
6402-02	Long	13.5" - 18.5"		
6405-01	0"-2" Extension		6.25" - 10.75"	12.75-15.25"
6409-01	1"-3" Extension		7.25-11.75"	13.75-16.25"
6406-01	2"-4" Extension		8.25-12.75"	14.75-17.25"
6410-01	3"-5" Extension		9.25-13.75"	15.75-18.25"
6407-01	4"-6" Extension		10.25-14.75"	16.75-19.25"
6411-01	5"-7" Extension		11.25-15.75"	17.75-20.25"
6508-01	6"-8" Extension		12.25-16.75"	18.75-21.25"

**Note:** Length from top of footrest plate to top of seat rail.  
**Note:** 3D HiPivot includes 4.5" of length adjustment, 3D includes 2.5" of length adjustment.

### Specialized Footrest Systems

**Note:** ALL legrest measurements are from the seat rail to footplate.  
**Note:** Not available with LEAF or SRO seat rail options.

#### Super Contracture Footrest Hanger (K0108)

**Note:** Super Contracture for frames 14" wide and greater.

**Contracture Short - 11" to 15"** PRO-30291-2 **\$3570**

**Note:** See Matrix for seat to rail footplate measurements.

#### Includes:

<b>Adj Contracture Upper</b>	PRO-6302-03	Incl.
<b>Footrest Adapter HD</b>	PRO-6220-01	Incl.
<b>Adj Plantar/Dorsi Ext. Short</b>	PRO-6401-02	Incl.

#### Footplate Options

<b>4" x 6" HD Alu Footplates</b>	PRO-6503	n/c
<b>5" x 7" HD Alu Footplates</b>	PRO-6504	n/c

**Note:** for frames 16" wide and greater only.

**Contracture Long - 15" to 19"** PRO-30292-2 **\$3570**

**Note:** See Matrix for seat to rail footplate measurements.

#### Includes:

<b>Adj Contracture Upper</b>	PRO-6302-03	Incl.
<b>Footrest Adapter HD</b>	PRO-6220-01	Incl.
<b>Adj Plantar/Dorsi Ext. Long</b>	PRO-6402-02	Incl.

#### Footplate Options

<b>4" x 6" HD Alu Footplates</b>	PRO-6503	n/c
<b>5" x 7" HD Alu Footplates</b>	PRO-6504	n/c

**Note:** Super Contracture is both swing-away and removable.

**Note:** Super Contracture has adjustable knee joints that have 48 positions that adjust in 7.5° increments.

**Note:** Frame and/or front caster choice may limit Super Contracture adjustment range and/or cause footplate/front caster interference.

**Note:** Super Contracture position may limit ankle joint/ footplate adjustment range.

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Specialized Footrest Systems - Continued

#### 3D Dynamic Footrest Hanger (K0108)

Note: 3D Dynamic for frames 14" wide and greater.

**3D Dynamic Footrest - Hi Pivot** PRO-30297-2 **\$6215**

**Includes:**

**HiPivot Dyn Adj Knee & Dyn Tel Up** PRO-6304-02 **Incl.**

**Footrest Adapter HD** PRO-6220-01 **Incl.**

#### Extension Options

Note: See Matrix for seat to rail footplate measurements.

**Dynamic Adj Plantar/Dorsi Ext, 0-2"** PRO-6405-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 1-3"** PRO-6409-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 2-4"** PRO-6406-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 3-5"** PRO-6410-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 4-6"** PRO-6407-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 5-7"** PRO-6411-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 6-8"** PRO-6408-01 **n/c**

#### Footplate Options

**4" x 6" HD Alu Footplates** PRO-6503 **n/c**

**5" x 7" HD Alu Footplates** PRO-6504 **n/c**

Note: for frames 16" wide and greater only.

**6" x 10" HD Alu Footplates** PRO-6505 **n/c**

Note: for frames 18" wide and greater only.

**3D Dynamic Footrest - STD** PRO-30298-2 **\$5775**

**Includes:**

**Dynamic Adj Knee & Dyn Tel Up** PRO-6304-01 **Incl.**

**Footrest Adapter HD** PRO-6220-01 **Incl.**

#### Extension Options

Note: See Matrix for seat to rail footplate measurements.

**Dynamic Adj Plantar/Dorsi Ext, 0-2"** PRO-6405-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 1-3"** PRO-6409-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 2-4"** PRO-6406-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 3-5"** PRO-6410-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 4-6"** PRO-6407-01 **n/c**

**Dynamic Adj Plantar/Dorsi Ext, 5-7"** PRO-6411-01 **n/c**

#### Footplate Options

**4" x 6" HD Alu Footplates** PRO-6503 **n/c**

**5" x 7" HD Alu Footplates** PRO-6504 **n/c**

Note: for frames 16" wide and greater only.

**6" x 10" HD Alu Footplates** PRO-6505 **n/c**

Note: for frames 18" wide and greater only.

Note: 3D Dynamic is both swing-away and removable.

Note: 3D Dynamic has adjustable dynamic knee joints that can be set at 40°, 50°, 70°, or 90° with 36° of dynamic elevation travel.

Note: 3D Dynamic has ankle joints with 18° of dynamic flexion for plantar/dorsiflexion ankle/foot positioning.

Note: Angle adjustable footplate included.

Note: Frame and/or front caster choice may limit 3D Dynamic adjustment range and/or cause footplate/front caster interference.

### Freedom Designs Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Single Post Fxd Hgt Short		8" - 10.5"	8" - 10.5"
Single Post Fxd Hgt Short		9" - 11.5"	9" - 11.5"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

Note: All heights from top of seat frame to top of armrest pad.

### Armrests and Armrest Pad Options

**Omit Armrests** PRO-38209

#### Back Post Mounted Armrests

Note: Flip-Back armrests length is measured from the front of the backpost tube to the end of the armrest with foam wrap pad installed.

**10.5" Long Flip-Back** PRO-38204 **n/c**

Note: Not available with full length arm pads

**13.5" Long Flip-Back** PRO-38202 **n/c**

**10.5" Long Flip-Back Locking** (E0973) PRO-38201 **\$ 245**

Note: Not available with full length arm pads

**13.5" Long Flip-Back Locking** (E0973) PRO-38200 **\$ 245**

**10.5" to 13.5" Long Flip-Back, Telescoping** (E0973) PRO-38208 **\$ 240**

Note: Full length pad on 13.5" or longer armrests only.

**10.5" to 13.5" Long Flip-Back Locking, Telescoping** (E0973) PRO-38210 **\$ 265**

Note: Full length pad on 13.5" or longer armrests only.

**17" to 20" Long Flip-Back Locking, Telescoping** (E0973) PRO-38215 **\$ 265**

**Foam Wraps - Standard** PRO-30216 **n/c**

**Desk Length Pads** PRO-38213 **n/c**

**Full Length Pads** PRO-38214 **n/c**

Note: Full length pad on 13.5" or longer armrests only.

#### Frame Mounted Armrests

**Single Post Fixed Height, Removable - Short** (K0015) PRO-38222 **\$ 225**

**Single Post Fixed Height, Removable - Tall** (K0015) PRO-38224 **\$ 225**

Note: Single Post Fixed Height arms include a removable side guard.

Note: Adjustable in height with tools only. See matrix for height ranges.

**Single Post Height Adjustable, Removable - Short** (E0973) PRO-38217 **\$ 265**

**Single Post Height Adjustable, Removable - Tall** (E0973) PRO-38203 **\$ 265**

**Dual Post Fixed Height Removable** PRO-38205 **n/c**

**Dual Post Flip-Back Height Adjustable, Removable** (E0973) PRO-38206 **\$ 275**

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

**Desk Length Pads** PRO-38213 **n/c**

**Full Length Pads** PRO-38214 **n/c**

Note: All dual post armrests not available with short frame.

Note: All dual post armrests not available with 22" or 24" wheels.

**Desk Length Pads** PRO-38213 **n/c**

**Full Length Pads** PRO-38214 **n/c**

Specify Seat Rail to Top of Armrest Measurement \_\_\_\_\_ "

**PRO-CG Rider Weight Capacity is 250 lbs.**  
All Prices subject to change without notice.

**PRO-CG Pediatric**  
Note: All Specifications and dimensions are approximate

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Caster Wheel				
	4"	5"	6"	7"	
12"	15"	15"	18.5"	18.5"	
	15.5"	15.5"	19"	19"	
	16"	16"			
	18.5"	18.5"			
	19"	19"			
16"	15"	15"	16"	17"	
	15.5"	15.5"	16.5"	17.5"	
	16"	16"	17"	20"	
	16.5"	16.5"	17.5"	20.5"	
	17"	17"	20"		
18"	15.5"	15.5"	16"	17"	
	16"	16"	16.5"	17.5"	
	16.5"	16.5"	17"	18"	
	17"	17"	17.5"	18.5"	
	17.5"	17.5"	18"	18.5"	
20"	17"	17"	17"	17"	
	17.5"	17.5"	17.5"	17.5"	
	18"	18"	18"	18"	
	18.5"	18.5"	18.5"	18.5"	
	19"	19"	19"	19"	
22"	18"	18"	18"	18"	
	18.5"	18.5"	18.5"	18.5"	
	19"	19"	19"	19"	
		19.5"	19.5"	20"	
			20"		
24"	18.5"	19"	19"	19"	
	19"	19.5"	19.5"	20"	
			20"	20.5"	

**Note:** Use chart to determine seat rail to floor height. Seat rail to floor height determined by frame type, rear wheels, and front caster size.

**Note:** Caster sizes vary slightly in height by type.

**Note:** Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4".

**Note:** Seat pan choice increases seat rail to floor height 3/8" to 3/4".

### Front Wheels - Caster Options

<b>Omit Caster Wheels</b>	PRO-34333	
<b>3" Polyurethane</b>	PRO-34888	n/c
<b>3" Polyurethane Light-Up</b>	PRO-34111	\$ 101
<b>4" Polyurethane</b>	PRO-40000	n/c
<b>4" Polyurethane Light-Up</b>	PRO-41000	\$ 101
<b>4" x 1.4" Poly, Aluminum Mag</b>	PRO-43000	\$ 175
<b>5" x 1" Poly, Plastic Mag</b>	PRO-34555	n/c
<b>5" Polyurethane Light-Up</b>	PRO-35050	n/c
<b>5" x 1.5" Poly, Plastic Mag</b>	PRO-34222	\$ 87
<b>5" x 1.4" Poly, Aluminum Mag</b>	PRO-42000	\$ 175

### Front Wheels - Caster Options - Continued

<b>6" x 1" Poly, Plastic Mag</b>	PRO-34666	n/c
<b>6" x 1.5" Poly, Plastic Mag</b>	PRO-34999	\$ 87
<b>6" x 1.4" Poly, Aluminum Mag</b>	PRO-44000	\$ 175
<b>6" x 2" Poly with Tread</b>	PRO-34444	\$ 92
<b>6" x 1" Pneumatic Plastic Mag</b>	(E2214/ E2215) PRO-34667	\$ 77
<b>7" x 1" Poly, Plastic Mag</b>	PRO-34777	n/c
<b>7" x 2" Polyurethane</b>	PRO-34200	\$ 118

### Rear Wheels - Tires - Accessories

**Note:** All wheel selections include Freedom Designs' Patented Double Locking Safety Axles.

<b>Omit Rear Wheels and Axles</b>	PRO-31524	
<b>Mag Wheels</b>		
<b>12" Low Profile Poly Mag</b>	(E2218) PRO-31526	\$ 36
<b>Note:</b> Not available with any foot/hand activated brake/wheel lock system.		
<b>12" Semi Pneumatic Plastic Mag</b>	PRO-31523	\$ 168
<b>Note:</b> Not available with airless inserts.		
<b>12" Pneumatic Plastic Mag</b>	PRO-31525	n/c
<b>12" Airless Inserts</b>	(E2213) PRO-31019	\$ 77
<b>16" Pneumatic Plastic Mag</b>	PRO-31530	\$ 41
<b>16" Airless Inserts</b>	(E2213) PRO-31020	\$ 77
<b>16" Aluminum Handrims</b>	PRO-31014	\$ 169
<b>18" Pneumatic Plastic Mag</b>	PRO-31533	\$ 41
<b>18" Airless Inserts</b>	(E2213) PRO-31021	\$ 77
<b>18" Aluminum Handrims</b>	PRO-31015	\$ 169
<b>20" Pneumatic Plastic Mag with Aluminum Handrim</b>	PRO-31527	n/c
<b>20" Airless Inserts</b>	(E2213) PRO-31022	\$ 77
<b>20" Plastic Coated Handrims</b>	PRO-31016	\$ 119
<b>Note:</b> Plastic Coated Handrims NOT available with Drum Braking Systems.		
<b>22" Pneumatic Plastic Mag with Aluminum Handrim</b>	PRO-31528	n/c
<b>22" Airless Inserts</b>	(E2213) PRO-31023	\$ 77
<b>22" Plastic Coated Handrims</b>	PRO-31017	\$ 119
<b>Note:</b> Plastic Coated Handrims NOT available with Drum Braking Systems.		
<b>24" Pneumatic Plastic Mag with Aluminum Handrim</b>	PRO-31529	n/c
<b>24" Airless Inserts</b>	(E2213) PRO-31024	\$ 77
<b>24" Plastic Coated Handrims</b>	PRO-31018	\$ 119
<b>Note:</b> Plastic Coated Handrims NOT available with Drum Braking Systems.		
<b>Omit Handrims from 20", 22" or 24" Wheels</b>	PRO-31006	

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Rear Wheels - Tires - Accessories - Cont.

#### Spoke Wheels

Note: All spoke wheels not available plastic coated handrims.

**16" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel** PRO-31534 n/c

Note: Rim, hub & spokes are anodized silver.

**16" Airless Inserts** (E2213) PRO-31020 \$ 77

**16" Silver Aluminum Handrim** PRO-31030 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**16" Black Pneumatic Tire & Black Spoke Aluminum Wheel** PRO-31535 n/c

Note: Rim, hub & spokes are anodized black.

**16" Airless Inserts** (E2213) PRO-31020 \$ 77

**16" Black Aluminum Handrim** PRO-31031 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**18" Grey Pneumatic Tire & Silver Spoke Aluminum Wheel** PRO-31537 n/c

Note: Rim, hub & spokes are anodized silver.

**18" Airless Inserts** (E2213) PRO-31021 \$ 77

**18" Silver Aluminum Handrim** PRO-31034 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**18" Black Pneumatic Tire & Black Spoke Aluminum Wheel** PRO-31538 n/c

Note: Rim, hub & spokes are anodized black.

**18" Airless Inserts** (E2213) PRO-31021 \$ 77

**18" Black Aluminum Handrim** PRO-31035 \$ 77

Note: 16" & 18" spoke wheels do not include handrims.

**20" Pneumatic Spoke Wheel with Aluminum Handrim** PRO-31540 \$ 204

**20" Airless Inserts** (E2213) PRO-31022 \$ 77

**22" Pneumatic Spoke Wheel with Aluminum Handrim** PRO-31541 \$ 204

**22" Airless Inserts** (E2213) PRO-31023 \$ 77

**24" Pneumatic Spoke Wheel with Aluminum Handrim** PRO-31542 \$ 204

**24" Airless Inserts** (E2213) PRO-31024 \$ 77

**Omit Handrims from 20", 22" or 24" Wheels** PRO-31006

### Wheel Lock Options

**Push to Lock - Standard** PRO-31501 n/c

**Pull to Lock** PRO-31506 n/c

**Extension Handles** (E0961) PRO-31801 \$ 61

**Foot Activated Wheel Lock** (E2228) PRO-31508 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Not available with ALL vent and battery trays.

**Hand Activated Brake and Wheel Lock System** (E2228) PRO-30275 n/c

Note: Includes manual wheel locks to be used with public transportation

Note: Not available with 12" mag wheels with low pro poly tires.

Note: Used to control direction, speed, brake and lock wheels from strolling position.

### Rear Anti-Tipper Options

**Curb Assist Anti-Tip Assembly Flip-Up, Adjustable** (E0971) PRO-31600 \$ 97

Note: Pair includes one anti-tip tube with curb assist feature.

**Spring Loaded Swing-Away Anti-Tip Assy - Foot Activated** (E0971) PRO-31605 \$ 180

### Seat Pan Options

Note: ALL seat pans include 2" x 2" notches for growth adjustment.

**Fixed Solid Seat Pan - Bolteron** (E2231) PRO-39010 \$ 170

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Note: ALL Bolteron seat pans are made to fit frame dimensions unless specified below.

Seat Width \_\_\_\_\_" X Seat Depth \_\_\_\_\_"

**Fixed Solid Seat Pan - Aluminum** (E2231) PRO-39030 \$ 200

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

**12" wide x 14" deep**                      **16" wide x 18" deep**  
**14" wide x 16" deep**                      **18" wide x 18" deep**

**Fixed, Width Adjustable Solid Seat Pan - Aluminum** (E2231) PRO-39040 \$ 290

Note: Fixed Seat Pan raises the effective seat to floor height by 1/2".

Note: Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with PRO frame.

Note: Not available with folding frame.

Note: For folding frame must select Removable Pan Hardware Kit below.

Available in the following sizes ONLY:

**12"-14" wide x 14" deep**                      **16"-18" wide x 18" deep**  
**14"-16" wide x 16" deep**

### Mounting Hardware Options

Note: Use to convert fixed seat pan to removable.

**Removable Pan Hardware Kit** (K0108) PRO-EZLITE \$ 198

Note: Removable Seat Pan Hardware raises the effective seat to floor height by 5/8".

Note: Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

**1" Drop Pan Hardware Kit** (K0108) PRO-DROP1 \$ 198

Note: Includes 4 each 1" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps.

Note: Not available with ANY Aluminum Solid Seat Pan.

**2" Drop Pan Hardware Kit** (K0108) PRO-DROP2 \$ 198

Note: Includes 4 each 2" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps.

Note: Not available with ANY Aluminum Solid Seat Pan.

Date: \_\_\_\_\_ PO# \_\_\_\_\_ Marked For: \_\_\_\_\_

Description	Part No.	Retail	Description	Part No.	Retail
-------------	----------	--------	-------------	----------	--------

### Accessories

<b>Tilt Range Limit Kit</b>	(K0108)	PRO-10920	<b>\$ 35</b>
<b>Note:</b> Pins install into outside of tilt channel to limit tilt travel.			
<b>Medical Necessities Pouch</b>	(K0108)	PRO-23207	<b>\$ 68</b>
<b>Note:</b> Black nylon.			
<b>Footrest Hanger Pads - Pair</b>	(K0108)	SFDI-9HP	<b>\$ 92</b>
<b>Note:</b> Wrap-around pads made from black neoprene.			
<b>Frame Pads - Pair</b>	(K0108)	10911-2	<b>\$ 82</b>
<b>Note:</b> Wrap-around pads made from black neoprene.			
<b>Note:</b> Typically used to pad exposed back post tubes.			
<b>Flip-Down Rigidizing Bar</b>	(K0108)	PRO-23210	<b>\$ 194</b>
<b>Note:</b> Mounts to back post tubes to add rigidity to frame.			
<b>Note:</b> Does NOT adjust in width to grow with frame.			
<b>Freedom Backpack</b>	(K0108)	PRO-10900	<b>\$ 56</b>
<b>Note:</b> Black Nylon only with Freedom Designs logo.			
<b>Note:</b> Monogramming is not available on Freedom Backpack.			
<b>Freedom Canopy - Flip-Back</b>	(K0108)	PRO-10905	<b>\$ 128</b>
<b>Note:</b> Canopy assembly includes removable mounting hardware.			
<b>Note:</b> Canopy is available in black only with Freedom Designs logo.			
<b>Note:</b> Canopy includes rear flap with zippered storage pocket.			
<b>Note:</b> Canopy includes clear-view window.			

**HCPCS CODES** are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

**Notes:**

### Anterior Support Surfaces – (Trays)

<b>3/8" Lexan Tray - Channel Bracket Hardware</b>	(E0950)	08-2000	<b>\$ 271</b>
<b>3/8" Lexan Tray - Adjustable Knob Hardware &amp; Ribs</b>	(E0950)	08-6000	<b>\$ 219</b>
<b>1/4" Lexan Tray - EZ Lite Channel Hardware &amp; Ribs</b>	(E0950)	08-13000	<b>\$ 183</b>
<b>1/4" Lexan Tray - Hook &amp; Loop Straps</b>	(E0950)	08-16000	<b>\$ 106</b>