



Invacare® Series **Silhouette®** Model #'s SCSC, SCBC
 For Simulation Systems: email to: elyriacustomseating@invacare.com
 fax to: # 800-834-4153
 For 3D Image Systems: E-mail Order Form and Digitized File to:
ElyriaCustomSeating@invacare.com

Company Name:	
Account #:	
Ship To Address:	
City/State/Zip:	
Client Reference:	
Custom Molded CLIP file name:	

Date:	
Contact Name:	
E-mail Address:	
Phone Number:	
Purchase Order #:	

*Note: CLIP file name= 1st Initial, 1st 3 of last name & Account #

Mobility Base Information (optional)			
Brand:		Model:	
Frame Width:		Frame Depth:	

On Chair Information	
<input type="checkbox"/>	Invacare to mount cushions to NEW Invacare Chair

When ordering email this order form and Invacare Chair order form to: elyriacustomseating@invacare.com or provide quote # below:

Invacare quote # : _____

Sensor Serial Numbers:	
four digit # located on rear of seat and back of back	
Seat: _____	Back: _____

<input type="checkbox"/>	Freedom Designs to mount cushions to NEW Freedom Designs Chair
--------------------------	--

If a Freedom Designs chair quote is to be converted please call 1-800-331-8551. Cushions will be drop shipped to Freedom Designs.

Freedom Designs quote / PO # : _____

<input type="checkbox"/>	Motion Concepts to mount cushions to NEW Motion Concepts Chair
--------------------------	--

If a Motion Concepts chair quote is to be converted please call 1-888-433-6818. Cushions will be drop shipped to Motion Concepts.

Motion Concepts quote / PO # : _____

Seat Specifications	
---------------------	--

HPCPS E2609 \$1068.00

Cushion Dimensions

Cushion Width _____" (variable from 8" to 24" wide*)
 *Seat cushion 21" or greater, additional cost applies \$75.00

Cushion Depth _____" (variable from 8" to 24" wide*)
 *Seat cushion 21" or greater, additional cost applies \$75.00

Measured from the back of the digitizer to the place on the client's thigh where you want the cushion to end, add 1"

SO90 Leg Length Discrepancy \$75.00
 Shorten left leg by _____"
 Shorten right leg by _____"
 *Longer Leg will be indicated by the overall cushion length

SO16 Growth Style \$75.00
 Length : _____" + Growth _____" = _____" Total Length

Foam Options

Soft Foam No Charge
 Firm Foam No Charge
 X-tra Firm Foam \$110.00

Mounting

TFNP New PinDot Pan \$215.00
 TFEP Fit Existing Pindot Pan _____" x _____" ... No Charge
 CM06T Wood Mount with t-nuts \$215.00
 Specify Non-std t-nut pattern: _____
 CM06 Wood Mount with out t-nuts \$150.00
 IDB/SM Infinity Drop Base 14 - 16" wide \$240.00
 IDB/MD Infinity Drop Base 17 - 19" wide \$240.00
 IDB/LG Infinity Drop Base 20 - 22" wide \$291.00

*If no mounting selected, cushion will have velcro on bottom

Seat Specifications

Cover Selection

- CMP Polartek
- CMS Lycra
- CMD Darlexx
- DT-2 Standard Dartex
- DT-2R Reversed Dartex
- L-05 Black Rubber Backed Lycra
- L-06 Gray Rubber Backed Lycra
- L-13 Light Blue Rubber Backed Lycra
- L-14 Red Rubber Backed Lycra
- L-15 Teal Rubber Backed Lycra
- L-18 Purple Rubber Backed Lycra
- L-19 Burgundy Rubber Backed Lycra
- L-20 Blue Rubber Backed Lycra
- L-23 Green Rubber Backed Lycra
- L-24 Spacer Mesh

When
opting for
a 2- tone
fabric
selection
indicate
"A" for
contact
surface,
"B" for
non-
contact
surfaces

Mounting Hardware

- HWPMK01 Pan Mounting Kit - 1" \$205.00
- HWPMK78 Pan Mounting Kit - 7/8" \$205.00
- HWWMK01 Pan Mounting Kit - 1" \$205.00
- HWWMK78 Pan Mounting Kit - 7/8" \$205.00
- OCMK01 J & L Brackets - 1" \$175.00
- OCMK01 J & L Brackets - 7/8" \$175.00
- FDI-538 Econo-Eze (includes Seat & Back) \$637.00
- OMIT Omit Hardware

Cover Style Selection

- Draw String No Charge
- 2 Piece Velcro No Charge
- Zipper with Velcro Bottom..... Standard

Additional Cover Selection

- SO25 Additional Cover \$100.00
Quantity _____ Fabric _____
- S040 Inner Water Resistant Cover (Zipper)..... \$210.00

Seat Modifications

- S075 Rail Cut - Railcut Width: _____ "..... \$75.00
- SO77 Undercut _____ "..... \$75.00
- CM20 Strap Slots 1/2" x 1 1/2" std -Qty: _____ ... \$150.00
- VH02 Vent Holes - Qty: \$150.00

Soft Spot / Overlay / Recess Option

- CM09 Softspot - Standard Foam \$165.00
Qty: _____ Depth: _____ (1/2" - 1")
- CM09V Softspot - Visco Foam \$275.00
(Qty: _____ Depth: _____ (1/2" - 1"))
- VF05 Visco Foam Overlay 1/2" \$275.00
- VF01 Visco Foam Overlay 1" \$330.00
- SF05 Standard Softspot Foam Overlay 1/2"..... \$165.00
- SF01 Standard Softspot Foam Overlay 1"..... \$220.00
- CM30 Recess only- Depth: _____ "..... \$105.00
- CM30G Recess with Silicone Pad..... \$275.00
- STAR Recess with Star Airflow Insert
 - 14" Wide \$404.50
 - 16" Wide \$404.50
 - 18" Wide \$539.50
- Recess with Action Gel Insert.....available after 6/15/18
 - ACTF Action Free Float * \$343.50 / \$575.00
 - ACTW Action Waffle * \$343.50 / \$575.00
 - ACTP Action Pebble..... * \$343.50 / \$575.00

*Size determines price >17" higher price applies

Back Specifications

HCPCS E2617 \$1068.00

Cushion Dimensions

- Cushion Width _____" (variable from 8" to 24" wide)
Seat cushion 21" or greater, additional cost applies ... \$75.00
- Cushion Height _____" (variable from 8" to 24" wide)
Seat cushion 21" or greater, additional cost applies \$75.00
Measured from the grid on Back Shape Sensor prior to filling gap.
If no gap, this will be your overall cushion height
- SO17 Fill Gap No Charge
Sensor Grid : _____ + Gap _____" = _____" Total Height

Foam Options

- Soft Foam No Charge
- Firm Foam No Charge
- X-tra Firm Foam \$110.00

Mounting

- SFK30 Kwik Fit Mounting shell w/ hardware..... \$400.00
Wheelchair Frame Width _____ inches (Variable 12" to 24")
KwikFit Mounting for Backs will fit between the wheelchair canes.
Cushion and Shell will be made 2 1/2" smaller than the wheelchair
width indicated. Standard KwikFit Mounted Back will have a 2" curved top.
- TFNP New PinDot Pan \$215.00
- TFEP Fit Existing Pindot Pan _____" x _____" ... No Charge
- CM06T Wood Mount with t-nuts \$215.00
Specify Non-std t-nut pattern: _____ Headrest: _____
- CM06 Wood Mount with out t-nuts \$150.00
- ELITE Matrix Elite Shell..... \$550.00
- ELITE TR Matrix Elite TR Shell..... \$550.00
- ELITE PB Matrix Elite PB Shell..... \$550.00
- TFCS To Fit Custom Shell..... \$75.00

Existing Shell Type: _____
*Elite shells over 20" wide - Additional charges apply

Back Specifications

Cover Selection

- CMP Polartek
- CMS Lycra
- CMD Darlexx
- DT-2 Standard Dartex
- DT-2R Reversed Dartex
- L-05 Black Rubber Backed Lycra
- L-06 Gray Rubber Backed Lycra
- L-13 Light Blue Rubber Backed Lycra
- L-14 Red Rubber Backed Lycra
- L-15 Teal Rubber Backed Lycra
- L-18 Purple Rubber Backed Lycra
- L-19 Burgundy Rubber Backed Lycra
- L-20 Blue Rubber Backed Lycra
- L-23 Green Rubber Backed Lycra
- L-24 Spacer Mesh


When
opting for
a 2-tone
fabric
selection
indicate
"A" for
contact
surface,
"B" for
non-
contact
surfaces

Mounting Hardware

- HWPMK01 Pan Mounting Kit - 1" \$205.00
- HWPMK78 Pan Mounting Kit - 7/8" \$205.00
- HWWMK01 Pan Mounting Kit - 1" \$205.00
- HWWMK78 Pan Mounting Kit - 7/8" \$205.00
- OCMK01 J & L Brackets - 1" \$175.00
- OCMK01 J & L Brackets - 7/8" \$175.00
- FDI-596 Econo-Eze (Back Hardware) \$596.00
- HWABK Adjustable Back Mounting Hardware ... \$260.00
- OMIT Omit Hardware

Back Options

- CM22 Rail Cut - Rail Cut Width: _____ " \$75.00
- Rail Cut to PlnDot Pan..... \$75.00
- Rail Cut to Kwik Fit \$75.00
- VH02 Vent Holes - Qty: _____ \$150.00
- CM30 Recess only - Depth: _____ \$105.00
- S076 I-Cut provide measurements below: \$75.00



Cover Style Selection

- Draw String No Charge
- 2 Piece Velcro No Charge
- Zipper with Velcro Bottom..... Standard

Additional Cover Selection

- SO25 Additional Cover \$100.00
Quantity _____ Fabric _____
- S040 Inner Water Resistant Cover (Zipper) \$210.00

Soft Spot / Overlay / Recess Option

- CM09 Softspot - Standard Foam \$165.00
Qty: _____ Depth: _____ (1/2" - 1")
- CM09V Softspot - Visco Foam \$275.00
(Qty: _____ Depth: _____ (1/2" - 1")
- VF05 Visco Foam Overlay 1/2" \$275.00
- VF01 Visco Foam Overlay 1" \$330.00
- SF05 Standard Softspot Foam Overlay 1/2"..... \$165.00
- SF01 Standard Softspot Foam Overlay 1"..... \$220.00
- CM30 Recess only- Depth: _____ " \$105.00
- CM30G Recess with Silicone Pad..... \$275.00
- Recess with Action Gel Insert.....available after 6/15/18
 - ACTF Action Free Float * \$343.50 / \$575.00
 - ACTW Action Waffle * \$343.50 / \$575.00
 - ACTP Action Pebble..... * \$343.50 / \$575.00

*Size determines price >17" higher price applies

- CM39 Laminated Reinforcement*..... \$105.00
*Available with Wood Mount w/ t-nuts only
- CM36 Extended Lateral Right* _____ " \$105.00
- CM36 Extended Lateral Left* _____ "..... \$105.00
*Cushions 7"-10"(10" Max) thick will incur a \$250 upcharge

Monogram

- \$65.00
- Style Block Cursive
- Thread Color: White Blue
- Pink Yellow
- Teal Red

Name: _____
(Name will be centered top 3/4" front of back cushion)

Accessories

- | | |
|--|---|
| <ul style="list-style-type: none"> <input type="checkbox"/> AR20 Desk Length Armrest (Pair) \$250.00 <input type="checkbox"/> AR24 Full Length Armrest (Pair) \$250.00 <input type="checkbox"/> LB07 Small Footrest (Pair) \$235.00 <input type="checkbox"/> LB08 Large Footrest (Pair) \$235.00 <input type="checkbox"/> LT06 Legrest Tops (Pair) \$235.00 <input type="checkbox"/> SB05 ContourU Pelvic Strap \$45.00 <input type="checkbox"/> PDSL M Swing Away Lateral Supports Medium..... \$290.00 <input type="checkbox"/> PDSL L Swing Away Lateral Supports Large..... \$290.00 | <ul style="list-style-type: none"> <input type="checkbox"/> TR40L Tray Receptacle - Left \$140.00 <input type="checkbox"/> TR40R Tray Receptacle - Right \$140.00 <input type="checkbox"/> GL01 Regular Links \$125.00 <input type="checkbox"/> GL03 Extra Strong Links..... \$140.00 <input type="checkbox"/> GL10 Dial Links (seat angle must be > 105)..... \$250.00 <input type="checkbox"/> Headrest Adapter Plate* \$45.00 <input type="checkbox"/> SLK Fixed Lateral Supports R/L for KwikFit shell.... \$125.00 |
|--|---|

* Must be ordered through Freedom Designs 800-331-8551