



Our kids are our passion.™

The **Easy Conforming Back**

Functionality and Comfort



"Easy Conforming"

Standard

The Easy Conforming Back has a contoured design with ultimate adjustability and quick release. This back hugs the contour of the body with curved panels to provide a custom fit.

The Split Wing modification (MOD2-QQ) has two separate panels that allow for easy individual adjustments at the lateral and hip areas to provide a truly custom fit. The Easy Conforming Back is available in standard and extended sizes with a choice of hardware.

HCPCS Code E2617



Extended



Split Wing Mod

See reverse side for
Part Numbers and
Measuring Guide

Freedom Designs, Inc.

2241 Madera Rd., Simi Valley, CA 93065

Phone: 800-331-8551 Fax: 888-582-1509


www.freedomdesigns.com




Our kids are our passion.™

The **Easy Conforming Back**

Part Number and Measuring Guide



Standard



Extended

Standard Sizes - Frame Widths: (5" narrower than chair width)

• 11/12	A=7	B=8	C=13	D=3	E=2	F=7
• 13/14	A=9.5	B=10.5	C=15	D=4	E=3	F=9
• 15/16	A=10.5	B=11.5	C=17	D=4	E=3	F=11
• 17/18	A=11	B=12.5	C=19	D=4	E=3	F=13
• 19/20	A=12	B=13.5	C=21	D=4	E=3.5	F=15
• 21/22*	A=13	B=14.5	C=23	D=4	E=3.5	F=17
• 23/24*	A=14	B=15.5	C=25	D=4	E=3.5	F=19
• 25/26*	A=14	B=15.5	C=27	D=4	E=3.5	F=21

Standard Sizes - Frame Widths: (5" narrower than chair width)

• 11/12	A=5	B=6	C=13	D=3	E=2	F=7	G=8
• 13/14	A=6	B=7	C=15	D=4	E=3	F=9	G=10.5
• 15/16	A=7	B=7.5	C=17	D=4	E=3	F=11	G=11.5
• 17/18	A=7.5	B=8	C=19	D=4	E=3	F=13	G=12.5
• 19/20	A=8	B=9	C=21	D=4	E=3.5	F=15	G=13.5
• 21/22*	A=9	B=10.5	C=23	D=4	E=3.5	F=17	G=14.5
• 23/24*	A=9	B=10.5	C=25	D=4	E=3.5	F=19	G=14.5
• 25/26*	A=9	B=10.5	C=27	D=4	E=3.5	F=21	G=14.5

Description:

- W/4 J & L's & 4 glide locks
- W/2 J & L's, glide lock, & seat/back hinge
- 21" wide W/4 J & L's & 4 glide locks
- 21" wide W/2 J & L's, glide locks & seat/back hinge
- W/4 J & L's, glide locks & 2 split collars
- 21" W/4 J & L's, 2 glide locks & 2 split collars
- 22-26" wide Easy Conforming heavy duty back
- Back with NO hardware
- W/4 Easy Clamps
- 21" wide W/4 Easy Clamps
- 22-26" wide W/4 Easy Clamps
- W/2 Easy Clamps & seat/back hinge
- 21" wide W/2 Easy Clamps & seat/back hinge
- 22-26" wide W/2 Easy Clamps & seat/back hinge
- Free Form: SOFT or MEDIUM
- Split outer wing - standard is cut in half
- All pads covered separately

Standard

Part Numbers:

102-12000
EZ CONF2
102-12000-21W
EZ CONF2-21W
102-17000
102-17000-21W
102-17000-22/26W
102-12000-NH
102-20000
102-20000-21W
102-20000-22/26W
EZCONF20
EZCONF20-21W
EZCONF20-22/26W
MOD-2S
MOD-2QQ
MOD-2SS

Extended

Part Numbers:

102-13000
EZ CONF-EXT2
102-13000-21W
EZCONF-EXT2-21W
102-18000
102-18000-21W
102-18000-22/26W
102-13000-NH
102-30000
102-30000-21W
102-30000-22/26W
EZCONF30
EZCONF30-21W
EZCONF30-22/26W
MOD-2S
MOD-2QQ
MOD-2SS

*Wide sizes only available for part numbers ending in "W" (eg. 102-13000-21W)

Freedom Designs, Inc.

2241 Madera Rd., Simi Valley, CA 93065

Phone: 800-331-8551 Fax: 888-582-1509

www.freedomdesigns.com



Our kids are our passion™

FSECB 22012