



INVACARE® MATRX® CUSHIONS
FREEDOM DESIGNS ON-CHAIR ORDER FORM
Effective November 1st, 2015



Date of Order: _____	Ship to Address: _____
Dealer Name: _____	_____
Dealer Account #: _____	City: _____
PO #: _____	State: _____ Zip: _____
Tag: _____	Comments: _____
Purchasing Contact: _____	_____
Phone: _____ Fax: _____	_____

Pediatric Cushions

► Page

Invacare® Matrx® Kid*ab*ra Cushion/Extra Outer Cover	2
Invacare® Matrx® Kid*ab*ra®Kid*ab*ra Drop Base/Extra Outer Cover	2
Invacare® Matrx® PS Cushion/Extra Outer Cover	5

Standard Cushions

► Page

Invacare® Matrx® Vi Cushion/Extra Outer Cover	3
Invacare® Matrx® Vi Drop Base/Extra Outer Cover	4
Invacare® Matrx® Vi Cushion Rigidizer	4
Invacare® Matrx® PS Cushion/Extra Outer Cover	6
Invacare® Matrx® Libra Cushion/Extra Outer Cover	7


Modifications & Accessories

► Page

Invacare® Matrx® Libra Fit Kit Components	8
Cushion Modifications	9
Cushion Accessories	9
Cushion Hardware	9

INVACARE® MATRX® CUSHIONS – FREEDOM DESIGNS ON-CHAIR ORDER FORM

MATRX® KID*AB*RA CUSHION: Pediatric – Skin Protection & Positioning Cushion



NOTE: Actual cushion width is 1/2" narrower than stated size and cushion depth is 1" longer than stated depth due to rear cushion radius.

- Cushion modifications are available and can be specified on Page 9

- Weight capacity up to 150 lbs
- Comes standard with inner cover and Kid*ab*ra outer covers



Kid*ab*ra
Outer Cover

Kid*ab*ra Cushion

► HCPCS Code: E2607



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ KA1212-Vi	11", 12"	12", 13"	\$425.00
☐ KA1214-Vi	11", 12"	14", 15"	\$425.00
☐ KA1412-Vi	13", 14"	12", 13"	\$425.00
☐ KA1414-Vi	13", 14"	14", 15"	\$425.00

Kid*ab*ra Extra Outer Cushion Cover (Refer to cushion sizes above for proper fit)

► HCPCS Code: E2619



MODEL	FITS WIDTH	FITS DEPTH	MSRP
KC	REFER TO CUSHION SIZES ABOVE		\$95.00
SPECIFY:	WIDTH _____	DEPTH _____	☐ QTY _____

MATRX® KID*AB*RA DROP BASE: Pediatric – Skin Protection & Positioning Cushion

NOTE: Base can only be ordered with corresponding Kid*ab*ra Cushion

- Must specify drop hook size: 1" diameter and 1" or 2" drop

- Weight capacity up to 150 lbs
- Includes dropbase and mounting hardware

Kid*ab*ra Drop Base Cushion

► HCPCS Code: E2231



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ KB1212-Vi	12"	12", 13"	\$200.00
☐ KB1214-Vi	12"	14", 15"	\$200.00
☐ KB1312-Vi	13"	12", 13"	\$200.00
☐ KB1314-Vi	13"	14", 15"	\$200.00
☐ KB1412-Vi	14"	12", 13"	\$200.00
☐ KB1414-Vi	14"	14", 15"	\$200.00

DROP HOOKS MOUNTING SPECIFICATIONS:

Diameter:	☐ 1"		
Front Seat:	☐ FLUSH	☐ 1"	☐ 2"
Rear Seat:	☐ FLUSH	☐ 1"	☐ 2"

Kid*ab*ra Drop Base Extra Outer Cover (Refer to drop base sizes above for proper fit)


► HCPCS Code: E2619



MODEL	FITS WIDTH	FITS DEPTH	MSRP
KBC	REFER TO CUSHION SIZES ABOVE		\$90.00
SPECIFY:	WIDTH _____	DEPTH _____	☐ QTY _____

INVACARE® MATRIX® CUSHIONS – FREEDOM DESIGNS ON-CHAIR ORDER FORM

MATRIX® Vi CUSHION: Standard – Skin Protection & Positioning Cushion



NOTE: Actual cushion width is 1/2" narrower than stated size and cushion depth is 1 1/2" longer than stated depth due to rear cushion radius.

- Cushion modifications are available and can be specified on Page 9

- Weight capacity up to 300 lbs
- Comes standard with inner cover and outer covers made of Meshtex material



Meshtex
Outer Cover

- Startex fabric or smooth fabric side up, Onyx or Fusion Covers are available as an upcharge

Vi Cushion

► HCPCS Code: E2607



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ MA1415-Vi	13", 14"	15", 16"	\$425.00
☐ MA1417-Vi	13", 14"	17", 18"	\$425.00
☐ MA1616-Vi	15", 16"	16", 17"	\$425.00
☐ MA1618-Vi	15", 16"	18", 19"	\$425.00
☐ MA1620-Vi	15", 16"	20", 21"	\$425.00
☐ MA1716-Vi	16", 17"	16", 17"	\$425.00
☐ MA1718-Vi	16", 17"	18", 19"	\$425.00
☐ MA1720-Vi	16", 17"	20", 21"	\$425.00
☐ MA1816-Vi	17", 18"	16", 17"	\$425.00
☐ MA1818-Vi	17", 18"	18", 19"	\$425.00
☐ MA1820-Vi	17", 18"	20", 21"	\$425.00
☐ MA2016-Vi	19", 20"	16", 17"	\$425.00
☐ MA2018-Vi	19", 20"	18", 19"	\$425.00
☐ MA2020-Vi	19", 20"	20", 21"	\$425.00

► HCPCS Code: E2608

☐ MA2218-Vi	21", 22"	18", 19"	\$475.00
☐ MA2220-Vi	21", 22"	20", 21"	\$475.00

COVER OPTIONS:

ADD

SPECIFY:	☐ STARTEX-F (Startex Cover - Fabric)	☐ STARTEX-S (Startex Cover-Smooth)	\$50.00
	☐ OCU-MA (Onyx Cover)	☐ FSU-MA(Fusion Cover)	\$100.00

Vi Extra Outer Cushion Cover (Refer to cushion sizes above for proper fit)

► HCPCS Code: E2619




MODEL	FITS WIDTH	FITS DEPTH	MSRP
C	REFER TO CUSHION SIZES ABOVE		\$100.00
SPECIFY:	WIDTH_____	DEPTH_____	☐ QTY_____

*C2218 & C2220 are \$110.00

INVACARE® MATRX® CUSHIONS – FREEDOM DESIGNS ON-CHAIR ORDER FORM

MATRX® Vi DROP BASE: Standard – Skin Protection & Positioning Cushion



NOTE: Base can only be ordered with corresponding Matrix Cushion

- Must specify drop hook size: 1" diameter and flush 1" or 2" drop

- Weight capacity up to 300 lbs
- Includes dropbase and mounting hardware

Vi Drop Base

► HCPCS Code: E2231



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ MB1415-Vi	14"	15", 16"	\$200.00
☐ MB1417-Vi	14"	17", 18"	\$200.00
☐ MB1516-Vi	15"	16", 17"	\$200.00
☐ MB1518-Vi	15"	18", 19"	\$200.00
☐ MB1616-Vi	16"	16", 17"	\$200.00
☐ MB1618-Vi	16"	18", 19"	\$200.00
☐ MB1716-Vi	17"	16", 17"	\$200.00
☐ MB1718-Vi	17"	18", 19"	\$200.00
☐ MB1816-Vi	18"	16", 17"	\$200.00
☐ MB1818-Vi	18"	18", 19"	\$200.00
☐ MB1916-Vi	19"	16", 17"	\$200.00
☐ MB1918-Vi	19"	18", 19"	\$200.00
☐ MB2016-Vi	20"	16", 17"	\$200.00
☐ MB2018-Vi	20"	18", 19"	\$200.00

DROP HOOKS MOUNTING SPECIFICATIONS:

Diameter: ☐ 1"

Front Seat: ☐ FLUSH ☐ 1" ☐ 2"

Rear Seat: ☐ FLUSH ☐ 1" ☐ 2"

Vi Drop Base Extra Outer Cover (Refer to drop base sizes above for proper fit)

► HCPCS Code: E2619



MODEL	FITS WIDTH	FITS DEPTH	MSRP
BC	REFER TO DROP BASE SIZES ABOVE		\$90.00*
SPECIFY:	WIDTH _____	DEPTH _____	☐ QTY _____

MATRX® Vi CUSHION RIGIDIZER

NOTE: Lightweight plastic seat insert fits between seat rails for level support and eliminates hammock effect of seat upholstery.

- Recommended for chairs with sling seat upholstery.

Vi Cushion Rigidizer


► HCPCS Code: E0992



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ CR14S	12"-14"	12"-17"	\$75.00
☐ CR16S	15", 16"	15", 16"	\$75.00
☐ CR16M	15", 16"	17", 18"	\$75.00
☐ CR16L	15", 16"	19", 20"	\$75.00
☐ CR18S	17", 18"	15", 16"	\$75.00
☐ CR18M	17", 18"	17", 18"	\$75.00
☐ CR18L	17", 18"	19", 20"	\$75.00
☐ CR20S	19", 20"	15", 16"	\$75.00
☐ CR20M	19", 20"	17", 18"	\$75.00
☐ CR20L	19", 20"	19", 20"	\$75.00
☐ CR22M	21", 22"	17", 18"	\$95.00
☐ CR22L	21", 22"	19", 20"	\$95.00

INVACARE® MATRIX® CUSHIONS – FREEDOM DESIGNS ON-CHAIR ORDER FORM

MATRIX® PS CUSHION: Pediatric Positioning Cushion



NOTE: Actual cushion width is 1/2" narrower than stated size and cushion depth is 1" longer than stated depth due to rear cushion radius.

- Cushion modifications are available and can be specified on Page 9 (cover will be non-reversible)

- Weight capacity up to 150 lbs
- Comes standard with inner cover and reversible Startex® outer covers



Startex
Outer Cover

- No charge outer cover options: Kid*ab*ra, Pink Lollipop or Blizzard Blue (all non-reversible)



Kid*ab*ra
Outer Cover



Pink Lollipop
Outer Cover



Blue Blizzard
Outer Cover

PS Cushion

► HCPCS Code: E2605



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ PS1010	9", 10"	10", 11"	\$345.00
☐ PS1012	9", 10"	12", 13"	\$345.00
☐ PS1014	9", 10"	14", 15"	\$345.00
☐ PS1212	11", 12"	12", 13"	\$345.00
☐ PS1214	11", 12"	14", 15"	\$345.00
☐ PS1216	11", 12"	16", 17"	\$345.00
☐ PS1412	13", 14"	12", 13"	\$345.00
☐ PS1414	13", 14"	14", 15"	\$345.00
☐ PS1416	13", 14"	16", 17"	\$345.00

COVER OPTIONS:

SPECIFY: ☐ PS (Startex) ☐ PS-K (Kid*abra) ☐ PS-P (Pink Lollipop) ☐ PS-B (Blizzard Blue)

PS Extra Outer Cover (Refer to cushion sizes above for proper fit)

► HCPCS Code: E2619



MODEL	FITS WIDTH	FITS DEPTH	MSRP
PSC	REFER TO CUSHION SIZES ABOVE		\$85.00*


SPECIFY: WIDTH _____ DEPTH _____ ☐ QTY _____

COVER OPTIONS:

SPECIFY: ☐ PSC (Startex) ☐ PSC-K (Kid*abra)
☐ PSC-P (Pink Lollipop) ☐ PSC-B (Blizzard Blue)

INVACARE® MATRIX® CUSHIONS – FREEDOM DESIGNS ON-CHAIR ORDER FORM

MATRIX® PS CUSHION: Standard – Positioning Cushion



NOTE: Actual cushion width is 1/2" narrower than stated size and cushion depth is 1" longer than stated depth due to rear cushion radius.

- Cushion modifications are available and can be specified on Page 9 (cover will be non-reversible)

- Weight capacity up to 300 lbs
- Comes standard with inner cover and reversible Startex® outer cover



Startex
Outer Cover

PS Cushion

► HCPCS Code: E2605



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ PS1616	15", 16"	16", 17"	\$345.00
☐ PS1618	15", 16"	18", 19"	\$345.00
☐ PS1620	15", 16"	20", 21"	\$345.00
☐ PS1816	17", 18"	16", 17"	\$345.00
☐ PS1818	17", 18"	18", 19"	\$345.00
☐ PS1820	17", 18"	20", 21"	\$345.00
☐ PS2016	19", 20"	16", 17"	\$345.00
☐ PS2018	19", 20"	18", 19"	\$345.00
☐ PS2020	19", 20"	20", 21"	\$345.00

PS Extra Outer Cover (Refer to cushion sizes above for proper fit)


► HCPCS Code: E2619



MODEL	FITS WIDTH	FITS DEPTH	MSRP
PSC	REFER TO CUSHION SIZES ABOVE		\$85.00*
SPECIFY:	WIDTH_____	DEPTH_____	☐ QTY_____

INVACARE® MATRX® CUSHIONS – FREEDOM DESIGNS ON-CHAIR ORDER FORM

MATRX® LIBRA CUSHION: Standard – Adjustable Skin Protection and Positioning Cushion



NOTE: Actual cushion width is as stated and cushion depth is 1½" longer than stated depth due to rear cushion radius.

- Weight capacity up to 300 lbs
- Comes standard with foam base, fluid sac, platilon inner liner and Startex® outer cover



Startex
Outer Cover

- Spacetex is a no charge outer cover option; Infection Control additional \$50 (all non-reversible)



Spacetex
Outer Cover



Infection Control
Outer Cover

Libra Cushion

► HCPCS Code: E2624



Startex Outer Cover



Spacetex Outer Cover



Infection Control Outer Cover



MODEL	FITS WIDTH	FITS DEPTH	MSRP
☐ LC1412	14"	12", 13"	\$485.00
☐ LC1414	14"	14", 15"	\$485.00
☐ LC1416	14"	16", 17"	\$485.00
☐ LC1516	15"	16", 17"	\$485.00
☐ LC1518	15"	18", 19"	\$485.00
☐ LC1520	15"	20", 21"	\$485.00
☐ LC1616	16"	16", 17"	\$485.00
☐ LC1618	16"	18", 19"	\$485.00
☐ LC1620	16"	20", 21"	\$485.00
☐ LC1622	16"	22", 23"	\$595.00
☐ LC1716	17"	16", 17"	\$485.00
☐ LC1718	17"	18", 19"	\$485.00
☐ LC1720	17"	20", 21"	\$485.00
☐ LC1722	17"	22", 23"	\$595.00
☐ LC1816	18"	16", 17"	\$485.00
☐ LC1818	18"	18", 19"	\$485.00
☐ LC1820	18"	20", 21"	\$485.00
☐ LC1822	18"	22", 23"	\$595.00
☐ LC1916	19"	16", 17"	\$485.00
☐ LC1918	19"	18", 19"	\$485.00
☐ LC1920	19"	20", 21"	\$485.00
☐ LC1922	19"	22", 23"	\$595.00
☐ LC2016	20"	16", 17"	\$485.00
☐ LC2018	20"	18", 19"	\$485.00
☐ LC2020	20"	20", 21"	\$485.00
☐ LC2022	20"	22", 23"	\$595.00

COVER OPTIONS:

SPECIFY: ☐ LC (Startex) ☐ LC-SP (Spacetex) ☐ LC-IC (Infection Control)*

*Additional \$50 for Infection Control Cover

Libra Extra Outer Cover (Refer to cushion sizes above for proper fit)


► HCPCS Code: E2619

MODEL	FITS WIDTH	FITS DEPTH	MSRP
LCC	REFER TO CUSHION SIZES ABOVE		\$100.00
LSC	REFER TO CUSHION SIZES ABOVE		\$125.00
LIC	REFER TO CUSHION SIZES ABOVE		\$150.00

SPECIFY: WIDTH _____ DEPTH _____ ☐ QTY _____


COVER OPTIONS:

SPECIFY: ☐ LCC (Startex) ☐ LSC (Spacetex) ☐ LIC (Infection Control)

MATRX® LIBRA FIT KIT COMPONENTS & OTHER COMPONENTS:

Libra Components

QUANTITY	MODEL	COMPONENT	MSRP
☐ QTY ____	LAB	Libra Abductor <i>1" Abductor built-up - Top mounted</i>	\$60.00
☐ QTY ____	FTSW	Anterior Seat Wedges <i>1" Anterior wedge increases pre-ischial contour - Underside mounted</i>	\$90.00
☐ QTY ____	FTHW	Hemi Wedges <i>1" Anterior thigh support - Underside mounted</i>	\$50.00
☐ QTY ____	FTTW-LFRR	Tapered Wedge LFRR <i>¾" Tapered wedge for users's Left front or Right rear - Underside mounted</i>	\$50.00
☐ QTY ____	FTTW-RFLR	Tapered Wedge LFLR <i>¾" Tapered wedge for users's Right front or Left rear - Underside mounted</i>	\$50.00
☐ QTY ____	LFS	Libra Fluid Supplement <i>May be used to accommodate mild pelvic obliquity/additional pressure reduction</i>	\$50.00
☐ QTY ____	FTLW	Lateral Wedge <i>1" Lateral wedge increases hip/thigh loading, - Underside mounted</i>	\$65.00
☐ QTY ____	CLCFP-M	Complete Libra Cushion Fluid Sac <i>18" Wde and below</i>	\$315.00
SPECIFY CUSHION SIZE:		WIDTH ____ DEPTH ____	
☐ QTY ____	CLCFP-L	Complete Libra Cushion Fluid Sac <i>19" Wde and above</i>	\$315.00
SPECIFY CUSHION SIZE:		WIDTH ____ DEPTH ____	
☐ QTY ____	LCPL-WWDD	Platilon Inner Liner	\$35.00
SPECIFY CUSHION SIZE:		WIDTH ____ DEPTH ____	


CUSHION MODIFICATIONS:



NOTE: Rear Cushion radius included on all PS, VI & Libra Cushions.

Notes:

Modifications (Additional one week lead time)



QUANTITY	MODEL	COMPONENT	MSRP
☐ QTY ____	SDR	Seat Depth Reduction Less: ☐ 1" ☐ 2"* ☐ 3"* ☐ 4"* Specify ____ <i>*Not available on Libra Cushions</i>	\$100.00
☐ QTY ____	SWR	Seat Width Reduction Less: ☐ 1" ☐ 2"* ☐ 3"* ☐ 4"* Specify ____ <i>*Not available on Libra Cushions</i>	\$100.00
☐ QTY ____	LLDR	Leg Length Discrepancy - Right Less: ☐ 1" ☐ 2" Specify ____	\$100.00
☐ QTY ____	LLDL	Leg Length Discrepancy - Left Less: ☐ 1" ☐ 2" Specify ____	\$100.00
☐ QTY ____	ECM	Outer Cover Modification <i>Accommodates LLDL/LLDR/SDR/SWR</i>	\$50.00
☐ QTY ____	MAC	Miscellaneous Modification <i>Please use a diagram above to indicate required modifications- Requires approval</i>	\$100.00

ACCESSORIES & HARDWARE:

QUANTITY	MODEL	COMPONENT	MSRP
☐ QTY ____	CIC	Extra Inner Cover <i>1" Pelvic Obliquity built-up- Top mounted</i>	\$75.00
☐ QTY ____	NSF	Non-Standard Fabric Selection <i>Available for Cushions/bases; excluding Libra</i>	\$75.00
☐ QTY ____	POK	Pelvic Obliquity Kit <i>Accommodates ½" and " Obliquities</i>	\$50.00
☐ QTY ____	BDH0	Base Drop Hooks (pair) <i>Flushmount (0) drop for 1" Diameter</i>	\$27.00
☐ QTY ____	BDH1-1	Base Drop Hooks (pair) <i>1" Drop for 1" Diameter</i>	\$27.00
☐ QTY ____	BDH2-1	Base Drop Hooks (pair) <i>1" Drop for 2" Diameter</i>	\$27.00
☐ QTY ____	BDC	Dura Clip (pair) <i>Top lock retractable hardware for use with Drop Hooks</i>	\$60.00
☐ QTY ____	BFC	FixClip (pair) <i>Top lock hardware for use with Drop Hooks</i>	\$25.00