



Freedom *Select* ANTERIOR TRUNK SUPPORTS



Standard Pull Straps

- 4 Vest Center Panel Styles to Choose From
- Push-Button Swivel Buckles for Easy Release



- Grommet Attachment & Side Release Buckle Upper Straps
- Cam Buckle Lower Strap Attachment




- All Mounting Fasteners Included
- Black Neoprene Standard
- Custom Sizes Available
- Dynamic Center Panel Modification Available



- All Mounting Fasteners Included
- Black Neoprene Standard
- Custom Sizes Available
- Dynamic Center Panel Modification Available



- Upper Straps are Over 34 Inches Long
- Upper Slide Buckle with Adjustable Side-Release Buckle for Lower Strap Adjustment
- Cam Buckle & Triangle Attaching Hardware for All 4 Straps Included



Dual-Adjustment Rear Pull Top Straps

- 4 Vest Center Panel Styles to Choose From
- Side-Release Buckles on Lower Straps

Freedom Select Anterior Trunk Supports

Clear Form

Email Form to
Customer Service



Freedom Designs, Incorporated

2241 N Madera Road, Simi Valley, CA 93065
phone: 800-331-8551 fax: 888-582-1509
e: customerservice@freedomdesigns.com



Yes, you can.

Quote

Order

Buyer: _____

Co./Account: _____

Location: _____

Contact #: _____


P.O. #: _____

Marked For: _____


Fax: _____

Email: _____

Select Standard



Select Narrow



Dimensions


Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SNZ-S	Select Standard Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-1-S	Select <i>Standard</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-M	Select <i>Standard</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-1-M	Select <i>Standard</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-L	Select <i>Standard</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-1-L	Select <i>Standard</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-XL	Select <i>Standard</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6SNZ-1-XL	Select <i>Standard</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

	SFDI-6NNZ-XS	Select <i>Narrow</i> X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-1-XS	Select <i>Narrow</i> Dual Rear-Pull X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-S	Select <i>Narrow</i> Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-1-S	Select <i>Narrow</i> Dual Rear-Pull Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-M	Select <i>Narrow</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-1-M	Select <i>Narrow</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-L	Select <i>Narrow</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-1-L	Select <i>Narrow</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-XL	Select <i>Narrow</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6NNZ-1-XL	Select <i>Narrow</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235


Note:

Freedom Select Anterior Trunk Supports (continued)

Select Standard Zipper



Select Narrow Zipper



Dimensions

Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SZ-S	Select <i>Standard Zipper</i> Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 167
	SFDI-6SZ-1-S	Select <i>Standard Zipper</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 167
	SFDI-6SZ-M	Select <i>Standard Zipper</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 167
	SFDI-6SZ-1-M	Select <i>Standard Zipper</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 167
	SFDI-6SZ-L	Select <i>Standard Zipper</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SZ-1-L	Select <i>Standard Zipper</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SZ-XL	Select <i>Standard Zipper</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 178
	SFDI-6SZ-1-XL	Select <i>Standard Zipper</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 178
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

	SFDI-6NZ-M	Select <i>Narrow Zipper</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 167
	SFDI-6NZ-1-M	Select <i>Narrow Zipper</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 167
	SFDI-6NZ-L	Select <i>Narrow Zipper</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 173
	SFDI-6NZ-1-L	Select <i>Narrow Zipper</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 173
	SFDI-6NZ-XL	Select <i>Narrow Zipper</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 178
	SFDI-6NZ-1-XL	Select <i>Narrow Zipper</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 178
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

Note: