

Simply Toddler Seating Order Form

Clear Form

Email Form to Customer Service



Freedom Designs, Incorporated
2241 N Madera Road, Simi Valley, CA 93065
phone: 800-331-8551 fax: 888-582-1509
e: customerservice@freedomdesigns.com



Quote

Order

Yes, you can:


Seating Specifications

- A Back Height: _____
B Lateral Height: _____
C Lateral Width: _____
D Hip Width (rear): _____
E Set Up Seat Depth: _____
F Hip Width (front): _____
G "E" + 3" for Growth: _____
H Seat Width: _____
I Back Width: _____

Important Note: To ensure quote accuracy, seat and back width, depth and height dimensions must be provided. Seats and backs over 18" in width, depth or height are subject to oversize charges.

Finished Dimensions

Standard cut is Boltaron with foam rolled over edge.



Wheelchair Specifications

- Model: _____
Width: _____
Depth: _____

Buyer: _____
Co./Account: _____
Location: _____
Contact #: _____
P.O.#: _____
Marked For: _____
Fax: _____
Email: _____

Freedom Designs Upholstery Selections

1. Choose Contact / Non-Contact
2. Choose Upholstery
3. Choose Component

	Seat Base	Back Base	Headrest Pad	Lateral Pads	Hip & Knee Pads
Body Contact Surfaces					
Primary Upholstery Option:					
Number: _____ Choose: _____					
Secondary Upholstery Option:					
Number: _____ Choose: _____					
Non-Contact Surfaces					
Primary Upholstery Option:					
Number: _____ Choose: _____					
Secondary Upholstery Option:					
Number: _____ Choose: _____					
Monogram Option	When ordering choose either Monogram or Embroidery				
Name: _____					
Script _____ Block _____					
Thread Color: _____					
Embroidery Option					
Design: _____					


Notes:

Simply Toddler Seating Order Form

Step 2. Choose Customization to EB1 Seat (continued)

EB1CS1CR Asymmetrical - Long Right


\$ 106



A _____
B _____
C _____
D _____

EB1CS1CL Asymmetrical - Long Left

\$ 106



A _____
B _____
C _____
D _____

EB1CS1MLB 1/2" Maple Laminate Base

N/C

EB1CS1WT Base with T-Nuts for Hardware
(no tracks included)

N/C

EB1CS1G Cushion Separate From Base

\$127

EB1CS1I Incontinence Resistance Liner

\$ 37

EB1CS1-
ROHOSEAT
EB1CS1L

Freedom Seat **ROHO** Insert
(see page 16 to specify order)
Grid Cut Foam Modification

\$924

\$403

EB1CS1FF

Removable Incontinence Slip Cover

\$247

Notes:

Step 3. Choose Foam

.5" 4# base w/.5" soft Sunmate top layer-Std.

OR

Custom Foam / Omit Standard Foam (see chart below)

Base layer to be: _____ Middle layer to be: _____ Top layer to be: _____

Seat Foam Chart

3# Foam	\$ 55 per 1"	Firm Sunmate	\$ 99 per 1"	Pudgee	\$252 per 1"
4# Foam	\$ 55 per 1"	X-Soft T-Foam	\$ 107 per 1"	Waffle	\$325 per .5"
5# Foam	\$ 57 per 1"	Soft T-Foam	\$ 107 per 1"	Waffle (over 18x16)	\$653 per .5"
Soft Sunmate	\$ 99 per 1"	Medium T-Foam	\$ 107 per 1"	Free Float	\$409 per .5"
Medium Sunmate	\$ 99 per 1"	Firm T-Foam	\$ 107 per 1"	FreeFloat (over 18x16)	\$815 per .5"

Notes:

Simply Toddler Seating Order Form

Solid Back Insert



\$1329.00
(E2617)

Toddler Simply Seating Back Construction Includes:

- 1/4" thick Boltaron base material up to 14" x 16" size ONLY • Cut between the rails unless otherwise specified • Base material available flat only • Full standard "T" nut pattern • Tracks for hardware and accessory attachment • Headrest adapt-a-plate • Anterior trunk support strap risers • 1/2" 4 lb. base foam with 1/2" soft Sunmate foam top layer is standard • 3 sided zipper closed cover •
- MAXIMUM WEIGHT CAPACITY 75 lbs. •




Step 1. Choose EB2-Kit (Back)

Note: Full length dual slot tracks not on backs narrower than 10 wide.

A	EB2-TDLR-Kit	2-Point Mounting Hdwe & Seat to Back Hinge	B	EB2-Knob-TDLR-Kit	2-Point Mounting Hdwe & Seat to Back Hinge
2 - 590A	Toddler Fast Clamp Adjustable Hardware		2 - 462G	Mini Econo Eze Adjustable Knob Hardware	
2 - SEZFC	Select EZ Frame Clamps		2 - 437E	Back Mount Bracket Mini Knob Hardware	
Note:	Choose 1 Strap Riser (pair) from C		Note:	Choose 1 Strap Riser (pair) from C	
C					
Strap Riser Choices (choose 1 pair)					
2 - MOD-6A 1-1/2" Wide Strap Opening, Flat					
2 - MOD-6A-1 1" Wide Strap Opening, Flat					
2 - MOD-6B 1-1/2" Wide Strap Opening, 25° Rear Bend					
2 - MOD-6B-1 1" Wide Strap Opening, 25° Rear Bend					

Notes:

Step 2. Choose Customization to EB2 Back


EB2CB2H	I-Cut Modification	\$155	EB2CB2IA	Curved Foam Modification	\$140
EB2CB2J	T-Cut Modification	\$155			
	A _____ B _____ C _____ D _____ E _____		E _____ F _____		

Notes:

Simply Toddler Seating Order Form

Step 2. Choose Customization to EB2 Back (continued)

EB2CB2L Biangular Cut Modification / Dual Cover \$366
EB2CB2LA Biangular Cut Modification / Single Cover \$366



E _____
F _____
G _____
H _____

EB2CB2MLB 1/2" Maple Laminate Base N/C
EB2CB2WT Base with T-Nuts for Hardware (no tracks included) N/C
EB2CB2K Grid Cut Foam Modification \$403
EB2CB2-ROHOBAG Freedom Back, **ROHO** Insert (see page 17 to specify) \$924
EB2CB2S Free Form Modification Soft/Med Sunmate - Box Kit (specify below) \$313
EB2CB2SBAG Free Form Modification Soft/Med Sunmate - Bag Kit (specify below) \$313
EB2CB2S1 Liquid Soft/Med Sunmate Only - Box Kit (specify below) \$293
EB2CB2S1-BAG Liquid Soft/Med Sunmate Only - Bag Kit (specify below) \$293
Soft Sunmate Medium Sunmate

Notes:

Step 3. Choose Foam

.5" 4# base w/.5" soft Sunmate top layer-Std

OR

Custom Foam / Omit Standard Foam (see chart below)

Base layer to be: _____ Middle layer to be: _____ Top layer to be: _____


Back Foam Chart

3# Foam	\$ 55 per 1"	Firm Sunmate	\$127 per 1"	Pudgee	\$273 per 1"
4# Foam	\$ 55 per 1"	X-Soft T-Foam	\$137 per 1"	Waffle	\$325 per .5"
5# Foam	\$ 75 per 1"	Soft T-Foam	\$137 per 1"	Waffle (over 18x16)	\$653 per .5"
Soft Sunmate	\$127 per 1"	Medium T-Foam	\$137 per 1"	Free Float	\$409 per .5"
Medium Sunmate	\$127 per 1"	Firm T-Foam	\$137 per 1"	FreeFloat (over 18x16)	\$815 per .5"

Notes:

Simply Toddler Seating Order Form

Medial Pelvic Support Package – Each

Part Number	Description	HCCPS	Price	Shape	Standard Pad Sizing	Custom Pad Size
EE3A-TDLR	Depth Adjustable - Removable Toddler Hdwe.– Choice of Pad	(E0957)	\$287	Wedge Oval Round	Small – 3L x 3H x 3-2 Wedge Medium – 4L x 3H x 3-2 Wedge Large – 5L x 4H x 4-2 Wedge X-Large – 6L x 5H x 5-3 Wedge	 A _____ B _____ C _____ D _____


Notes:

Medial Pelvic Support Components – Each

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
MOD-3M	Protective Boot For Hardware	(K0108)	\$ 72 ea.	438T	Removable, Adjustable Toddler Medial Hdwe	(E0128)	\$ 80 ea.
MOD-3E	Medial Pelvic Support Pad Only	(E0957)	\$206 ea.				


Notes:

Lateral Pelvic Support Package – Pair – Mounted at Hips

Part Number	Description	HCCPS	Price	Location	Toddler Pad Sizing	Custom Pad Size
EE4A-TDLR	Fixed, Flush Toddler Hdwe & Pads	(E0956)	\$287	Hip	Set Up Seat Depth x 3H Pad Profile Standard (1" thick) Thin (1/2" thick)	 A _____ B _____ C _____ D _____

Notes:

Lateral Pelvic Support Package – Pair – Mounted at Knees

Part Number	Description	HCCPS	Price	Location	Toddler Pad Sizing	Custom Pad Size
EE4A-TDLR	Fixed, Flush Toddler Hdwe & Pads	(E0956)	\$287	Knee	Set Up Seat Depth x 3H Pad Profile Standard (1" thick) Thin (1/2" thick)	 A _____ B _____ C _____ D _____


Notes:

Lateral Pelvic Support Components – Pair

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
MOD-4B	Lateral Pelvic Pads Only	(E0956)	\$197 pr.	450T-2	Fixed, Flush Toddler Lateral Hdwe	(K0108)	\$ 89 pr.
MOD-4D	Extended Lateral Pelvic Hardware	(K0108)	\$102 pr.				
MOD-4A	Contour Pad Modification	(K0108)	\$110 pr.				

Simply Toddler Seating Order Form

Lateral Thoracic Support Package – Pair

Part Number	Description	HCCPS	Price	Shape	Pad Profile	Std Pad Size	Custom Pad Size
EE5A-TDLR	Fixed, Flush Toddler Hdwe & Pads	(E0956)	\$287	Flat	Standard (1" thick)	4L x 3H	
EE5E-TDLR	Lift & Swing-Away, Flush Toddler Hdwe & Pads (Flush Only)	(E0956 & E1028)	\$573	Curved Reverse Curve	Thin (1/2" thick)	5L x 4H 6L x 5H	

Notes:

Lateral Thoracic Support Components – Pair

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
EE5I	Tracks For Lateral Hdwe Mounting	(K0108)	\$ 64 pr.				
MOD-5G	Lateral Thoracic Pads Only	(E0956)	\$207 pr.	MOD-5F	Protective Boots for Hardware	(K0108)	\$ 63 pr.
MOD-5J	Extended Hardware	(K0108)	\$ 182 pr.				
450T-2	Fixed, Flush Toddler Lateral Hdwe	(K0108)	\$ 89 pr.	455T-2	Lift, Swing-Away Toddler Lateral Hdwe	(E1028)	\$367 pr.

Notes:

Headrest Support Hardware – Each

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
7733	Adaptaplate - Headrest Hardware Interface Plate for Freedom Designs & Other Manufacturers	(E1028)	INCL.	7700-TDLR	Height, Depth Adjustable - Removable Toddler Hdwe	(E1028)	\$174ea.
				7714-TDLR	Height, Depth, & Angle Adjustable - Removable Toddler Hdwe	(E1028)	\$222 ea.

Notes:



Headrest Support – Each

Part Number	Description	HCCPS	Price	Standard Sizes		
7756	Chubbie Headrest Pad Only	(E0955)	\$227 ea.	X-Sm – 8W x 3H	Sm – 10W x 4H	Med – 12W x 5H
7756-ADJ	Chubbie 3-Piece Adjustable Headrest Pad	(E0955 & E1028)	\$432 ea.		Sm – 10W x 4H	Med – 12W x 5H
7757	Low-Profile Headrest Pad	(E0955)	\$180 ea.		Sm – 8.5W x 4H	Med – 9W x 4H
7753	Hinge-Loc 3-Piece Adjustable Headrest Pad	E0955 & E1028	\$380 ea.		Sm – 4W x 4H All Pads	Med – 5W x 4H All Pads
7760	Skinny Headrest Pad Only	(E0955)	\$227 ea.	X-Sm – 6W x 3.5H	Sm – 8W x 4.5H	Med – 10W x 4.5H

Notes:


Simply Toddler Seating Order Form

Lower Extremity Supports – Pair

Part Number	Description	HCCPS	Price	Standard Sizes			Closure Type			
EE9B	2-Point Hi-Tops STATIC Ankle Supports	(K0108)	\$135 pr.	X-Small	Small	Medium	Hook & Loop Buckle			
EE9B-D	2-Point Hi-Tops DYNAMIC Ankle Supports	(K0108)	\$135pr.	X-Small	Small	Medium				
EE9F	3-Point Hi-Tops STATIC Ankle Supports	(K0108)	\$204 pr.	X-Small	Small	Medium				
EE9F-D	3-Point Hi-Tops DYNAMIC Ankle Supports	(K0108)	\$204 pr.	X-Small	Small	Medium				
				Ankle Circumference Range:	5.5" - 8.5"	7.5" - 10.5"	9.5" - 12.5"			
EE9H	Foot & Ankle Positioners (Ankle Positioner Measures Same as Above)	(K0108)	\$255 pr.	X-Small	Sm / Med					
				Mounting Plate Measures:	4.5"W x 7.5"L	5"W x 9"L				
EE9	Shoe Holders – Hook & Loop Closure	(K0108)	\$201 pr.	Small	Small/Med	Medium				
EE9C	Shoe Holders – Buckle Closure	(K0108)	\$201 pr.	Small	Small/Med	Medium				
EE9G	Shoe Holders – Buckle Closure Ankle, Hook & Loop Toe Straps	(K0108)	\$201 pr.	Small	Small/Med	Medium				
				Shoe Holder Size (inches):	2.75W x 5L	3W x 5.75L	3.25W x 6.5L			
08-10000	Individual Leg Cradles	(K0108)	\$728 pr.	 						
08-11000	One-Piece, Double Leg Footbox	(K0108)	\$623 ea.							
MOD-8B	1/4" Black ABS (Bumper)	(K0108)	\$ 58 ea.							
Hanger Mount Footplate / Footboard Mount W _____ or Width to Fit Frame D _____ H _____										
Notes:										

Simply Toddler Seating Order Form




Anterior Support Tray – Packages



Standard Tray Sizes (Approximate)				
	Size	Chair	Tray	Body Cut
	X-Small	10-11"	21 x 18"	10 x 7"
	Small	12-13"	21 x 18"	12 x 7"
	Medium	14-15"	24 x 22"	14 x 9"
	Large	16-17"	26 x 24"	16 x 10"
	X-Large	18-19"	28 x 24"	18 x 10"


Wheelchair Tray Packages – Including Hardware			
Part Number	Description	HCPCS	Price
08-16000	1/4" Lexan Tray - Hook & Loop Straps	(E0950)	\$128
08-12000	1/4" Lexan Tray - EZ Lite Channel Hardware	(E0950)	\$206
08-13000	1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs	(E0950)	\$221
08-4000	1/2" Birch Tray - Hook & Loop Straps	(E0950)	\$234
08-1000	3/8" Lexan Tray - Hook & Loop Straps	(E0950)	\$239
08-6000	3/8" Lexan Tray - Adjustable Knob Hardware & Ribs	(E0950)	\$265
08-6500	1/2" Birch Tray - Adjustable Knob Hardware & Ribs	(E0950)	\$265
08-2500	1/2" Birch Tray - Channel Bracket Hardware	(E0950)	\$279
08-7500	1/2" Birch Tray - Toggle Hardware	(E0950)	\$279
08-2000	3/8" Lexan Tray - Channel Bracket Hardware	(E0950)	\$328
08-7000	3/8" Lexan Tray - Toggle Hardware	(E0950)	\$328

Notes:



Simply Toddler Seating Order Form

Anterior Support Tray – Build to Order



A _____ B _____

C _____ D _____

Standard Tray Sizes (Approximate)				
	Size	Chair	Tray	Body Cut
	X-Small	10-11"	21 x 18"	10 x 7"
	Small	12-13"	21 x 18"	12 x 7"
	Medium	14-15"	24 x 22"	14 x 9"
	Large	16-17"	26 x 24"	16 x 10"
	X-Large	18-19"	28 x 24"	18 x 10"

Wheelchair Trays - No Hardware			
Part Number	Description	HCPSCS	Price
08-17000	1/4" Lexan Tray - No Hardware	(E0950)	\$115
08-17100	1/4" Bolteron Tray - No Hardware (Black)	(E0950)	\$163
08-8000	3/8" Lexan Tray - No Hardware	(E0950)	\$193
08-8100	3/8" Bolteron Tray - No Hardware (Black)	(E0950)	\$232
08-8500	1/2" Birch Tray - No Hardware	(E0950)	\$193
Notes:			

Wheelchair Tray Hardware			
Part Number	Description	HCPCS	Price
454A	Adjustable Knob Release Hardware	(E1028)	\$126 pr.
454B	Adjustable Toggle Clamp Hardware	(E1028)	\$298 pr.
454C	Hook & Loop Strap Hardware	(E1028)	\$ 98 pr.
461	Channel Bracket Hardware	(E1028)	\$160 pr.
559	EZ Lite Adjustable Channel Bracket	(E1028)	\$160 pr.

Notes:

Simply Toddler Seating Order Form

Custom Wheelchair Tray Modifications

Part Number	Description	HCPCS	Price	Quantity
MOD-8A	3/8" Lexan Material Upgrade For Trays (Substitutes Std Tray Material with 3/8" Lexan)	(K0108)	\$ 78	Each
MOD-8A1	1/4" Bolteron Material Upgrade For Trays(Substitutes Std Tray Material with 1/4" Bolteron)	(K0108)	\$ 59	Each
MOD-8A2	3/8" Bolteron Material Upgrage For Trays (Substitutes Std Tray Material with 3/8" Bolteron)	(K0108)	\$117	Each
MOD-8C	Joystick Cut Out Right Left _____ L x _____ W	(K0108)	\$ 106	Each
MOD-8D	Channel Tray Stopper Hardware Kit	(K0108)	\$ 50	Kit
MOD-8E	Custom Slots for Hardware – Must Specify # & Location	(K0108)	\$ 106	Pair
MOD-8F	Round Body Cut-Out		N/C	Each
MOD-8G	Aluminum Rim - 1/16" x 1"H Rims Sides and Front Edge	(K0108)	\$ 66	Kit
MOD-8H	Plastic Rim Guard - 3/8" Attaches Directly To Tray Edge	(K0108)	\$ 43	Each
MOD-8L	Pad Tray - Tray Top Only _____ Upholstery	(K0108)	\$364	Each
MOD-8L-1	Pad Tray - Top & Bottom _____ Upholstery	(K0108)	\$605	Each
MOD-8M	Pad Tray - Half Tray Pad _____ Upholstery	(K0108)	\$306	Each
MOD-8N	Ribs - Wood Lexan Aluminum	(K0108)	\$ 58	Set
MOD-8P	Overlay - 1/8" Lexan	(K0108)	\$140	Each
MOD-8Q	Underlay - 1/8" Lexan	(K0108)	\$140	Each
MOD-8R	Easel - Removable, Adjustable, 3/8" Lexan w/ Pencil Ledge	(K0108)	\$269	Each
MOD-8S	Tray Trim - 1/16" For Aluminum Rim	(K0108)	\$ 43	Each
MOD-8T	Slip Cover For Tray Padding _____ Upholstery	(K0108)	\$271	Each

Notes:

Custom Wheelchair Tray Accessories


Part Number	Description	HCPCS	Price	Quantity
MOD-8J	Rear Stability Belt – 1" Webbing W/ Side Release Buckle	(K0108)	\$ 48	Each
MOD-8K	Elbow Blocks w/ L Bracket _____ L x _____ H	(K0108)	\$287	Pair
MOD-8O	Elbow Pads - Hook & Loop Mounting _____ L x _____ H	(K0108)	\$134	Pair
MOD-8OW	Elbow Pads - Wood Base _____ L x _____ H	(K0108)	\$287	Pair
	Upholstery _____			
MOD-8V	Wrist Straps - Hook & Loop	(K0108)	\$ 66	Pair

Notes / Diagrams:


Simply Toddler Seating Order Form

Freedom Select Anterior Trunk Supports

Select Standard



Select Narrow



Dimensions

Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SNZ-S	Select Standard Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-1-S	Select <i>Standard</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-M	Select <i>Standard</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-1-M	Select <i>Standard</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-L	Select <i>Standard</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-1-L	Select <i>Standard</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-XL	Select <i>Standard</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6SNZ-1-XL	Select <i>Standard</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235


	SFDI-6NNZ-XS	Select <i>Narrow</i> X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-1-XS	Select <i>Narrow</i> Dual Rear-Pull X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-S	Select <i>Narrow</i> Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-1-S	Select <i>Narrow</i> Dual Rear-Pull Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-M	Select <i>Narrow</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-1-M	Select <i>Narrow</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-L	Select <i>Narrow</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-1-L	Select <i>Narrow</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-XL	Select <i>Narrow</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6NNZ-1-XL	Select <i>Narrow</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

Note:


Simply Toddler Seating Order Form

Freedom Select Anterior Trunk Supports (continued)

Select Standard Zipper



Select Narrow Zipper



Dimensions

Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SZ-S	Select <i>Standard Zipper</i> Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 167
	SFDI-6SZ-1-S	Select <i>Standard Zipper</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 167
	SFDI-6SZ-M	Select <i>Standard Zipper</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 167
	SFDI-6SZ-1-M	Select <i>Standard Zipper</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 167
	SFDI-6SZ-L	Select <i>Standard Zipper</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SZ-1-L	Select <i>Standard Zipper</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SZ-XL	Select <i>Standard Zipper</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6SZ-1-XL	Select <i>Standard Zipper</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 179
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235


	SFDI-6NZ-M	Select <i>Narrow Zipper</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 167
	SFDI-6NZ-1-M	Select <i>Narrow Zipper</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 167
	SFDI-6NZ-L	Select <i>Narrow Zipper</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 173
	SFDI-6NZ-1-L	Select <i>Narrow Zipper</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 173
	SFDI-6NZ-XL	Select <i>Narrow Zipper</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 179
	SFDI-6NZ-1-XL	Select <i>Narrow Zipper</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 179
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

Note:

Simply Toddler Seating Order Form

Freedom Select Anterior Pelvic Supports

Select 2-Point



2-POINT – PADDED – EZ RELEASE PUSH BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2MP-XS	XS	MICRO 1"	1	50	1.75	6	\$ 94
	SFDI-3A2MP-S	S	MICRO 1"	1	55	2.25	7	\$ 94
	SFDI-3A2P-S	S	STD 1.5"	1.5	55	2.25	7	\$ 94
	SFDI-3A2P-M	M	STD 1.5"	1.5	60	2.5	9	\$ 94
	SFDI-3A2P-L	L	STD 1.5"	1.5	65	3	11	\$ 94
	SFDI-3A2LP-L	L	LG 2"	2	65	3	11	\$ 94
	SFDI-3A2LP-XL	XL	LG 2"	2	70	3.25	13	\$ 94
	SFDI-3A2P-CUSTOM		QUOTE					\$ 118

2-POINT – PADDED – SIDE RELEASE BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2MS-XS	XS	MICRO 1"	1	50	1.75	6	\$ 68
	SFDI-3A2MS-S	S	MICRO 1"	1	55	2.25	7	\$ 68
	SFDI-3A2S-S	S	STD 1.5"	1.5	55	2.25	7	\$ 68
	SFDI-3A2S-M	M	STD 1.5"	1.5	60	2.5	9	\$ 68
	SFDI-3A2S-L	L	STD 1.5"	1.5	65	3	11	\$ 68
	SFDI-3A2LS-L	L	LG 2"	2	65	3	11	\$ 68
	SFDI-3A2LS-XL	XL	LG 2"	2	70	3.25	13	\$ 68
	SFDI-3A2S-CUSTOM		QUOTE					\$ 94

Note:

2-POINT – UNPADDED – EZ RELEASE PUSH BUTTON BUCKLE

Note: Pads not included, must be ordered separately.

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2UMP-XS	XS	MICRO 1"	1	50	1.75	6	\$ 45
	SFDI-3A2UMP-S	S	MICRO 1"	1	55	2.25	7	\$ 45
	SFDI-3A2UP-S	S	STD 1.5"	1.5	55	2.25	7	\$ 45
	SFDI-3A2UP-M	M	STD 1.5"	1.5	60	2.5	9	\$ 45
	SFDI-3A2UP-L	L	STD 1.5"	1.5	65	3	11	\$ 45
	SFDI-3A2ULP-L	L	LG 2"	2	65	3	11	\$ 45
	SFDI-3A2ULP-XL	XL	LG 2"	2	70	3.25	13	\$ 45
	SFDI-3A2UP-CUSTOM		QUOTE					\$ 57

Note:

2-POINT ANTERIOR PELVIC SUPPORT PADS

Qty.	PRICE EA.	Qty.	PRICE PR.
SFDI-3HBP-XS	\$ 32	SFDI-3HBP-XS	\$ 64
SFDI-3HBP-S	\$ 32	SFDI-3HBP-S	\$ 64
SFDI-3HBP-M	\$ 32	SFDI-3HBP-M	\$ 64
SFDI-3HBP-L	\$ 32	SFDI-3HBP-L	\$ 64
SFDI-3HBP-XL	\$ 32	SFDI-3HBP-XL	\$ 64
SFDI-3HBP-CUSTOM	\$ 44	SFDI-3HBP-CUSTOM	\$ 88


Simply Toddler Seating Order Form

Freedom Select Anterior Pelvic Supports (continued)

2-POINT – UNPADDED – SIDE RELEASE BUTTON BUCKLE

Note: Pads not included, must be ordered separately.

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price	2-POINT ANTERIOR PELVIC SUPPORT PADS					
									Qty.	PRICE EA.	Qty.	PRICE PR.		
	SFDI-3A2UMS-XS	XS	MICRO 1"	1	50	1.75	6	\$ 40	SFDI-3HBP-XS		\$ 32		\$ 64	
	SFDI-3A2UMS-S	S	MICRO 1"	1	55	2.25	7	\$ 40	SFDI-3HBP-S		\$ 32		\$ 64	
	SFDI-3A2US-S	S	STD 1.5"	1.5	55	2.25	7	\$ 40	SFDI-3HBP-M		\$ 32		\$ 64	
	SFDI-3A2US-M	M	STD 1.5"	1.5	60	2.5	9	\$ 40	SFDI-3HBP-L		\$ 32		\$ 64	
	SFDI-3A2US-L	L	STD 1.5"	1.5	65	3	11	\$ 40	SFDI-3HBP-XL		\$ 32		\$ 64	
	SFDI-3A2ULS-L	L	LG 2"	2	65	3	11	\$ 40	SFDI-3HBP-CUSTOM		\$ 44		\$ 88	
	SFDI-3A2ULS-XL	XL	LG 2"	2	70	3.25	13	\$ 40						
	SFDI-3A2US-CUSTOM		QUOTE					\$ 51						



4-POINT – EZ RELEASE PUSH BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	E	F	Price
	SFDI-3A4MP-XS	XS	MICRO 1"	1	1	50	14	4.25	5.75	\$203
	SFDI-3A4MP-S	S	MICRO 1"	1	1	55	16	5	6.5	\$203
	SFDI-3A4P-S	S	STD 1.5"	1.5	1	55	16	5	6.5	\$203
	SFDI-3A4P-M	M	STD 1.5"	1.5	1	60	18	5.5	7.5	\$203
	SFDI-3A4P-L	L	STD 1.5"	1.5	1.5	65	20	6	8.75	\$203
	SFDI-3A4P-XL	XL	STD 1.5"	1.5	1.5	70	22	6	8.75	\$203
	SFDI-3A4LP-XL	XL	LG 2"	2	1.5	70	22	6	8.75	\$203
	SFDI-3A4P-CUSTOM		QUOTE							\$235

4-POINT – SIDE RELEASE BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	E	F	Price
	SFDI-3A4MS-XS	XS	MICRO 1"	1	1	50	14	4.25	5.75	\$202
	SFDI-3A4MS-S	S	MICRO 1"	1	1	55	16	5	6.5	\$202
	SFDI-3A4S-S	S	STD 1.5"	1.5	1	55	16	5	6.5	\$202
	SFDI-3A4S-M	M	STD 1.5"	1.5	1	60	18	5.5	7.5	\$22
	SFDI-3A4S-L	L	STD 1.5"	1.5	1.5	65	20	6	8.75	\$202
	SFDI-3A4S-XL	XL	STD 1.5"	1.5	1.5	70	22	6	8.75	\$202
	SFDI-3A4LS-XL	XL	LG 2"	2	1.5	70	22	6	8.75	\$202
	SFDI-3A4S-CUSTOM		QUOTE							\$235

Note:

Simply Toddler Seating Order Form



PO# _____

Marked for: _____

Cell
Size

1"

2"

3"

4"

Deleted Cells
Base Intact

X

Cut Out Cells
& Base

C

Valve Hose
Location
Corner & Cell
Height

V1

V2

V3

V4

Select Seat
Choice

Standard

☐

Asymmetrical
Left

☐

Asymmetrical
Right

☐

Other

☐

MOD-1ROHO SEAT

Freedom Seat ROHO Insert

\$924.00

Note: Retail pricing includes any combination of specification option.

FRONT OF INSERT

COLUMNS

INCHES		3	4.75	6.5	8	9.75	11.5	13.25	15	16.5	18.25	20	21.75
CELLS		1	2	3	4	5	6	7	8	9	10	11	12
		A	B	C	D	E	F	G	H	I	J	K	L
3	1	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
4.75	2	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
6.5	3	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
8	4	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
9.75	5	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
11.5	6	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
13.25	7	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
15	8	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
16.5	9	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
18.25	10	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
20	11	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
21.75	12	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

R
O
W
S

INCHES ARE FROM BASE EDGE TO EDGE OF INSERT

1. Start diagram at cell 1A which represents front left of cushion from client's seated position.
2. Specify "Cell Height" (1", 2", 3" or 4") by placing a number in each cell box.
3. Specify request to "Delete Cells with Base Intact" by typing an "X" in those cell boxes.
4. Specify request to "Cut Out Cells & Base" by typing a "C" in those cell boxes.
5. Specify "Valve Hose Location Corner & Cell Height" by typing "V1", "V2", "V3" or "V4" in that corner box.

Note: For multiple valve or custom valve locations and custom air chamber divisions call Freedom Designs Customer Service for assistance at 800-331-8551

PLEASE ALLOW 3-4 WEEKS FOR DELIVERY

Simply Toddler Seating Order Form



PO# _____

Marked for: _____

MOD-2ROHO BACK

Freedom Back ROHO Insert

\$924.00

Note: Retail pricing includes any combination of specification option.

TOP OF INSERT

COLUMNS

INCHES → 3 4.75 6.5 8 9.75 11.5 13.25 15 16.5 18.25 20 21.75

CELLS → 1 2 3 4 5 6 7 8 9 10 11 12

A B C D E F G H I J K L

R O W S	3	1													
	4.75	2													
	6.5	3													
	8	4													
	9.75	5													
	11.5	6													
	13.25	7													
	15	8													
	16.5	9													
	18.25	10													
	20	11													
	21.75	12													

INCHES ARE FROM BASE EDGE TO EDGE OF INSERT

1. Start diagram at cell 1A which represents top right of back from client's seated position.
2. Specify "Cell Height" (1", 2", or 3") by placing a number in each cell box.
3. Specify request to "Deleted Cells with Base Intact" by typing an "X" in those cell boxes.
4. Specify request to "Cut Out Cells & Base" by typing a "C" in those cell boxes.
5. Specify "Valve Hose Location Corner & Cell Height" by typing "V1", "V2", or "V3" in that corner box.

Note: For multiple valve or custom valve locations and custom air chamber divisions call Freedom Designs Customer Service for assistance at 800-331-8551

PLEASE ALLOW 3-4 WEEKS FOR DELIVERY

Cell Size

1"

2"

3"

Deleted Cells Base Intact

X

Cut Out Cells & Base

C

Valve Hose Location Corner & Cell Height

V1

V2

V3

Select Back Choice

Standard

"I" Back

"T" Back