

Simply Seating Order Form

Clear Form

Email Form to Customer Service



Freedom Designs, Incorporated
2241 N Madera Road, Simi Valley, CA 93065
phone: 800-331-8551 fax: 888-582-1509
e: customerservice@freedomdesigns.com



Yes, you can.

Quote

Order

Buyer: _____
Co./Account: _____
Location: _____
Contact #: _____
P.O.#: _____
Marked For: _____
Fax: _____
Email: _____


Seating Specifications

A Back Height: _____
B Lateral Height: _____
C Lateral Width: _____
D Hip Width (rear): _____
E Set Up Seat Depth: _____
F Hip Width (front): _____
G "E" + 3" for Growth: _____
H Seat Width: _____
I Back Width: _____

Important Note: To ensure quote accuracy, seat and back width, depth and height dimensions must be provided. Seats and backs over 18" in width, depth or height are subject to oversize charges.

Finished Dimensions

Standard cut is Boltaron with foam rolled over edge.



Wheelchair Specifications

Model: _____
Width: _____
Depth: _____

Freedom Designs Upholstery Selections

1. Choose Contact / Non-Contact
2. Choose Upholstery
3. Choose Component

	Seat Base	Back Base	Headrest Pad	Lateral Pads	Hip & Knee Pads
Body Contact Surfaces					
Primary Upholstery Option:					
Number: _____ Choose: _____					
Secondary Upholstery Option:					
Number: _____ Choose: _____					
Non-Contact Surfaces					
Primary Upholstery Option:					
Number: _____ Choose: _____					
Secondary Upholstery Option:					
Number: _____ Choose: _____					
Monogram Option	When ordering choose either Monogram or Embroidery				
Name: _____					
Script _____ Block _____					
Thread Color: _____					
Embroidery Option					
Design: _____					

Notes:

Simply Seating Order Form


Solid Seat Insert

**\$1260.00
(E2609)**

Simply Seating EB1 Seat Construction Includes:



Boltaron base material 1/4" thick up to 14" x 14", 3/8" thick for sizes larger than 14" • Cut full width • 2" x 2" back post notches standard • Boltaron "T" nutted for tracks, hinge and medial thigh support • Base sizes larger than 18" x 18" subject to oversize charge • Tracks for hardware and accessory attachment • Base attachment hardware • Seat to back hinge • 1 1/2" 4 lb. foam is standard • 3-sided zipper cover •

Step 1. Choose EB1-Kit (Seat)

A	EB1-Kit	2-Point Adjustable Mounting Hdwe & Seat to Back Hinge (Crash Test Approved)	C	EB1-Knob-Kit	2-Point Adjustable Mounting Hdwe & Seat to Back Hinge (Crash Test Approved)
	2 - 585(X)	Aluminum Fast Clamp Hdwe Assy		2 - 462A	Knob Mounting Hdwe Assy
	2 - 568	Frame Stopper Clamps 3/4"-1"		2 - 439	Seat Mount Bracket
	2 - 436	Aluminum 30 Degree Brackets		2 - 436	Aluminum 30 Degree Brackets
	Note:	Choose 1 Seat to Back Hinge Option From D		Note:	Choose 1 Seat to Back Hinge Option From D
A1	EB1-EXT-Kit	2-Point Mounting Hdwe & Seat to Back Hinge	C1	EB1-Knob-Kit Options:	ADD
	2 - 504B	EZ Lite Extended Mounting Hdwe	C2	EB1-Knob-4	Upgrade EB1-Knob-Kit to 4-Point Knob Hdwe (replaces hinge) \$419
	2 - 568	Frame Stopper Clamps 3/4"-1"	C3	EB1-Knob-Pin2	Upgrade EB1-Knob-Kit to 2-Point Pin-Lock Knob Hdwe \$102
	2 - 436	Aluminum 30 Degree Brackets		EB1-Knob-Pin4	Upgrade EB1-Knob-Kit to 4-Point Pin-Lock Knob Hdwe (replaces hinge) \$527
	Note:	Choose 1 Seat to Back Hinge Option From D	D	Seat to Back Hinge Options	
B	EB1-4-Kit	4-Point Adjustable Mounting Hdwe NO Seat to Back Hinge (Crash Test Approved)		1 - 430A	Seat to Back Hinge Flush (pair)
	4 - 585(X)	Aluminum Fast Clamp Hdwe Assy		1 - 430	Seat to Back Hinge Offset (pair)
	2 - 568	Frame Stopper Clamps 3/4"-1"		1 - 430B	Flip-And-Fold Seat to Back Hinge (pair)

Notes:

Step 2. Choose Customization to EB1 Seat

EB1CS1K	Anti-Thrust Foam Modification	\$255	EB1CS1Y	Contour Foam Modification - ADD, ABD & Anti-Thrust	\$394
	• Height of Anti-Thrust shelf is 1"			• Adductor height is 1"	
	• C= First 1/3 of set-up seat depth (see "E" page 1)			• Height of Anti-Thrust shelf is 1"	
				• C= First 1/3 of set-up seat depth (see "E" page 1)	
EB1CS1V	Contour Foam Modification - ABD, Anti-Thrust Only	\$342	EB1CS1YA	Contour Foam Modification - ADD, ABD, NO Anti-Thrust	\$394
EB1CS1VA	Contour Foam Modification - ABD Only	\$342	EB1CS1W	Contour Foam Modification - ADD, Anti-Thrust Only	\$342
			EB1CS1WA	Contour Foam Modification - ADD Only	\$342


Notes:

Simply Seating Order Form

Step 2. Choose Customization to EB1 Seat (continued)

EB1CS1CR Asymmetrical - Long Right


\$ 106



A _____
B _____
C _____
D _____

EB1CS1CL Asymmetrical - Long Left

\$ 106



A _____
B _____
C _____
D _____

EB1CS1WW Flared Seat Cut Modification for LEAF Frame (specify shape to fit LEAF frame)

\$493

EB1CS1MLB 1/2" Maple Laminate Base

N/C

EB1CS1WT Base with T-Nuts for Hardware (no tracks included)

N/C

EB1CS1EE Base Only Cut Between Rails

\$198

EB1CS1G Cushion Separate From Base

\$127

EB1CS1I Incontinence Resistance Liner

\$ 37

EB1CS1FF Removable Incontinence Slip Cover

\$247

EB1CS1-ROHOSEAT

Freedom Seat **ROHO** Insert (see page 16 to specify order)

\$924

EB1CS1DD

Modular Contour Foam Modification

\$523

EB1S1UU

Cut, Carve, Contour Foam Mod.

\$557

EB1CS1L

Grid Cut Foam Modification

\$403

Oversize Seat Upcharge

EB1CS1OVER18 Seat Insert Over 18W OR 18D

\$ 81

EB1CS1OVER22 Seat Insert Over 22W OR 22D

\$ 110

EB1CS1OVER24 Seat Insert Over 24W OR 24D

\$136

Notes:

Step 3. Choose Foam

1-1/2" Std 4# foam Only

OR

Custom Foam / Omit Standard Foam (see chart below)

Base layer to be: _____ Middle layer to be: _____ Top layer to be: _____

Seat Foam Chart

3# Foam	\$ 55 per 1"	Firm Sunmate	\$ 99 per 1"	Pudgee	\$252 per 1"
4# Foam	\$ 55 per 1"	X-Soft T-Foam	\$ 107 per 1"	Waffle	\$325 per 1/2"
5# Foam	\$ 57 per 1"	Soft T-Foam	\$ 107 per 1"	Waffle (over 18x16)	\$653 per 1/2"
Soft Sunmate	\$ 99 per 1"	Medium T-Foam	\$ 107 per 1"	Free Float	\$409 per 1/2"
Medium Sunmate	\$ 99 per 1"	Firm T-Foam	\$ 107 per 1"	FreeFloat (over 18x16)	\$815 per 1/2"

Notes:

Simply Seating Order Form

Solid Back Insert

**\$1329.00
(E2617)**

Simply Seating EB2 Back Construction Includes:

Boltaron planar / flat base $\frac{1}{4}$ " thick up to 18" x 18" size, $\frac{3}{8}$ " thick for sizes larger than 18" • Cut between rails unless otherwise specified • Boltaron "T" nutted for tracks, hinge and headrest plate • Base sizes larger than 18" x 18" subject to oversize charge • Tracks for hardware and accessory attachment • Base attachment hardware • Headrest adapt-a-plate • Anterior trunk support strap risers • 1" 4 lb. foam standard • 3-sided zipper cover •

Step 1. Choose EB2-Kit (Back)


A	EB2-Kit	2-Point Attachment & Seat to Back Hinge (Crash Test Approved)	C	EB2-Knob-Kit	2-Point Attachment & Seat to Back Hinge (Crash Test Approved)	
	2 - 585(X)	Plastic Fast Clamp Hdwe Assy		2 - 462A	Knob Mounting Hardware Assy	
	2 - 568	Frame Stopper Clamps $\frac{3}{4}$ "-1"		2 - 437A	Back Mount Bracket	
	Note:	Choose 1 Strap Riser (pair) from D		Note:	Choose 1 Strap Riser (pair) from D	
						
				EB2-Knob-Kit Options:		ADD
			C1	EB2-Knob-4	Upgrade EB2-Knob-Kit to 4-Point Knob Hdwe (replaces hinge)	\$490
			C2	EB2-Knob-Pin2	Upgrade EB2-Knob-Kit to 2-Point Pin-Lock Knob Hdwe	\$106
			C3	EB2-Knob-Pin4	Upgrade EB2-Knob-Kit to 4-Point Pin-Lock Knob Hdwe (replaces hinge)	\$603
B	EB2-4-Kit	4-Point Attachment - NO Seat to Back Hinge (Crash Test Approved)	D	Strap Riser Choices (choose 1 pair)		INCL.
	4 - 585(X)	Plastic Fast Clamp Hdwe Assy		2 - MOD-6A	1- $\frac{1}{2}$ " Wide Strap Opening, Flat	
	2 - 568	Frame Stopper Clamps $\frac{3}{4}$ "-1"		2 - MOD-6A-1	1" Wide Strap Opening, Flat	
	Note:	Choose 1 Strap Riser (pair) from D		2 - MOD-6B	1- $\frac{1}{2}$ " Wide Strap Opening, 25° Rear Bend	
				2 - MOD-6B-1	1" Wide Strap Opening, 25° Rear Bend	

Notes:


Step 2. Choose Customization to EB2 Back

EB2CB2H I-Cut Modification **\$155**

EB2CB2J T-Cut Modification **\$155**




A _____
B _____
C _____
D _____
E _____



EB2CB2I Curved Base Modification **\$163**

EB2CB2IA Curved Foam Modification **\$140**

E _____
F _____




Notes:

Simply Seating Order Form

Step 2. Choose Customization to EB2 Back (continued)


EB2CB2L Biangular Cut Modification / Dual Cover \$366

EB2CB2LA Biangular Cut Modification / Single Cover \$366



E _____
F _____
G _____
H _____

EB2CB2DD Adjustable Lumbar Modification (n/a with Toddler Back) \$362



E _____
F _____
J _____

EB2CB2MLB 1/2" Maple Laminate Base N/C

EB2CB2WT Base with T-Nuts for Hardware (no tracks included) N/C

EB2CB2GG Ergo Adjustable Modification \$510

EB2CB2K Grid Cut Foam Modification \$403

Oversize Back Upcharge

EB2CB2OVER22 Back Insert Over 22W OR 22D \$110

EB2CB2OVER24 Back Insert Over 24W OR 24D \$136

EB2CB2OVER-CUSTOM Back Insert Over 28", Every 3" (in addition to Over24 charge) \$ 74

EB2CB2-ROHOBACK Freedom Back, **ROHO** Insert (see page 17 to specify) \$924

EB2CB2S Free Form Modification Soft/Med Sunmate - Box Kit (specify below) \$313

EB2CB2SBAG Free Form Modification Soft/Med Sunmate - Bag Kit (specify below) \$313

EB2CB2S1 Liquid Soft/Med Sunmate Only - Box Kit (specify below) \$293

EB2CB2S1BAG Liquid Soft/Med Sunmate Only - Bag Kit (specify below) \$293

Soft Sunmate Medium Sunmate

Notes:

Step 3. Choose Foam

1" Std 4# foam Only

OR

Custom Foam / Omit Standard Foam (see chart below)

Base layer to be: _____ Middle layer to be: _____ Top layer to be: _____


Back Foam Chart

3# Foam	\$ 55 per 1"	Firm Sunmate	\$127 per 1"	Pudgee	\$273 per 1"
4# Foam	\$ 55 per 1"	X-Soft T-Foam	\$137 per 1"	Waffle	\$325 per 1/2"
5# Foam	\$ 75 per 1"	Soft T-Foam	\$137 per 1"	Waffle (over 18x16)	\$653 per 1/2"
Soft Sunmate	\$127 per 1"	Medium T-Foam	\$137 per 1"	Free Float	\$409 per 1/2"
Medium Sunmate	\$127 per 1"	Firm T-Foam	\$137 per 1"	FreeFloat (over 18x16)	\$815 per 1/2"

Notes:

Simply Seating Order Form

Medial Pelvic Support Package – Each

Part Number	Description	HCCPS	Price	Shape	Standard Pad Sizing	Custom Pad Size
EE3B-2	Lever Flip-Down – Choice of Pad	(E0957 & 1028)	\$536	Wedge	Small – 3L x 3H x 3-2 Wedge	
EE3S	Push Button, Flip-Down – Choice of Pad	(E0957 & 1028)	\$509	Oval	Medium – 4L x 3H x 3-2 Wedge	
EE3P	Hide Away - Lever, Flip-Down – Choice of Pad	(E0957 & 1028)	\$510	Round	Large – 5L x 4H x 4-2 Wedge	
					X-Large – 6L x 5H x 5-3 Wedge	


Notes:

Medial Pelvic Support Components – Each

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
MOD-3M	Protective Boot For Hardware	(K0108)	\$ 72 ea.	MOD-3E	Medial Pelvic Support Pad Only	(E0957)	\$206 ea.
MOD-3F	Extended Flip Bracket for Abductor	(K0108)	\$ 34 ea.				
424-2	Lever Flip-Down Hdwe (EE3B-2) W/ STD Top Plate	(E1028)	\$299 ea.	581	Push Button Flip-Down Hdwe (EE3S) W/ STD Top Plate	(E1028)	\$281 ea.
543	Hide Away, Lever Flip-Down Hdwe (EE3P) W/ STD Top Plate	(E1028)	\$299 ea.				


Notes:

Lateral Pelvic Support Package – Pair – Mounted at Hips

Part Number	Description	HCCPS	Price	Offset	Standard Pad Sizing	Custom Pad Size
EE4A	Fixed, Flush STD Hdwe & Pads	(E0956)	\$287	1"	Set Up Seat Depth x 4H	
EE4B	Fixed, Offset STD Hdwe & Pads (Specify Offset)	(E0956)	\$287	2"	Pad Profile	
EE4C	Quick-Release, Flush or Offset STD Hdwe & Pads (Flush & 1" Only)	(E0956 & 1028)	\$573	Location	Standard (1" thick)	
				Hip	Thin (1/2" thick)	

Notes:

Lateral Pelvic Support Package – Pair – Mounted at Knees

Part Number	Description	HCCPS	Price	Offset	Standard Pad Sizing	Custom Pad Size
EE4A	Fixed, Flush STD Hdwe & Pads	(E0956)	\$287	1"	Set Up Seat Depth x 4H	
EE4B	Fixed, Offset STD Hdwe & Pads (Specify Offset)	(E0956)	\$287	2"	Pad Profile	
EE4C	Quick-Release, Flush or Offset STD Hdwe & Pads (Flush & 1" Only)	(E0956 & 1028)	\$573	Location	Standard (1" thick)	
				Knee	Thin (1/2" thick)	

Notes:

(Lateral Pelvic Components continued next page)


Simply Seating Order Form

Lateral Pelvic Support Components - Pair

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
MOD-4B	Lateral Pelvic Pads Only	(E0956)	\$ 197 pr.	MOD-4A	Contour Pad Modification	(K0108)	\$ 110 pr.
MOD-4D	Extended Lateral Pelvic Hardware	(K0108)	\$ 102 pr.				
450A	Lateral Bracket, Fixed Flush STD	(K0108)	\$ 68 pr.	442	Lateral Bracket, Fixed " Offset STD	(K0108)	\$ 74 pr.
443	Lateral Bracket, Fixed 2" Offset STD	(K0108)	\$ 74 pr.				
456A	Lateral Bracket, Quick-Release Flush STD	(E1028)	\$342 pr.	456B	Lateral Bracket, Quick-Release 1" Offset STD	(E1028)	\$342 pr.

Notes:

Lateral Thoracic Support Package – Pair

Part Number	Description	HCCPS	Price	Shape	Pad Profile	Std Pad Size	Custom Pad Size
EE5A	Fixed, Flush STD Hdwe & Pads	(E0956)	\$287	Flat	Standard (1" thick)	4L x 3H	
EE5B	Fixed, Offset STD Hdwe & Pads (Specify Offset)	(E0956)	\$287	Curved	Thin (1/2" thick)	5L x 4H	
EE5E	Lift, Swing-Away STD Hdwe & Pads (Specify Offset)	(E0956 & 1028)	\$573	Reverse Curve		6L x 5H	
EE5W	Push-Button, Swing-Away STD Hdwe & Pads (Specify Offset)	(E0956 & 1028)	\$573				
EE5P	Trigger, Swing-Away STD Hdwe & Pads (Specify Offset)	(E0956 & 1028)	\$643				
EE5K	Lift, Swing-Away Summer/Winter STD Hdwe & Pads (Specify Offset)	(E0956 & 1028)	\$765				
EE5Q	Trigger, Swing-Away Summer/Winter STD Hdwe & Pads (Specify Offset)	(E0956 & 1028)	\$808				

Notes:

Lateral Thoracic Support Components – Pair

Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
EE5I	Tracks For Lateral Hdwe Mounting	(K0108)	\$ 64 pr.				
MOD-5G	Lateral Thoracic Pads Only	(E0956)	\$207 pr.	MOD-5F	Protective Boots for Hardware	(K0108)	\$ 63 pr.
MOD-5J	Extended Hardware	(K0108)	\$102 pr.				
450A	Lateral Bracket, Fixed Flush STD	(K0108)	\$ 68 pr.	442	Lateral Bracket, Fixed 1" Offset STD	(K0108)	\$ 74 pr.
443	Lateral Bracket, Fixed 2" Offset STD	(K0108)	\$ 74 pr.				
455A	Lateral Bracket, Lift Swing-Away Flush STD	(E1028)	\$356 pr.	455C	Lateral Bracket, Lift Swing-Away 1" Offset STD	(E1028)	\$356 pr.
455D	Lateral Bracket, Lift Swing-Away 2" Offset STD	(E1028)	\$356 pr.				
589A	Lateral Bracket, Push-Button Swing-Away Flush STD	(E1028)	\$356 pr.	587A	Lateral Bracket, Push-Button Swing-Away 1" Offset STD	(E1028)	\$356 pr.
588A	Lateral Bracket, Push-Button Swing-Away 2" Offset STD	(E1028)	\$356 pr.				
510A	Lateral Bracket, Trigger Swing-Away Flush STD	(E1028)	\$406 pr.	510B	Lateral Bracket, Trigger Swing-Away 1" Offset STD	(E1028)	\$406 pr.
510C	Lateral Bracket, Trigger Swing-Away 2" Offset STD	(E1028)	\$406 pr.				
497A	Lateral Bracket, Lift Swing-Away, Summer/Winter Flush STD	(E1028)	\$524 pr.	497C	Lateral Bracket, Lift Swing-Away, Summer/Winter 1" Offset STD	(E1028)	\$524 pr.
497D	Lateral Bracket, Lift Swing-Away, Summer/Winter 2" Offset STD	(E1028)	\$524 pr.				
511A	Lateral Bracket, Trigger Swing-Away, Summer/Winter Flush STD	(E1028)	\$566 pr.	511B	Lateral Bracket, Trigger Swing-Away, Summer/Winter 1" Offset STD	(E1028)	\$566 pr.
511C	Lateral Bracket, Trigger Swing-Away, Summer/Winter 2" Offset STD	(E1028)	\$566 pr.				

Simply Seating Order Form

Headrest Support Hardware – Each



Part Number	Description	HCCPS	Price	Part Number	Description	HCCPS	Price
7733	Adaptaplate - Headrest Hardware Interface Plate for Freedom Designs & Other Manufacturers	(E1028)	INCLUDED WITH ALL EB-2-KITS	7709	Rock-N-Roll Flip-Back, Height & Depth Adjustable Hdwe (non-removable)	(E1028)	\$486 ea.
7700	Height, Depth Adjustable - Removable Hdwe	(E1028)	\$174 ea.	7715-STD	Multi-Axis Headrest Hdwe - Offset 9.5 Horizontal Bar	(E1028)	\$247 ea.
7711	Adjust-Eze Multit-Position Adjustable Hdwe	(E1028)	\$279 ea.	7715-LONG	Multi-Axis Headrest Hdwe - Offset 11-1/2 Horizontal Bar	(E1028)	\$247 ea.
7702	Flip-Back, Height & Depth Adjustable (non-removable)	(E1028)	\$341 ea.				
7708	Rock-N-Roll Height & Depth Adjustable Hdwe (Removable)	(E1028)	\$475 ea.				

Headrest Support – Each

Part Number	Description	HCCPS	Price	Standard Sizes			
7756	Chubbie Headrest Pad Only	(E0955)	\$227 ea.	X-Sm – 8W x 3H	Sm – 10W x 4H	Med – 12W x 5H	Lg – 13-1/2W x 4-1/2H
7756-ADJ	Chubbie 3-Piece Adjustable Headrest Pad	(E0955 & E1028)	\$432 ea.		Sm – 10W x 4H	Med – 12W x 5H	Lg – 13-1/2W x 4-1/2H
7757	Low-Profile Headrest Pad	(E0955)	\$180 ea.		Sm – 8-1/2W x 4H	Med – 9W x 4H	Lg – 13W x 4H
7753	Hinge-Loc 3-Piece Adjustable Headrest Pad	E0955 & E1028	\$380 ea.		Sm – 4W x 4H All Pads	Med – 5W x 4H All Pads	Lg – 5W x 5H All Pads
7760	Skinny Headrest Pad Only	(E0955)	\$227 ea.	X-Sm – 6W x 3-1/2H	Sm – 8W x 4-1/2H	Med – 10W x 4-1/2H	Lg – 12W x 4-1/2H


Notes:

Lower Extremity Supports – Pair

Part Number	Description	HCCPS Price		Standard Sizes					
EE9B	2-Point Hi-Tops STATIC Ankle Supports	(K0108)	\$135 pr.	X-Small	Small	Medium	Large	Closure Type	
EE9B-D	2-Point Hi-Tops DYNAMIC Ankle Supports	(K0108)	\$135 pr.	X-Small	Small	Medium	Large	Hook & Loop	Buckle
EE9F	3-Point Hi-Tops STATIC Ankle Supports	(K0108)	\$204 pr.	X-Small	Small	Medium	Large		
EE9F-D	3-Point Hi-Tops DYNAMIC Ankle Supports	(K0108)	\$204 pr.	X-Small	Small	Medium	Large		
		Ankle Circumference Range:		5-½" - 8-½"	7-½" - 10-½"	9-½" - 12-½"	11-½" - 14-½"		
EE9H	Foot & Ankle Positioners	(K0108)	\$255 pr.	X-Small	Small/Med	Large			
EE9	Shoe Holders – Hook & Loop Closure	(K0108)	\$201 pr.	Small	Small/Med	Medium	Med/Lg	Large	X-Large
EE9C	Shoe Holders – Buckle Closure	(K0108)	\$201 pr.	Small	Small/Med	Medium	Med/Lg	Large	X-Large
EE9G	Shoe Holders – Buckle Closure Ankle, Hook & Loop Toe Straps	(K0108)	\$201 pr.	Small	Small/Med	Medium	Med/Lg	Large	X-Large
		Shoe Holder Size (inches):		2-¾W x 5L	3W x 5-¾L	3-¼W x 6-½L	3-½W x 7-¾L	4-¼W x 8-½L	4-½W x 10L
08-10000	Individual Leg Cradles	(K0108)	\$728 pr.						
08-11000	One-Piece, Double Leg Footbox	(K0108)	\$623 ea.						
MOD-8B	¼" Black ABS (Bumper)	(K0108)	\$ 58 ea.						
Hanger Mount									
Footplate / Footboard Mount									
W _____ or Width to Fit Frame									
D _____									
H _____									

Simply Seating Order Form



Anterior Support Tray – Packages



Standard Tray Sizes (Approximate)				
	Size	Chair	Tray	Body Cut
	X-Small	10-11"	21 x 18"	10 x 7"
	Small	12-13"	21 x 18"	12 x 7"
	Medium	14-15"	24 x 22"	14 x 9"
	Large	16-17"	26 x 24"	16 x 10"
	X-Large	18-19"	28 x 24"	18 x 10"


Wheelchair Tray Packages – Including Hardware			
Part Number	Description	HCPCS	Price
08-16000	1/4" Lexan Tray - Hook & Loop Straps	(E0950)	\$128
08-12000	1/4" Lexan Tray - EZ Lite Channel Hardware	(E0950)	\$206
08-13000	1/4" Lexan Tray - EZ Lite Channel Hardware & Ribs	(E0950)	\$221
08-4000	1/2" Birch Tray - Hook & Loop Straps	(E0950)	\$234
08-1000	3/8" Lexan Tray - Hook & Loop Straps	(E0950)	\$239
08-6000	3/8" Lexan Tray - Adjustable Knob Hardware & Ribs	(E0950)	\$265
08-6500	1/2" Birch Tray - Adjustable Knob Hardware & Ribs	(E0950)	\$265
08-2500	1/2" Birch Tray - Channel Bracket Hardware	(E0950)	\$279
08-7500	1/2" Birch Tray - Toggle Hardware	(E0950)	\$279
08-2000	3/8" Lexan Tray - Channel Bracket Hardware	(E0950)	\$328
08-7000	3/8" Lexan Tray - Toggle Hardware	(E0950)	\$328

Notes:



Simply Seating Order Form

Anterior Support Tray – Build to Order



A _____ B _____

C _____ D _____

Standard Tray Sizes (Approximate)



	Size	Chair	Tray	Body Cut
	X-Small	10-11"	21 x 18"	10 x 7"
	Small	12-13"	21 x 18"	12 x 7"
	Medium	14-15"	24 x 22"	14 x 9"
	Large	16-17"	26 x 24"	16 x 10"
	X-Large	18-19"	28 x 24"	18 x 10"

Wheelchair Trays - No Hardware

Part Number	Description	HCPCS	Price
08-17000	1/4" Lexan Tray - No Hardware	(E0950)	\$115
08-17100	1/4" Bolteron Tray - No Hardware (Black)	(E0950)	\$163
08-8000	3/8" Lexan Tray - No Hardware	(E0950)	\$193
08-8100	3/8" Bolteron Tray - No Hardware (Black)	(E0950)	\$232
08-8500	1/2" Birch Tray - No Hardware	(E0950)	\$193

Notes:

Wheelchair Tray Hardware

Part Number	Description	HCPCS	Price
454A	Adjustable Knob Release Hardware 	(E1028)	\$126 pr.
454B	Adjustable Toggle Clamp Hardware 	(E1028)	\$298 pr.
454C	Hook & Loop Strap Hardware 	(E1028)	\$ 98 pr.
461	Channel Bracket Hardware 	(E1028)	\$160 pr.
559	EZ Lite Adjustable Channel Bracket 	(E1028)	\$160 pr.

Notes:

Simply Seating Order Form

Custom Wheelchair Tray Modifications

Part Number	Description	HCPCS	Price	Quantity
MOD-8A	3/8" Lexan Material Upgrade For Trays (Substitutes Std Tray Material with 3/8" Lexan)	(K0108)	\$ 78	Each
MOD-8A1	1/4" Bolteron Material Upgrade For Trays (Substitutes Std Tray Material with 1/4" Bolteron)	(K0108)	\$ 59	Each
MOD-8A2	3/8" Bolteron Material Upgrade For Trays (Substitutes Std Tray Material with 3/8" Bolteron)	(K0108)	\$117	Each
MOD-8C	Joystick Cut Out Right Left _____ L x _____ W	(K0108)	\$106	Each
MOD-8D	Channel Tray Stopper Hardware Kit	(K0108)	\$ 50	Kit
MOD-8E	Custom Slots for Hardware – Must Specify # & Location	(K0108)	\$ 106	Pair
MOD-8F	Round Body Cut-Out		N/C	Each
MOD-8G	Aluminum Rim - 1/16" x 1"H Rims Sides and Front Edge	(K0108)	\$ 66	Kit
MOD-8H	Plastic Rim Guard - 3/8" Attaches Directly To Tray Edge	(K0108)	\$ 43	Each
MOD-8L	Pad Tray - Tray Top Only _____ Upholstery	(K0108)	\$364	Each
MOD-8L-1	Pad Tray - Top & Bottom _____ Upholstery	(K0108)	\$605	Each
MOD-8M	Pad Tray - Half Tray Pad _____ Upholstery	(K0108)	\$306	Each
MOD-8N	Ribs - Wood Lexan Aluminum	(K0108)	\$ 58	Set
MOD-8P	Overlay - 1/8" Lexan	(K0108)	\$140	Each
MOD-8Q	Underlay - 1/8" Lexan	(K0108)	\$140	Each
MOD-8R	Easel - Removable, Adjustable, 3/8" Lexan w/ Pencil Ledge	(K0108)	\$268	Each
MOD-8S	Tray Trim - 1/16" For Aluminum Rim	(K0108)	\$ 43	Each
MOD-8T	Slip Cover For Tray Padding _____ Upholstery	(K0108)	\$271	Each

Notes:

Custom Wheelchair Tray Accessories


Part Number	Description	HCPCS	Price	Quantity
MOD-8J	Rear Stability Belt – 1" Webbing W/ Side Release Buckle	(K0108)	\$ 48	Each
MOD-8K	Elbow Blocks w/ L Bracket _____ L x _____ H	(K0108)	\$287	Pair
MOD-8O	Elbow Pads - Hook & Loop Mounting _____ L x _____ H	(K0108)	\$134	Pair
MOD-8OW	Elbow Pads - Wood Base _____ L x _____ H	(K0108)	\$287	Pair
	Upholstery _____			
MOD-8V	Wrist Straps - Hook & Loop	(K0108)	\$ 66	Pair

Notes / Diagrams:


Simply Seating Order Form

Freedom Select Anterior Trunk Supports

Select Standard



Select Narrow



Dimensions

Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SNZ-S	Select Standard Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-1-S	Select <i>Standard</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-M	Select <i>Standard</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-1-M	Select <i>Standard</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-L	Select <i>Standard</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-1-L	Select <i>Standard</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-XL	Select <i>Standard</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6SNZ-1-XL	Select <i>Standard</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235


	SFDI-6NNZ-XS	Select <i>Narrow</i> X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-1-XS	Select <i>Narrow</i> Dual Rear-Pull X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-S	Select <i>Narrow</i> Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-1-S	Select <i>Narrow</i> Dual Rear-Pull Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-M	Select <i>Narrow</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-1-M	Select <i>Narrow</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-L	Select <i>Narrow</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-1-L	Select <i>Narrow</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-XL	Select <i>Narrow</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6NNZ-1-XL	Select <i>Narrow</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

Note:


Simply Seating Order Form

Freedom Select Anterior Trunk Supports (continued)

Select Standard Zipper



Select Narrow Zipper



Dimensions

Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SZ-S	Select <i>Standard Zipper</i> Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 167
	SFDI-6SZ-1-S	Select <i>Standard Zipper</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 167
	SFDI-6SZ-M	Select <i>Standard Zipper</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 167
	SFDI-6SZ-1-M	Select <i>Standard Zipper</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 167
	SFDI-6SZ-L	Select <i>Standard Zipper</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SZ-1-L	Select <i>Standard Zipper</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SZ-XL	Select <i>Standard Zipper</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6SZ-1-XL	Select <i>Standard Zipper</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 178
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235


	SFDI-6NZ-M	Select <i>Narrow Zipper</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 167
	SFDI-6NZ-1-M	Select <i>Narrow Zipper</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 167
	SFDI-6NZ-L	Select <i>Narrow Zipper</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 173
	SFDI-6NZ-1-L	Select <i>Narrow Zipper</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 173
	SFDI-6NZ-XL	Select <i>Narrow Zipper</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 179
	SFDI-6NZ-1-XL	Select <i>Narrow Zipper</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 179
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

Note:

Simply Seating Order Form

Freedom Select Anterior Pelvic Supports

Select 2-Point



2-POINT – PADDED – EZ RELEASE PUSH BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2MP-XS	XS	MICRO 1"	1	50	1.75	6	\$ 94
	SFDI-3A2MP-S	S	MICRO 1"	1	55	2.25	7	\$ 94
	SFDI-3A2P-S	S	STD 1.5"	1.5	55	2.25	7	\$ 94
	SFDI-3A2P-M	M	STD 1.5"	1.5	60	2.5	9	\$ 94
	SFDI-3A2P-L	L	STD 1.5"	1.5	65	3	11	\$ 94
	SFDI-3A2LP-L	L	LG 2"	2	65	3	11	\$ 94
	SFDI-3A2LP-XL	XL	LG 2"	2	70	3.25	13	\$ 94
	SFDI-3A2P-CUSTOM		QUOTE					\$ 118

2-POINT – PADDED – SIDE RELEASE BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2MS-XS	XS	MICRO 1"	1	50	1.75	6	\$ 68
	SFDI-3A2MS-S	S	MICRO 1"	1	55	2.25	7	\$ 68
	SFDI-3A2S-S	S	STD 1.5"	1.5	55	2.25	7	\$ 68
	SFDI-3A2S-M	M	STD 1.5"	1.5	60	2.5	9	\$ 68
	SFDI-3A2S-L	L	STD 1.5"	1.5	65	3	11	\$ 68
	SFDI-3A2LS-L	L	LG 2"	2	65	3	11	\$ 68
	SFDI-3A2LS-XL	XL	LG 2"	2	70	3.25	13	\$ 68
	SFDI-3A2S-CUSTOM		QUOTE					\$ 94

Note:

2-POINT – UNPADDED – EZ RELEASE PUSH BUTTON BUCKLE

Note: Pads not included, must be ordered separately.

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2UMP-XS	XS	MICRO 1"	1	50	1.75	6	\$ 45
	SFDI-3A2UMP-S	S	MICRO 1"	1	55	2.25	7	\$ 45
	SFDI-3A2UP-S	S	STD 1.5"	1.5	55	2.25	7	\$ 45
	SFDI-3A2UP-M	M	STD 1.5"	1.5	60	2.5	9	\$ 45
	SFDI-3A2UP-L	L	STD 1.5"	1.5	65	3	11	\$ 45
	SFDI-3A2ULP-L	L	LG 2"	2	65	3	11	\$ 45
	SFDI-3A2ULP-XL	XL	LG 2"	2	70	3.25	13	\$ 45
	SFDI-3A2UP-CUSTOM		QUOTE					\$ 57

Note:

2-POINT ANTERIOR PELVIC SUPPORT PADS

Qty.	PRICE EA.	Qty.	PRICE PR.
SFDI-3HBP-XS	\$ 32		\$ 49
SFDI-3HBP-S	\$ 32		\$ 49
SFDI-3HBP-M	\$ 32		\$ 49
SFDI-3HBP-L	\$ 32		\$ 49
SFDI-3HBP-XL	\$ 32		\$ 49
SFDI-3HBP-CUSTOM	\$ 44		\$ 73


Simply Seating Order Form

Freedom Select Anterior Pelvic Supports (continued)

2-POINT – UNPADDED – SIDE RELEASE BUTTON BUCKLE

Note: Pads not included, must be ordered separately.

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price	2-POINT ANTERIOR PELVIC SUPPORT PADS					
									Qty.	PRICE EA.	Qty.	PRICE PR.		
	SFDI-3A2UMS-XS	XS	MICRO 1"	1	50	1.75	6	\$ 40						
	SFDI-3A2UMS-S	S	MICRO 1"	1	55	2.25	7	\$ 40	SFDI-3HBP-XS		\$ 32			\$ 49
	SFDI-3A2US-S	S	STD 1.5"	1.5	55	2.25	7	\$ 40	SFDI-3HBP-S		\$ 32			\$ 49
	SFDI-3A2US-M	M	STD 1.5"	1.5	60	2.5	9	\$ 40	SFDI-3HBP-M		\$ 32			\$ 49
	SFDI-3A2US-L	L	STD 1.5"	1.5	65	3	11	\$ 40	SFDI-3HBP-L		\$ 32			\$ 49
	SFDI-3A2ULS-L	L	LG 2"	2	65	3	11	\$ 40	SFDI-3HBP-XL		\$ 32			\$ 49
	SFDI-3A2ULS-XL	XL	LG 2"	2	70	3.25	13	\$ 40	SFDI-3HBP-CUSTOM		\$ 44			\$ 73
	SFDI-3A2US-CUSTOM		QUOTE					\$ 51						



4-POINT – EZ RELEASE PUSH BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	E	F	Price
	SFDI-3A4MP-XS	XS	MICRO 1"	1	1	50	14	4.25	5.75	\$203
	SFDI-3A4MP-S	S	MICRO 1"	1	1	55	16	5	6.5	\$203
	SFDI-3A4P-S	S	STD 1.5"	1.5	1	55	16	5	6.5	\$203
	SFDI-3A4P-M	M	STD 1.5"	1.5	1	60	18	5.5	7.5	\$203
	SFDI-3A4P-L	L	STD 1.5"	1.5	1.5	65	20	6	8.75	\$203
	SFDI-3A4P-XL	XL	STD 1.5"	1.5	1.5	70	22	6	8.75	\$203
	SFDI-3A4LP-XL	XL	LG 2"	2	1.5	70	22	6	8.75	\$203
	SFDI-3A4P-CUSTOM		QUOTE							\$235

4-POINT – SIDE RELEASE BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	E	F	Price
	SFDI-3A4MS-XS	XS	MICRO 1"	1	1	50	14	4.25	5.75	\$202
	SFDI-3A4MS-S	S	MICRO 1"	1	1	55	16	5	6.5	\$202
	SFDI-3A4S-S	S	STD 1.5"	1.5	1	55	16	5	6.5	\$202
	SFDI-3A4S-M	M	STD 1.5"	1.5	1	60	18	5.5	7.5	\$202
	SFDI-3A4S-L	L	STD 1.5"	1.5	1.5	65	20	6	8.75	\$202
	SFDI-3A4S-XL	XL	STD 1.5"	1.5	1.5	70	22	6	8.75	\$202
	SFDI-3A4LS-XL	XL	LG 2"	2	1.5	70	22	6	8.75	\$202
	SFDI-3A4S-CUSTOM		QUOTE							\$235

Note:

Simply Seating Order Form



Cell Size

1"

2"

3"

4"

Deleted Cells
Base Intact

X

Cut Out Cells
& Base

C

Valve Hose
Location
Corner & Cell
Height

V1

V2

V3

V4

Select Seat
Choice

Standard

☐

Asymmetrical
Left

☐

Asymmetrical
Right

☐

Other

☐

MOD-1ROHO SEAT

Freedom Seat ROHO Insert

\$924.00

Note: Retail pricing includes any combination of specification option.

FRONT OF INSERT COLUMNS

INCHES	3	4.75	6.5	8	9.75	11.5	13.25	15	16.5	18.25	20	21.75
CELLS	1	2	3	4	5	6	7	8	9	10	11	12
	A	B	C	D	E	F	G	H	I	J	K	L

3	1	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
4.75	2	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
6.5	3	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
8	4	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
9.75	5	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
11.5	6	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
13.25	7	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
15	8	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
16.5	9	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
18.25	10	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
20	11	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐
21.75	12	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐	☐

INCHES ARE FROM BASE EDGE TO EDGE OF INSERT

1. Start diagram at cell 1A which represents front left of cushion from client's seated position.
2. Specify "Cell Height" (1", 2", 3" or 4") by placing a number in each cell box.
3. Specify request to "Delete Cells with Base Intact" by typing an "X" in those cell boxes.
4. Specify request to "Cut Out Cells & Base" by typing a "C" in those cell boxes.
5. Specify "Valve Hose Location Corner & Cell Height" by typing "V1", "V2", "V3" or "V4" in that corner box.

Note: For multiple valve or custom valve locations and custom air chamber divisions call Freedom Designs Customer Service for assistance at 800-331-8551

PLEASE ALLOW 3-4 WEEKS FOR DELIVERY

Simply Seating Order Form



MOD-2ROHO BACK

Freedom Back ROHO Insert

\$924.00

Note: Retail pricing includes any combination of specification option.

Cell Size

☐ 1"

☐ 2"

☐ 3"

Deleted Cells Base Intact

☐ X

Cut Out Cells & Base

☐ C

Valve Hose Location Corner & Cell Height

☐ V1

☐ V2

☐ V3

Select Back Choice

☐ Standard

☐ "I" Back

☐ "T" Back

		TOP OF INSERT COLUMNS											
INCHES		3	4.75	6.5	8	9.75	11.5	13.25	15	16.5	18.25	20	21.75
CELLS		1	2	3	4	5	6	7	8	9	10	11	12
		A	B	C	D	E	F	G	H	I	J	K	L
R O W S	3	1											
	4.75	2											
	6.5	3											
	8	4											
	9.75	5											
	11.5	6											
	13.25	7											
	15	8											
	16.5	9											
	18.25	10											
	20	11											
	21.75	12											

INCHES ARE FROM BASE EDGE TO EDGE OF INSERT

1. Start diagram at cell 1A which represents top right of back from client's seated position.
2. Specify "Cell Height" (1", 2", or 3") by placing a number in each cell box.
3. Specify request to "Deleted Cells with Base Intact" by typing an "X" in those cell boxes.
4. Specify request to "Cut Out Cells & Base" by typing a "C" in those cell boxes.
5. Specify "Valve Hose Location Corner & Cell Height" by typing "V1", "V2", or "V3" in that corner box.

Note: For multiple valve or custom valve locations and custom air chamber divisions call Freedom Designs Customer Service for assistance at 800-331-8551

PLEASE ALLOW 3-4 WEEKS FOR DELIVERY