



Yes, you can.®



PRO-CG™ Adult Manual Wheelchair w/Matrx

Thank you for using Freedom Designs / Invacare eForms

Invacare is pleased to offer you an upgraded solution to your ordering process. Our enhanced order forms allow you to fill out a form electronically, print and fax the form, save and email* it to Customer Service, or maintain the business practices that work for you today. The format has been revised to reveal a cleaner look with electronic selection and input functions.

For additional information please visit www.freedomdesigns.com



Adobe Acrobat Reader DC

* Interactive functions of our new forms work best with the latest version of Adobe Acrobat Reader DC visit <https://get.adobe.com/reader/> to download and install on your PC or Mac or visit Google Play/ iTunes to download the Adobe Acrobat Reader DC app for your device.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-331-8551



Account Information

Request Type: Quote

Order

Date: _____

Account #: _____

Company: _____

SHIP TO Name: _____

Address: _____

City: _____ State: _____ Zip: _____

Country: _____

Purchase Order #: _____

CONTACT Name: _____

Back Up Contact: _____

Phone: _____

Email: _____

Comments / Marked For: _____

MF: Please do not provide with a customer name

Base Options

A seating system must be purchased separately. No seating system is provided with the PRO-CG frame. No seat or back upholstery is available as an option.

TRANSPORT OPTIONAL

PRO-30310 PRO-CG Adult Rigid\$3,215 (E1161)

Transportation Option selected separately below. (Includes Future Fit Program)

PRO-99998 Transport Brackets.....\$323

Must be installed by Freedom Designs manufacturing. Freedom's combination of customized seating and wheelchair frames are WC19 crash test approved in accordance with RESNA WC-4: 2012 Section 19, Annexes A, B, C & E as well as ISO 7176-19 (2008) Annexes A, C & D. Most Freedom seating, hardware, and positioning accessories have been WC20 tested and approved for transit.

TRANSPORT INCLUDED

PRO-30311 PRO-CG Adult Rigid.....\$3,548 (E1161)

Frame Dimensions

PRO-30240 16" Base Width

Specify Seat Width Setting

PRO-14W 14" n/c

PRO-15W 15" n/c

PRO-16W 16" n/c
Seat frame adjusts 14" to 16" wide in 1" increments.

Specify Seat Depth Setting

PRO-30506 14" n/c

PRO-30521 15" n/c

PRO-30522 16" n/c

PRO-30523 17" n/c

PRO-30524 18" n/c

PRO-30527 19" n/c
2" extension requiredPRO-30528 20" \$592 (E2203)
2" extension requiredPRO-30529 21" \$592 (E2203)
4" extension requiredPRO-30535 22" \$1,065 (E2204)
4" extension required

PRO-30260 18" Base Width

Specify Seat Width Setting

PRO-16W 16" n/c

PRO-17W 17" n/c

PRO-18W 18" n/c
Seat frame adjusts 16" to 18" wide in 1" increments.

Specify Seat Depth Setting

PRO-30522 16" n/c

PRO-30523 17" n/c

PRO-30524 18" n/c

PRO-30527 19" n/c

PRO-30528 20" \$592 (E2203)

PRO-30529 21" \$592 (E2203)
2" extension requiredPRO-30535 22" \$1,065 (E2204)
4" extension requiredPRO-30536 23" \$1,065 (E2204)
4" extension requiredPRO-30537 24" \$1,065 (E2204)
4" extension required

**PRO-30261 20" Base Width****Specify Seat Width Setting**

PRO-18W 18"	n/c
PRO-19W 19"	n/c
PRO-20W 20"	\$538 (E2201)
Seat frame adjusts 18" to 20" wide in 1" increments.	

Specify Seat Depth Setting

PRO-30524 18"	n/c
PRO-30527 19"	\$590
PRO-30528 20"	\$592 (E2203)
PRO-30529 21"	\$592 (E2203)
PRO-30535 22"	\$1,065 (E2204)
PRO-30536 23"	\$1,065 (E2204)
2" extension required	
PRO-30537 24"	\$1,065 (E2204)
2" extension required	
PRO-30538 25"	\$1,065 (E2204)
4" extension required	
PRO-30539 26"	\$1,065 (E2204)
4" extension required	

SEAT RAIL EXTENSIONS

PRO-30525 2" Extension - Pair.....	\$45 (K0108)
PRO-30531 4" Extension - Pair.....	\$45 (K0108)

Specify Back Post Location (if other than specified above) _____ "

Seat Rail Options**LEAF - LOWER EXTREMITY
ACCOMMODATION FRAME**

Not available with ALL Contracture & 3D Dynamic footrest systems.

PRO-30324 LEAF Seat Frame- Pair.....	\$834 (K0108)
--------------------------------------	---------------

Front frame angle adjusts 22° from midline for abduction, adduction, and windswept lower extremity positioning. Seat frame(s) can be individually repositioned for up to a 2" leg length discrepancy.

SRO - SEAT RAIL ONLY

PRO-30325 SRO Seat Frame- Pair.....	\$350 (K0108)
-------------------------------------	---------------

SRO is manufactured without standard footrest hanger attachments.

Vent Options

PRO-30335 Vent Tray Base W/ Tray.....	\$634 (E1029)
---------------------------------------	---------------

Extended wheelchair base including vent tray. For LP10 style vent, base must be 16" or greater. Tray measures full width of the base frame and 12" deep

11-1100 Back Post Mounted Vent Bracket.....	\$633 (E1029)
---	---------------

For LTV vents ONLY. Does NOT include dovetail bracket. Weight limit 15 lbs.

Auxillary Device Mounts

11-3100 O2 Holder - Aluminum w/ QR Hdwe....	\$631 (K0108)
---	---------------

11-4000 IV Pole w/ Mounting Hardware.....	\$311 (K0108)
---	---------------

Tilt Activations

PRO-30272 Dual Lever Hand Activated Tilt.....	n/c
---	-----

PRO-30273 Foot Activated Tilt.....	n/c
------------------------------------	-----


Not available with vent tray.

Frame Colors (Includes Seat Rail, Caster Housing, and Hangers)

Cherry Flame..... Glossy	Midnight Teal..... Glossy
Classic Black..... G / S / M	Minty Fresh.....G / S / M
Electric Pink.....G / S / M	Pure Snow.....G / S
Forest Green..... G / S	Radiant Ruby..... G / S
Granite Grey..... G / S / M	Sunburst Yellow.....G / S / M
Lavender Mist.....G / S / M	True Blue.....G / S / M

FINISH - NOTE-Default finish is Glossy

Glossy Matte Sparkle



Base Colors

PRO-BBASE Black Base.....n/c
PRO-CBASE Matching Base and Frame Color.....\$100
PRO-CACCENT Accent Base Color.....\$150

Back Posts

BACK POST OPTIONS

PRO-30545 Fold-Down Back Posts.....n/c
PRO-30562 Reduced Shear Recline Back.....\$1,500 (E1226)
PRO-30542 Dynamic Rocker Back - Seating Dynamics.....\$2,200 (E2398) Angle adjustable from 72° to 108° in 6° increments. Can be locked-out when needed. Includes a set of 4 elastomers in different durometers to allow spring adjustment in the field.

PRO-30564 Reduced Shear Recline Back with Angle-
Adjustable Stroller Handle Extension.....\$1,712

Recline range is 65° from 97° to 162°. Pivot point is 2.5" above the seat rail. Solid Back Insert installation requires 4-points of attachment hardware. Seat-to-Back hinges are not compatible with any Super Reclining Reduced Shear Back.

PUSH HANDLE OPTIONS

PRO-30518 90° Bend Back Post Tubes.....n/c

PRO-30544 Straight Back Post Tubesn/c

Requires selection of Push Button Angle Adjustable Stroller Handle (PRO-30571 or PRO-30576). Not available with Reduced Shear Recline Back.

PUSH HANDLE HEIGHT

PRO-30515 Short - 20" High.....n/c

PRO-30514 Medium - 23" High.....n/c

PRO-30516 Tall - 26" High.....n/c

Push Handle Extension Options



-30561



-30590



-30571



-30576



-30517

All push handle extensions except PRO-30561 add approximately 6.5" to push handle height. All push handle extensions except Individual include a black neoprene wrap around pad.

PRO-30561 Push Handle Extension - Width Adjustable..n/c

Not available with Straight Back Posts. Not available with Recline.

PRO-30571 Push Button, Angle-Adjustable Stroller Handle Extension.....\$121

Available on wheelchair frames 15W & smaller only. Plastic hinge up to 200° range of adjustment. PRO-30571 or PRO-30576 required with Straight Back Post Tubes (PRO-30544)

PRO-30517 Individual Stroller Handle Extensions Removable - Pair.....\$119

Not available with Straight Back Posts. Not available with Recline.

PRO-30590 Stroller Handle Extension - Width Adjustable, Removable\$103

Not available with Straight Back Posts. Not available with Recline.

PRO-30576 Heavy Duty Push Button, Angle-Adjustable Stroller Handle Extension\$193

Available on any size frame. Metal hinge up to 180° range of adjustment. PRO-30571 or PRO-30576 required with Straight Back Post Tubes (PRO-30544)

Footrests

Footrest / Extension Lengths


Extension/ Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	6"-8"	5.5"-8.5"	4"-8"	4"-8"	6"-12"
High Mount	N/A	N/A	1.5"-3"	N/A	N/A
4.5" H D Ext.	9"-12"	9.5"-12.5"	9"-12"	9"-12"	13"-15"
6.5" H D Ext.	9"-15"	9.5"-15.5"	9"-15"	9"-15"	13"-17"
4" Adjustable	10"-13"	9.5"-12.5"	11"-13"	10"-13"	13"-15"
6" Adjustable	12"-15"	10.5"-14.5"	13"-15"	12"-15"	13"-18"
8" Adjustable	14"-17"	12.5"-16.5"	N/A	N/A	14"-20"
10" Adjustable	16"-19"	14.5"-18.5"	N/A	N/A	16"-22"



MAPS Footrest / Extension Lengths

Extension/ Hanger	75° Pop	75° Lch	90° Pop	90° Lch	ELR (90°)
Hanger Only	5.5"-7.5"	5"-8"	5.5"-7.5"	4.5"-7.5"	4.5"-8.5" & 11"
High Mount	N/A	N/A	N/A	N/A	N/A
4.5" H D Ext.	8.5"-11.5"	9"-12"	8.5"-11.5"	8.5"-11.5"	12.5"-14.5"
6.5" H D Ext.	8.5"-13.5"	9"-14"	8.5"-13.5"	8.5"-13.5"	12.5"-16.5"
4" Adjustable	8.5"-11.5"	9"-12"	9.5"-11.5"	8.5"-11.5"	12.5"-14.5"
6" Adjustable	10.5"-13.5"	10"-14"	11.5"-13.5"	8.5"-13.5"	12.5"-17.5"
8" Adjustable	12.5"-15.5"	12"-16"	N/A	N/A	13.5"-19.5"
10" Adjustable	14.5"-17.5"	14"-18"	N/A	N/A	15.5"-21.5"

Length is measured from top of footrest plate to top of seat rail. Without extensions bottom hanger hole not accessible for hardware. Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy duty extension tubes.





HANGERS

PRO-30280 75° Pop-Up Swing Away Hangers.....n/c

PRO-30281 90° Pop-Up Swing Away Hangers.....n/c
Not available with 8" or 10" Telescoping Extensions.

PRO-30279 90° Pop-Up Swing Away Elevating
Legrests (ELRs).....\$318 (E099)
Not available with all Telescoping extensions.

PRO-30282 75° Latch Style Swing Away Hangers.....n/c

PRO-30283 90° Latch Style Swing Away Hangers.....n/c
Not available with 8" or 10" Telescoping Extensions.

PRO-33000 Omit Legrest Hangers & Footrests.....n/c

FOOTPLATES & FOOTBOARDS

PRO-36033 Aluminum Angle Adj. Flip-up
Footplates.....\$205 (K0040)
Does not include Heel Loops, select below.

PRO-17070 Composite Flip-Up.....n/c

PRO-36029 Aluminum Angle
Adjustable Footboard Removable\$205 (K0040)

PRO-36404 Omit Footrest Plates.....n/c

FOOTPLATE OPTIONS

PRO-36040 MAPS (Multi Axis Positioning System)\$350
Requires selection of Aluminum Angle Adjustable Footplates. Not compatible
with
Composite footplates or Aluminum Footboard. Not compatible with Seating
Dynamics hangers.

PRO-36022 Heel Loops.....\$49 (E0951)

PRO-36024 Calf Strap - Webbing.....\$49 (K0038)

HANGER EXTENSION OPTIONS

PRO-30285 4" Telescoping Extensions.....n/c

PRO-30286 6" Telescoping Extensions.....n/c

PRO-30287 8" Telescoping Extensions.....n/c
Not available with all 90° hangers.

PRO-30288 10" Telescoping Extensions.....n/c
Not available with all 90° hangers.

PRO-36402 4.5" Heavy Duty Extensions.....\$63

PRO-36403 6.5" Heavy Duty Extensions.....\$74

Heavy duty extensions do not telescope adjust. Footrests adjust vertically in length on heavy
duty extension tubes.

PRO-36410 High Mount Clamp & Tubes\$110

Clear Selection

SPECIALIZED FOOTRESTS

ALL legrest measurements are from the seat rail to footplate. Not available with LEAF or SRO seat rail options. Frame and/or front caster choice may limit adjustment range and/or cause interference.

Specialized Footrest Seat Rail to Footplate Measurements

Length measured from top of footrest plate to top of seat rail.

MUST BE COMPLETED FOR CORRECT SELECTION OF COMPONENTS

PRO - LL Back of the Knee (to Heel): _____

PRO - CT Cushion Thickness: _____

CONTRACTURE FOOTREST HANGERS

ALL legrest measurements are from the seat rail to footplate. Not available with LEAF or SRO seat rail options. Frame and/or front caster choice may limit

PRO-30291-2 Contracture Short - 11" to 15" \$3838

Includes short hanger extension.

PRO-30292-2 Contracture Long - 15" to 19" \$3838

Includes long hanger extension.



3D DYNAMIC FOOTREST HANGERS

PRO-30299-2 3D Dynamic Footrest - \$7800

SPECIALIZED FOOTPLATES

Must select a footplate and extension

PRO-6503 4" x 6" HD Aluminum Footplates n/c

PRO-6504 5" x 8" HD Aluminum Footplates n/c

For frames 16" wide and greater only.

PRO-6505 6" x 10" HD Aluminum Footplates n/c

For frames 18" wide and greater only. Not available with Contracture Hangers



Armrests

BACK POST MOUNTED ARMRESTS

Omit Armrests

Foam wrap pad STD on all back post mounted armrests.

PRO-38204 10.5" Flip-Backn/c

PRO-38202 13.5" Flip-Backn/c

PRO-38210 10.5" to 13.5" Flip-Back
Locking, Telescoping.....\$307 (E0973)

Full length pad on 13.5" or longer armrests only.

PRO-38201 10.5" Flip-Back Locking.....\$285 (E0973)

PRO-38200 13.5" Flip-Back Locking\$285 (E0973)

PRO-38215 17" to 20" Flip-Back
Locking, Telescoping.....\$307 (E0973)

Armrest Heights

Armrest and Pad	Foam	Desk	Full
Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Telescoping Flip Back	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Flip Back Locking Telescoping	6.5" - 11.5"	7.5" - 12.5"	7.5" - 12.5"
Single Post Height Adj. Short		8" - 14"	8" - 14"
Single Post Height Adj. Tall		10" - 16"	10" - 16"
Dual Post Fixed Height		9.5"	9.5"
Dual Post Flip Back Height Adj.		9" - 13"	9" - 13"

All heights from top of seat frame to top of armrest pad.



FRAME MOUNTED ARMRESTS

PRO-38217 Single Post Height Adjustable,
Removable - Low\$307 (E0973)

PRO-38203 Single Post Height Adjustable,
Removable - Standard.....\$307 (E0973)

PRO-38206 Dual Post Flip-Back Height
Adjustable, Removable\$318 (E0973)

PRO-38205 Dual Post Fixed Height Removablen/c

Not available with short frame. 22" & 24" Wheels not available with dual post armrests

ARM PAD OPTIONS

PRO-38213 Desk Length Padsn/c

PRO-38214 Full Length Padsn/c
13.5" or longer armrests only.

Casters

PRO-40000 4" x 1" Polyurethane.....n/c

PRO-43000 4" x 1.4" Poly, Aluminum Mag.....\$180

PRO-34555 5" x 1" Polyurethane.....n/c

PRO-34222 5" x 1.5" Polyurethane.....\$94

PRO-35050 5" Polyurethane Light-Up.....\$110

PRO-42000 5" x 1.4" Poly, Aluminum Mag.....\$180

PRO-34666 6" x 1" Polyurethane.....n/c

PRO-34999 6" x 1.5" Polyurethane.....\$94

PRO-44000 6" x 1.4" Poly, Aluminum Mag\$180

PRO-34444 6" x 2" Polyurethane, Treaded.....\$105

PRO-34777 7" x 1" Polyurethane.....n/c

PRO-34200 7" x 2" Foam Filled.....\$215

Wheel Options	12"	16"	18"	20"	22"	24"
MAGS FOR POLYURETHANE TIRES - FULL PROFILE	X	X	X	X	X	X
MAGS FOR PNEUMATIC TIRES - <i>Compatible with all pneumatic tire options</i>	X	X	X	X	X	X
Spoke Wheel - Black - <i>Compatible with all pneumatic tire options</i>		X	X			
Polyurethane Tires - Black, Full Profile	X	X	X	X	X	X
Pneumatic Tires - Black	X	X	X	X	X	X
Pneumatic Tires - Gray	X	X	X	X	X	X
Airless Inserts - <i>Compatible with all pneumatic tire options</i>	X	X	X	X	X	X

Rear Wheels

MAGS FOR **PNEUMATIC** TIRES

PRO-31525 12" Mag Wheel.....	n/c
12" mag wheels not available with 7" casters.	
PRO-31530 16" Mag Wheel.....	\$45
PRO-31533 18" Mag Wheel.....	\$45
PRO-31527 20" Mag Wheel - Aluminum Handrim	n/c
PRO-31528 22" Mag Wheel - Aluminum Handrim.....	n/c
PRO-31529 24" Mag Wheel - Aluminum Handrim.....	n/c

MAGS FOR **POLYURETHANE** TIRES - FULL PROFILE

PRO-31520 12" Wheel for Full Profile Poly.....	n/c
PRO-31522 16" Wheel for Full Profile Poly.....	n/c
PRO-31521 18" Wheel for Full Profile Poly.....	n/c
PRO-31519 20" Wheel for Full Profile Poly.....	n/c
PRO-31518 22" Wheel for Full Profile Poly.....	n/c
PRO-31517 24" Wheel for Full Profile Poly.....	n/c

SPOKE

PRO-31535 16" Spoke Wheel - Black w/ Black Tire.....	n/c
PRO-31538 18" Spoke Wheel - Black w/ Black Tire.....	n/c

HANDRIM OPTIONS

PRO-31005 Aluminum Coated Handrims.....	\$130
---	-------

PRO-31006 Omit Handrims
20", 22" & 24" Mag wheels only

PNEUMATIC TIRE OPTIONS

PRO-31543 Black Pneumatic Tire.....	n/c
PRO-31544 Black Pneumatic Tire with Airless Inserts.....	\$92 (E2213)
PRO-31040 Gray Pneumatic Tire.....	n/c (E2211, E2212)
PRO-31000 Gray Pneumatic Tire with Airless Inserts.....	\$92 (E2213)

POLYURETHANE TIRE OPTIONS

PRO-31516 Full Poly Black - 12", 18", 20", 22", 24"	n/c
PRO-31515 Full Poly Black - 16"	\$80

Wheel/Caster Seat Rail to Floor Heights

Rear Wheel	Caster Wheel					
	4"	5"	6"	6 x 2"	7"	7 x 2"
12"	15"	15"	18.5"	18.5"	18.5"	18.5"
	15.5"	15.5"	19"	19"	19"	19"
	16"	16"				
	18.5"	18.5"				
	19"	19"				
16"	15"	15"	16"	16"	17"	17"
	15.5"	15.5"	16.5"	16.5"	17.5"	17.5"
	16"	16"	17"	17"	20"	20"
	16.5"	16.5"	17.5"	17.5"	20.5"	20.5"
	17"	17"	20"	20"		
18"	17.5"	17.5"				
	15"	15.5"	16"	16"	17"	17"
	15.5"	16"	16.5"	16.5"	17.5"	17.5"
	16"	16.5"	17"	17"	18"	18"
	16.5"	17"	17.5"	17.5"	18.5"	18.5"
20"	17"	17.5"	18"	18"		
	17.5"	18"	18.5"	18.5"		
	18"	18.5"	19"	19"		
	18.5"	19"	19.5"	19.5"		
	19"	19.5"	20"	20"		
22"						
	18"	18"	18"	18"	18"	18"
	18.5"	18.5"	18.5"	18.5"	18.5"	18.5"
	19"	19"	19"	19"	19"	19"
		19.5"	19.5"	19.5"	20"	20"
24"						
	18.5"	18.5"	19"	19"	19"	19"
	19"	19"	19.5"	19.5"	20"	20"
			20"	20"	20.5"	20.5"
Clear Selection						

Use chart to determine seat rail to floor height. Caster sizes may vary slightly in height by type. Seat rail to floor heights are approximate, and are impacted by wheel type, caster size and braking choice. Heights can vary +/- 1/4". Seat pan choice increases seat rail to floor height 3/8" to 3/4".



Wheel Locks

STANDARD

PRO-31501 Push to Lock.....n/c

PRO-31506 Pull to Lock.....n/c

WHEEL LOCK OPTIONS

PRO-31801 Extension Handles..... \$81 (E0961)

DRUM LOCK

Includes manual wheel locks to be used with public transportation. Not available with 12" mag wheels with low pro poly tires. Not Available with 20", 22" & 24" Spoke Wheels.

PRO-31508 Foot Activated Wheel Lock.....n/c
Not available with ALL vent and battery trays.

PRO-30275 Hand Activated Brake
and Wheel Lock System \$1,000 (E2228)

Anti-Tippers

PRO-31600 Curb Assist Anti-Tip
Assembly Flip-Up, Adjustable.....\$135 (E0971)
Pair includes one anti-tip tube with curb assist feature.

PRO-31605 Spring Loaded Swing-Away
Anti-Tip Assy - Foot Activated..... \$194 (E0971)

Seat Pans

Only choose one Seat Pan option

ALL seat pans include 2" x 2" notches for growth adjustment.

NON-ADJUSTABLE Fixed Seat Pan raises the effective seat to floor height by 1/2". ALL Composite seat pans are made to fit frame dimensions unless specified below.

PRO-39010 Mounted Solid Seat Pan - Composite..\$183 (E2231)

Seat Width _____" X Seat Depth _____"

PRO-39030 Mounted Solid Seat Pan - Aluminum.... \$215 (E2231)

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Available in the following sizes ONLY:

10" wide x 12" deep

12" wide x 14" deep

14" wide x 16" deep

16" wide x 18" deep

18" wide x 18" deep

ADJUSTABLE

PRO-39040 Fixed, Width Adjustable

Solid Seat Pan - Aluminum..... \$312 (E2231)

Fixed Seat Pan raises the effective seat to floor height by 1/2".

Width Adjustable Solid Seat Pan adjusts 2" in width (in 1" increments) by assembly hardware to grow with PRO frame.

Available in the following sizes ONLY:

10"- 12" wide x 12" deep

12" - 14" wide x 14" deep

14" - 16" wide x 16" deep

16" - 18" wide x 18" deep

MOUNTING HARDWARE OPTIONS

Use to convert fixed seat pan to removable.

PRO-DROP1 1" Drop Pan Hardware Kit..... \$223 (K0108)

PRO-EZLITE Removable Pan Hardware Kit..... \$223 (K0108)

Removable Seat Pan Hardware raises the effective seat to floor height by 5/8". Includes 4-points of EZ-Lite mounting hardware, 2 each frame stopper clamps.

PRO-DROP2 2" Drop Pan Hardware Kit..... \$223 (K0108)

Includes 4 each 2" drop hooks, 4 each 408 clamps, & 2 each Glide Lock clamps. Not available with ANY Aluminum Solid Seat Pan.

MATRIX SEAT CUSHIONS




Matrix PS Seat Cushion 14" - 20" wide/reversable cover

HCPCS code E2605.....

\$345

- | | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ PS1414 | ☐ PS1416 | ☐ PS1616 | ☐ PS1618 | ☐ PS1620 | ☐ PS1716 |
| ☐ PS1718 | ☐ PS1720 | ☐ PS1816 | ☐ PS1818 | ☐ PS1820 | ☐ PS2016 |
| ☐ PS2018 | ☐ PS2020 | | | | |



Matrix PSP Seat Cushion 14" - 20" wide/reversable cover

HCPCS code E2607.....

\$395

- | | | | | | |
|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|----------------------------------|
| ☐ PSP1414 | ☐ PSP1416 | ☐ PSP1616 | ☐ PSP1618 | ☐ PSP1620 | ☐ PSP1716 |
| ☐ PSP1718 | ☐ PSP1720 | ☐ PSP1816 | ☐ PSP1818 | ☐ PSP1820 | ☐ PSP2016 |
| ☐ PSP2018 | ☐ PSP2020 | | | | |



Matrix Vi Seat Cushion 14" - 20" wide

HCPCS code E2607.....

\$425

- | | | | | | |
|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|------------------------------------|
| ☐ KA1414-VI | ☐ MA1415-VI | ☐ MA1417-VI | ☐ MA1616-VI | ☐ MA1618-VI | ☐ MA1620-VI |
| ☐ MA1716-VI | ☐ MA1718-VI | ☐ MA1720-VI | ☐ MA1816-VI | ☐ MA1818-VI | ☐ MA1820-VI |
| ☐ MA2016-VI | ☐ MA2018-VI | ☐ MA2020-VI | | | |



Matrix Libra Seat Cushion 14" - 20" wide

HCPCS code E2624.....

\$485

- | | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC1414 | ☐ LC1416 | ☐ LC1516 | ☐ LC1518 | ☐ LC1520 | ☐ LC1616 |
| ☐ LC1618 | ☐ LC1620 | ☐ LC1716 | ☐ LC1718 | ☐ LC1720 | ☐ LC1816 |
| ☐ LC1818 | ☐ LC1820 | ☐ LC1916 | ☐ LC1918 | ☐ LC1920 | ☐ LC2016 |
| ☐ LC2018 | ☐ LC2020 | | | | |

Matrix Libra Seat Cushion 22" deep

HCPCS code E2624.....

\$595

- | | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ LC1622 | ☐ LC1722 | ☐ LC1822 | ☐ LC1922 | ☐ LC2022 | ☐ LC2122 |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|



Matrix Libra Adjust Seat Cushion 14" -20" wide

HCPCS code E2624.....

\$625

- | | | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC1414-ADJ | ☐ LC1416-ADJ | ☐ LC1516-ADJ | ☐ LC1518-ADJ | ☐ LC1520-ADJ | ☐ LC1616-ADJ |
| ☐ LC1618-ADJ | ☐ LC1620-ADJ | ☐ LC1716-ADJ | ☐ LC1718-ADJ | ☐ LC1720-ADJ | ☐ LC1816-ADJ |
| ☐ LC1818-ADJ | ☐ LC1820-ADJ | ☐ LC1916-ADJ | ☐ LC1918-ADJ | ☐ LC1920-ADJ | ☐ LC2016-ADJ |
| ☐ LC2018-ADJ | ☐ LC2020-ADJ | | | | |



Matrix Libra Adjust Seat Cushion 22" deep

HCPCS code E2624.....

\$735

- | | | | | | |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|
| ☐ LC1622-ADJ | ☐ LC1722-ADJ | ☐ LC1822-ADJ | ☐ LC1922-ADJ | ☐ LC2022-ADJ | ☐ LC2122-ADJ |
|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|-------------------------------------|



Matrix Multi Cushion 14" - 20" wide/reversable cover

HCPCS code E2607.....

\$445


- | | | | | | |
|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|---------------------------------|
| ☐ MM1414 | ☐ MM1416 | ☐ MM1616 | ☐ MM1618 | ☐ MM1620 | ☐ MM1716 |
| ☐ MM1718 | ☐ MM1720 | ☐ MM1816 | ☐ MM1818 | ☐ MM1820 | ☐ MM2016 |
| ☐ MM2018 | ☐ MM2020 | | | | |

FABRIC OPTIONS FOR CUSHIONS

Cushion	Meshtex	Startex	Infection	Spacetex
Libra/Adjust	N/A	STD	\$50	N/C
PS/PSP	N/A	STD Revers-ible	STD Revers-ible	N/A
Vi	STD	\$50	\$50	N/A

- ☐ SFLF Spacetex Fabric Option for Libra Cushion.....N/C
- ☐ ICFL Infection Control Fabric Upgrade for Libra Cushion\$50
- ☐ SFVS Smooth side Startex Fabric Upgrade for Matrix Vi Cushion\$50
- ☐ SFVF Fabric side Startex Fabric Upgrade for Matrix Vi Cushion\$50

LIBRA CUSHION COMPONENTS



- ☐ LAB Libra Abductor.....\$60.00
- ☐ FT Anterior Seat Wedges.....\$125.00
- ☐ FTHW Hemi Wedge.....\$65.00
- ☐ FTTW-LFRR Tapered Wedge LFRR.....\$50.00
- ☐ FTTW-RFLR Tapered Wedge RFLR.....\$50.00
- ☐ LFS Libra Fluid Supplement (each).....\$75.00
- ☐ LLW Lateral Wedge\$75.00
- ☐ LFK Complete Libra Fit Kit.....\$435.00
*Components:
 LAB, FTSW, FTHW, FTTW-LFRR, FTTW-RFLR, LFS x 2*

PELVIC OBLIQUITY KIT

- ☐ POK Pelvic Obliquity Kit for PS / PSP / Multi / Vi cushions\$50



Elite E2 **Standard Back:** 3" Contour



BACK WIDTH AVAILABLE

BACK LENGTH

SIZE	14"	15"	16"	17"	18"	19"	20"	MSRP	HCPCS
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"	\$625.00	E2615
14"	☐ E2S1414	☐ E2S1514	☐ E2S1614	☐ E2S1714	☐ E2S1814	☐ E2S1914	☐ E2S2014		
16"	☐ E2S1416	☐ E2S1516	☐ E2S1616	☐ E2S1716	☐ E2S1816	☐ E2S1916	☐ E2S2016		
18"	☐ E2S1418	☐ E2S1518	☐ E2S1618	☐ E2S1718	☐ E2S1818	☐ E2S1918	☐ E2S2018	\$645.00	E2615
20"	☐ E2S1420	☐ E2S1520	☐ E2S1620	☐ E2S1720	☐ E2S1820	☐ E2S1920	☐ E2S2020		

Elite E2 **LB:** 3" Contour



BACK WIDTH AVAILABLE

BACK LENGTH

SIZE	14"	15"	16"	17"	18"	19"	20"	MSRP	HCPCS
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"	\$725.00	E2620
14"	☐ E2LB1414	☐ E2LB1514	☐ E2LB1614	☐ E2LB1714	☐ E2LB1814	☐ E2LB1914	☐ E2LB2014		
16"	☐ E2LB1416	☐ E2LB1516	☐ E2LB1616	☐ E2LB1716	☐ E2LB1816	☐ E2LB1916	☐ E2LB2016		
18"	☐ E2LB1418	☐ E2LB1518	☐ E2LB1618	☐ E2LB1718	☐ E2LB1818	☐ E2LB1918	☐ E2LB2018	\$745.00	E2620
20"	☐ E2LB1420	☐ E2LB1520	☐ E2LB1620	☐ E2LB1720	☐ E2LB1820	☐ E2LB1920	☐ E2LB2020		

Elite E2 **Deep Back:** 6" Contour



BACK WIDTH AVAILABLE

BACK LENGTH

SIZE	14"	15"	16"	17"	18"	19"	20"	MSRP	HCPCS
FITS W/C	13" - 16"	14" - 17"	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"	\$755.00	E2620
14"	☐ E2D1414	☐ E2D1514	☐ E2D1614	☐ E2D1714	☐ E2D1814	☐ E2D1914	☐ E2D2014		
16"	☐ E2D1416	☐ E2D1516	☐ E2D1616	☐ E2D1716	☐ E2D1816	☐ E2D1916	☐ E2D2016		
18"	☐ E2D1418	☐ E2D1518	☐ E2D1618	☐ E2D1718	☐ E2D1818	☐ E2D1918	☐ E2D2018	\$775.00	E2620
20"	☐ E2D1420	☐ E2D1520	☐ E2D1620	☐ E2D1720	☐ E2D1820	☐ E2D1920	☐ E2D2020		

MAC Dual Layer Foam Back Model



BACK WIDTHS

BACK LENGTH	SIZE	16"	17"	18"	19"	20"	MSRP	HCPCS
	FITS W/C	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"	\$955.00	E2620
	16"	☐ MAC1616	☐ MAC1716	☐ MAC1816	☐ MAC1916	☐ MAC2016		
	18"	☐ MAC1618	☐ MAC1718	☐ MAC1818	☐ MAC1918	☐ MAC2018		
	20"	☐ MAC1620	☐ MAC1720	☐ MAC1820	☐ MAC1920	☐ MAC2020		
	22"	☐ MAC1622	☐ MAC1722	☐ MAC1822	☐ MAC1922	☐ MAC2022		

Lateral Contour Extensions (Included - Please select the size)

CONTOUR	LEFT LATERAL	RIGHT LATERAL
3.5" Pad	☐ MLCS	☐ MRCS
6" Pad	☐ MLCD	☐ MRCD
9" Pad	☐ MLCXD	☐ MRCXD



MAC HUG StaminaFibre® Back Model



BACK WIDTHS

BACK LENGTH	SIZE	16"	17"	18"	19"	20"	MSRP	HCPCS
	FITS W/C	15" - 18"	16" - 19"	17" - 20"	18" - 21"	19" - 22"	\$1055.00	E2620
	16"	☐ MAC1616-HUG	☐ MAC1716-HUG	☐ MAC1816-HUG	☐ MAC1916-HUG	☐ MAC2016-HUG		
	18"	☐ MAC1618-HUG	☐ MAC1718-HUG	☐ MAC1818-HUG	☐ MAC1918-HUG	☐ MAC2018-HUG		
	20"	☐ MAC1620-HUG	☐ MAC1720-HUG	☐ MAC1820-HUG	☐ MAC1920-HUG	☐ MAC2020-HUG		
	22"	☐ MAC1622-HUG	☐ MAC1722-HUG	☐ MAC1822-HUG	☐ MAC1922-HUG	☐ MAC2022-HUG		

Lateral Contour Extensions (Included - Please select the size)

CONTOUR	LEFT LATERAL	RIGHT LATERAL
3.5" Pad	☐ MLCS	☐ MRCS
6" Pad	☐ MLCD	☐ MRCD
9" Pad	☐ MLCXD	☐ MRCXD



Posture Back: 3" Contour



BACK WIDTH AVAILABLE

BACK LENGTH	SIZE	14"	15"	16"	17"	18"	19"	20"	MSRP	HCPCS
	FITS W/C	14" - 15"	15" - 16"	16" - 17"	17" - 18"	18" - 19"	19" - 20"	20" - 21"	\$625.00	E2615
	12"	☐ MPB1412	☐ MPB1512	☐ MPB1612	☐ MPB1712	☐ MPB1812	☐ MPB1912	☐ MPB2012		
	14"	☐ MPB1414	☐ MPB1514	☐ MPB1614	☐ MPB1714	☐ MPB1814	☐ MPB1914	☐ MPB2014		
	16"	☐ MPB1416	☐ MPB1516	☐ MPB1616	☐ MPB1716	☐ MPB1816	☐ MPB1916	☐ MPB2016		
	18"	☐ MPB1418	☐ MPB1518	☐ MPB1618	☐ MPB1718	☐ MPB1818	☐ MPB1918	☐ MPB2018		
	20"	☐ MPB1420	☐ MPB1520	☐ MPB1620	☐ MPB1720	☐ MPB1820	☐ MPB1920	☐ MPB2020		

Posture Deep Back: 3" Contour



BACK WIDTH AVAILABLE

BACK LENGTH	SIZE	14"	15"	16"	17"	18"	19"	20"	MSRP	HCPCS
	FITS W/C	14" - 15"	15" - 16"	16" - 17"	17" - 18"	18" - 19"	19" - 20"	20" - 21"	\$645.00	E2615
	12"	☐ PBD1412	☐ PBD1512	☐ PBD1612	☐ PBD1712	☐ PBD1812	☐ PBD1912	☐ PBD2012		
	14"	☐ PBD1414	☐ PBD1514	☐ PBD1614	☐ PBD1714	☐ PBD1814	☐ PBD1914	☐ PBD2014		
	16"	☐ PBD1416	☐ PBD1516	☐ PBD1616	☐ PBD1716	☐ PBD1816	☐ PBD1916	☐ PBD2016		
	18"	☐ PBD1418	☐ PBD1518	☐ PBD1618	☐ PBD1718	☐ PBD1818	☐ PBD1918	☐ PBD2018		
	20"	☐ PBD1420	☐ PBD1520	☐ PBD1620	☐ PBD1720	☐ PBD1820	☐ PBD1920	☐ PBD2020		

E2S EASYFIX LATERALS HCPCS CODE E0956

- ☐ E2FL Fixed Lateral Supports
Pad 6.5" long x 5" deep (pair) \$290



E2S SWINGAWAY LATERAL HARDWARE

HCPCS code E0128

- ☐ SAHEB Swing-Away Lateral Hardware Only
Universal (each) \$225



E2S OFFSET, SWINGAWAY LATERAL HARDWARE

HCPCS CODE E0128

- ☐ OSAH-L Offset Swing-Away Lateral Hardware
User Left (each) \$225

- ☐ OSAH-R Offset Swing-Away Lateral Hardware
User Right (each) \$255



E2S SWINGAWAY LATERAL PADS HCPCS CODE E0956

- ☐ EPAD-XS Height 3.25", depth 4.25" (each) \$110
☐ EPAD-S Height 3.5", depth 5.25" (each) \$110
☐ EPAD-M Height 3.75", depth 5.75" (each) \$110
☐ EPAD-L Height 4.25", depth 6.5" (each) \$110



E2S/E2D HEADREST INTERFACE

- ☐ E2HI E2 Headrest Interface \$15
Molded lightweight, low profile mounting
plate for headrest attachment;
compatible and tested for Matrix headrest
attachment

- ☐ E2HDI E2 Heavy Duty Headrest Interface \$25
Extra strong steel headrest mount for extreme tone
or heavy use applications; Recommended for use
with non Matrix headrests



E2S/E2D CHEST HARNESS INTERFACE

- ☐ E2CHI-S Fits E2 backs 14, 15" wide \$75
☐ E2CHI-M Fits E2 backs 16-20" wide \$75



E2S/E2D BACKBONE PSIS SUPPORT HCPCS CODE E0128

- ☐ E2BKB Contoured PSIS/Sacrolumbar
Support \$55



POSTUREFLEX ABDOMINAL SUPPORT

- ☐ PFLEX-M Length 6.5", length 16" - 23" \$245
☐ PFLEX-L Length 7", length 20" - 32 \$245



- ☐ QRTAB Quad release Tab \$25
For easy independent application
and removal of PostureFlex



PRIVACY FLAPS

- ☐ PFLAP-XS X-Small Privacy Flap \$40
10" W x 16" L - fits 14" wide backs
- ☐ PFLAP-S Small Privacy Flap \$40
13" W x 16" L - fits 15"-16" wide backs
- ☐ PFLAP-M Medium Privacy Flap \$40
15" W - 16" L - fits 17"-18" wide backs
- ☐ PFLAP-L Large Privacy Flap \$40
17" W - 16" L - fits 19"-20" wide backs



MATRIX HEADREST PADS [HCPCS CODE E0955](#)

Standard Pads - Startex Fabric Cover

Elan HUG Headrest



Elan Headrest



☐ ESP6	Matrix Elan Standard Pad 6"	\$235
☐ ESP10	Matrix Elan Standard Pad 10"	\$235
☐ ESP14	Matrix Elan Standard Pad 14"	\$235
☐ ESP6-HUG	Matrix Elan Hug Standard Pad 6"	\$235
☐ ESP10-HUG	Matrix Elan Hug Standard Pad 10"	\$235
☐ ESP14-HUG	Matrix Elan Hug Standard Pad 14"	\$235

Standard Pads - Wipeable Infection Control Cover

☐ ESP6-IC	Matrix Elan Standard Pad 6" with Infection Control.....	\$235
☐ ESP10-IC	Matrix Elan Standard Pad 10" with Infection Control	\$235
☐ ESP14-IC	Matrix Elan Standard Pad 14" with Infection Control	\$235
☐ ESP6-IC-HUG	Matrix Elan Hug Standard Pad 6" with Infection Control.....	\$235
☐ ESP10-IC-HUG	Matrix Elan Hug Standard Pad 10" with Infection Control	\$235
☐ ESP14-IC-HUG	Matrix Elan Hug Standard Pad 14" with Infection Control	\$235

Occipital Pads - Startex Fabric Cover



☐ EOP9	Matrix Elan Occipital Pad 9"	\$235
☐ EOP12	Matrix Elan Occipital Pad 12"	\$235

Occipital Pads - Wipeable Infection Control Cover

☐ EOP9-IC	Matrix Elan Occipital Pad 9" with Infection Control	\$235
☐ EOP12-IC	Matrix Elan Occipital Pad 12" with Infection Control.....	\$235

MATRIX HEADREST PADS CONT. [HCPCS CODE E0955](#)

4-Point Pads - Startex Fabric Cover



☐ E4PTSM	Matrix Elan 4-Point Pad 9"W x 7"H	\$325
☐ E4POINT	Matrix Elan 4-Point Pad 11"W x 10"H	\$325
☐ E4PTSHB	4 Point (Small) Pad with Startex..... fabric cover with Elan Headband	\$450
☐ E4PTSHBIC	4 Point (Small) Pad with Infection Control fabric cover with Elan Headband	\$450

4-Point Pads - Wipeable Infection Control Cover

☐ E4PTSM-IC	Matrix Elan 4-Point Pad 9"W x 7"H with Infection Control	\$325
☐ E4POINT-IC	Matrix Elan 4-Point Pad 11"W x 10"H with Infection Control	\$325
☐ E4PTLHB	4 Point (Large) Pad with Startex fabric on cover Elan Headband	\$450
☐ E4PTLHBIC	4 Point (Large) Pad with Infection Control ... fabric on Cover and Elan Headband	\$450

MATRIX HEADREST HARDWARE [HCPCS CODE E1028](#)

Elan Hardware



FLIPP Hardware




LOXX Hardware



☐ LXSHW	Loxx Standard Hardware 15"	\$215
☐ LXMMH	Loxx Mini Hardware 11"	\$215
☐ MEHW	Elan Hardware 15".....	\$220
☐ MEMH	Elan Mini Hardware 11"	\$220
☐ MEHW-EXT	Elan Extended Hardware	\$295
☐ FPLXSHW	FLIPP Loxx Standard Hardware 15".....	\$515
☐ FPLXSHW-EX	FLIPP Loxx Extended Hardware	\$525
☐ FPLXMMH	FLIPP Loxx Mini Hardware 11"	\$515
☐ FPLXMMH-EX	FLIPP Loxx Mini Extended Hardware	\$525
☐ FPMEHW	FLIPP Elan Hardware 15"	\$525
☐ FPMMH	FLIPP Elan Mini Hardware 11"	\$525
☐ FPMEHW-EXT	FLIPP Elan Extended Hardware.....	\$605

Anterior Support Tray – Build to Order



Standard Tray Sizes (Approximate)				
	Size	Chair	Tray	Body Cut
	X-Small	10-11"	21 x 18"	10 x 7"
	Small	12-13"	21 x 18"	12 x 7"
	Medium	14-15"	24 x 22"	14 x 9"
	Large	16-17"	26 x 24"	16 x 10"
	X-Large	18-19"	28 x 24"	18 x 10"

Wheelchair Trays - No Hardware			
Part Number	Description	HCPCS	Price
08-17000	1/4" Lexan Tray - No Hardware	(E0950)	\$115
08-17100	1/4" Bolteron Tray - No Hardware (Black)	(E0950)	\$163
08-8000	3/8" Lexan Tray - No Hardware	(E0950)	\$193
08-8100	3/8" Bolteron Tray - No Hardware (Black)	(E0950)	\$232
08-8500	1/2" Birch Tray - No Hardware	(E0950)	\$193
Notes:			

Wheelchair Tray Hardware				
Part Number	Description	HCPCS	Price	
454A	Adjustable Knob Release Hardware	(E1028)	\$126 pr.	
454B	Adjustable Toggle Clamp Hardware	(E1028)	\$298 pr.	
454C	Hook & Loop Strap Hardware	(E1028)	\$ 98 pr.	
461	Channel Bracket Hardware	(E1028)	\$160 pr.	
559	EZ Lite Adjustable Channel Bracket	(E1028)	\$160 pr.	

Notes:

Custom Wheelchair Tray Modifications

Part Number	Description	HCPCS	Price	Quantity
MOD-8A	3/8" Lexan Material Upgrade For Trays (Substitutes Std Tray Material with 3/8" Lexan)	(K0108)	\$ 78	Each
MOD-8A1	1/4" Bolteron Material Upgrade For Trays (Substitutes Std Tray Material with 1/4" Bolteron)	(K0108)	\$ 59	Each
MOD-8A2	3/8" Bolteron Material Upgrade For Trays (Substitutes Std Tray Material with 3/8" Bolteron)	(K0108)	\$117	Each
MOD-8C	Joystick Cut Out Right Left _____ L x _____ W	(K0108)	\$106	Each
MOD-8D	Channel Tray Stopper Hardware Kit	(K0108)	\$ 50	Kit
MOD-8E	Custom Slots for Hardware – Must Specify # & Location	(K0108)	\$ 106	Pair
MOD-8F	Round Body Cut-Out		N/C	Each
MOD-8G	Aluminum Rim - 1/16" x 1"H Rims Sides and Front Edge	(K0108)	\$ 66	Kit
MOD-8H	Plastic Rim Guard - 3/8" Attaches Directly To Tray Edge	(K0108)	\$ 43	Each
MOD-8L	Pad Tray - Tray Top Only _____ Upholstery	(K0108)	\$364	Each
MOD-8L-1	Pad Tray - Top & Bottom _____ Upholstery	(K0108)	\$605	Each
MOD-8M	Pad Tray - Half Tray Pad _____ Upholstery	(K0108)	\$306	Each
MOD-8N	Ribs - Wood Lexan Aluminum	(K0108)	\$ 58	Set
MOD-8P	Overlay - 1/8" Lexan	(K0108)	\$140	Each
MOD-8Q	Underlay - 1/8" Lexan	(K0108)	\$140	Each
MOD-8R	Easel - Removable, Adjustable, 3/8" Lexan w/ Pencil Ledge	(K0108)	\$268	Each
MOD-8S	Tray Trim - 1/16" For Aluminum Rim	(K0108)	\$ 43	Each
MOD-8T	Slip Cover For Tray Padding _____ Upholstery	(K0108)	\$271	Each

Notes:


Custom Wheelchair Tray Accessories

Part Number	Description	HCPCS	Price	Quantity
MOD-8J	Rear Stability Belt – 1" Webbing W/ Side Release Buckle	(K0108)	\$ 48	Each
MOD-8K	Elbow Blocks w/ L Bracket _____ L x _____ H	(K0108)	\$287	Pair
MOD-8O	Elbow Pads - Hook & Loop Mounting _____ L x _____ H	(K0108)	\$134	Pair
MOD-8OW	Elbow Pads - Wood Base _____ L x _____ H	(K0108)	\$287	Pair
	Upholstery _____			
MOD-8V	Wrist Straps - Hook & Loop	(K0108)	\$ 66	Pair


Notes / Diagrams:

Freedom Select Anterior Trunk Supports

Select Standard



Select Narrow



Dimensions


Qty.	Part Number	Description	A	B	C	D	E	F	G	Price
	SFDI-6SNZ-S	Select Standard Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-1-S	Select <i>Standard</i> Dual Rear-Pull Small	7.5	4.5	3.25	6	11	5.25	13.75	\$ 161
	SFDI-6SNZ-M	Select <i>Standard</i> Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-1-M	Select <i>Standard</i> Dual Rear-Pull Medium	9.5	6.5	4	6.5	14.75	7.25	17	\$ 161
	SFDI-6SNZ-L	Select <i>Standard</i> Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-1-L	Select <i>Standard</i> Dual Rear-Pull Large	10.5	7.5	4.5	6.75	18	9	19.75	\$ 173
	SFDI-6SNZ-XL	Select <i>Standard</i> Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6SNZ-1-XL	Select <i>Standard</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	7.75	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

	SFDI-6NNZ-XS	Select <i>Narrow</i> X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-1-XS	Select <i>Narrow</i> Dual Rear-Pull X-Small	6	3	2	2	6.5	4.25	10.5	\$ 161
	SFDI-6NNZ-S	Select <i>Narrow</i> Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-1-S	Select <i>Narrow</i> Dual Rear-Pull Small	7.5	4	2.75	3	11	5.25	13.75	\$ 161
	SFDI-6NNZ-M	Select <i>Narrow</i> Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-1-M	Select <i>Narrow</i> Dual Rear-Pull Medium	9.5	6.5	4	3	14.75	7.25	17	\$ 161
	SFDI-6NNZ-L	Select <i>Narrow</i> Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-1-L	Select <i>Narrow</i> Dual Rear-Pull Large	10.5	7.5	4.5	3.5	18	8	19.75	\$ 167
	SFDI-6NNZ-XL	Select <i>Narrow</i> Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6NNZ-1-XL	Select <i>Narrow</i> Dual Rear-Pull Extra Large	10.5	7.5	4.5	4.25	21.75	10	22	\$ 173
	SFDI-6CUSTOM	Select Anterior Support Custom Sizing (specify)								\$ 235

Note:

Freedom Select Anterior Pelvic Supports

Select 2-Point



2-POINT – PADDED – EZ RELEASE PUSH BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2MP-XS	XS	MICRO 1"	1	50	1.75	6	\$ 94
	SFDI-3A2MP-S	S	MICRO 1"	1	55	2.25	7	\$ 94
	SFDI-3A2P-S	S	STD 1.5"	1.5	55	2.25	7	\$ 94
	SFDI-3A2P-M	M	STD 1.5"	1.5	60	2.5	9	\$ 94
	SFDI-3A2P-L	L	STD 1.5"	1.5	65	3	11	\$ 94
	SFDI-3A2LP-L	L	LG 2"	2	65	3	11	\$ 94
	SFDI-3A2LP-XL	XL	LG 2"	2	70	3.25	13	\$ 94
	SFDI-3A2P-CUSTOM		QUOTE					\$ 118

2-POINT – PADDED – SIDE RELEASE BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2MS-XS	XS	MICRO 1"	1	50	1.75	6	\$ 68
	SFDI-3A2MS-S	S	MICRO 1"	1	55	2.25	7	\$ 68
	SFDI-3A2S-S	S	STD 1.5"	1.5	55	2.25	7	\$ 68
	SFDI-3A2S-M	M	STD 1.5"	1.5	60	2.5	9	\$ 68
	SFDI-3A2S-L	L	STD 1.5"	1.5	65	3	11	\$ 68
	SFDI-3A2LS-L	L	LG 2"	2	65	3	11	\$ 68
	SFDI-3A2LS-XL	XL	LG 2"	2	70	3.25	13	\$ 68
	SFDI-3A2S-CUSTOM		QUOTE					\$ 94

Note:

2-POINT – UNPADDED – EZ RELEASE PUSH BUTTON BUCKLE

Note: Pads not included, must be ordered separately.

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2UMP-XS	XS	MICRO 1"	1	50	1.75	6	\$ 45
	SFDI-3A2UMP-S	S	MICRO 1"	1	55	2.25	7	\$ 45
	SFDI-3A2UP-S	S	STD 1.5"	1.5	55	2.25	7	\$ 45
	SFDI-3A2UP-M	M	STD 1.5"	1.5	60	2.5	9	\$ 45
	SFDI-3A2UP-L	L	STD 1.5"	1.5	65	3	11	\$ 45
	SFDI-3A2ULP-L	L	LG 2"	2	65	3	11	\$ 45
	SFDI-3A2ULP-XL	XL	LG 2"	2	70	3.25	13	\$ 45
	SFDI-3A2UP-CUSTOM		QUOTE					\$ 57

Note:

2-POINT ANTERIOR PELVIC SUPPORT PADS

Qty.	PRICE EA.	Qty.	PRICE PR.
SFDI-3HBP-XS	\$ 32		\$ 49
SFDI-3HBP-S	\$ 32		\$ 49
SFDI-3HBP-M	\$ 32		\$ 49
SFDI-3HBP-L	\$ 32		\$ 49
SFDI-3HBP-XL	\$ 32		\$ 49
SFDI-3HBP-CUSTOM	\$ 44		\$ 73

Freedom Select Anterior Pelvic Supports (continued)


2-POINT – UNPADDED – SIDE RELEASE BUTTON BUCKLE

Note: Pads not included, must be ordered separately.

Qty.	Part Number	Size	Buckle Type	A	B	C	D	Price
	SFDI-3A2UMS-XS	XS	MICRO 1"	1	50	1.75	6	\$ 40
	SFDI-3A2UMS-S	S	MICRO 1"	1	55	2.25	7	\$ 40
	SFDI-3A2US-S	S	STD 1.5"	1.5	55	2.25	7	\$ 40
	SFDI-3A2US-M	M	STD 1.5"	1.5	60	2.5	9	\$ 40
	SFDI-3A2US-L	L	STD 1.5"	1.5	65	3	11	\$ 40
	SFDI-3A2ULS-L	L	LG 2"	2	65	3	11	\$ 40
	SFDI-3A2ULS-XL	XL	LG 2"	2	70	3.25	13	\$ 40
	SFDI-3A2US-CUSTOM		QUOTE					\$ 51

2-POINT ANTERIOR PELVIC SUPPORT PADS

Qty.	PRICE EA.	Qty.	PRICE PR.
SFDI-3HBP-XS	\$ 32		\$ 49
SFDI-3HBP-S	\$ 32		\$ 49
SFDI-3HBP-M	\$ 32		\$ 49
SFDI-3HBP-L	\$ 32		\$ 49
SFDI-3HBP-XL	\$ 32		\$ 49
SFDI-3HBP-CUSTOM	\$ 44		\$ 73



4-POINT – EZ RELEASE PUSH BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	E	F	Price
	SFDI-3A4MP-XS	XS	MICRO 1"	1	1	50	14	4.25	5.75	\$203
	SFDI-3A4MP-S	S	MICRO 1"	1	1	55	16	5	6.5	\$203
	SFDI-3A4P-S	S	STD 1.5"	1.5	1	55	16	5	6.5	\$203
	SFDI-3A4P-M	M	STD 1.5"	1.5	1	60	18	5.5	7.5	\$203
	SFDI-3A4P-L	L	STD 1.5"	1.5	1.5	65	20	6	8.75	\$203
	SFDI-3A4P-XL	XL	STD 1.5"	1.5	1.5	70	22	6	8.75	\$203
	SFDI-3A4LP-XL	XL	LG 2"	2	1.5	70	22	6	8.75	\$203
	SFDI-3A4P-CUSTOM		QUOTE							\$235

4-POINT – SIDE RELEASE BUTTON BUCKLE

Qty.	Part Number	Size	Buckle Type	A	B	C	D	E	F	Price
	SFDI-3A4MS-XS	XS	MICRO 1"	1	1	50	14	4.25	5.75	\$202
	SFDI-3A4MS-S	S	MICRO 1"	1	1	55	16	5	6.5	\$202
	SFDI-3A4S-S	S	STD 1.5"	1.5	1	55	16	5	6.5	\$202
	SFDI-3A4S-M	M	STD 1.5"	1.5	1	60	18	5.5	7.5	\$202
	SFDI-3A4S-L	L	STD 1.5"	1.5	1.5	65	20	6	8.75	\$202
	SFDI-3A4S-XL	XL	STD 1.5"	1.5	1.5	70	22	6	8.75	\$202
	SFDI-3A4LS-XL	XL	LG 2"	2	1.5	70	22	6	8.75	\$202
	SFDI-3A4S-CUSTOM		QUOTE							\$235

Note:



Accessories

PRO-10920 Tilt Range Limit Kit.....\$38 (K0108)

Pins install into outside of tilt channel to limit tilt travel.

SFDI-9HP Footrest Hanger Pads - Pair.....\$112 (K0108)

Wrap-around pads made from black neoprene.

PRO-23210 Flip-Down Rigidizing Bar.....\$189 (K0108)

Mounts to back post tubes to add rigidity to frame. Does NOT adjust in width to grow with frame.

PRO-10905 Freedom Canopy - Flip-Back.....\$275* (K0108)

Canopy assembly includes removable mounting hardware. Available in black only with Freedom Designs logo. Includes rear flap with zippered storage pocket and clear-view window.

*included at no charge when ordered with a Freedom Solid Seat & Solid Back

Specifications

Frame Style	Tilt-In-Space
Base Frame Widths	16", 18", 20"
Overall Width	Base width + 7.5" (all wheels without drum lock) Base width + 8.5" (with drum lock except 16" & 18" mag wheel) Base width + 9.5" (with drum brake and 16" or 18" mag wheel)
Seat Widths	14" - 20"
Seat Depths	14" - 26"
Back Heights	20", 23", 26"
Adjustable Angle Back	90°, 100°, 110°, 120°
Tilt Range	0° to 55°
Recline Back	65° from 97° to 172°
Weight Limit	250 LBS.
Transport Weight	34 lbs.
Warranty	5 year warranty on side frames and quick release axles

Special Configurations Through Customer Service

Anterior Tilt	As much as -5° depending on configuration. Reduces posterior tilt angle by same amount.
Asymmetrical Seat Depth	Seat rails may be configured up to 4" offset

HCPCS Codes

are provided based on verification through PDAC or our interpretation of the code definitions. The use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Save

Adobe Acrobat Reader DC allows you to save this form with your content - to complete later or use as a starting point for your next form. Please note that content must be added and saved in Acrobat - saving content from completed forms in the browser may not be possible.

Submit

Adobe Acrobat Reader DC allows you to submit this form electronically via your email client.

Simply click the submit button below and step through the simple process.

Print

If you do not have access to Adobe Acrobat Reader DC simply print this form and complete it by hand and fax it to our Customer Service Department at:

800-331-8551